## ÁCIDO ÚRICO MONOREAGENTE

REF K139

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#### FINALIDADE

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### INSTRUCTIONS FOR USE


#### FUNCTION

Method for determination of Uric Acid in biologicals samples of serum, plasma (Heparin or EDTA) and urine. Enzymatic colorimetric test, only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Enzymatic Colorimetric UOD-PAP (Trinder)

The enzymatic determination of Uric Acid is made according to the following reactions:

The intensity of the cherry color formed is directly proportional to the concentration of Uric Acid in the sample.

#### REAGENTS

**Number 1 - Enzyme Reagent** - Store between 2 and 8°C. Do not freeze. Contains: DHBS < 10 mmol/L, 4-Aminoantipyrine < 10 mmol/L, Peroxidase < 10 KU/L, Uricase < 10 KU/L, buffer, surfactant, stabilizers, and preservative.

**Number 2 - Standard** - Store between 2 and 8°C. Contains: Uric Acid 6,0 mg/dL, buffer, stabilizer and preservative.

#### PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	100 mL	3 mL
2	2 x 100 mL	3 mL
3	4 x 100 mL	3 mL
4	5 x 20 mL	3 mL
5	2 x 40 mL	3 mL
6	3 x 40 mL	3 mL
7	4 x 40 mL	3 mL
8	5 x 40 mL	3 mL
9	6 x 40 mL	3 mL
10	2 x 60 mL	3 mL
11	3 x 60 mL	3 mL
12	4 x 60 mL	3 mL
13	8 x 60 mL	3 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water-bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For *in vitro* diagnostic use only.**

**2- Strictly follow the methodology proposed to obtain exact results.**

**3- Water used in material cleaning must be recent and free of contaminants.**

**4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.**

**5- The Sodium Hypochlorite is a contaminating agent that can significantly alter the results, so the materials used to perform the tests must be properly washed and free from this type of residue.**

**6- Water level in water bath should be above the level of reagents in test tubes.**

**7- Handle with care Reagent N° 1, containing Sodium Azide, irritating to skin and mucous.**

**8- Slightly pink color development in the Reagent N° 1 does not damage their quality, as long as the corresponding white and standard dosages are used.**

**9- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance Blank current legislation.**

**10- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.**

**11- Do not use the product in case of damaged packaging.**

**12- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.**

#### SAMPLES

Serum obtained free from hemolysis, Plasma (collected with Heparin or EDTA) and Urine.

Serum or plasma: Stable for 3 days at 2 to 8°C and 7 days at -20°C.

Urine: Stable for 3 days at 2 to 8°C. **Do not freeze** <sup>5</sup>.

#### PROCESS DESCRIPTION

The stability of the Uric Acid Single Reagent product after opening (In Use / On Board – Product installed in equipment) is 30 days. This stability can vary depending on the test conditions. Therefore, it is suggested to monitor the performance of the product using control serum.

#### TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark three test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	---	---	25 µL
Reagent N° 2	---	25 µL	---
Reagent N° 1	1.0 mL	1.0 mL	1.0 mL

Homogenize well and put in a water-bath at 37°C for 5 minutes. Read the absorbance of Sample and Standard at 505 nm (490 - 540 nm), hitting the zero with the Blank. The color is stable for 30 minutes.

**24 - Hour urine** - Preparation of sample:

Take note of the volume (V) of 24-hour urine in mL, homogenize, separate a rate of 10 mL and adjust the pH between 7 and 9 with 5% NaOH. If the sample is cloudy, warm it in water bath for 10 minutes at 56°C, to dissolve the crystals. Dilute urine 1:10 (100 µL of urine + 900 µL of distilled or deionized water).

Multiply the result by 10. Proceed the dosage as described for serum.

#### CALCULATIONS

$$\text{Uric Acid (mg/dL)} = \frac{\text{Sample Absorbance} \times 6.0}{\text{Standard Absorbance}}$$

As the reaction follows the Beer-Lambert Law, the Calibration Factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (6.0 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{Uric Acid (mg/dL)} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

#### Urine

$$\text{Uric Acid (mg/24h)} = \frac{\text{mg/dL} \times \text{Urinary Volume (mL)}}{100}$$

Results are expressed in mg/dL.

#### INTERFERENCES

Use plasma samples containing only EDTA or Heparin. Other anticoagulants interfere with the reaction, with Fluoride Uricase inhibitor. The citrate can cause falsely decreased results.

For serum lipids or Bilirubin levels (> 10 mg/dL), mix 2.0 mL of 0.85% NaCl and 0.05 mL of serum. Measure absorbance at 505 nm, hitting the zero with distilled water. Decrease the absorbance thus obtained in the absorbance of the sample, and calculate the concentration in mg/dL.

The use of highly reducing drugs, such as Ascorbic Acid (Vitamin C), interferes with the reaction, because they compete with the consumption of H<sub>2</sub>O<sub>2</sub>, by providing falsely decreased values. For this reason, you should discontinue its use at least 12 hours before collection of sample.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 913A.

#### REFERENCE VALUES

The reference values, in mg/dL, for this method were obtained through the determination of Uric Acid in healthy populations of male and female.

Serum	Adults	Male	2.5 to 7.0 mg/dL
		Female	1.5 to 6.0 mg/dL
	Children	Male	1.5 to 6.0 mg/dL
		Female	0.5 to 5.0 mg/dL
Urine	250 to 750 mg/24 hours		

To convert the values from mg/dL to mmol/L (IS), multiply by 0.05948.

These values should be used as guidelines, each laboratory should create their own range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

## PRODUCT PERFORMANCE

### ACCURACY

#### Comparison of Methods

The Uric Acid Single Reagent kit was compared with other method for determination of Uric Acid dosage commercially available. 42 analyzes were performed and the results were evaluated. The linear equation obtained was  $Y = 0.9718X + 0.0996$  and correlation coefficient 0.9965. With these results we can conclude that the kit shows good methodological specificity.

### PRECISION

#### Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	11.86	6.27	4.18
<b>Standard Deviation (mg/dL)</b>	0.08	0.05	0.06
<b>Coefficient of Variation (%)</b>	0.69	0.88	1.39

#### Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
<b>Average Concentration (mg/dL)</b>	11.51	6.17	4.11
<b>Standard Deviation (mg/dL)</b>	0.34	0.09	0.06
<b>Coefficient of Variation (%)</b>	2.96	1.38	1.49

### SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free of Uric Acid. The average found was 0.042 mg/dL with standard deviation of 0.029 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 0.129 mg/dL.

### LINEARITY

The reaction is linear up to 20.0 mg/dL. For samples with values above 20 mg/dL or optical density above 0.8, it is recommended diluting the sample with distilled water, repeat the dosage and multiply the result by the dilution factor.

### DIAGNOSTIC SIGNIFICANCE

Uric Acid is the end product of metabolism of Nucleic Acids and Purines. Its concentration in body fluids depends on the balance between production and delimitation through the kidneys, and generated about 400 mg daily.

The biochemical characteristics and the basic condition for the diagnosis of gout is hyperuricemia. From the epidemiological point of view, Uric Acid levels above 7 mg/dL may indicate high risk of gouty arthritis or nephrolithiasis.

The concentration of plasma Uric Acid is elevated in several other clinical conditions such as renal failure, congestive heart failure, toxemia of pregnancy (eclampsia), leukemias and lymphomas, ketoacidosis, Down's syndrome. Often the observation of hyperuricemia associated with factors such as obesity, hypertriglyceridemia, hypertension and mellitus diabetes.

Reduced levels of Uric Acid are observed in Fanconi syndrome, Wilson's disease and in patients treated with Allopurinol, Corticosteroids, ACTH and Phenylbutazone.

### BIBLIOGRAPHIC REFERENCES

- 1 - TRIVEDI, R. C.; REBAR, L.; BERKA, E.; STRONG, L., New Enzymatic Method for Serum Uric Acid at 500 nm. Clin. Chem., 1978, 24, 1908.
- 2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.
- 3 - TRINDER, P., Ann. Clin. Biochem., 1969, 6, 24.
- 4 - TIETZ Textbook of Clinical Chemistry, 2nd ed., 1994.
- 5 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte/MG - Brasil  
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Uric Acid Single Reagent kit: 10269360228

Review: May/2023

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE  
(last day of the month)

TEMPERATURE LIMIT  
(store)

CONTENT IS SUFFICIENT  
FOR <N> TEST

SEE INSTRUCTIONS  
FOR USE

IN VITRO DIAGNOSTIC  
PRODUCT

KEEP AWAY  
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF  
PACKAGE IS  
DAMAGED

PRODUCT  
STERILIZED

DANGER