ALBUMINA MONORREAGENTE

REF K304

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FINALIDADE

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VETERINARY USE

REF K304

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

INSTRUCTIONS FOR USE

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for determination of Albumin in biologicals samples of serum or plasma (Heparin or EDTA). Colorimetric test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Bromocresol Green (BCG)

The measurement uses what is called "protein error of indicators". In the presence of Albumin, Bromocresol Green forms a colored complex, which exhibits an absorption spectrum different from the dye in its free state, thus allowing the Albumin dosage.

REAGENTS

Number 1 - Color Reagent - Store between 2 and 8 °C. Contains: Verde de Bromocresol < 1.0 mmol/L, Citrato buffer < 100 mmol/L (pH 4.0), surfactant and preservative.
Number 2 - Standard - Store between 2 and 8 °C. Shake before use. Contains: Albumin 3.8 g/dL, stabilizer and preservative.

Presentation	Reagent N° 1	Reagent N° 2
K304-1-VET	250 mL	2 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, watch or stopwatch, pipettes, test tubes, Bioclin's Bicontrol N and Biocontrol P. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8 °C. The transport at temperatures up to 30 °C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Reagent N° 2 contains Sodium Azide: handle with care.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (Heparin or EDTA) obtained free of hemolysis stable for 21 days between 2 and 8°C and 4 months at -20 °C⁷.

PROCESS DESCRIPTION

The stability of the Albumin Monoreagent product after opening (In Use / On Board – Product installed in equipment) is 30 days. This stability may vary according to test conditions. Therefore, it is suggested to monitor the performance of the product using control serum.

WORKING REAGENT PREPARATION

The reagent is ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 03 test tubes: B (Blank), P (Standard), A (Sample), and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 µL
Standard	--	10 µL	--
Color Reagent	2.5 mL	2.5 mL	2.5 mL

Homogenize well and leave it at rest for 5 minutes. Read the absorbance of Sample and Standard at 630 nm (600 - 640 nm), hitting the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

$$\text{Albumin (g/dL)} = \frac{\text{Sample Absorbance} \times 3.8}{\text{Standard Absorbance}}$$

As the reaction follows the Beer-Lambert Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (3.8 g/dL)}}{\text{Standard Absorbance}}$$

$$\text{Albumin (g/dL)} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

Results are expressed in g/dL.

PROCEDURE LIMITATIONS

Use an automated pipette and good quality tips for quality pipetting from Standard and Sample, to minimize problems with volume imprecision. Use of tourniquet for longer than 3 minutes increases albumin levels.

INTERFERENT

No interference was observed by Bilirubin until 38 mg/dL, Hemoglobin until 180 mg/dL and Triglycerides until 250 mg/dL. Corticosteroid use can cause falsely increased albumin results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D .

REFERENCE VALUES

The reference values were determined in serum samples of healthy dogs, cats, cattle and horses, males and females, adults and without distinction of race.

Species	Value
Dogs	2.4 - 4.2 g/dL
Cats	2.1 - 3.6 g/dL
Cattle	2.0 - 3.7 g/dL
Horses	2.1 - 3.6 g/dL

These values should be used as a guideline and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Albumin Monoreagent kit was compared with another method commercially available to measure Albumin. 42 analysis were performed and the results were evaluated. The linear equation obtained was $Y = 1.0311X - 0.0731$ and the correlation coefficient 0.9957. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (g/L)	2.49	3.84	3.22
Standard Deviation (g/L)	0.01	0.02	0.01
Coefficient of Variation (%)	0.50	0.47	0.43

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (g/L)	2.51	3.87	3.22
Standard Deviation (g/L)	0.01	0.03	0.01
Coefficient of Variation (%)	0.43	0.75	0.21

SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free of Albumin. The average found, which indicates the detection limit, was of 0.034 g/dL with standard deviation of 0.002 g/dL. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 0.039 g/dL.

LINEARITY

The reaction is linear up to 6.0 g/dL. For samples with values above 6.0 g/dL, it is recommended to dilute it with Saline Solution at 0.85%, repeat the dosage and multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Albumin is the most abundant protein in plasma and its functions are to maintain colloid osmotic pressure and transport substances in the body.

The main causes of the increase in albumin concentration is dehydration or loss of body fluids. Other less frequent causes for increased plasma albumin are adrenal gland dysfunction and hepatocellular carcinoma.

The collection of plasma in a tube containing heparin and the use of corticosteroids are causes of false albumin

increase.

Decreased serum albumin may occur in situations of dietary deficit, reduced albumin production (due to liver failure, recent proinflammatory response, and increased oncotic pressure), loss of albumin (caused by enteropathies, glomerulopathies, hemorrhage, or exudative skin diseases), severe sequestration (in case of effusions) or catabolism.

BIBLIOGRAPHIC REFERENCES

- 1 - Gornall AG, Bardawill CJ, David MMJ. Biol Chem, 1977, 751.
- 2 - Weichselbaum TE, Amer J. Clin Pathol. 1946, 16,40.
- 3 - Slater L, Carter PM, Hobbs JR. Ann Clin Biochem. 1975, 12,333.
- 4 - Batsakis JG, Aroushon RS, Walker WA, Barnes B, Amer J. Clin. Pathol., 1.976, 66,238.
- 5 - Hoel PG. Estatística Elementar, Ed. Fundo de Cultura S/A, 1969.
- 6 - Tonks DB. Clin Chem. 1983, 9,217.
- 7 - Who. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 8 - ECLINPATH. Cornell University College of Veterinary Medicine. 2020.
- 9 - Kaneko JJ, Harvey JW, Bruss ML. Clinical Biochemistry of Domestic Animals, 6 ed. Editora Elsevier, 2008.
- 10 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | Email: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER