ALBUMINA MONORREAGENTE

REF K304

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Method for determination of Albumin in biologicals samples of serum or plasma (Heparin or EDTA). Colorimetric test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Bromocresol Green (BCG)

The measurement uses what is called "protein error of indicators". In the presence of Albumin, Bromocresol Green forms a colored complex, which exhibits an absorption spectrum different from the dye in its free state, thus allowing the Albumin dosage.

REAGENTS

Number 1 - Color Reagent - Store between 2 and 8°C. Contains: Verde de Bromocresol < 1.0 mmol/L, Citrato buffer < 100 mmol/L (pH 4.0), surfactant and preservative.
Number 2 - Standard - Store between 2 and 8°C. Shake before use. Contains: Albumin 3.8 g/dL, stabilizer and preservative.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2*
1	250 mL	2 mL
2	500 mL	2 mL
3	5 x 40 mL	2 mL
4	10 x 40 mL	2 mL

*The R2 Standard is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, watch or stopwatch, pipettes, test tubes, Bioclin's Bicontrol N and Biocontrol P. They can be found at markets specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must to be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Reagent Nº 2 contains Sodium Azide: handle with care.

- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (Heparin or EDTA) obtained free of hemolysis stable for 21 days between 2 and 8°C and 4 months at -20°C⁷.

PROCESS DESCRIPTION

The stability of the Albumin Monoreagent product after opening (In Use / On Board – Product installed in equipment) is 30 days. This stability may vary according to test conditions. Therefore, it is suggested to monitor the performance of the product using control serum.

WORKING REAGENT PREPARATION

The reagent is ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 03 test tubes: B (Blank), P (Standard), A (Sample), and proceed as follows:

	Blank	Standard	Sample
Sample	--	--	10 µL
Standard	--	10 µL	--
Color Reagent	2.5 mL	2.5 mL	2.5 mL

Homogenize well and leave it at rest for 5 minutes. Read the absorbance of Sample and Standard at 630 nm (600 - 640 nm), hitting the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

Albumin (g/dL) =
$$\frac{\text{Sample Absorbance} \times 3.8}{\text{Standard Absorbance}}$$

As the reaction follows the Beer-Lambert Law, the calibration factor can be used.

Calibration Factor =
$$\frac{\text{Standard Concentration (3.8 g/dL)}}{\text{Standard Absorbance}}$$

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

Albumin (g/dL) = Sample Absorbance x Calibration Factor

Results are expressed in g/dL.

PROCEDURE LIMITATIONS

Use an automated pipette and good quality tips for quality pipetting from Standard and Sample, to minimize problems with volume imprecision. Use of tourniquet for longer than 3 minutes increases albumin levels.

INTERFERENT

No interference was observed by Bilirubin until 38 mg/dL, Hemoglobin until 180 mg/dL, Triglycerides until 250 mg/dL, Ascorbic Acid until 60 mg/dL and Rheumatoid Factor until 300 UI/mL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D .

REFERENCE VALUES

The reference values in g/dL, for this method, were obtained through the determination of Albumin in healthy populations.

Children and Teenagers	1 - 30 days	2.6 - 4.3 g/dL
	31 - 182 days	2.8 - 4.6 g/dL
	183 - 365 days	2.8 - 4.8 g/dL
	1 - 18 years	2.9 - 4.7 g/dL
Adults		3.5 - 5.5 g/dL

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY

Comparison of Methods

The Albumin Monoreagent kit was compared with another method commercially available to measure Albumin. 42 analysis were performed and the results were evaluated. The linear equation obtained was $Y = 1.0311X - 0.0731$ and the correlation coefficient 0.9957. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (g/L)	2.49	3.84	3.22
Standard Deviation (g/L)	0.01	0.02	0.01
Coefficient of Variation (%)	0.50	0.47	0.43

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (g/L)	2.51	3.87	3.22
Standard Deviation (g/L)	0.01	0.03	0.01
Coefficient of Variation (%)	0.43	0.75	0.21

SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free of Albumin. The average found, which indicates the detection limit, was of 0.034 g/dL with standard deviation of 0.002 g/dL. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 0.039 g/dL.

LINEARITY

The reaction is linear up to 6.0 g/dL. For samples with values above 6.0 g/dL, it is recommended to dilute it with Saline Solution at 0.85%, repeat the dosage and multiply the result by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The Albumin constitutes the main serum protein. Synthesized almost entirely by the liver, has a half-life of approximately two weeks. An increase of albumin may be observed in dehydration, shock and hemoconcentration. Lower values occur in malnutrition, nephrotic syndrome, hepatic failure, glomerulonephritis, multiple myeloma, anemia serious, pregnancy, severe infections and protracted.

BIBLIOGRAPHIC REFERENCES

- 1 - GORNALL, A. G.; BARDAWILL, C. J.; DAVID, M. M. J., Biol. Chem., 1977, 751.
- 2 - WEICHELBAUM, T. E.; AMER, J., Clin. Pathol., 1946, 16,40.
- 3 - SLATER, L.; CARTER, P. M.; HOBBS, J. R., Ann. Clin. Biochem., 1975, 12,333.
- 4 - BATSAKIS, J. G.; AROUSOHN, R. S.; WALKER, W. A.; BARNES, B.; AMER, J., Clin. Pathol., 1.976, 66,238.
- 5 - HOEL. P. G., em Estatística Elementar, Ed. Fundo de Cultura S/A, 1969.
- 6 - TONKS, D. B., Clin. Chem., 1983, 9,217.
- 7 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | Email: sac@bioclin.com.br

ANVISA registration for Albumin Monoreagent kit:
10269360376

Review: April/2025

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER