# Bioclin

## ALFA 1 GLICOPROTEÍNA ÁCIDA

REF K060

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REF K060

### INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

#### FUNCTION

Method for the determination of Alpha 1 Acid Glycoprotein in biological human samples of serum or plasma (heparin). Immunoturbidimetric test, for *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Immunoturbidimetry

The reagent allows the quantitative determination of Alpha 1 Acid Glycoprotein in human serum by immunoprecipitation, in the presence of an activator polymer that increases the sensitivity and immunoturbidimetric test speed. The Alpha1 Glycoprotein Acid to forms along with the specific antiserum an insoluble complex, producing turbidity, whose absorbance is proportional to the concentration of Alpha 1 Acid Glycoprotein in the sample.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: PEG 6.000 < 15%, buffer and preservative.

**Number 2 - Antiserum** - Store between 2 and 8°C. Contains: Goat antiserum for Alpha 1 Acid Glycoprotein < 2 g/L, buffer and preservative.

#### PRESENTATION

Presentation	Reagent Nº 1	Reagente Nº 2
1	25 mL	2 mL
2	2 x 25 mL	2 x 2 mL
3	3 x 25 mL	3 x 2 mL
4	4 x 25 mL	4 x 2 mL
5	1 x 40 mL	1 x 10 mL
7	20 mL	5 mL
8	2 x 20 mL	2 x 5 mL
9	3 x 20 mL	3 x 5 mL
10	4 x 20 mL	4 x 5 mL

Calibrator not included in kit.

#### EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Multical kit and Multicontrol Bioclin.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must to be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Handle with care as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Use fresh serum or plasma (heparin). The analyte is stable for 2 days at 2 to 8°C and 30 days at -20°C. Avoid repeated freezing and thawing. For samples thawed, homogenized prior to use. Do not use samples with signs of contamination.

#### PROCESS DESCRIPTION

The stability of the Alpha 1 Acid Glycoprotein product after opening (In Use / On Board – Product installed in equipment) is 15 days. This stability may vary according to test conditions. Therefore, it is suggested to monitor the performance of the product using control serum.

#### WORKING REAGENT PREPARATION

Reagents are ready to use.

#### TECHNIQUE

For calibration and control of the reaction, utilize the Multical kits and Multicontrol Bioclin.

This kit is for use in automated biochemical analyzers only.

**Check the application sheet for this equipment on the website [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC 0800 031 5454.**

#### PROCEDURE LIMITATIONS

Hemolysis, icterus and lipemia (low or moderate) do not not interfere in performance of the tests, using the recommended dilution for sample. Gross lipemia and turbid samples should be processed with a higher dilution.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### REFERENCE VALUES

Serum/Plasma (heparin)	Adults	30 - 120 mg/dL
------------------------	--------	----------------

These values should be used as a guideline, and each laboratory should establish its range of reference values, according with the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

#### PRODUCT PERFORMANCE

##### ACCURACY

##### Recovery

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 98 and 106%.

##### Comparison of Methods

Alpha 1 Acid Glycoprotein kit was compared with other method for determination of Alpha 1 Acid Glycoprotein dosage commercially available. 40 tests were performed and the results were evaluated. The linear equation obtained was  $Y = 0.946X + 5.980$  and correlation coefficient 0.989. With these results we can conclude that the kit shows good methodological specificity.

## PRECISION

### Repeatability

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2
Concentration (mg/dL)	62.87	40.52
Standard Deviation (mg/dL)	0.24	0.42
Coefficient of Variation (%)	0.38	1.05

### Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2
Concentration (mg/dL)	62.98	40.80
Standard Deviation (mg/dL)	0.25	0.38
Coefficient of Variation (%)	0.39	0.93

## SENSITIVITY

Sensitivity was calculated from 10 determinations of a sample free of the Alpha 1 Acid Glycoprotein. The average, which indicates the method detection limit, found was 0.53 mg/dL, with standard deviation of 0.05 mg/dL. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 0.67 mg/dL.

## LINEARITY

The reaction is linear up to the level of the highest calibrator concentration.

## DIAGNOSTIC SIGNIFICANCE

Main component of the Wintzler Mucoproteins, rises considerably in acute inflammatory processes. It is an important index of rheumatic activity, because it remains high when other tests normalize.

The increase in values also happens in erythematosus lupus, neoplasm metastasis and administration of estrogen. Decrease of values occurs in liver failure, adrenal insufficiency, malnutrition and pituitary insufficiency

## BIBLIOGRAPHIC REFERENCES

- 1 - KILLINGSWORTH, L. M; SAVORY, J.: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971; 17, 936.
- 2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry, C. V. Mosby Company, 1987.
- 3 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.
- 4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.
- 5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26, 1459.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

## QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte/MG - Brasil  
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

## CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | Email: sac@bioclin.com.br

ANVISA registration for Alpha 1 Acid Glycoprotein kit:  
10269360061

**Review:** November/2025

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER