# ANTI ESTREPTOLISINA O

REF **K057**

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## INSTRUCTIONS FOR USE

### FUNCTION

Method of quantitative determination of Anti-Streptolysin O (ASO) quantity in biological sample of serum. Immunoturbidimetric test, for *in vitro* diagnostic use only.

### PRINCIPLE OF ACTION

#### Methodology: Immunoturbidimetry

The polystyrene particles coated with Streptolysin O coalesce with samples containing Anti-Streptolysin O. The intensity of scattered light is proportional to the concentration of Anti-Streptolysin O, so that, by comparison with the calibrator concentration known, one can determine the amount of ASO in the sample.

#### Reagents

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: Buffer < 200 mmol/L, stabilizers, surfactant and preservative.

**Number 2 - ASO Latex** - Store between 2 and 8°C. Contains: Suspension of Latex particles sensitized with Streptolysin O < 50 mg/dL, stabilizers and preservative.

**Number 3 - Calibrator** - Store between 2 and 8°C. Contains: Anti Streptolysin O, stabilizers and preservative. **Potentially infectious.**

**Warning: The concentration of Anti Streptolysin O varies according to the lot - see the bottle's label.**

### PRESENTATION

Presentations	Reagent N° 1	Reagent N° 2	Reagent N° 3
1	18 mL	2 mL	1.0 mL
2	45 mL	5 mL	1.0 mL
3	2 x 45 mL	2 x 5 mL	2 x 1.0 mL
4	20 mL	5 mL	1.0 mL
5	2 x 20 mL	2 X 5 mL	1.0 mL
6	3 x 20 mL	3 x 5 mL	2 x 1.0 mL
7	4 x 20 mL	4 x 5 mL	2 x 1.0 mL
8	40 mL	10 mL	1.0 mL
9	2 x 40 mL	2 x 10 mL	2 x 1.0 mL
10	3 x 40 mL	3 x 10 mL	3 x 1.0 mL
11	4 x 40 mL	4 x 10 mL	4 x 1.0 mL
12	60 mL	15 mL	2 x 1.0 mL
13	2 x 60 mL	2 x 15 mL	3 x 1.0 mL

### EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Multical and Multicontrol Bioclin. They can be found at markets specialized in Laboratories of Clinical Analysis.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

### SPECIAL CARE

**1- For *in vitro* diagnostic use only.**

**2- Strictly follow the methodology proposed to obtain exact results.**

**3- Water used in material cleaning must be recent and free of contaminants.**

**4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.**

**5- Reagent N° 3 should be handled cautiously since it is subject to biological contamination.**

**6- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.**

**7- The reagent has been tested for HIV, HCV and HBsAg using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially infectious.****

**8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.**

**9- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.**

**10- Do not use the product in case of damaged packaging.**

**11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.**

### SAMPLES

Serum obtained free of hemolysis and intense lipemia. Samples are stable in temperatures between 2 and 8°C for up to 2 days and 3 months at -20°C.

### PROCESS DESCRIPTION

The calibration stability of the Anti Streptolysin O kit installed on refrigerated equipment is up to 22 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

### TECHNIQUE

The reagents are ready to use. Thoroughly mix the Reagent N° 2 before use.

For calibration of the reaction, utilize preferentially the kit's internal calibrator or the Bioclin Multical kit. For control use the Bioclin Multicontrol Kit.

This kit is for use in automated biochemical analyzers.

**Check the application sheet for this equipment on [www.bioclin.com.br](http://www.bioclin.com.br) or through SAC.**

### PROCEDURE LIMITATIONS

Hemolysis, icterus and lipemia (low or moderate) do not interfere with assay performance. Gross lipemia and turbid samples should be processed with a higher dilution and read with a dynamic range. Values below 20 UI/mL lead to poor reproductive results.

### INTERFERENT

No interference was observed by Bilirubin until 60 mg/dL, Hemoglobin until 750 mg/dL, Triglycerides until 900 mg/dL and Rheumatoid Factor 600 UI/L.

### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

### TRACEABILITY

The kit calibrator and the Multical are traceable to the reference material NIBSC 97/662 (Anti-Streptolysin-O, Human, NIBSC Reagent).

### REFERENCE VALUES

The reference values for this method were obtained by determining anti-streptolysin O in healthy populations of male and female.

Less than 200 UI/mL

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

### PRODUCT PERFORMANCE

#### ACCURACY

##### Recovery

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 96 and 105%.

## Comparison of Methods

The Anti Streptolysin O kit was compared with other commercially available methods for determination of Anti-Streptolysin O dosage. 41 tests were performed and the results were evaluated. The linear equation obtained was  $Y = 0.947X + 17.127$  and correlation coefficient 0.9779. With these results we can conclude that the kit shows good methodological specificity.

## PRECISION

### Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Concentration (UI/mL)	96.28	201.59	390.12
Standard Deviation (UI/mL)	1.97	4.35	6.53
Coefficient of Variation (%)	2.04	2.16	1.67

### Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Concentration (UI/mL)	101.44	206.18	399.72
Standard Deviation (UI/mL)	4.42	6.87	11.39
Coefficient of Variation (%)	4.35	3.33	2.85

## SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free of Anti-Streptolysin O. The average found, which indicates the method detection limit, was 2.276 UI/mL, with standard deviation of 1.584 UI/mL. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 7.027 UI/mL.

Values below 20 UI/mL lead to poor reproductive results.

## LINEARITY

The reaction is linear up to 1200 UI/mL. For samples with higher values than 1200 UI/mL, it is recommended to dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the result by the factor dilution.

## HIGH DOSE HOOK EFFECT

Hook effect wasn't observed with a high concentration of Antistreptolysin O up to 1800 IU/mL.

## DIAGNOSTIC SIGNIFICANCE

The Streptolysin O is a toxic immunogenic exoenzymes produced by group A  $\beta$ -hemolytic streptococci. Given the large amount of exoenzymes released *in vivo*, determination of antibody response against the presence of Streptolysin O has become a routine procedure for the diagnosis and treatment of rheumatic fever, acute glomerulonephritis, scarlet, tonsillitis, erysipelas, sepsis puerperal and other streptococcal infections in group A.

## BIBLIOGRAPHIC REFERENCES

- 1 - NASAHIKO, Oukumi, Antistreptolysin O Japan Clinical, 1985, 43, 387.
- 2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry, C. V. Mosby Company, 1987.
- 3 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.
- 4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.
- 5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980, 26, 1459.
- 6 - GALVIN, J. P. et., al Clin. Lab. Assays, 73, 4th, 1983.
- 7 - PASSING, H.; BABLOCK, W., J. Clin. Chem. Clin. Biochem 21, 709; 1983.
- 8 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

## QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte/MG - Brasil  
Phone.: +55 31 3439.5454 | E-mail: [bioclin@bioclin.com.br](mailto:bioclin@bioclin.com.br)  
CNPJ: 19.400.787/0001-07 - Made in Brazil

## CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

ANVISA registration number for Anti Streptolysin O kit:  
10269360067

**Review:** January/2026

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER