BABESIA IgG VET FAST

REF VET014

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FINALIDADE

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- 1- BOOZER AL, MACINTIRE DK. Canine babesiosis. Vet Clin North Am Small Anim Pract, v. 33, n. 4, p. 885-905, viii, 2003.
- 2- KOSTER L, LOBETTI R, KELLY P. Canine babesiosis: a perspective on clinical complications, biomarkers, and treatment. Vet Med (Auckl), v.6, p.119-128, 2015.
- 3- LEISEWITZ AL, MRLJAK V, DEAR JD, BIRKENHEUER A. The Diverse Pathogenicity of Various Babesia Parasite Species That Infect Dogs. Pathogens, v. 12, n.12, p. 1437, 2023.
- 4- PANTI-MAY JA, RODRÍGUEZ-VIVAS RI. Canine babesiosis: A literature review of prevalence, distribution, and diagnosis in Latin America and the Caribbean. Vet Parasitol Reg Stud Reports, v. 21, p. 100417, 2020.
- 5- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1- BOOZER AL, MACINTIRE DK. Canine babesiosis. Vet Clin North Am Small Anim Pract, v. 33, n. 4, p. 885-905, viii, 2003.
- 2- KOSTER L, LOBETTI R, KELLY P. Canine babesiosis: a perspective on clinical complications, biomarkers, and treatment. Vet Med (Auckl), v.6, p.119-128, 2015.
- 3- LEISEWITZ AL, MRLJAK V, DEAR JD, BIRKENHEUER A. The Diverse Pathogenicity of Various Babesia Parasite Species That Infect Dogs. Pathogens, v. 12, n.12, p. 1437, 2023.
- 4- PANTI-MAY JA, RODRÍGUEZ-VIVAS RI. Canine babesiosis: A literature review of prevalence, distribution, and diagnosis in Latin America and the Caribbean. Vet Parasitol Reg Stud Reports, v. 21, p. 100417, 2020.
- 5- QUIBASA: Datos del Departamento de Investigación y Desarrollo.

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INSTRUCTIONS FOR USE

FUNCTION

Rapid immunochromatographic test for the qualitative determination of anti-*Babesia vogeli* antibodies in serum, plasma, or whole blood samples.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
Babesia IgG VET FAST is an *in vitro* diagnostic test for the rapid and qualitative detection of anti-*Babesia vogeli* antibodies in serum, plasma, or whole blood samples. The cassette contains a nitrocellulose membrane impregnated at the test line with recombinant *Babesia vogeli* antigen. The conjugate consists of a mouse anti-dog IgG antibody conjugated to colloidal gold particles. When the sample is added to the cassette, it migrates by capillarity until it reaches the conjugate. In this region, anti-*Babesia vogeli* IgG antibodies present in the sample interact with the conjugate, forming the complex: Dog IgG – Conjugate. The sample continues to migrate along the immunochromatographic strip until it reaches the test line, where the dog IgG antibody also interacts with the immobilized *Babesia vogeli* antigen at the test line, forming the complex: *Babesia vogeli* Antigen – Dog IgG – Conjugate. This results in the formation of a pink-red line at the test line, indicating a reactive result. Excess sample continues migrating along the test strip until it reaches the control line, which is composed of immobilized anti-conjugate antibody. This results in the formation of a second line (control line), indicating that the test is valid.

REAGENTS

Number 1 (R1) Cassette - Store between 15 and 30 °C. **Do not freeze.** Contains:

A) 01 plastic cassette with 01 test strip containing:
- Control Region (C): Goat anti-mouse antibody immobilized.
- Test Region (T): Immobilized *Babesia vogeli* antigen.
- Conjugate: Mouse anti-dog IgG antibody conjugated to colloidal gold particles.

B) 01 sachet of silica

Number 2 (R2) Diluent - Store between 15 and 30 °C. **Do not freeze.** Contains: Buffered solution and preservative.

Accessories: Disposable plastic pipette.

PRESENTATION

Reagent	Presentation	Number of Cassettes per Pack	Diluent
Cassette + Diluent	VET014-1	10	1 x 2.5 mL
	VET014-2	25	1 x 5.0 mL
	VET014-3	10	10 x 0.5 mL
	VET014-4	25	25 x 0.5 mL
	VET014-5	5	1 x 1.25 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table
- Instructions for use (manual)

Required materials not contained in the kit:

Clock or stopwatch, PPE (gloves, mask, protective glasses) and trash for disposal. All listed inputs can be found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage and transport temperature of the product should be between 15 and 30 °C. Keep protected from direct sunlight. Do not freeze.

The cassette is sensitive to air humidity. Therefore, it must always be kept inside its primary packaging (aluminum pouch) until the moment of use.

SPECIAL CARES

- Product for *in vitro* diagnosis, for veterinary use only.**
- Strictly follow the proposed methodology to ensure accurate results.
- The cassettes (test devices) do not contain contaminating substances, but since the samples are biological, strict safety protocols must be followed. The use of disposable gloves and other personal protective equipment is mandatory.
- Use suitable containers for storing biological samples.
- Open the aluminum pouch to remove the cassette (test device) only at the time of use.
- The kit must not be frozen or exposed to high temperatures, as this will cause deterioration of its components.
- Reagents and samples must reach room temperature before starting the test.
- Do not use reagents from different batches or kits.
- Do not use the kit past its expiration date. The expiration date is indicated on the product label.
- All components of the kit that come into contact with the sample are disposable and must not be reused.
- Components that have been in contact with the patient's sample and the diluent included in the kit must be discarded as biological waste.
- Secondary packaging (box), Instructions for Use, aluminum pouch, and silica gel packet can be disposed of in regular waste.
- Follow local, state, and federal environmental protection regulations to dispose of reagents and biological material in compliance with current legislation.
- For biosafety information or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or request it via Quibasa's Customer Support Service.
- Do not use the product if the packaging is damaged.

SAMPLES

Whole Blood

Collect the blood sample into a tube containing EDTA, Heparin, or Citrate, following standard laboratory procedures. Whole blood can be stored between 2 and 8 °C and tested within 24 hours after collection. **Do not freeze.**

Serum and Plasma

Separate from the blood interface after centrifugation. For serum, allow the blood to clot before centrifugation. For plasma, collect using EDTA, Heparin, or Citrate as an anticoagulant. Plasma or serum can be stored for up to 3 days between 2 and 8 °C or for 6 months at -20 °C.

PROCESS DESCRIPTION

- The sample must be at a temperature between 15 and 30 °C before starting the test.
- Remove the cassette (test device) from the protective packaging, place it on a clean, level surface, and label it appropriately.
- Collect the sample using the disposable plastic pipette provided in the kit.
- Squeeze the pipette bulb and touch the sample. Release the bulb pressure to draw the sample. Draw up to the line to obtain a volume of 20 µL.
- To dispense the sample, squeeze the bulb again, touching the lower tip of the pipette to the test membrane at the sample application well.


Note: Alternatively, an automatic pipette can be used to collect and add the sample.

Note: Before applying the Diluent (R2), wait for the sample to be absorbed into the membrane at the sample application well.

6- Hold the Diluent (R2) bottle vertically and apply 2 drops into the sample application well (70 µL).

7- Wait for the lines to form. Interpret the results between 20 and 30 minutes. **Do not interpret after 30 minutes.**

Note: For samples that do not show visible color development on the test line at 20 minutes, it is recommended to wait until the end of the reading time, totaling 30 minutes.

INTERPRETATION OF THE RESULT

Reactive Test: Formation of a red line in the test region (T) and another line in the control region (C) within the first 20 to 30 minutes. Do not interpret after 30 minutes.	
Non-Reactive Test: Formation of a red line in the control region (C) and complete absence of a red line in the test region (T). Do not interpret after 30 minutes.	
Invalid Result: Complete absence of a line in the control region (C), with or without a red line in the test region (T). Test the sample again.	

Reactive: A reactive result in the Babesia IgG VET FAST test indicates the presence of IgG antibodies specific to *Babesia vogeli*, suggesting current infection or previous exposure to the pathogen. However, the test does not differentiate between active and past infection. Therefore, it is essential to correlate the result with the patient's history, clinical signs, and complete blood count. If necessary, perform complementary tests to confirm active infection.

Non-reactive: A non-reactive result in the Babesia IgG VET FAST test indicates the absence of IgG antibodies specific to *Babesia vogeli*, suggesting that there has been no infection or previous exposure to the pathogen. However, if recent infection is suspected, the organism may not yet have produced antibodies in sufficient quantity to be detected. It is imperative to correlate the laboratory result with the clinical presentation, the patient's history (especially tick exposure), and, if necessary, repeat the test or perform other complementary examinations for a more accurate interpretation.

Notes:

- 1- The line in the control region may appear before the 20 minutes of incubation; however, this does not mean results can be interpreted earlier than the stipulated time.
- 2- Results must not be interpreted after 30 minutes to avoid incorrect interpretations.
- 3- Results provided by this kit must be interpreted by a licensed veterinary professional and should not be the sole criterion for diagnosing or treating the patient.

TEST LIMITATIONS

- 1- Using the correct sample volume is crucial, as smaller or larger volumes may produce erroneous results.
- 2- A non-reactive result may occur if the amount of analyte in the sample is below the detection limit of the test. This may depend on the disease stage, the type of sample used, and its storage conditions.
- 3- For accurate diagnosis, the professional must consider not only the test result but also clinical signs, the animal's health condition, and any additional diagnostic results.

INTERNAL QUALITY CONTROL

The clinical laboratory should implement an internal quality control program with clearly established procedures, standards, limits, and tolerances for variations. It is essential to recognize that all measurement systems have inherent analytical variability, which should be monitored by the laboratories themselves. For this purpose, the use of controls is recommended to assess precision and accuracy.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations using 2 different positive samples. The following results were obtained:

Repeat-ability	N° of Repeti-tions	Expected outcome	Result Found
Sample 1	10	Reactive	Reactive
Sample 2	10	Reactive	Reactive

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days using 2 different positive samples. The following results were obtained:

Reproduc-ibility	N° of Days	N° of Rep-etitions	Expected outcome	Result Found
Sample 1	3	10	Reactive	Reactive
Sample 2	3	10	Reactive	Reactive

CLINICAL SENSITIVITY AND SPECIFICITY

This product was tested against other methods. A total of 98 samples characterized for canine babesiosis were analyzed, with 38 samples diagnosed as reactive and 60 as non-reactive. The results show that the clinical sensitivity of the product is 97.36% and the clinical specificity is > 99%.

METHOD		REFERENCE		Total
BABESIA IgG VET FAST	Result	Reactive	Non-reactive	
	Reactive	37	0	37
	Non-reactive	1	60	61
Result Total		38	60	98

Clinical Sensitivity: 97.36% (37/38) - 95% CI: 92.34% - 100%
Clinical Specificity: >99% (60/60) - 95% CI: 95% - 100%
Accuracy: 99% [(37+60)/(38+60)]

CLINICAL SIGNIFICANCE

Canine babesiosis is a tick-borne infection prevalent in tropical and subtropical regions, including Brazil, where *Babesia vogeli* is the predominant species. The disease can be subclinical or mild in adult animals, but in puppies, it can cause severe anemia and other serious complications. The incubation period ranges from 10 to 28 days, and antibody production begins between 8 and 10 days after infection. Infection by *Babesia* causes erythrocyte destruction, leading to anemia, thrombocytopenia, coagulopathies, and, in severe cases, organ failure. Clinical manifestations vary depending on the degree of hemolysis and the body's inflammatory response, and they may include renal, hepatic, and pancreatic lesions, as well as pulmonary edema. Cerebral babesiosis is a rare and severe form characterized by vascular lesions that increase the risk of mortality. Diagnosing bbesiosis can be challenging, as the parasite is not always detectable in blood smears due to low parasitemia. PCR is one of the most sensitive diagnostic methods for detecting *Babesia* species in dogs. However, false-negative PCR results have been reported in subclinical and chronic cases of babesiosis, attributed to undetectable levels or the elimination of the parasite from circulating blood. In such cases, serology may be the only method to reveal these infection profiles. Serological tests are widely used for screening and monitoring due to their broad diagnostic window, as antibodies may persist for months or years after infection. For better interpretation, serological test results should be associated with clinical signs, history of vector exposure, and suggestive hematological changes.

BIBLIOGRAPHIC REFERENCE

- 1- BOOZER AL, MACINTIRE DK. Canine babesiosis. Vet Clin North Am Small Anim Pract, v. 33, n. 4, p. 885-905, viii, 2003.
- 2- KOSTER L, LOBETTI R, KELLY P. Canine babesiosis: a perspective on clinical complications, biomarkers, and treatment. Vet Med (Auckl), v.6, p.119-128, 2015.
- 3- LEISEWITZ AL, MRLJAK V, DEAR JD, BIRKENHEUER A. The Diverse Pathogenicity of Various Babesia Parasite Species That Infect Dogs. Pathogens, v. 12, n.12, p. 1437, 2023.
- 4- PANTI-MAY JA, RODRIGUEZ-VIVAS RI. Canine babesiosis: A literature review of prevalence, distribution, and diagnosis in Latin America and the Caribbean. Vet Parasitol Reg Stud Reports, v. 21, p. 100417, 2020.
- 5- QUIBASA: Data from the Research and Development Department.

QUALITY GUARANTEE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product Licensed at the Ministry of Agriculture, Livestock and Supply since 12/12/2019 under number 10.291/2019.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611).

Revision: January/2025

REF

CATALOG NUMBER

LOT

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

STERILE R

PRODUCT
STERILIZED

DANGER