## VETLISA BABESIA IgG

**REF VET038**

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#### FINALIDADE

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- 1- TABOADA J, MERCHANT SR. Babesiosis of Companion Animals and Man. Veterinary Clinics of North America: Small Animal Practice. n.1. v. 21, p. 103-123. Janeiro, 1991.
- 2- COMAZZI S, PALTRINIERI S, MANFREDI MT, AGNES F. Diagnosis of canine babesiosis by Percoll gradient separation of parasitized erythrocytes. Journal of Veterinary Diagnostic and Investigation. v.11, p. 102-104. 1999.
- 3- BOOZER AL, MACINTIRE DK. CANINE BABESIOSIS. Vet Clin Small Anim, v.33, p.885-904, 2003.
- 4- BRANDÃO LP, HAGIWARA MK, MYIASHIRO SI. Humoral immunity and reinfection resistance in dogs experimentally inoculated with *Babesia canis* and either treated or untreated with imidocarb dipropionate. Vet Parasitol., v.114, p.253-265, 2003.
- 5- VERCAMMEN F, DE DEKEN R, MAES L. Duration of protective immunity in experimental canine babesiosis after homologous and heterologous challenge. Veterinary Parasitology, v.68, p. 51-55, 1997.
- 6- Dados do Departamento de Pesquisa e Desenvolvimento.

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**1-** TABOADA J, MERCHANT SR. Babesiosis of Companion Animals and Man. Veterinary Clinics of North America: Small Animal Practice. n.1. v. 21, p. 103-123. Janeiro, 1991.

**2-** COMAZZI S, PALTRINIERI S, MANFREDI MT, AGNES F. Diagnosis of canine babesiosis by Percoll gradient separation of parasitized erythrocytes. Journal of Veterinary Diagnostic and Investigation. v.11, p. 102-104. 1999.

**3-** BOOZER AL, MACINTIRE DK. CANINE BABESIOSIS. Vet Clin Small Anim, v.33, p.885-904, 2003.

**4-** BRANDÃO LP, HAGIWARA MK, MYIASHIRO SI. Humoral immunity and reinfection resistance in dogs experimentally inoculated with *Babesia canis* and either treated or untreated with imidocarb dipropionate. Vet Parasitol., v.114, p.253-265, 2003.

**5-** VERCAMMEN F, DE DEKEN R, MAES L. Duration of protective immunity in experimental canine babesiosis after homologous and heterologous challenge. Veterinary Parasitology, v.68, p. 51-55, 1997.

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### INSTRUCTIONS FOR USE

#### FUNCTION

Test for qualitative determination of IgG antibodies against *Babesia vogeli* using dog serum or plasma, by solid phase enzyme immunoassay, in microplate. For *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Enzyme immunoassay or immunoenzymatic  
The VETLISA BABESIA IgG is an immunoenzymatic assay in solid phase, based on the principle of immunocapture for the qualitative detection of IgG antibodies against *Babesia vogeli* in serum or plasma. Antibodies against *Babesia vogeli* present in the sample, bind to the recombinant antigen of *Babesia vogeli* coated on the microplate, forming immobilized complexes: *Babesia vogeli* antigen – anti-*Babesia vogeli* antibodies. After the initial incubation, the microplate is washed to eliminate unbound materials. The conjugate, formed by the anti-dog IgG antibody conjugated to peroxidase, binds the immobilized antigen-antibody complex on the plate. After the second incubation, the microplate is washed and incubated with substrate, and the blue color results from the interaction with the antibodies against *Babesia vogeli* present in the sample. The stop solution is added to stop the reaction, promoting a color change to yellow, measured in a microplate reader.

#### REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8°C. Contains: Polystyrene plate, divided into 12 strips of 8 wells each, impregnated with *Babesia vogeli* antigen.
- 2- Concentrated Conjugate (100x)** - Store between 2 and 8°C. Contains: Peroxidase-linked anti-dog IgG antibody solution.
- 3- Concentrated Wash (20x)** - Store between 2 and 8°C. Contains: Phosphate Buffer Solution, surfactant and preservative.
- 4- Diluent** - Store between 2 and 8°C. Contains: Buffer solution, surfactant, stabilizer and preservative.
- 5- Substrate TMB** - Store between 2 and 8°C. Contains: Solution containing Tetramethylbenzidine (TMB < 1.0 mg / mL), Citric Acid Solution <5% and Urea Peroxide <1%.
- 6- Stop Solution** - Store between 2 and 8°C. Contains: 1N Sulfuric Acid.
- 7- Negative Control** - Store between 2 and 8°C. Contains: Serum containing negative serum for IgG antibodies against *Babesia vogeli* with preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8°C. Contains: Serum containing positive serum for IgG antibodies against *Babesia vogeli* with preservative. **Potentially infectious.**
- 9- Plate Sealers**

#### PRESENTATIONS

COMPONENTS	PRESENTATION		
	1	2	3
	96 cavities	192 cavities	480 cavities
<b>1- Sensitized Plate</b>	1 unit	2 units	5 units
<b>2- Concentrated Conjugate (100x)</b>	1 x 300 µL	2 x 300 µL	5 x 300 µL
<b>3- Concentrated Wash (20x)</b>	1 x 50 mL	2 x 50 mL	5 x 50 mL
<b>4- Diluent</b>	1 x 60 mL	2 x 60 mL	5 x 60 mL
<b>5- Substrate TMB</b>	1 x 15 mL	2 x 15 mL	5 x 15 mL
<b>6- Stop Solution</b>	1 x 15 mL	2 x 15 mL	5 x 15 mL
<b>7- Negative Control</b>	1 x 1 mL	2 x 1 mL	5 x 1 mL
<b>8- Positive Control</b>	1 x 1 mL	2 x 1 mL	5 x 1 mL
<b>9- Plate Sealers</b>	3 units	6 units	15 units

#### EQUIPMENTS AND OPERATIONAL INPUTS

##### Materials contained in the kit:

- Reagents described in the previous table.
- Instructions for use (manual).

##### Materials needed, but not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 5 to 500 µL with a variation coefficient of less than 1.5%.
- 2- Repipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5- Absorbent paper to dry the wells.
- 6- Stopwatch or watch.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator at 37 °C ± 2 °C.

#### STORAGE AND TRANSPORT CONDITIONS

The product storage temperature should be 2 to 8 °C. Transport can be carried out at room temperature (30°C) for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

#### SPECIAL CARE

- 1- For professional veterinary *in vitro* diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate must be opened only after reaching room temperature. Replace strips from unused microwells in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be fresh and free of contaminants.
- 5- Saturated deionizer columns release alkaline water, different ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Sulfuric Acid which is a strong acid. So

handle it with due care.

7- The handling of any product that contains serum is potentially capable of transmitting diseases. Therefore, it is necessary to take proper biosafety precautions when handling these products.

8- Pipette the reagents always in the same order to minimize the difference in reaction time between the microwells.

9- For protection measure, the plate must be covered during the reaction.

10- Ensure that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow wells to dry during the assay.

11- Do not expose reagents, especially Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.

12- We recommend applying local, state and federal regulations for environmental protection so that the disposal of reagents and biological material is done in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Service of Customer Service) by Quibasa.

14- Do not use the product in case of damage to the packaging.

15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

##### Serum or Plasma (EDTA or Heparin)

Hemolysis, lipemia and bilirubin (jaundice) can interfere with test result. The result obtained from samples in these conditions must be assessed with caution. The samples can be kept refrigerated, between 2 and 8 °C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at a temperature of -20 °C.

#### PROCESS DESCRIPTION

##### Stability After Open

The results of the stability test show that the VETLISA BABESIA IgG kit is stable after being opened up to 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using internal controls of the kit and the criteria for validating the technique.

**ATTENTION: The Positive and Negative Controls are ready to use.**

#### PREPARATION OF WORK REAGENTS

##### Washing Solution

Dilute the contents of vial N° 3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for up to 30 days. It can be stored at room temperature. If crystallization occurs, heat to 37 °C until dissolved.

##### Substrate

The Substrate is ready for use.

##### Conjugate Solution

Dilute the Concentrated Conjugate (Reagent N° 2) in the ratio of 1:101 in Diluent (Reagent N° 4). Prepare the solution at the moment of performing the test.

To perform an assay using all cavities in the kit, mix 110 µL of

## VETERINARY USE

Concentrated Conjugate in 11 mL of Diluent.

To perform an assay using 8 cavities (1 strip), mix 10 µL of the Concentrated Conjugate in 1 mL of Diluent.

**IMPORTANT:** The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

##### Sample Dilution

In a test tube, dilute 15 µL of the sample in 300 µL of Diluent (Reagent N° 4), if performing the assay in duplicate. Cap tube and vortex gently or mix manually by inversion. Dilutions cannot be stored.

#### TECHNIQUE

##### For use in automatic equipment, consult the Customer Service Department (SAC).

Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 - 30°C) for at least 30 minutes. Return the unused strips to the original sealed packaging.

1- Separate the wells to be used considering: Controls and Samples ((it is recommended to test in duplicate).

2- Separate the first well for the Blank (OPTIONAL).

3- Pipette 100 µL of the Negative Control and the Positive Control into the previously determined cavities.

**Note: The controls are ready to use, there is no need to dilute them.**

4- Pipette 100 µL of the previously diluted samples into the previously determined cavities. In the Blank cavity, if you have chosen, pipette only 100 µL of the diluent

5- Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.

6- Incubate for 30 minutes at 37 °C ± 2 °C.

7- Discard the content of the cavities by aspiration (washing machine) or by decantation (manual). Use approximately 300µL/well of previously diluted Wash Solution, for a total of five (5) wash cycles. Shake for three seconds after each wash cycle. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.

**IMPORTANT:** Poor washing/drying can cause inadequate results.

8- Pipette 100 µL of the Conjugate previously diluted into each well, even in the Blank cavity.

9- Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.

10- Incubate for 30 minutes at a 37 °C ± 2 °C.

11- Repeat item 7.

12- Add 100 µL of Substrate TMB into each well.

13- Homogenize gently for ± 3 seconds. Cover the cavities with plate sealer.

14- Incubate for 10 minutes at 37°C, protected from light.

15- Remove the plate sealer from the cavities.

16- Pipette 100 µL of Stop Solution into all wells.

17- Homogenize gently for ± 3 seconds

18- Read the absorbance in an ELISA reader on a dual filter 450 nm (primary filter) / 630nm (secondary filter) within a maximum of 10 minutes after adding the Stop Solution.

## TECHNIQUE VERIFICATION

Verify if the results obtained from reading the Controls and the Blank are compatible with the specification above:

ITEM	ABSORBANCE (DUAL FILTER)
Blank	< 0.050
Negative Control	0.100 a 0.300
Positive Control	> 0.800

If the values are out the expected ones, the test should be repeated.

## Calculations

Calculate Cut Off according to the following formula:

Cut Off = Average Negative Control Absorbance + 0.300

Example:

ITEM	ABSORBANCIA
Negative Control	A1 = 0.208 A2 = 0.212
Cut Off = Negative Control Average Absorbance + 0.300	(0.208 + 0.212) / 2 + 0.300 = 0.510

Calculate the index by dividing the absorbance of the sample by the cut off value.

Example:

ITEM	ABSORBANCE
Sample	0.992
Negative Control	0.210
Cut Off Value	0.510
Index: Sample Abs. / Cut Off Value	0.992 / 0.510 = 1.95

## INTERPRETATION OF THE RESULT

After calculating the sample index, consider the indices below to determine the results.

RESULTS	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	≥ 0.9 y ≤ 1.1

**Note:** An indeterminate result indicates an assay response within the indeterminate range (gray zone between the cutoff points) and, on its own, neither confirms nor excludes infection. It may reflect the immunological window, onset of seroconversion, declining antibodies, or pre-analytical interferences. The interpretation should consider the history and clinical presentation, with repeat testing (preferably on a new sample), resampling in 7–14 days, and, depending on the suspicion and stage of the disease, confirmatory testing. The diagnostic conclusion rests with the veterinarian, integrating clinical signs and complementary tests.

## PROCEDURE LIMITATIONS

The presence of IgG antibodies against *Babesia vogeli* should be evaluated along with clinical information and other data from complementary tests. IgG production usually begins between 15 and 20 days after infection and persists for up to 2 years. The presence of IgG indicates that the dog has demonstrated an immune response to exposure to *Babesia vogeli*; however, the absence of antibodies does not rule out infection. Factors such as the stage of the disease, the individual's immune response, the amount of antibodies available in the sample, and the presence of interfering factors should be considered when interpreting the result.

## INTERNAL QUALIT CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement

systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

## PRODUCT PERFORMANCE

### ACCURACY

#### Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following results:

REPEATABILITY	ABSORBANCE		
	1	2	3
Mean	0.780	0.980	0.090
Standard Deviation	0.037	0.077	0.002
Coefficient of Variation (%)	4.730	7.800	2.604

#### Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following results:

REPRODUCTIBILITY	ABSORBANCE		
	1	2	3
Mean	0.780	1.000	0.093
Standard Deviation	0.050	0.049	0.005
Coefficient of Variation (%)	6.414	4.879	5.072

## CLINICAL SENSITIVITY AND SPECIFICITY

The kit VETLISA BABESIA IgG analyzed clinical samples compared to other ELISA method. The results show that the kit VETLISA BABESIA IgG presents a clinical sensitivity of 97.4% and a clinical specificity of 97.8%.

Method	Reference		Total	
	Result	Positive		Negative
VETLISA BABESIA IgG	Positive	88	1	89
	Negative	2	37	39
Total Result		90	38	128

**Clinical Sensitivity:** 97.4% (37/38) (IC 95%; 92.34 a 100%)

**Clinical Specificity:** 97.8% (88/90) (IC 95%; 94.77 a 100%)

**Precision:** 97.7% [(88 + 37) / (90 + 38)]

## CLINICAL SIGNIFICANCE

Babesiosis is a disease caused by protozoa of the genus *Babesia*, which can affect domestic dogs, wild animals, and humans. Transmission occurs primarily through the bite of infected ticks during a blood meal. The incubation period between infection and the onset of clinical signs ranges from 7 to 30 days, and the disease can manifest as peracute, acute, chronic, or asymptomatic. The most frequently observed clinical signs include anemia, thrombocytopenia, hemoglobinuria, fever, pale mucous membranes, tachycardia, tachypnea, depression, anorexia, jaundice, and splenomegaly. In chronic cases, progressive weight loss and persistent anorexia are common. Due to the clinical similarities between different hemoparasitoses and the possibility of coinfections, since they share the same vector, the use of complementary diagnostic tools, such as serological tests, is essential. These tests should be used to support clinical evaluation, never in isolation. Laboratory results, combined with the history and clinical examination, should be interpreted carefully, with the diagnostic conclusion always being the veterinarian's responsibility.

## BIBLIOGRAPHIC REFERENCES

1- TABOADA J, MERCHANT SR. Babesiosis of Companion Animals and Man. Veterinary Clinics of North America: Small Animal Practice. n.1. v. 21, p. 103-123. Janeiro, 1991.

2- COMAZZI S, PALTRINIERI S, MANFREDI MT, AGNES F. Diagnosis of canine babesiosis by Percoll gradient separation of parasitized erythrocytes. Journal of Veterinary Diagnostic and Investigation. v.11, p. 102-104. 1999.

3- BOOZER AL, MACINTIRE DK. CANINE BABESIOSIS. Vet Clin Small Anim, v.33, p.885-904, 2003.

4- BRANDÃO LP, HAGIWARA MK, MYIASHIRO SI. Humoral immunity and reinfection resistance in dogs experimentally inoculated with *Babesia canis* and either treated or untreated with imidocarb dipropionate. Vet Parasitol., v.114, p.253-265, 2003.

5- VERCAMMEN F, DE DEKEN R, MAES L. Duration of protective immunity in experimental canine babesiosis after homologous and heterologous challenge. Veterinary Parasitology, v.68, p. 51-55, 1997.

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## QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda.**

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Industria Brasileira

## CONSUMER SERVICE

Customer Support Service

Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

**Product licensed in the Ministry of Agriculture, Livestock and Supply since 10/04/2019 with the number 10.235/2019.**

**Technical manager:** Dra. Camila Eckstein (CRMV/MG 20611)

**Review:** January/2026

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER