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REF K006

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- 1 - Representatives of the American Academy of Pediatrics, the College of American Pathologists, The American Association of Clinical Chemists, and the National Institutes of Health: Recommendation on a Uniform Bilirubin Standard. Standard Methods in Clinical Chemistry. Edited by S. Meites. New York, Academic Press, 1965, Vol. 5, p. 75.
- 2 - Recommendations of the College of American Pathologists Standards Committee: Aunifor Bilirubin Standard. Am. J. Clin. Pathol. 39.90, (1963).
- 3 - WINKELMAN, J., CANNON, D. C., JACOBS, S. L.: Liver Function Tests, including Bile Pigments, In Clinical Chemistry – Principles and Technics, 2nd. RJ Henry, DC Cannon, JW Winkelman, Editors Haper & How, Hagerstown, MMD, 1974, p. 1053 – 1062.
- 4 - DOUMAS, B. T., PERR Y, B. W. SASSE, E. A., STRAUMFJORD Jr., J. V., Clin. Chem. 19: 984, 1973.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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REFERENCIAS BIBLIOGRAFICAS
1 - Representatives of the American Academy of Pediatrics, the College of American Pathologists, The American Association of Clinical Chemists, and the National Institutes of Health: Recomendation on a Uniform Bilirrubin Standard. Standard Methods in Clinical Chemistry. Edited by S. Meites. New York, Academic Pres, 1965, vol. 5, p. 75.
2 - Recommendations of the college of american Pathologists Standards Commitee: Aunifor Bilirrubin Standard. Am. J. Clin. Pathol. 39.90, (1963).
3 - WINKELMAN, J., CANNON, D. C., JACOBS, S. L.: Liver Function Tests, including Bile Pigments, In Clinical Chemistry- Principales and Technics, 2nd. RJ Henry, DC Cannon, JW Winkelman, Editors Haper & How, Hagerstown, MMD, 1974, p. 1053-1062.
4 - DOUMAS, B. T., PERRY, B. W. SASSE, E. A., STRAUMFJORD J., J. V., Clin. Chem. 19: 984, 1973.
5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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REF K006

INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Reagent for Bilirubin calibration. Colorimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Bioclin uses a Standard following rigorously the specifications of the “Joint Bilirubin Committee” and presents a high level of purity, being dosed by attachment with the Sulphanilic Acid diazotate, in the presence of Caffeine and Sodium Benzoate.

REAGENTS

Number 1 - Standard - Store between 15 and 30°C. Contains: Bilirubin 0.3 mg. After the reconstitution the concentration equals to 10 mg/dL. Do not remove the paper seal that protects the Standard from the action of light.
Number 2 - Diluent - Store between 15 and 30°C. Contains: Dimethyl Sulfoxide P.A.

PRESENTATION

Reagent	Volume
Nº 1 (Desiccated)	0.30 mg
Nº 2	4 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, watch and stopwatch, pipettes, test tubes, and Bioclin Bilirubin kit. They can be found at markets specialized in laboratories of clinical analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be between 15 and 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers.
- 5- The calibration shall be repeated periodically to verify any change in the colorimeter response or spectrophotometer.
- 6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

- 7- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 8- Do not use the product in case of damaged packaging.
- 9- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

PREPARATION OF STANDARD (10 mg/dL)

The Bilirubin is high photosensitivity, **store the bottle protected from light**, during the standard preparation. Pipette 3 mL of the Diluent (Reagent Nº 2) and add the bottle of the Standard (Reagent Nº 1). Agitate well and leave it resting for 5 minutes. Agitate again to homogenize. This Standard is stable for 12 hours at environment temperature and protected from light. At 20°C negative is stable for 2 months. The standard can only be freeze and defrosted once. The final concentration of the Standard solution is 10 mg/dL.

Note: During the production process, upon solvent evaporation, the bilirubin may be dispersed and impregnated over the entire surface of the vial. Therefore, the amount of powder apparent in the vial prior to reconstitution may be imperceptible. In this way, the reconstitution must be performed with care, following the guidelines described above.

TECHNIQUE

Carry out technique along with the **Bioclin Bilirubin*** kit, using the Reagent Nº 2 and Reagent Nº 3 and the Diazo Reagent for the technique execution.

MACRO TECHNIQUE - Carry out the Standard in triplicate.

	Blank	Standard (Triplicate)
Reagent Nº 3*	4.5 mL	4.5 mL
Reagent Nº 2*	500 µL	---
Diazo Reagent*	---	500 µL
Standard	300 µL	300 µL

Homogenize well, let it rest for 5 minutes. Determine the absorbancies of Standard at 525 nm (500 - 550 nm), adjusting the zero with the Blank.

MICRO TECHNIQUE - Carry out the standard in triplicate.

	Blank	Standard (Triplicate)
Reagent Nº 3*	1.8 mL	1.8 mL
Reagent Nº 2*	150 µL	---
Diazo Reagent*	---	150 µL
Standard	50 µL	50 µL

Homogenize well and let it rest for 5 minutes. Determine the absorbancies of Standard at 525 nm (500 - 550 nm), adjusting the zero with the Blank.

CALCULATIONS

It is important to dose the Standard in triplicate, using the average absorbance to calculate the factor. Determine a factor for the Micro technique and another for the Macro technique, because the dilution is different in each technique.

Calibration Factor = $\frac{\text{Standard Concentration (10.0 mg/dL)}}{\text{Standard Absorbance (average)}}$

mg/dL = Sample Absorbance x Calibration Factor

Example:

Standard average Absorbance = 0.480
Sample Absorbance = 0.045

Calibration Factor = $\frac{10 \text{ mg/dL}}{0.480}$ = 20.83

mg/dL = 0.045 x 20.83 = 0.93

PROCEDURE LIMITATIONS

Dose the Standard in triplicate, using the absorbancies average to calculate the factor. To use this average, the differences between the absorbancies shall be inferior to 5%.

INTERFERENTES

Light can interfere in the Bilirubin Standard dosage.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

Bioclin recommends the usage of the Bioclin Bilirubin Calibration kit, which is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 916A.

PRODUCT PERFORMANCE

ACCURACY

Recovery

The recovery analysis was made with 5 Standard determinations. The accuracies were calculated and if they are found in a good concordance with the reference values, obtaining a recovery between 95 and 105%

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations using Bilirubin Standard, obtaining the following results:

Repeatability	Standard
Average Concentration(mg/dL)	10.285
Standard Deviation (mg/dL)	0.075
Coefficient of Variation (%)	0.725

Reproducibility

The reproducibility was calculated from 20 successive determinations using Bilirubin Standard, for 3 consecutive days, obtaining the following results:

Reproducibility	Standard
Average Concentration(mg/dL)	10.290
Standard Deviation (mg/dL)	0.005
Coefficient of Variation (%)	0.049

DIAGNOSTIC SIGNIFICANCE

The Bilirubin lift on plasma in the presence of hepatic injuries, bile obstruction or when the speed of destruction of red globule is increased.
The increase of the Indirect Bilirubin is observed on the hemolytic syndrome, neonatal jaundice, Cligler-Najjar syndrome and Gilbert disease.
The Direct Bilirubin is increased on acute and chronicle hepatitis, the toxic reactions to many drugs (Chlorpromazine, Arsenicals Organics, Methyltestosterone), in bile obstruction treatment.

BIBLIOGRAPHIC REFERENCES

- 1 - Representatives of the American Academy of Pediatrics, the College of American Pathologists, The American Association of Clinical Chemists, and the National Institutes of Health: Recomendation on a Uniform Bilirrubin Standard. Standard Methods in Clinical Chemistry. Edited by S. Meites. New York, Academic Press, 1965, Vol. 5, p. 75.
- 2 - Recommendations of the College of American Pathologists Standards Commitee: Aunifor Bilirrubin Standard. Am. J. Clin. Pathol. 39.90, (1963).
- 3 - WINKELMAN, J., CANNON, D. C., JACOBS, S. L.: Liver Function Tests, including Bile Pigments, In Clinical Chemistry – Principles and Technics, 2nd. RJ Henry, DC Cannon, JW Winkelman, Editors Haper & How, Hagerstown, MMD, 1974, p. 1053 – 1062.
- 4 - DOUMAS, B. T., PERRY, B. W. SASSE, E. A., STRAUMFJORD Jr., J. V., Clin. Chem. 19: 984, 1973.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte/MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for the Bilirubin Calibration kit:
10269360090

Review: August/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER