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REF K107

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FINALIDADE

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REF K107

INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for determination of Direct Bilirubin in biological samples serum or plasma (EDTA or Heparin). Colorimetric test, for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Colorimetric Test - DPD

Direct Bilirubin, through the coupling reaction with Dichlorophenyl Diazonium Tetrafluoroborate forms an azocompound, a red colored complex with maximum absorption at 546 nm. The intensity of color formed is directly proportional to the concentration of Direct Bilirubin in the sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Sulfamic Acid < 200 mmol/L and stabilizers.

Number 2 - Color Reagent - Store between 2 and 8°C. Contains: Dichlorophenyl Diazonium Tetrafluoroborate < 5 mmol/L and stabilizer.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2
1	1 x 40 mL	1 x 10 mL
2	2 x 40 mL	2 x 10 mL
3	3 x 40 mL	3 x 10 mL
4	4 x 40 mL	4 x 10 mL
5	5 x 40 mL	5 x 10 mL
6	2 x 40 mL	1 x 20 mL
7	4 x 40 mL	2 x 20 mL
8	1 x 24 mL	1 x 6 mL
9	2 x 24 mL	2 x 6 mL
10	3 x 24 mL	3 x 6 mL
11	4 x 24 mL	4 x 6 mL
12	5 x 24 mL	5 x 6 mL
13	6 x 24 mL	6 x 6 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Biocal kits and Biocontrol N and P Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- The reagent Nº 2 must be kept protected from light.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Calibration must be repeated periodically to verify any alterations in colorimetric responses or spectrophotometer.

6- Hemolysis, even if discreet, interfere in dosage.

7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

9- Do not use the product in case of damaged packaging.

10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma (EDTA or Heparin) free of hemolysis.

Sample must be protected from light.

The samples is stable for 3 days between 2 and 8°C and 3 months at -20°C.

Sample for therapeutic control should always be collected at the same hour.

PROCESS DESCRIPTION

The calibration stability of the Bilirubin Direct Automated kit installed on refrigerated equipment is at least 30 days.

This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

PREPARATION OF REAGENTS

Reagents are ready for use.

TECHNIQUE

For the use of the kit, use as calibrator the Biocal Bioclin kit and as control serum the Biocontrol N and P Bioclin kit.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

INTERFERENCES

Hemolysis, even if slight, interfere with measurement of Bilirubin. No interference was observed for Ascorbic Acid up to 30 mg/dL and lipemia up to 1000 mg/L from Triglycerides. According to bibliographic data, some drugs, as Anabolic Steroids, Ascorbic Acid, Salicylates and Vitamin A, can raise bilirubin values. Caffeine and Penicillins may lead to falsely decreased.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to

variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitored by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 916A.

REFERENCE VALUES

Reference values were obtained through determinations of Direct Bilirubin in healthy populations of male and female.

For all ages:

Direct Bilirubin: up to 0.4 mg/dL

To convert the values to mg/dL to mmol/L (SI), multiply the results obtained by 0.0171.

These values should be used only as a guide. Each laboratory should establish its own reference range, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

Bilirubin Direct Automated kit was compared with other methods for determination of Direct Bilirubin dosage commercially available measurement. 42 analyzes were performed and the results were evaluated. The equation linear obtained was $Y = 1.0228X + 0.0034$ and the correlation coefficient 0.9916. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0.42	0.69	0.58
Standard Deviation (mg/dL)	0.01	0.01	0.01
Coefficient of Variation (%)	1.50	1.13	2.14

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0.43	0.70	0.60
Standard Deviation (mg/dL)	0.01	0.03	0.02
Coefficient of Variation (%)	1.51	4.75	3.43

Sensitivity

The sensitivity was calculated from 40 determinations of a sample free of Direct Bilirubin. The average found was 0.055 mg/dL, with deviation standard of 0.005 mg/dL. The sensitivity, which indicates the method limit detection, corresponds the average plus 3 times the standard deviation, and is equal to 0.070 mg/dL.

Linearity

The reaction is linear to the concentration of 12 mg/dL. For higher values dilute the sample with NaCl 0.85% and repeat the dose. Multiply the resulting number by the dilution factor used.

Diagnostic Significance

Bilirubin is a breakdown product of hemoglobin in reticuloendothelial system. It is conjugated in the liver to the then be excreted in bile. The level of Bilirubin is useful for differential diagnosis of hepatobiliary disease and other causes of jaundice which becomes clinically evident when the concentration of Total Bilirubin greater than 2.5 mg/dL. The Direct Bilirubin rises in plasma in the presence of hereditary liver diseases such as disease Dubin-Johnson and Rotor, hepatocyte injury (viral, toxic or alcoholic), biliary obstruction (stones, neoplasm), acute or chronic hepatitis and toxic reactions to various drugs (such as Chlorpromazine, and organic arsenic compounds, Methyltestosterone, among others). Levels of Direct Bilirubin over 50% of the Total Bilirubin concentration indicate the existence of post-liver. The use of drugs that activate the hepatic microsomal system can reduce the concentration of Total and Direct Bilirubin in plasma.

Bibliographic References

1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998:192-202.
2. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors, Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1125-77.
3. Rand RN, di Pasqua A. A new diazo method for the determination of bilirubin. Clin Chem 1962; 6:570-8.
4. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

Quality Assurance

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

Customer Service

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Bilirubin Direct Automated kit: 10269360176

Review: August/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER