BIOCLIN FAST COVID-19 Ag

REF K244

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FINALIDADE

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- REFERÊNCIAS BIBLIOGRÁFICAS**
- 1 - Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
 - 2 - CENTERS FOR DISEASE CONTROL AND PREVENTION et al. Interim guidelines for collecting, handling, and testing clinical specimens for COVID-19. Retrieved June, v. 19, p. 2020, 2020. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13. A.Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
 - 3 - Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
 - 4 - Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb15;395(10223):497-506.
 - 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1 - Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- 2 - CENTERS FOR DISEASE CONTROL AND PREVENTION et al. Interim guidelines for collecting, handling, and testing clinical specimens for COVID-19. Retrieved June, v. 19, p. 2020, 2020. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13. A.Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- 3 - Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- 4 - Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb15;395(10223):497-506.
- 5- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Test for the qualitative determination of antigens for SARS-CoV-2, virus which causes COVID-19 (severe acute respiratory syndrome coronavirus 2) in biological **nasopharyngeal swab and nasal swab** samples. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The BIOCLIN FAST COVID-19 Ag kit is an immunochromatographic test that detects antigens of the SARS-CoV-2 Virus employing mouse monoclonal antibodies anti N protein (nucleocapsid) and RBD (receptor binding domain) present in the sample. The test line contains Anti-N and Anti-RBD monoclonal antibodies. After adding the sample, colloidal gold particles conjugated Anti-N and Anti-RBD bind to viral antigens and migrate by capillarity across the membrane.

The colloidal gold complex and viral antigen are captured in the test line (T), leading to the formation of a reddish line indicating the presence of SARS-CoV-2 in the sample. The non-formation of a reddish line in the T region indicates the absence of a viral particle in the sample or an amount of Virus below the detection limit. The sample-antigen mixture continues running through the membrane to the control line, where a reddish line is formed, both for positive and negative samples. The formation of the reddish line in the Control region (C) demonstrates that the test is suitable for use.

REAGENTS

Number 1 (R1) Cassette – Store at 2 to 30°C. **Do not freeze.** Contains:

A) 01 plastic device with 01 test strip containing:
 - Control region (C): Goat antibody Anti-mouse (> 0.06 µg).
 - Test region (T): Anti-N (> 0.3 µg), Anti-RBD (> 0.3 µg) - Conjugate:
 Colloidal gold particles conjugated with mouse monoclonal antibody Anti-N and Anti-RBD.

B) 01 silica sachet

Number 2 (R2) Collection Tube - Store at 2 to 30°C. **Do not freeze.** Contains: Buffer < 100 mmol/L pH 7.4, and preservatives.

Accessories for Sample Collection:

- Dropper cap (attached to the collection tube)
- Disposable Swab.

PRESENTATION

Reagent	N° of Cassettes per packaging	N° of Collection tube per packaging
Cassette + Collection Tube (Included dropper cap and disposable swab) *	05	05
	10	10
	25	25
	50	50
	100	100

Reagent	N° of tests per packaging
Aluminized pouch, contains: 1 Cassette 1 Collecting Tube 1 Dropper Cap 1 Disposable Swab	05
	10
	25
	50
	100

* The number of accessories is proportional to the number of cassettes and collecting tubes.

INSTRUMENTS AND OPERACIONAL INPUTS

Clock or stopwatch, PPE (glove, mask, goggles), and garbage for disposal. All the inputs listed are in the specialized market for articles for Clinical Analysis Laboratories.

CONDITIONS OF STORAGE AND TRANSPORTATION

The storage and transport temperature should be 2 to 30°C. Avoid moisture. **Do not freeze.**

WARNINGS

- 1 - For *in vitro* diagnostic use only.
- 2 - Strictly follow the methodology proposed to obtain exact results.
- 3 - The used device has potential infectious risk. The safety routine must be followed strictly, and the use of disposable gloves is indispensable.
- 4 - Use clean and dry recipients to settle the samples.
- 5 - Open the aluminum envelope to take the cassette off, only for its immediate use.
- 6 - The kit must not be frozen or exposed to temperatures above 45°C. That exposure may cause its deterioration.
- 7 - The reagents and the sample must reach the environment temperature before the test begins.
- 8 - Do not mix reagents between kits with the same batch or with different batch numbers.
- 9 - Do not use the kit after the expiration date. The expiry date is indicated on the product label.
- 10 - The cassette (test device), collecting tube, dropper cap, and swab are disposable and cannot be reused.
- 11 - The kit materials that had contact with the patient's biological material must have a special destination. The collecting tubes and swabs and Cassettes (test devices) must be discarded in biological waste.
- 12 - Boxes, instructions for use, and aluminized pouches can be disposed of in regular trash.
- 13 - We recommend applying local, state, and federal environmental protection standards so that the disposal of reagents and biological material is done following current legislation.
- 14 - To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by Quibasa's Customer Service Department.
- 15 - Do not use the product in case of damage to the packaging.

SAMPLES

Nasopharynx swab or nasal swab. The sample must be tested preferably after collection or stored between 2-8°C for 24 hours.

DESCRIPTION OF THE TECHNICAL PROCESS

The BIOCLIN FAST COVID-19 Ag cassette's stability (after opening the aluminized pouch) is 2 hours. It is stability may vary according to the environmental conditions after opening the aluminized pouch.


TECHNIQUE

A) Segregate the collection material: collecting tube (dropper cap attached), swab.

B) Sampling

- Use the swab included in the kit to collect the nasopharynx or nasal samples.
- The sample should be tested preferably after collection.


Nasopharynx Swab


- 1- Ask the patient to tilt his head slightly.
- 2- Carefully insert the swab into the patient's nostril until the swab head meets resistance against the nasopharynx surface.
- 3- Rotate the swab 5 times, in contact with the posterior nasopharynx, remove the swab.
- 4- Repeat steps 1, 2, and 3 in the second nostril.

Nasal Swab

- 1- Ask the patient to blow his nose.
- 2- Insert approximately 2 cm from the tip of the swab into one nostril.
- 3- Using circular movements, scrape the swab head 5 times into the patient's nostril.
- 4- Repeat steps 2 and 3 for the other nostril.


- Ensure that the collection tube's buffer is at the bottom of the tube and remove the dropper cap.
- After collection, insert the swab into the collection tube.
- Mix for 10 seconds to ensure good extraction through the buffer.
- Remove excess buffer at the tip of the swab by pressing the walls of the collection tube against the head of the swab.
- Insert the dropper cap into the collection tube.
- Add three drops of the buffer (90-100 µL) to the sample well. Read the test between 15 and 20 minutes. Do not interpret after 20 minutes.


INTERPRETATION OF THE RESULT


Reagent Test: Formation of a reddish line in the test region (T) and another line in the control region (C). Do not interpret after 20 minutes

Non-Reagent Test: Formation of a reddish line in the control region (C) and complete absence of a red line in the test region (T). Do not interpret after 20 minutes.

Invalid Result: Complete absence of a line in the region control (C) with or without a red line in the test region (T). Test the sample again.

- Notes:**
- 1 - Results should only be read between 15 and 20 minutes, and not before or after that time.
 - 2 - The line in the control region can appear before the test is incubated for 10 minutes, which does not mean that negative results can be interpreted ahead of time.
 - 3 - The results should not be interpreted after 20 minutes to avoid incorrect results.
 - 4 - SARS-CoV-2 virus particles are most commonly present during the acute phase of the disease. It may be present even after the appearance of antibodies against the virus. As it is a new disease, there is no consensus on the exact moment of viral load reduction with the formation of an immune response.
 - 5 - The responsible medical professional must interpret the results provided by this test, and it is not the only criterion for determining the diagnosis or treatment of the patient.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have internal quality control, in which all procedures, rules, limits, and tolerance to variations must be established clearly. It is essential to mention that all measurement systems present an analytical variety, and the laboratory must monitor it. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

LIMITATIONS OF THE PROCESS

- 1 - The BIOCLIN FAST COVID-19 Ag kit is a qualitative test and should only be used for in vitro diagnostic use.
- 2 - The BIOCLIN FAST COVID-19 Ag Kit determines only the presence of SARS-CoV-2 antigens and should not be used in isolation as a diagnostic criterion. Reagent results should be evaluated with other clinical information available. The isolated results of the test cannot be used for a definitive diagnosis.
- 3 - **A non-reactive result does not exclude the possibility of infection with the SARS-CoV-2 Virus. Viral antigen standards may be below the kit's analytical sensitivity, 2 x 10^{3.1} TCID50/mL (Tissue Culture Infectious Dose 50%/mL).**
- 4 - Like any other diagnostic test, the results must be interpreted together with other patient clinical data.
- 5 - If the non-reactive result persists and the symptoms carry on, an additional test should be performed using another method. A non-reactive result does not exclude, at any time, the possibility of infection by the SARS-CoV-2 virus.
- 6 - **The accuracy of the test is closely linked to correct sample collection. False-negative results may be caused by inadequate collection or improper storage of the sample.**
- 7 - The buffer has a colorless to slightly yellowish color with a clear to slightly cloudy appearance, and may present suspended particles inherent to the raw materials used, without compromising the performance of the product.

CROSS REACTION

Cross-reactivity for BIOCLIN FAST COVID-19 Ag was evaluated using samples containing antigens for pathogens: Human Adenovirus 1, Parainfluenza, Rhinovirus, Influenza A H1N1, Influenza A H3N2, Mycoplasma pneumonia, and Streptococcus pneumoniae. None of the evaluated samples presented cross-reactivity to BIOCLIN FAST COVID-19 Ag. Despite the results found, it cannot be ruled out entirely the possibility of other reactivities crossed. The final diagnosis should consider the data patient's clinical data along with other data laboratory tests.

INTERFERENT

For the assessment of interferents, substances were added to 2 samples of known results. No divergent results were found up to concentrations of: acetylsalicylic acid 65 mg/dL, bilirubin 30 mg/dL, hemoglobin 500 mg/dL, triglycerides 1000 mg/dL and mucin 100 µg/mL and rheumatoid factor 1000 IU/mL.

PRODUCT PERFORMANCE

ACCURACY
Comparison of Methods, Methodology Specificity and Sensitivity
To establish the sensitivity and specificity of the kit BIOCLIN FAST COVID-19 Ag were analyzed 292 known clinical samples. Clinical sensitivity and clinical specificity found were:

Sensitivity: (135/142) = 95.1% (IC 95%: 91.6% – 98.6%)
Specificity: (146/150) = 97.3% (IC 95%: 94.8% – 99.9%)

To establish the sensitivity and specificity for the BIOCLIN FAST COVID-19 Ag kit for nasal swab samples, 100 known clinical samples were analyzed. The clinical sensitivity and clinical specificity found for nasal swab samples were:

Sensitivity: (48/50) = 96.0% (IC 95%: 90.67% – 100%)
Specificity: (49/50) = 98.0% (IC 95%: 94.2% – 100%)

PRECISION

Repeatability
The repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Samples	Repetition N°	Expected Results	Results Found
01	10	100% Reagent	100% Reagent
02	10	100% Non-Reagent	100% Non-Reagent

Reproducibility

The reproducibility was calculated from 5 successive determinations for 3 consecutive days, using 2 different samples, obtaining the following results:

Samples	N° of days	Repetition N°	Expected Results	Results Found
01	03	05	100% Reagent	100% Reagent
02	03	05	100% Non-Reagent	100% Non-Reagent

ANALYTICAL SENSITIVITY

Sensitivity 2 x 10^{3.1} TCID50/mL (Tissue Culture Infectious Dose 50%/mL).

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are fever, cough, dyspnoea, myalgia, or fatigue. It can cause diseases ranging from the common cold to more serious illnesses, such as the Middle East Respiratory Syndrome (MERS-CoV) and the Severe Acute Respiratory Syndrome (SARS-CoV). In December 2019, a series of pneumonia cases of unknown cause appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples it was shown to be Coronavirus infection, which was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents, cases of infection have already been reported.

BIBLIOGRAPHIC REFERENCES

- 1 - Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- 2 - CENTERS FOR DISEASE CONTROL AND PREVENTION et al. Interim guidelines for collecting, handling, and testing clinical specimens for COVID-19. Retrieved June, v. 19, p. 2020, 2020. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation, and storage of specimens for molecular methods; approved guideline. CLSI document MM13. A.Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- 3 - Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- 4 - Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb15;395(10223):497-506.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until the expiration date stated on the presentation packaging when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte/MG - Brasil
Tel.: +55 31 3439 5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Indústria Brasileira

CUSTOMER SERVICE

Customer Advisory Service
Phone: +55 31 3439 5454 | E-mail: customerservice@bioclin.com.br

ANVISA registration for Bioclin Fast COVID-19 Ag kit: 10269360335

Review: September/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER