BIOCLIN FAST DENGUE NS1

REF K130

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FINALIDADE

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.
- Dussart P., et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
- Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.
- Dussart P., et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
- Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento..

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INSTRUCTIONS FOR USE

FUNCTION

Test for the qualitative determination of the NS1 antigen of the Dengue Virus, in biological samples (serum, plasma or whole blood) through an immunochromatographic test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The Bioclin Fast Dengue NS1 kit is an immunochromatographic test for the detection of the NS1 antigen of the Dengue virus. The kit has a "T" test line with immobilized anti-NS1 antibodies and a "C" control line with chicken Anti-IgY antibodies. When a sample is added, it interacts with anti-NS1 antibodies conjugated to colloidal gold particles. Then, the sample migrates through the membrane by capillarity, and interacts with the anti-NS1 antibodies present in the "T" line. If the sample contains NS1 antigens, a colored band will appear in the "T" test region, indicating a reactive result. Failure to form a line in the "T" region indicates absence of Dengue NS1 antigen in the analyzed sample. The sample continues migrating through the membrane until it reaches the control "C" region, where the conjugate containing chicken IgY interacts with anti-chicken IgY antibodies, generating a colored line in the "C" region that indicates the correct functioning of the product.

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 30 °C. **Do not freeze.** Each sachet contains:

- A) 01 plastic cassette (device) with 01 test strip containing:
- Control Region (C): Chicken Anti-IgY.
 - Test Region (T): NS1 Anti-Dengue Antibody.
 - Conjugate: Anti-Dengue NS1 Antibody and chicken IgY Antibody conjugated with colloidal gold particles.

B) 01 sachet of silica.

Number 2 (R2) Diluent - Store between 2 and 30 °C. **Do not freeze.** Contains: Buffer, surfactant and preservative.

Collection Accessories:

- Disposable Plastic Pipette.
- Disposable Lancet*.

*The shape and color of the lancet may vary and do not affect the performance of the accessory and kit.

PRESENTATION

Reagent	Nº of Cassettes per Package	Diluent
Cassette + Diluent (Includes Pipettes)*	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
Cassette + Diluent (Includes Lancetes and Pipettes)**	50	6.0 mL
	100	2 x 6.0 mL
	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
30	4.0 mL	
40	5.0 mL	
50	6.0 mL	
100	2 x 6.0 mL	

* These presentations are accompanied by disposable plastic pipettes, proportionally to the number of cassettes.

** These presentations are accompanied by disposable plastic lancets and pipettes, proportionally to the number of cassettes.

OPERATIONAL EQUIPMENT AND SUPPLIES

Clock or stopwatch. They can be found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORTATION CONDITIONS

The storage temperature should be 2 to 30 °C. Transport at temperatures up to 45 °C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARES

- 1- Only for *in vitro* diagnostic use.
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples used are biological, the safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Open the aluminum envelope to remove the cassette, only when using it.
- 5- The kit must not be frozen or exposed to high temperatures (above 45°C), as this will cause deterioration.
- 6- The reagents, as well as the samples, must reach room temperature before starting the test.
- 7- Do not mix reagents from kits with the same batch or with different batch numbers.
- 8- Do not use the kit beyond its expiration date. The expiry date is indicated on the product label.
- 9- The cassette (test device), the plastic pipette and the lancet are disposable and cannot be reused.
- 10- The kit materials that came into contact with the patient's sample must have a special destination. Lancets must be disposed of in sharps waste. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biowaste.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

11- Secondary packaging (box), Instructions for Use, aluminum envelope and silica sachet can be disposed of in regular trash.

12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Safety Service) Customer Advice) by Quibasa.

14- Do not use the product if the packaging is damaged.

15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

It is suggested that the sample be collected between the 1st and 4th day after the onset of symptoms to perform the test.

Serum or Plasma: Samples can be stored under refrigeration, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer). To obtain the plasma, collect using EDTA, Heparin or Citrate as an anticoagulant.

Whole Blood: It can be collected by venipuncture using the appropriate anticoagulant (EDTA, Heparin or Citrate) or by digital puncture with the aid of a lancet. Whole blood collected by lancet should be used for immediate testing. Whole blood collected by venipuncture can be stored under refrigeration, between 2 and 8°C, for a maximum period of 5 days. **Do not freeze.***



* At the time of use, samples must be stabilized at room temperature.

PROCESS DESCRIPTION

The stability of the BIOCLIN FAST DENGUE NS1 cassette (after opening the aluminized envelope) is 2 hours. This stability may vary depending on the environmental conditions after opening the aluminized envelope. However, for safety reasons, we recommend that the BIOCLIN FAST DENGUE NS1 kit be used immediately after removing the cassette from the sachet.

TECHNIQUE - VENOPUNCIÓN

- 1- The sample must be at room temperature before starting the test.
- 2- Remove the cassette (test unit) from its protective packaging, place it on a flat and clean surface and label it properly.
- 3- **For Serum, Plasma or Total Blood:** Transfer 20 µL of sodium, plasma or total blood for the little sample.
- 4- Secure the bottle in an upright position and apply 2 drops (approximately 80 µL) of the diluent (Reagent N° 2) to the little sample.

5- Wait for the lines to form after resting for 25 to 30 minutes. **Do not interpret after 30 minutes.**

TECHNIQUE - DIGITAL PUNCTURE / LANCET

- 1- Remove the cassette from the protective packaging, place it on a clean and level surface and identify it properly.
- 2- Carry out asepsis on the finger tip to be perforated with 70% (p/p) alcohol.
- 3- Press the tip of the finger that will be pierced by the lancet so that there is an accumulation of blood in the region.
- 4- Remove the lancet protection cap.






5- Position and press the lancet firmly on the area to be lancet. Immediately, the blood will come out on the pierced area.

6- Collect the blood with the help of the disposable plastic pipette that accompanies the kit. Open the plastic pipette above the marked line and edge to the open cavity of the plastic pipette in the drop of blood. Release the pressure in the plastic pipette so that the blood is aspirated. Aspirate the blood at the marked line into a disposable plastic pipette (equivalent to 20 µL of blood).

7- Do not sample a little, dispense or bleed by pressing the pipette again.

8- Secure the bottle vertically and apply 2 drops (approximately 80 µL) of diluent (Reagent N° 2) to the little bit of the sample.

9- Wait for the lines to form after or rest for 25 to 30 minutes. **Do not interpret after 30 minutes.**


RESULT INTERPRETATION



Reagent Test for Dengue NS1: Formation of a red line in the test region (T) and another line in the control region (C) in the first 25 to 30 minutes. Do not interpret after 30 minutes.

Test Non-Reactive: Formation of a red line in the control region (C) and complete absence of a red line in the test region (T). Do not interpret after 30 minutes.

Invalid Result: Complete absence of line in the control region (C) with or without red line in the test region (T). Test the sample again.

**Notes:**

- 1- Non-reactive results should only be read within 30 minutes. Do not interpret after 30 minutes.
- 2- The line in the control region may appear before 25 minutes of test incubation, and this does not mean that non-reactive results can be interpreted ahead of time.
- 3- The results must not be interpreted after 30 minutes to avoid incorrect results.
- 4- The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.
- 5- Any intensity of color in the test line, even if very weak, characterizes the reagent result.

PROCESS LIMITATIONS

- 1- The Bioclin Fast Dengue NS1 kit is a qualitative test and should only be used for *in vitro* diagnostic use.
- 2- The Bioclin Fast Dengue NS1 kit is a screening test to characterize the presence of NS1 antigens of the Dengue virus.
- 3- Like all diagnostic tests, the results must be interpreted together with other clinical data of the patient.
- 4- Isolated test results cannot be used for definitive diagnosis.
- 5- If the non-reactive result remains and the clinical symptoms persist, an additional test must be performed using another method. A non-reactive result does not exclude, at any time, the possibility of infection with the Dengue virus.
- 6- Reactive results must be confirmed with another test and available clinical information.

INTERFERING

No interference was observed by Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS REACTIVITY

A cross-reactivity study was performed, evaluating 149 serum, plasma and whole blood samples negative for Dengue NS1 but positive for other infections. Among them 5 positive samples for HIV, 4 positive samples for HTLV, 5 positive samples for HBsAg, 15 positive samples for anti-HBs, 10 positive samples for Syphilis, 10 positive samples for Yellow Fever, 10 positive samples for Zika, 5 positive samples for Chikungunya, 10 positive

samples for Chagas Disease, 10 positive samples for CMV, 10 positive samples for Rubella, 10 positive samples for Toxoplasmosis, 7 positive samples for H1N1, 28 positive samples for Sars-CoV-2 and 10 positive samples for HCV. No cross-reactivity was observed with positive samples for HIV, HTLV, HBsAg, anti-HBs, Syphilis, Yellow Fever, Zika, Chikungunya, Chagas Disease, CMV, Rubella, Toxoplasmosis, H1N1, Sars-Cov-2 and HCV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE**ACCURACY****Method Comparison**

The Method Comparison study was performed using 365 biological samples, 93 of which were reactive for Dengue NS1 and 272 were non-reactive for Dengue NS1. The agreement found between the reference kit and the Bioclin Fast Dengue NS1 was 98.1%. With these results, it can be concluded that the kit presents good methodological specificity.

	Reference Kit	Bioclin Fast Dengue NS1
Reagent	93	89
Non-Reagent	272	269
Total	365	358

Agreement: 98.1% (358/365)

PRECISION**Repeatability**

Repeatability was calculated from 10 successive determinations, using 4 different samples, obtaining the following results:

SAMPLE	N° OF REPETITIONS	EXPECTED RESULT	RESULT FOUND
Plasma	10	100% Reagent	100% Reagent
Serum	10	100% Reagent	100% Reagent
Whole Blood	10	100% Reagent	100% Reagent
Serum	10	100% Non-Reagent	100% Non-Reagent

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 4 different samples, obtaining the following results:

SAMPLE	N° DE DAYS	N° OF REPETITIONS	EXPECTED RESULT	RESULT FOUND
Plasma	03	10	100% Reagent	100% Reagent
Serum	03	10	100% Reagent	100% Reagent
Whole Blood	03	10	100% Reagent	100% Reagent
Serum	03	10	100% Non-Reagent	100% Non-Reagent

SENSITIVITY AND CLINICAL SPECIFICITY

The BIOCLIN FAST DENGUE NS1 kit was used for the analysis of clinical samples previously characterized and confirmed using the reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOCLIN FAST DENGUE NS1 kit is 95.6% and the clinical specificity is 98%.

		EXPECTED RESULT		
		Reagent	Non-Reagent	Total
BIOCLIN FAST DENGUE NS1	Reagent	65	2	67
	Non-Reagent	3	98	101
	Total	68	100	168

Clinical Sensitivity: 95.6% (65/68) - 95% CI: 87.6% - 99.1%

Clinical Specificity: 98% (98/100) - 93% CI: 87.3% - 99.8%

DIAGNOSTIC SIGNIFICANCE

Dengue virus is a flavivirus transmitted by the *Aedes aegypti* mosquito. It is distributed throughout tropical and subtropical areas of the world, causing more than 100 million infections annually. The classic Dengue infection is characterized by high fever (usually between 38 and 40°C) with an abrupt onset, malaise, anorexia (little appetite), headache, muscle and eye pain. Cases of Hemorrhagic Dengue can cause gingivorrhagia and epistaxis, internal bleeding and disseminated intravascular coagulation, with damage to various organs, which can cause death. NS1 is a highly conserved glycoprotein that is present in high concentrations in the sera of dengue-infected patients. The NS1 antigen is found from the 1st to the 9th day after the onset of fever in samples from patients with primary and secondary dengue virus infection. However, its concentration decreases over the days of infection, considerably reducing the chances of detecting the antigen. Articles report a sensitivity of less than 50% for detection after the 4th day of infection by various diagnostic methods. Primary dengue infection causes an increase in IgM antibodies 3 to 5 days after the onset of fever. IgM antibodies normally remain in the circulation for 30 to 90 days. Patients from endemic regions may experience secondary infections, which result in elevated levels of IgG antibodies, either alone or concurrently with an IgM response. In case of secondary infection, the NS1 antigen may not be detected due to the possibility of lower viremia.

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR.
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
4. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
5. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
6. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
7. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
8. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
9. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
10. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Microbiol Immunol. 1994; 38(12): 951-5.
11. Dussart P., et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
12. Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
13. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | Email: sac@bioclin.com.br

ANVISA registration for BIOCLIN FAST DENGUE NS1 kit:
10269360203

Review: September/2025

UNIVERSAL SYMBOLOLOGY