BIOCLIN FAST HCV

REF K255

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FINALIDADE

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Sensibilidade: (93/93) = 99,99% (IC 95%: 98,0 – 100%)
Especificidade: (80/80) = 99,99% (IC 95%: 98,0 – 100%)

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1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Fourantigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Fourantigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
5. QUIBASA: Datos do Departamento de Pesquisa e Desenvolvimento.

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USE INSTRUCTIONS

FUNCTION

Rapid immunochromatographic test for the qualitative determination of total anti-HCV antibodies in fingerstick samples, whole blood (EDTA, Heparin or Citrate), serum or plasma (EDTA, Heparin or Citrate). For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
BIOCLIN FAST HCV is an immunochromatographic test intended for the detection of total antibodies against the hepatitis C virus (HCV) in whole blood, serum or plasma. Anti-HCV antibodies bind to the HCV antigens conjugated to colloidal gold particles (Core, NS3, NS4, and NS5 proteins) and form complexes that will migrate across the membrane toward immobilized HCV antigens in the test region (T), giving rise to the formation of a pink band that indicates the presence of anti-HCV antibodies. The absence of a pink band in the T region indicates the absence of anti-HCV antibodies. The colloidal gold particles continue to travel through the membrane to the control line, where a pink line forms for both reactive and non-reactive samples. The formation of the pink line in the control region (C) demonstrates that the kit is suitable for use.

REAGENTS

Number 1: Cassette - Store between 2 and 30 °C. **Do not freeze.** Contains:

- A)** 01 plastic cassette (device) with 01 test strip containing:
- Control region (C): goat anti-mouse IgG < 1000 ng.
 - Test region (T): HCV antigens (Core, NS3, NS4 and NS5) < 1000 ng.
 - Conjugate: conjugated colloidal gold; Mouse IgG < 1000 ng. HCV antigens (Core, NS3, NS4 and NS5) < 1000 ng.

B) 01 silica envelope

Number 2: Diluent - Store between 2 and 30 °C. **Do not freeze.** Contains: 100 mmol/L buffer pH 7.4 and preservative.

Collection Accessories:

- Disposable plastic pipette
- Disposable lancet (sterile)

PRESENTATIONS

Reagent	Nº of Cassetter per Packaging	Diluent
Cassette + Diluent	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL
	100	4 x 2.5 mL
Cassette + Diluent (pipette included)*	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL
	100	4 x 2.5 mL
Cassette + Diluent (lancet and pipette included)*	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL
	100	4 x 2.5 mL

*These presentations come with lancets and/or disposable plastic pipettes, proportionally to the number of cassettes.

EQUIPMENT AND OPERATIONAL SUPPLIES

Clock or stopwatch. They can be found in the specialized market of articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage and transport temperature should be between 2 and 30°C. Keep away from light and avoid humidity.

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- Strictly follow the proposed methodology to obtain accurate results.
- The cassettes (test devices) do not contain contaminating substances, but since the samples are biological, safety routines must be strictly followed and the use of disposable gloves is essential.
- Use clean and dry containers to store the samples.
- Treat all samples as potentially infectious.
- Reagent No. 2 contains Sodium Azide, irritant for skin and mucous membranes. Handle the reagent with care.
- The kit must not be frozen or exposed to high temperatures (above 45 °C), as this will cause it to deteriorate.

- The reagents, as well as the sample, must reach room temperature before starting the test.
- Open the aluminum sachet to remove the cassette (test device) only at the time of use.
- It is interesting to carry out the test immediately after the sample is collected.
- The kit should not be used after the expiration date.
- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request from the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damage to the packaging.
- It is essential that the instruments and equipment used are properly calibrated and subject to regular maintenance.

SAMPLES

Serum or Plasma: Samples can be stored refrigerated, between 2 and 8 °C, for a maximum period of 3 days. If samples cannot be tested within 3 days, they can be stored up to 30 days at -20 °C (freezer). To obtain plasma, collect the sample using EDTA, Heparin or Citrate as anticoagulant.



Whole Blood: It can be collected by venipuncture using the appropriate anticoagulant (EDTA, Heparin or Citrate) or by fingerstick with the help of a lancet. Whole blood collected by fingerstick should be used for immediate testing. Whole blood collected by anticoagulant venipuncture may be stored at 2-8 °C for up to 3 days. **Do not freeze.**

PROCESS DESCRIPTION



The stability of the Bioclin FAST HCV cassette (after opening the aluminized sachet) is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

TECHNIQUE - VENOPUNCTION

- The sample must be at room temperature before starting the test.
- Remove the cassette (test device) from the protective packaging, place it on a clean, level surface and properly identify it.
- For Serum or Plasma:** Transfer 25 µL of serum or plasma to the sample well. Hold the bottle of Diluent (Reagent Nº 2) vertically and apply 2 drops (80 µL) to the sample well.









- For Whole Blood:** Transfer 50 µL of whole blood to the sample well. Hold the bottle of Diluent (Reagent Nº 2) vertically and apply 2 drops (80 µL) to the sample well

- Wait for the lines to form. Interpret the results between 10 and 20 minutes. **Do not interpret after 20 minutes.**


TECHNIQUE - DIGITAL PUNCH / LANCET

- Remove the cassette (test device) from the protective packaging, place it on a clean, level surface and properly identify it.
- Wash and dry the patient's hands. Clean the tip of the finger with alcohol at 70% (p/p) for asepsis of the area used.
- Press the tip of the finger that will be punched by the lancet for blood accumulation in this region.
- Remove the protective cap from the lancet.
- Place and firmly press the lancet at the area to be punched. Blood will then come out of the punched area.
- Collect the blood with the disposable plastic pipette that comes with the kit. Place the pipette over the drop of blood and allow the blood to flow by capillarity, without pressing the bulb, up to the line marked on the disposable plastic pipette (equivalent to 50 µL of blood).
- In the sample well, dispense the blood by pressing the bulb of the pipette.
- Hold the bottle of Diluent (Reagent Nº 2) upright and apply 2 drops (80 µL) of Diluent to the sample well.
- Wait for the lines to form. Interpret the results between 10 and 20 minutes. **Do not interpret after 20 minutes.**


INTERPRETATION OF RESULTS



Reactive Test: Formation of two red lines, one in the Control region (C) and one in the Test region (T), after 10 minutes. Do not interpret after 20 minutes.

Non-Reactive Test: Formation of a red line in the control region (C) and total absence of a red line in the test region (T). Do not interpret after 20 minutes.

Invalid Result: Total absence of a red line in the control region (C), regardless of whether or not there is a red line in the test region (T). Test the sample again.

PROCEDURE LIMITATIONS

- 1- The line in the control region will appear before 10 minutes of test incubation, and this does not mean that the results can be interpreted before the expected time.
- 2- It is important to use the correct volume of sample, since lower or higher volumes can determine erroneous results.
- 3- The BIOCLIN FAST HCV kit detects antibodies against the HCV virus, the results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria to determine the diagnosis and/or treatment of the patient.

INTERFERENTS

No cross reaction was found with samples positive for: HAMA (human anti-murine antibody), rheumatoid factor, HBsAg, HBsAb, HBeAg, HBeAb, HBcAb, Syphilis, HIV, Mononucleosis, Cytomegalovirus, Rubella and Toxoplasmosis. Samples with concentrations up to 1000 mg/dL of Hemoglobin, Bilirubin 50 mg/dL and Triglycerides 1000 mg/dL do not interfere with the test.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where the procedures, standards, limits and tolerance to variations are clearly established. It is important to emphasize that every measurement system has a characteristic analytical variability, which must be monitored by the laboratories themselves. Therefore, the use of controls is recommended, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE
ACCURAY

Comparison of Methods

The BIOCLIN FAST HCV kit was compared to another commercially available HCV test kit. 50 tests were performed and the results were evaluated. The concordance found was > 99.9%. With these results, it can be concluded that the kit has good methodological specificity.

	Reference Kit	HCV
Reactive Sample	30	30
Non-Reactive Sample	20	20
Total of Samples	50	

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Sample	Nº of Repetitions	Expected Result	Found Result
01	10	Reactive	Reactive
02	10	Non-Reactive	Non-Reactive

Reproducibility

Reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 different samples, yielding the following results:

Sample	Nº of Days	Nº of Repetitions per Day	Expected Result	Found Result
01	03	10	Reactive	Reactive
02	03	10	Non-Reactive	Non-Reactive

SENSITIVITY AND CLINICAL SPECIFICITY

BIOCLIN FAST HCV kit has been tested with clinically known samples. For this test, 173 samples were analyzed, being 93 HCV reactive and 80 HCV non-reactive. The following results were found:

	Expected Result	Found Result
Reactive	93	93
Non-Reactive	80	80
Total of Samples	173	

Sensitivity: (93/93) = 99.99% (IC 95%: 98.0 – 100%)
Specificity: (80/80) = 99.99% (IC 95%: 98.0 – 100%)

DIAGNOSTIC SIGNIFICANCE

Hepatitis C virus has a small envelope and it is the main cause of parenteral transmission of non-A and non-B hepatitis. HCV infection causes a wide variety of chronic liver diseases, cirrhosis, and liver cancer. The main route of transmission of the virus is through blood and blood product transfusions, organ transplants, and sharing contaminated needles and syringes. HCV antibodies are found in more than 80% of patients with non-A and non-B hepatitis. Cloning of the viral genome made possible the development of serological tests using recombinant antigens.

BIBLIOGRAPHIC REFERENCES

1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. *Science*. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. *Science*. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Fourantigen Recombinant Immunoblot Assay. *Lancet*. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. *J. Clinical Immunoassay*. 1993;16:204.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOCLIN FAST HCV kit: 10269360350

Review: January/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER