BIOCLIN FAST HCV

REF K255

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FINALIDADE

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1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.

2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.

3. Van der Poel, C. L., H.T.M. Cuyppers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Fourantigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.

4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.

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1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuyppers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Fourantigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
5. QUIBASA: Datos do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Rapid immunochromatographic test for the qualitative determination of total anti-HCV antibodies in fingerstick samples, whole blood (EDTA, Heparin or Citrate), serum or plasma (EDTA, Heparin or Citrate). For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Immunochromatography

The Bioclin Fast HCV kit is an immunochromatographic test for detecting total antibodies against the Hepatitis C Virus (HCV) in whole blood, serum or plasma samples. The conjugate contains HCV antigens (Core proteins, NS3, NS4 and NS5) and rabbit IgG conjugated with colloidal gold. Total anti-HCV antibodies bind to HCV antigens conjugated with colloidal gold particles and form complexes that will migrate through the nitrocellulose membrane towards the HCV antigens immobilized in the "T" test region, leading to the formation of a red band indicating a reagent result. The absence of a red band in the "T" region indicates a non-reactive result. The control line "C" consists of immobilized anti-rabbit IgG antibody. In this region, the excess conjugate is captured by this antibody, regardless of whether the sample is reactive or non-reactive, forming the complex: anti-rabbit IgG antibodies – conjugate. The formation of a red band in the control region "C" indicates that the test is working perfectly.

REAGENTS

Number 1 (R1): Cassette - Store between 2 and 30 °C. **Do not freeze.** Each sachet contains:

- A)** 01 plastic cassette (device) with 01 test strip containing:
 - Control Region (C): Rabbit Anti-IgG Antibody: < 1 µg.
 - Test Region (T): HCV antigens (Core, NS3, NS4 and NS5) < 1 µg.
 - Conjugate: HCV antigens (Core, NS3, NS4 and NS5): < 1 µg and rabbit IgG < 1 µg conjugated with colloidal gold

B) 01 sachet of Silica

Number 2 (R2): Diluent - Store between 2 and 30 °C. **Do not freeze.** Contains: Buffer < 100 mmol/L pH 7.4, and preservatives.

Collection Accessories:

- Disposable plastic pipette (this item is included in the kit according to the commercial presentation)
- Lancet* (this item is included in the kit according to the commercial presentation)

*The shape and color of the lancet may vary and do not affect the performance of the accessory and kit.

PRESENTATIONS

Reagent	Nº of Cassetter per Packaging	Diluent
Cassette + Diluent	5	1.0 mL
	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL
Cassette + Diluent (pipette included)*	5	1.0 mL
	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL
Cassette + Diluent (lancet and pipette included)*	5	1.0 mL
	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL
Cassette + Diluent (lancet and pipette included)*	5	1.0 mL
	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL
Cassette + Diluent (lancet and pipette included)*	5	1.0 mL
	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL
Cassette + Diluent (lancet and pipette included)*	5	1.0 mL
	10	1.0 mL
	20	2.5 mL
	25	2.5 mL
	30	2 x 2.5 mL
	40	2 x 2.5 mL
	50	2 x 2.5 mL

*These presentations come with lancets and/or disposable plastic pipettes, proportionally to the number of cassettes.

EQUIPMENT AND OPERATIONAL SUPPLIES

Clock or stopwatch. They can be found in the specialized market of articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 30 °C. Transport at temperatures up to 45 °C should not exceed 5 days. Avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The cassettes (test devices) do not contain contaminating substances, but since the samples are biological, safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Use clean and dry containers to store the samples.
- 5- Treat all samples as potentially infectious.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

6- Reagent No. 2 contains Sodium Azide, irritant for skin and mucous membranes. Handle the reagent with care.

7- The kit must not be frozen or exposed to high temperatures (above 45 °C), as this will cause it to deteriorate.

8- The reagents, as well as the sample, must reach room temperature before starting the test.

9- Open the aluminum sachet to remove the cassette (test device) only at the time of use.

10- It is interesting to carry out the test immediately after the sample is collected.

11- The kit should not be used after the expiration date.

12- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request from the SAC (Customer Advisory Service) of Quibasa.

14- Do not use the product in case of damage to the packaging.

15- It is essential that the instruments and equipment used are properly calibrated and subject to regular maintenance.

SAMPLES

Serum or Plasma: Samples can be stored refrigerated, between 2 and 8 °C, for a maximum period of 3 days. If samples cannot be tested within 3 days, they can be stored up to 30 days at -20 °C (freezer). To obtain plasma, collect the sample using EDTA, Heparin or Citrate as anticoagulant.


Whole Blood: It can be collected by venipuncture using the appropriate anticoagulant (EDTA, Heparin or Citrate) or by fingerstick with the help of a lancet. Whole blood collected by fingerstick should be used for immediate testing. Whole blood collected by anticoagulant venipuncture may be stored at 2-8 °C for up to 3 days. **Do not freeze.**

PROCESS DESCRIPTION


The stability of the BIOCLIN FAST HCV cassette (after opening the aluminized sachet) is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

TECHNIQUE - VENOPUNCTION

- 1- The sample must be at room temperature before starting the test.
- 2- Remove the cassette (test device) from the protective packaging, place it on a clean, level surface and properly identify it.
- 3- **For Serum or Plasma:** Transfer 30 µL of serum or plasma to the sample well. Hold the bottle of Diluent (Reagent N° 2) vertically and apply 1 drop (about 40 µL) to the sample well.


- 4- **For Whole Blood:** Transfer 50 µL of whole blood to the sample well. Hold the bottle of Diluent (Reagent N° 2) vertically and apply 1 drop (about 40 µL) to the sample well.

5- Wait for the lines to form. Interpret the results between 15 and 20 minutes. **Do not interpret after 20 minutes.**


TECHNIQUE - DIGITAL PUNCH / LANCET

- 1- Remove the cassette (test device) from the protective packaging, place it on a clean, level surface and properly identify it.
- 2- Wash and dry the patient's hands. Clean the tip of the finger with alcohol at 70% (p/p) for asepsis of the area used.
- 3- Press the tip of the finger that will be punched by the lancet for blood accumulation in this region.
- 4- Remove the protective cap from the lancet.
- 5- Place and firmly press the lancet at the area to be punched. Blood will then come out of the punched area.
- 6- Collect the blood with the disposable plastic pipette that comes with the kit. Place the pipette over the drop of blood and allow the blood to flow by capillarity, without pressing the bulb, up to the line marked on the disposable plastic pipette (equivalent to 50 µL of blood).
- 7- In the sample well, dispense the blood by pressing the bulb of the disposable plastic pipette.
- 8- Hold the bottle of Diluent (Reagent N° 2) upright and apply 1 drop (about 40 µL) of Diluent to the sample well.
- 9- Wait for the lines to form. Interpret the results between 15 and 20 minutes. **Do not interpret after 20 minutes.**


INTERPRETATION OF THE RESULT


Reagent Test: Formation of two lines, one in the Control region (C) and one in the Test region (T), between 15 and 20 minutes. **Do not interpret after 20 minutes.** Any intensity of color on the test line, even if very weak, characterizes a reactive result.

Non-Reactive Test: Formation of a red line in the Control region (C) and complete absence of a red line in the Test region (T). **Do not interpret after 20 minutes.**

Invalid Result: Complete absence of a red line in the Control (C) region, regardless of whether or not there is a red line in the Test (T) region, indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

PROCESS LIMITATIONS

- 1- The line in the control region will appear before 10 minutes of test incubation, and this does not mean that the results can be interpreted before the expected time.
- 2- It is important to use the correct sample volume, as lower or higher volumes may cause erroneous results.
- 3- The BIOCLIN FAST HCV kit detects antibodies against the HCV virus, the results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the diagnosis and/or treatment of the patient.

INTERFERENCES

Samples with concentrations of up to 1000 mg/dL of Hemoglobin, Bilirubin 50 mg/dL and Triglycerides 1000 mg/dL do not interfere with the test.

CROSS REACTIVITY

No cross-reaction was observed with positive samples for: HCG, Syphilis, HIV, HBV, HTLV, HBsAg, Rheumatoid Factor, and PSA.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where the procedures, standards, limits and tolerance to variations are clearly established. It is important to emphasize that every measurement system has a characteristic analytical variability, which must be monitored by the laboratories themselves. Therefore, the use of controls is recommended, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE ACCURACY

Comparison of Methods

The BIOCLIN FAST HCV kit was compared to another commercially available HCV test kit. 80 tests were performed and the results were evaluated. The concordance found was > 99.9%. With these results, it can be concluded that the kit has good methodological specificity.

	Reference Kit	HCV
Reactive Sample	40	40
Non-Reactive Sample	40	40
Total of Samples	80	

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 different samples, obtaining the following results:

Sample	Nº of Repetitions	Expected Result	Found Result
01	10	Reactive	Reactive
02	10	Non-Reactive	Non-Reactive
03	10	Non-Reactive	Non-Reactive

Reproducibility

Reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 different samples, yielding the following results:

Sample	Nº of Days	Nº of Repetitions per Day	Expected Result	Found Result
01	03	10	Reactive	Reactive
02	03	10	Non-Reactive	Non-Reactive
03	03	10	Reactive	Reactive

SENSITIVITY AND CLINICAL SPECIFICITY

BIOCLIN FAST HCV kit has been tested with clinically known samples. For this test, 530 samples were analyzed, being 250 HCV reactive and 280 HCV non-reactive. The following results were found:

	Expected Result	Found Result
Reactive	250	249
Non-Reactive	280	279
Total of Samples	530	

Sensitivity: $(249/250) = 99,60\%$ (IC 95%: 98,82 – 100%)

Specificity: $(279/280) = 99,64\%$ (IC 95%: 98,95 – 100%)

DIAGNOSTIC SIGNIFICANCE

Hepatitis C virus has a small envelope and it is the main cause of parenteral transmission of non-A and non-B hepatitis. HCV infection causes a wide variety of chronic liver diseases, cirrhosis, and liver cancer. The main route of transmission of the virus is through blood and blood product transfusions, organ transplants, and sharing contaminated needles and syringes. HCV antibodies are found in more than 80% of patients with non-A and non-B hepatitis. Cloning of the viral genome made possible the development of serological tests using recombinant antigens.

BIBLIOGRAPHIC REFERENCES

1. Choo, Q.L., G. Kuo, A.J. Weiner, L.R. Overby, D.W. Bradley, and M. Houghton. Isolation of a cDNA Clone Derived from a Blood-borne Non-A, Non-B Viral Hepatitis Genome. Science. 1989;244:359.
2. Kuo, G., Q.L. Choo, H.J. Alter, and M. Houghton. An Assay for Circulating Antibodies to a Major Etiologic Virus of Human Non-A, Non-B Hepatitis. Science. 1989;244:362.
3. Van der Poel, C. L., H.T.M. Cuyppers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Fourantigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
4. Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOCLIN FAST HCV kit: 10269360350

Review: January/2025

UNIVERSAL SYMBOLOGY