Bioclin

BIOCLIN FAST ML FLOW HANSENÍASE

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FINALIDADE

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- 1) Associação entre reação hanseníca após alta e a carga bacilar avaliada utilizando sorologia anti-PGL-I e baciloscopia. Universidade Federal de Pernambuco (UFPE), 2008.
- 2) Association between anti-PGL-I IgM and clinical and demographic parameters in leprosy. *Leprosy Review* 77: 343-355, 2006.
- 3) Description of Leprosy Classification at Baseline among Patients Enrolled at the Uniform Multidrug Therapy Clinical Trial for Leprosy Patients in Brazil; *Am. J. Trop. Med. Hyg.*, 92(6), 2015, pp. 1280–1284.
- 4) Evaluation of a rapid serological test for leprosy classification using human serum albumin as the antigen carrier; *Journal of Immunological Methods* 412 (2014) 35–41.
- 5) IgM serum antibodies to phenolic glycolipid-I and clinical leprosy: two years' observation in a community with hyperendemic leprosy. *International Journal of Leprosy and Other Mycobacteriology Diseases* 58: 25-30, 1990.
- 6) Risk factors for type-1 reactions in borderline leprosy patients. *Lancet* 338: 654-657, 1991.
- 7) Simple and Fast Lateral Flow Test for Classification of Leprosy Patients and Identification of Contacts with High Risk of Developing Leprosy; *JOURNAL OF CLINICAL MICROBIOLOGY*, May 2003, p. 1991–1995.
- 8) The ML Flow test as a point of care test for leprosy control programmes: potential effects on classification of leprosy patients; *Lepr Rev* (2007) 78, 70–79.
- 9) The use of whole blood in a dipstick assay for detection of antibodies to *Mycobacterium leprae*: a field evaluation. *FEMS Immunology and Medical Microbiology* 21: 197-201, 1998.

- 10) Use of synthetic glycoconjugates containing the *Mycobacterium leprae* specific and immunodominant epitope of phenolic glycolipid I in the serology of leprosy. *Clinical and Experimental Immunology* 64: 476-483, 1986.
- 11) The use of serology as an additional tool to support diagnosis of difficult multibacillary leprosy cases: lessons from clinical care. *Revista da Sociedade Brasileira de Medicina Tropical* 41(Suplemento II):27-33, 2008)
- 12) QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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The Bioclin Fast ML Flow Hansen’s Disease kit was tested using samples from patients clinically positive for Hansen’s Disease and negative (without clinical signs) for the disease. The samples came from different regions of Brazil and had no distinction of age, sex or ethnic race.

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- Associação entre reação hansênica após alta e a carga bacilar avaliada utilizando sorologia anti PGL-I e baciloscopia. Universidade Federal de Pernambuco (UFPE), 2008.
- Association between anti-PGL-I IgM and clinical and demographic parameters in leprosy. Leprosy Review 77: 343-355, 2006.
- Description of Leprosy Classification at Baseline among Patients Enrolled at the Uniform Multidrug Therapy Clinical Trial for Leprosy Patients in Brazil; Am. J. Trop. Med. Hyg., 92(6), 2015, pp. 1280–1284.
- Evaluation of a rapid serological test for leprosy classification using human serum albumin as the antigen carrier; Journal of Immunological Methods 412 (2014) 35–41.
- IgM serum antibodies to phenolic glycolipid-I and clinical leprosy: two years' observation in a community with hyperendemic leprosy. International Journal of Leprosy and Other Mycobacteriology Diseases 58: 25-30, 1990.
- Risk factors for type-1 reactions in borderline leprosy patients. Lancet 338: 654-657, 1991.
- Simple and Fast Lateral Flow Test for Classification of Leprosy Patients and Identification of Contacts with High Risk of Developing Leprosy; JOURNAL OF CLINICAL MICROBIOLOGY, May 2003, p. 1991–1995.
- The ML Flow test as a point of care test for leprosy control programmes: potential effects on classification of leprosy patients; Lepr Rev (2007) 78, 70–79.
- The use of whole blood in a dipstick assay for detection of antibodies to Mycobacterium leprae: a field evaluation. FEMS Immunology and Medical Microbiology 21: 197-201, 1998.

- Use of synthetic glycoconjugates containing the Mycobacterium leprae specific and immunodominant epitope of phenolic glycolipid I in the serology of leprosy. Clinical and Experimental Immunology 64: 476-483, 1986.
- The use of serology as an additional tool to support diagnosis of difficult multibacillary leprosy cases: lessons from clinical care. Revista da Sociedade Brasileira de Medicina Tropical 41(Suplemento II):27-33, (2008)
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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BIOCLIN FAST ML FLOW HANSEN’S DISEASE

REF K254

INSTRUCTIONS FOR USE

FUNCTION

Rapid immunochromatographic test for the qualitative determination of anti-*Mycobacterium leprae* IgM antibodies in biological samples (serum, plasma or whole blood). Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Immunochromatography
The BIOCLIN FAST ML FLOW HANSEN’S DISEASE kit is an immunochromatographic test that detects anti-*Mycobacterium leprae* IgM antibodies in serum, plasma and whole blood (venipuncture or finger puncture) samples, employing a specific antigen from *M. Leprae*: fenolic glicolipid I (PGL-I). The test line contains the PGL-I immobilized. The conjugate is formed by colloidal gold nanoparticles linked with human anti-IgM monoclonal antibodies. When a positive sample is added to the sample well, it migrates through the membranes by capillary action. When it achieves the conjugate, the IgM antibodies present on the sample are captured by the human anti-IgM monoclonal antibodies, conjugated to colloidal gold, forming the complex: human IgM – conjugate. The sample continues to migrate by the membrane until reaches the test line (T), where the IgM antibodies positive for PGL-I are captured by the PGL-I immobilized, forming the complex: PGL-I – human IgM – conjugate. This reaction leads to the formation of a reddish line, indicating the presence of positive IgM antibodies to PGL-I. Negative samples do not have IgM antibodies reactive to PGL-I. In this way there is no complex formation at the test Line. The non-formation of a reddish line in the T region indicates the absence of IgM antibodies to PGL-I in the sample. The sample continues the migration through the membrane until reach the control line (C), formed by protein A. In this line it is observed the formation of a reddish line, due to the capture of the excess of conjugate by the protein A. The signal on the control region should be present in positive and negative samples. The formation of the reddish line in the Control region demonstrates that the result is valid.

REAGENTS

Number 1 (R1) Cassette – Store at 2 to 30°C. **Do not freeze.** Contains:
A) 01 plastic device with 01 test strip containing:
- Control region (C): Protein A (> 0.06 µg).
- Test region (T): PGL-I antigen (> 0.3 µg).
- Conjugate: Anti-human IgM antibody conjugated to colloidal gold nanoparticles.

B) 01 silica sachet

Number 2 (R2) Sample Diluent - Diluent - Store between 2 and 30°C. Do not freeze. Contains: Buffer Solution pH 7.2, surfactant and preservative.

Accessories for Sample Collection:

- Plastic capillary tube
- Disposable lancet* +
- Alcohol sachet*

*Available in accord to the commercial presentation.
+ The format and color of the lancet may vary, but, without interfering in the performance of the product.

PRESENTATIONS

Reagent	Number of Cassettes of Packaging	Diluent
Cassette + Diluent (Includes Pipette)*	5	1.0 mL
	10	1.5 mL
	15	2.5 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
	50	6.0 mL
Cassette + Diluent (Includes Lancet and Pipette)**	100	2 x 6.0 mL
	5	1.0 mL
	10	1.5 mL
	15	2.5 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
Cassette + Diluent (Includes Lancet and Pipette and Alcoholic Sachet)***	50	6.0 mL
	100	2 x 6.0 mL
	5	5 x 0.25 mL
	10	10 x 0.25 mL
	15	15 x 0.25 mL
	20	20 x 0.25 mL
	25	25 x 0.25 mL
	30	30 x 0.25 mL
	40	40 x 0.25 mL
	50	50 x 0.25 mL
	100	100 x 0.25 mL

*These presentations are accompanied by disposable plastic pipettes, proportionally to the number of cassettes.
**These presentations are accompanied by disposable plastic lancets and pipettes, proportionally to the number of cassettes.
***These presentations are accompanied by lancets, disposable plastic pipettes and alcoholic sachets, proportionally to the number of cassettes.

INSTRUMENTS AND OPERACIONAL INPUTS

Pipettes, tips, watch or stopwatch, 70% alcohol (w/w), cotton and garbage for disposal. These items are found in the specialized market for articles for Clinical Analysis Laboratories. The kit allows visible reading without the aid of equipment.

CONDITIONS OF STORAGE AND TRANSPORTATION

The storage and transport temperature should be 2 to 30°C. Avoid moisture. **Do not freeze.**

WARNINGS

1 - For *in vitro* diagnostic use only.
2 - Strictly follow the methodology proposed to obtain exact results.
3 - The cassettes (test devices) do not contain contaminating substances, but as the samples are biological, safety routines must be strictly followed and the use of disposable gloves is essential.
4 - Treat all samples as potentially infectious.
5 - Use clean and dry containers to store the samples.
6 - Open the aluminum envelope, to remove the cassette (test device), only at the moment of use.
7 - The kit must not be frozen or exposed to high temperatures, as this will cause deterioration.
8 - The reagents, as well as the sample, must reach room temperature before starting the test.
9 - Do not mix the reagents (components) of kits with the same lot or different lot numbers.
10 - Do not use the kit after the expiration date. The expiry date is indicated on the product label.
11 - The cassette (test device), pipette and lancet are disposable and cannot be reused.
12 - The kit materials that had contact with the patient’s biological material must have a special destination. Disposal of lancets and pipettes should be done in sharps. Cassettes (test devices) and diluent vials should be disposed of in biowaste.
13 - Boxes, Instructions for Use, aluminum envelope and silica sachet can be disposed of in regular trash.
14 - We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
15 - To obtain information related to biosafety or in case of accidents with the product, consult the FISPQ (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Safety Service Client Advisory) by Quibasa.
16 - Do not use the product in case of damage to the packaging.
17 - It is essential that the instruments and equipment used are properly calibrated and subject to periodic maintenance.

SAMPLES

Serum or Plasma: Samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 3 days. If samples cannot be analyzed within 3 days, they can be stored for up to 30 days at -20 °C (freezer). To obtain the plasma, perform collection using EDTA, Heparin or Citrate as an anticoagulant.

Whole Blood: Collect whole blood by venipuncture using the appropriate anticoagulant (EDTA, Heparin or Citrate) or by digital puncture with the aid of a lancet. Whole blood collected by lancet should be used for immediate testing. Whole blood collected by venipuncture with anticoagulant can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 3 days. **Do not freeze.**


PROCESS DESCRIPTION

Bioclin FAST ML FLOW HANSEN’S DISEASE reactive strip stability (after opening the aluminized sachet) is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

TECHNIQUE - VENOPUNCTURE

1 - The sample must be at room temperature before starting the test.
2 - Remove the cassette from the protective packaging, place it on a clean and leveled surface. Clearly identify it.
3 - For Whole Blood, serum or plasma: Transfer 10 µL of the sample to the sample well


4 - Hold the bottle vertically and apply 2 drops (approximately 70µL) of Diluent (Reagent N° 2) in the sample well.

5 - Wait for the lines to form. Interpret the results between 15 and 20 minutes. Do not interpret after 20 minutes.

TECHNIQUE - DIGITAL PUNCH / LANCET

1 - Remove the cassette from the protective packaging, place it on a clean, level surface and clearly identify it.
2 - Wash and dry the patient’s hands. Apply 70% alcohol (w/w) on the tip of the finger for asepsis of the area used.
3 - Press the fingertip that will be pierced by the lancet to accumulation of blood in this region.
4 - Remove the protective cap from the lancet.
5 - Position and press the lancet firmly under the area to be lancet. Then blood will come out of the punctured area.
6 - Press the capillary above the trace and touch the sample. Relieve the pressure for which the sample is aspirated. Aspirate up to the dash to get the volume of 10 µL
7 - In the sample well, dispense the blood by pressing again the pipette.
8 - Hold the Diluent bottle (Reagent N° 2) vertically and apply 2 drops (approximately 70 µL) of Diluent into the sample well.
9 - Wait for the lines to form. Interpret the results between 15 and 20 minutes. Do not interpret after 20 minutes.

INTERPRETATION OF THE RESULT


Reagent Test: Formation of a red line in the test region (T) and another line in the control region (C) in the first 15 to 20 minutes. Do not interpret after 20 minutes.

Non-Reagent Test: Formation of a reddish line in the control region (C) and complete absence of a red line in the test region (T). Do not interpret after 20 minutes.

Invalid Result: Complete absence of a line in the region control (C) with or without a red line at test region (T). Test the sample again.

Notes:

- Results should only be read between 15 and 20 minutes. Not before or after that time.
- The line in the control region can appear before 15 minutes. However, the reading of the results must be made between 15 and 20 minutes.
- The results should not be interpreted after 20 minutes to avoid incorrect results.
- The results provided by this kit must be interpreted by the responsible healthcare professional, who must take into account the clinical signs and the test result for the determination of the diagnosis and/or treatment of the patient.

LIMITATIONS OF THE PROCESS

- It is important to use the correct volume of sample, because lower volumes or higher may lead to erroneous results.
- The results provided by this kit must be interpreted by the responsible healthcare professional, who must take into account the clinical signs and the test result for the determination of the diagnosis and/or treatment of the patient.
- The presence of a red line in the test region (T) at any intensity, weak or strong, when the control line is also present, identifies a reagent result.
- The intensity of the control line (C) may vary according to the sample used in the test and should not be used as a criterion for analysis of the test line.
- If the test shows a non-reactive result, even if the clinical symptoms persist, the investigation of the diagnosis must be maintained following the criteria of the Ministry of Health.
- Do not use *post-mortem* samples.

INTERFERENT

Hemolysis with hemoglobin values up to 1000 mg/dL do not interfere with the dosage of IgM anti-M antibodies. leprae. No significant interference was found in samples with up to 20 mg/dL of acetylsalicylic acid, 2 g/dL of ascorbic acid, 200 mg/dL of creatine, 1 g/dL of bilirubin, 2 g/dL of albumin, 60 mg/dL of oxalic acid. Interferences above these values can cause false reagent or false non-reactive results. There was no evidence of interference with Rheumatoid Factor up to 1000 IU/mL, C-Reactive Protein up to 50 mg/dL and Anti Streptolysin O up to 1000 IU/mL.

HIGH DOSE PROZONE EFFECT

Although no prozone effect was observed in laboratory tests, field tests show that blood, serum or plasma samples from patients with clinical signs of Virchowian leprosy may show a false negative result in the test or test band with weaker intensity. In this case, carrying out serial dilutions of the sample in saline solution can be an alternative to show a reactive result in the immunochromatographic test.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality assurance program, in which all procedures and rules are previously established and clearly shared for all the organization and where every limits and tolerance values to variations are clearly established. It is also recommended that every lab have its own set of controls, in order to access the precision and accuracy of the dosages methods, considering that every dosage method has its own analytical variety.

PRODUCT PERFORMANCE ACCURACY

Methods comparison

The BIOCLIN FAST ML FLOW HANSEN'S DISEASE kit was compared with ELISA for identification of IgM antibodies anti-*M. leprae*. 116 different samples were analyzed in both methods. The performed and results were evaluated. The total agreement found was > 99.1%, showing that the lateral flow test has good specificity methodological.

	Kit Reference	BIOCLIN FAST ML FLOW HANSEN'S DISEASE
Reactive Samples	51	52
Non-Reactive Samples	65	64
TOTAL	116	

Total agreement found: > 99.1% (115/116)

PRECISION

Repeatability

The repeatability was calculated from 10 successive determinations, using 12 different samples. On this test it was observed 100% of agreement between reference and found results, obtaining the following results:

Amostras	Nº de Repetições	Resultado Esperado	Resultado Encontrado
Serum	10	Reagent	Reagent
Plasma EDTA	10	Reagent	Reagent
Whole Blood EDTA	10	Reagent	Reagent
Whole Blood Heparin	10	Reagent	Reagent
Plasma Heparin	10	Reagent	Reagent
Digital Puncture	10	Reagent	Reagent
Serum	10	Non- Reagent	Non- Reagent
Plasma EDTA	10	Non- Reagent	Non- Reagent
Whole Blood EDTA	10	Non- Reagent	Non- Reagent
Whole Blood Heparin	10	Non- Reagent	Non- Reagent
Plasma Heparin	10	Non- Reagent	Non- Reagent
Digital Puncture	10	Non- Reagent	Non- Reagent

Reproducibility

The reproducibility was calculated from 5 successive determinations for 3 consecutive days, using 12 different samples. On this test it was observed 100% of agreement between reference and found results, obtaining the following results:

Amostra	Nº de Dias	Nº de Repetições	Resultado Esperado	Resultado Encontrado
Serum	3	5	Reagent	Reagent
Plasma EDTA	3	5	Reagent	Reagent
Whole Blood EDTA	3	5	Reagent	Reagent
Whole Blood Heparin	3	5	Reagent	Reagent
Plasma Heparin	3	5	Reagent	Reagent
Digital Puncture	3	5	Reagent	Reagent
Serum	3	5	Non- Reagent	Non- Reagent
Plasma EDTA	3	5	Non- Reagent	Non- Reagent
Whole Blood EDTA	3	5	Non- Reagent	Non- Reagent
Whole Blood Heparin	3	5	Non- Reagent	Non- Reagent

Plasma Heparin	3	5	Non- Reagent	Non- Reagent
Digital Puncture	3	5	Non- Reagent	Non- Reagent

CLINICAL SENSITIVITY AND SPECIFICITY

The Bioclin Fast ML Flow Leprosy kit has been tested using samples from patients clinically positive for leprosy and negative (without clinical signs) for the disease. The samples came from different regions of Brazil and had no distinction of age, sex or ethnic race. The results are presented below:

	Kit Reference	BIOCLIN FAST ML FLOW HANSEN'S DISEASE
Reactive Samples	48	44
Non-Reactive Samples	60	54
TOTAL	108	

The results found suggest the following values:

Clinical Sensitivity: 91.7% (IC 84.15% - 99.18%)

Clinical Specificity: 90.0% (IC 82.76% - 97.24%)

DIAGNOSTIC SIGNIFICANCE

Leprosy is an infection caused by *Mycobacterium leprae*. Bioclin FAST ML FLOW HANSEN'S DISEASE is a simple test for the detection of specific IgM antibodies against the phenolic glycolipid-I (PGL-I) from *Mycobacterium leprae*, in human serum, plasma or whole blood (venipuncture or finger puncture).

The test is useful to aid in diagnosis when results are analyzed together with the patient's clinical information^{2, 5, 9, 10}. The presence of antibodies suggests the presence of a multibacillary infection (MB) and therefore the test can be used to discriminate multibacillary patients of the paucibacillary ones. Monitoring the effectiveness of therapy should be accompanied by a decrease in the titers of circulating antibodies obtained by ELISA. In leprosy reactions, PGL-I also proved to be useful as predictive test^{1, 6}. The BAAR detection method in intradermic smear indicates active disease for Leprosy. Differently of this, the BIOCLIN FAST ML FLOW HANSEN'S DISEASE serological test detects antibodies against *Mycobacterium leprae*, which indicates infection that may progress or not into disease state. The test can be used to help identify differential diagnoses, reinforce clinical suspicion of cases asymptomatic MB. The BIOCLIN FAST ML FLOW HANSEN'S DISEASE test should be used as auxiliary tool to assist in the diagnosis of complex MB cases (Virchowian and Borderline Virchowiano), mainly in cases where there are difficulties with the absence of cardinal signs and symptoms, or when they are not recognized in the health services¹¹.

The BIOCLIN FAST ML FLOW HANSEN'S DISEASE test can be used to investigate erythema knotty leprosy. The presence of erythematous and painful nodules may not be Leprosy. In these cases, a negative result on the BIOCLIN FAST ML FLOW HANSEN'S DISEASE should be used to direct investigation to other infectious diseases that could explain the clinical symptoms, such as Hepatitis C (HCV)¹¹. In endemic areas, the presence of a low specificity antibodies may correspond to a previous exposure, and therefore should be viewed as a negative result.

BIBLIOGRAPHIC REFERENCES

- Associação entre reação hansênica após alta e a carga bacilar avaliada utilizando sorologia anti PGL-I e baciloscopia. Universidade Federal de Pernambuco (UFPE), 2008.
- Association between anti-PGL-I IgM and clinical and demographic parameters in leprosy. Leprosy Review 77: 343-355, 2006.
- Description of Leprosy Classification at Baseline among Patients Enrolled at the Uniform Multidrug Therapy Clinical Trial for Leprosy Patients in Brazil; Am. J. Trop. Med. Hyg., 92(6), 2015, pp. 1280–1284.
- Evaluation of a rapid serological test for leprosy classification using human serum albumin as the antigen carrier; Journal of Immunological Methods 412 (2014) 35–41.
- IgM serum antibodies to phenolic glycolipid-I and clinical leprosy: two years' observation in a community with hyperendemic leprosy. International Journal of Leprosy and Other Microbacteriology Diseases 58: 25-30, 1990.
- Risk factors for type-1 reactions in borderline leprosy patients. Lancet 338: 654-657, 1991.
- Simple and Fast Lateral Flow Test for Classification of Leprosy Patients and Identification of Contacts with High Risk of Developing Leprosy; JOURNAL OF CLINICAL MICROBIOLOGY, May 2003, p. 1991–1995.
- The ML Flow test as a point of care test for leprosy control programmes: potential effects on classification of leprosy patients; Lepr Rev (2007) 78, 70–79.
- The use of whole blood in a dipstick assay for detection of antibodies to *Mycobacterium leprae*: a field evaluation. FEMS Immunology and Medical Microbiology 21: 197-201, 1998.

10) Use of synthetic glycoconjugates containing the *Mycobacterium leprae* specific and immunodominant epitope of phenolic glycolipid I in the serology of leprosy. Clinical and Experimental Immunology 64: 476-483, 1986.

11) The use of serology as an additional tool to support diagnosis of difficult multibacillary leprosy cases: lessons from clinical care. Revista da Sociedade Brasileira de Medicina Tropical 41(Suplemento II):27-33, 2008)

12) QUIBASA: Development data from the Research and Development department.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until the expiration date stated on the presentation packaging when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte/MG - Brasil

Tel.: +55 31 3439 5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Brazilian Industry

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Bioclin Fast ML Flow Hansen's Disease kit:

10269360345

Review: May/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER