## BIOCLIN FAST SANGUE OCULTO

REF K351

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## FINALIDADE

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1- NaKama H. et al., A comparative study of immunochemical fecal testes for detection of colorectal adenomatous polyps. Hepato-Gastroenterology. 2000, 47:386-389.

2- Uchida K. et al., Immunochemical detection of human blood in feces. Clinica Chimica Acta. 1990, 189:267-274.

3- Tommaso CL, Salzeider K, Arif M, et al. Serial myoglobin Gorzynski TJ, Krcó CJ, David CS., An idiotope expresse don a monoclonal antibody specific for human haemoglobin chain na naturally occurring immunoglobulin(s). J Immunogenet. 1985 Dec;12(6):281-91.

4- Gorzynski TJ, Krcó CJ, David CS., The Hb-2d cross-reactive idiotyp. A common idiotype expressed by monoclonal and polyclonal antibodies to human haemoglobin. J Immunogenet. 1985 Dec;12(6):267-79

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### INSTRUCTIONS FOR USE

#### FUNCTION

Immunochromatographic method for rapid qualitative determination of Occult Blood on feces. For *in vitro* diagnostic use only.

#### ACTION PRINCIPLE

##### Methodology: Immunochromatography

The BIOCLIN FAST OCCULT BLOOD Kit is an immunochromatographic test for the qualitative detection of human hemoglobin (Hb) in stool samples. The method uses membranes pre-coated with antibodies to selectively identify human hemoglobin. If present in the sample, human hemoglobin binds to the anti-hemoglobin antibody conjugated with colloidal gold nanoparticles and forms complexes that will migrate through the nitrocellulose membrane by capillarity toward the second antibody immobilized in the test line region "T," leading to the formation of a red band indicating a reactive result. The absence of a red band in the "T" region indicates a non-reactive result. The sample continues to migrate along the membrane until it reaches the control region "C," where the conjugate containing rabbit IgG antibody interacts with the anti-rabbit IgG antibody, generating a colored line in the "C" region that indicates the proper functioning of the product.

#### REAGENTS

**Number 1 (R1): Cassette** - Store between 2 and 30°C. **Do not freeze.** Each sachet contains:

A) 01 plastic device with 01 test strip containing:

**Control Region (C):** Anti-rabbit IgG antibody: < 1 µg.

**Test Region (T):** Anti-human hemoglobin antibody: < 1 µg.

**Conjugate:** Anti-human hemoglobin antibody < 1 µg and rabbit IgG antibody < 1 µg conjugated with colloidal gold nanoparticles.

B) 01 sachet of silica.

**Number 2 (R2) Diluent** - Contained in the tube. Store between 2 and 30°C. **Do not freeze.** Contains: Buffer < 100 mmol, surfactant and preservatives.

#### PRESENTATION

N° of Cassettes per package	Diluent
1	1 x 2 mL
5	5 x 2 mL
10	10 x 2 mL
20	20 x 2 mL
25	25 x 2 mL
30	30 x 2 mL
40	40 x 2 mL
50	50 x 2 mL
100	100 x 2 mL

#### OPERATIONAL EQUIPMENT AND SUPPLIES

Clock or stopwatch. They can be found in the specialized market for articles for Clinical Analysis Laboratories.

#### STORAGE AND TRANSPORTATION CONDITIONS

The storage and transport temperature should be between 2 and 30°C. Avoid moisture. **Do not freeze**

#### SPECIAL CARES

**1- Only for *in vitro* diagnostic use.**

**2- Strictly follow the proposed methodology to obtain accurate results.**

**3- The cassettes (test devices) do not contain contaminating substances, but since the samples are biological, safety routines must be strictly followed and the use of disposable gloves is essential.**

**4- Use clean and dry containers to store the samples.**

**5- Treat all samples as if they were potentially infectious.**

**6- The kit should not be frozen or exposed to high temperatures (above 45 °C), as this will cause deterioration.**

**7- The reagents, as well as the sample, must reach room temperature before starting the test.**

**8- Open the aluminum envelope to remove the cassette (test device) only at the time of use.**

**9- It is advisable to perform the test immediately after collecting the sample.**

**10- The kit should not be used after the expiration date.**

**11- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.**

**12- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or by requesting it through Quibasa's SAC (Customer Advisory Service).**

**13- Do not use the product if the packaging is damaged.**

**14- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.**

#### SAMPLES

Only for testing with human feces samples.

Stool samples should be collected in clean, dry containers. Avoid contamination by water and other substances from the bathroom. Samples should be kept in a cool place and under controlled temperature conditions.

Samples should be tested within 6 hours of collection. To store them for longer, they must be frozen at -20°C for 1 week or 3 days at 2 - 8°C.

Dilute the sample only in the diluent of the tube (R2) of the kit.

Samples should always be at room temperature for testing.

Samples should not be collected during the menstrual period. Samples from patients suffering from hemorrhoids, hemorrhages or who have blood in their urine may generate inconsistent results.

The use of medications that irritate the gastric mucosa (anti-inflammatories, corticosteroids or iron and vitamin C-based medications) should be discontinued at least 2 days before collecting the sample.

Excessive alcohol consumption can also cause gastrointestinal irritation. Therefore, these substances should not be used for at least 72 hours before performing the test. No dietary restrictions are necessary to perform the test.

#### PROCESS DESCRIPTION

##### Stability after opening

The BIOCLIN FAST OCCULT BLOOD cassette is stable (after opening the foil sachet) for 2 hours. This stability may vary depending on environmental conditions after opening the foil sachet. However, for safety reasons, we recommend that the BIOCLIN FAST OCCULT BLOOD kit be used immediately after removing the cassette from the sachet.

#### TECHNIQUE

The sample must be at a temperature between 15 and 30°C (room temperature) before starting the test.

**1- Open the analysis tube containing the diluent solution and insert its measuring loop into 5-6 different points of the sample (A).**





**2- Close the analysis tube with the measuring loop impregnated with feces, and shake vigorously for 3 to 5 minutes to emulsify the feces and ensure dispersion and dissolution (a vortex can be used) (B).**

**3- Remove the cassette from the aluminum sachet and identify it appropriately.**

**4- Open the upper end of the tube containing the diluent and sample so that it functions as a dropper (C).**

**5- Place 2 drops in the sample hole of the cassette, allowing the first drop to be absorbed so that the next drop can be placed (C).**


**6- Read the result between 5 and 10 minutes. Do not read the result after 10 minutes (D).**


#### RESULT INTERPRETATION

**Reagent Test:** Formation of two lines, one in the Control region (C) and one in the Test region (T), between 5 and 10 minutes.

**Do not interpret after 10 minutes.** Any intensity of color on the test line, even if very weak, characterizes a reactive result.


**Non-Reactive Test:** Formation of only one red line, between 5 and 10 minutes, in the Control region. **Do not interpret after 10 minutes.**

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

**Invalid Result:** Complete absence of a red line in the Control (C) region, regardless of whether or not there is a red line in the Test (T) region, indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

#### PROCESS LIMITATIONS

**1-** The line in the control region will appear before the 5 minute incubation of the test, and this does not mean that the results can be interpreted before the expected time.

**2-** The test is qualitative, although there is a proportionality in the intensity of the color obtained in the test line.

**3-** It is important to use the correct volume of sample, as lower or higher volumes may determine erroneous results.

**4-** The reaction reading time must be followed according to the established technique, in order to avoid false interpretations of the results.

**5-** The results provided by this kit must be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient, and must be interpreted in conjunction with the available clinical evaluation information of the patient and other diagnostic procedures.

#### INTERFERING

The BIOCLIN FAST OCCULT BLOOD kit was evaluated for possible interference. No interference was observed for Peroxidase (up to 20,000 µg/mL), Beetroot, Red Radish, Turnip, Cauliflower, Broccoli, Parsnip, Cantaloupe, Ascorbic Acid (up to 20 mg/dL) and Bilirubin (up to 20 mg/dL). Inconsistent results may be present in the presence of hemorrhoids, menstrual bleeding or blood in the urine, due to possible contamination of the feces.

#### CROSS REACTIVITY

The BIOCLIN FAST OCCULT BLOOD kit was evaluated for possible cross-reactions and no interference was observed with hemoglobin from swine, chicken, cattle, duck, fish, and goat.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. For this purpose, it is recommended to use controls, which allow the evaluation of the precision and accuracy of the dosages.

**PRODUCT PERFORMANCE****ACCURACY****Method Comparison**

The BIOCLIN FAST OCCULT BLOOD kit was compared with a reference kit for detecting occult blood in feces. 192 analyses were performed and the results were evaluated. The agreement found was > 99.47%. With these results, it can be concluded that the kit presents good methodological specificity.

		Expected Result		
		Reactive	Non Reactive	Total
<b>BIOCLIN FAST SANGUE OCULTO</b>	<b>Reactive</b>	76	1	77
	<b>Non Reactive</b>	0	115	115
	<b>Total</b>	76	116	192

**PRECISION****Repeatability**

Repeatability was calculated from 20 successive determinations, using 3 different samples, obtaining the following results:

Sample	N° of Repetitions	Expected Results	Found Results
01	20	Reactive	Reactive
02	20	Reactive	Reactive
03	20	Reactive	Reactive
04	20	Non-Reactive	Non-Reactive
05	20	Non-Reactive	Non-Reactive

**Reproducibility**

Reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 different samples, yielding the following results:

Sample	N° of days	N° of repetitions	Expected Results	Found Results
01	03	20	Reactive	Reactive
02	03	20	Reactive	Reactive
03	03	20	Reactive	Reactive
04	03	20	Non-Reactive	Non-Reactive
05	03	20	Non-Reactive	Non-Reactive

**SENSITIVITY AND CLINICAL SPECIFICITY**

The BIOCLIN FAST OCCULT BLOOD kit was tested using clinically known samples. For this test, 192 samples were analyzed, 76 positive and 116 negative. The following results were found:

		Expected Result		
		Reactive	Non Reactive	Total
<b>BIOCLIN FAST SANGUE OCULTO</b>	<b>Reactive</b>	76	1	77
	<b>Non Reactive</b>	0	115	115
	<b>Total</b>	76	116	192

Sensitivity: (76/76) > 99.99% (95% CI: 95.26% - 100%)

Specificity: (115/116) = 99.14% (95% CI: 95.29% - 99.98%)

**ANALYTICAL SENSITIVITY**

The BIOCLIN FAST OCCULT BLOOD kit was able to detect human hemoglobin concentrations equal to or greater than 50 ng/mL.

**DIAGNOSTIC SIGNIFICANCE**

The presence of Occult Blood in feces is directly associated with gastrointestinal disorders such as Hemorrhoids, Anal Fissure, Diverticulitis, Polyps, Crohn's Disease, Colorectal Cancer,

Gastric or Duodenal Ulcers, Intestinal Infection, among others.

Small amounts of blood in the stool or bleeding that can only be detected after cleaning the anus with toilet paper are the most common forms of rectal bleeding. In 90% of cases, the etiology is benign and mainly corresponds to hemorrhoids and anal fissures. When the amount of blood in the stool is moderate to large, or when there is melena (stool with digested blood), the origin of the bleeding is usually more internal, usually colon or stomach. Colorectal cancer is one of the main causes.

Immunological tests used as screening and developed to detect human hemoglobin are more accurate and do not require special diets for patients.

Although the screening test is not specific to determine which disease is causing the bleeding, it serves as an important alert for the beginning of a clinical trial.

Early diagnosis and immediate treatment have shown a significant reduction in complications and mortality from colorectal cancer.

**BIBLIOGRAPHIC REFERENCES**

1- Nakama H. et al., A comparative study of immunochemical fecal testes for detection of colorectal adenomatous polyps. Hepato-Gastroenterology. 2000, 47:386-389.

2- Uchida K. et al., Immunochemical detection of human blood in feces. Clinica Chimica Acta. 1990, 189:267-274.

3- Tommaso CL, Salzeider K, Arif M, et al. Serial myoglobin Gorzynski TJ, Krco CJ, David CS., An idiotope expresse don a monoclonal antibody specific for human haemoglobin chain na naturally occurring immunoglobulin(s). J Immunogenet. 1985 Dec;12(6):281-91.

4- Gorzynski TJ, Krco CJ, David CS., The Hb-2d cross-reactive idiootype. A common idiootype expressed by monoclonal and polyclonal antibodies to human haemoglobin. J Immunogenet. 1985 Dec;12(6):267-79

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**QUALITY ASSURANCE**

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date mentioned on the presentation package, provided they are stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

**CUSTOMER SERVICE**

Customer Advisory Service  
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOCLIN FAST OCCULT BLOOD kit:  
10269360484

Review: March/2026

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER