BIOCLIN FAST TROPONINA I

REF K277

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FINALIDADE

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1. Mair J., Genser N., Morandell D., et al. Cardiac troponin-I in the diagnosis of myocardial injury and infarction. Clin. Chem. Acta 1996;245:19-38;
2. Vallins W.J., Brand N.J., Dabhade N., et al. Molecular cloning of human cardiac troponin-I using polymerase chain reaction. FEBS Lett 1990;270:57-61;
3. Bodor G.S., Porter S., Landt Y., et al. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in suspected cases of myocardial infarction. Clin. Chem. 1992;38:2203-14;
4. Zaninotto M., Altinier S., Lachin M., et al. Fluoroenzymometric method to measure cardiac troponin-I in the sera of patients with myocardial infarction. Clin. Chem. 1996;42:1460-6;
5. OH S., Foster K., Datta P., et al. Use of a dual monoclonal solid phase and polyclonal detector to create an immunoassay for the detection of human cardiac troponin-I. Clinical Biochemistry 2000;33(4);255-262.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.
7. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

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INSTRUCTIONS FOR USE

FUNCTION

Test for rapid and qualitative determination of cTnI (Human Cardiac Troponin I) in samples of whole blood, serum or plasma (Heparin or EDTA) through immunochromatographic test. Only for *in vitro* diagnostic use.

ACTION PRINCIPLE

Methodology: Immunochromatography

The Bioclin Fast Troponin I kit is an immunochromatographic test for the detection of Human Cardiac Troponin I (cTnI). When a sample is added, interaction occurs with anti-cTnI monoclonal antibodies conjugated to colloidal gold particles. Then, the sample migrates through the membrane by capillarity and interacts with the Anti-cTnI monoclonal antibodies present in the "T" line. If the sample contains Troponin I (cTnI), a colored band will appear in the "T" test region, indicating a reactive result. Failure to form a line in the "T" region indicates absence of Troponin I in the analyzed sample. The control line consists of immobilized anti-mouse IgG. In this region, excess conjugate is captured by this antibody, regardless of whether the sample is positive or negative, forming the complex: anti-mouse IgG antibodies – conjugate. As a result, a red band appears, indicating the perfect functioning of the test.

REAGENTS

Number 1 (R1): Cassette - Store between 2 and 30 °C. **Do not freeze.** Each sachet contains:

- A)** 01 plastic cassette (device) with 01 test strip containing:
 - Control Region (C): Anti-mouse IgG polyclonal antibody (> 0.05 µg).
 - Test Region (T): Anti-cTnI monoclonal antibodies (> 0.05 µg).
 - Conjugate: Anti-cTnI monoclonal antibodies (> 0.05 µg) conjugated with colloidal gold particles.

- B)** 01 sachet of silica.

Number 2 (R2) Diluent - Store between 2 and 30°C. **Do not freeze.** Contains: Buffer, surfactant and preservative.

Collection Accessories:

- Disposable Plastic Pipette.
- Disposable lancet.

PRESENTATION

Reagent	Nº of Cassettes per Package	Diluent
Cassette + Diluent	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
	50	6.0 mL
Cassette + Diluent (Includes Pipettes)*	100	2 x 6.0 mL
	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
Cassette + Diluent (Includes Lancetes and Pipettes)**	50	6.0 mL
	100	2 x 6.0 mL
	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
	50	6.0 mL
	100	2 x 6.0 mL

* These presentations are accompanied by disposable plastic pipettes, proportionally to the number of cassettes.
 ** These presentations are accompanied by disposable plastic lancets and pipettes, proportionally to the number of cassettes.

OPERATIONAL EQUIPMENT AND SUPPLIES

Clock or stopwatch. They can be found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORTATION CONDITIONS

The storage and transport temperature should be 2 at 30°C. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- **Only for *in vitro* diagnostic use.**
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples used are biological, the safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Open the aluminum envelope to remove the cassette, only when using it.
- 5- The kit must not be frozen or exposed to high temperatures (above 45°C), as this will cause deterioration.
- 6- The reagents, as well as the samples, must reach room temperature before starting the test.
- 7- Do not mix reagents from kits with the same batch or with different batch numbers.
- 8- Do not use the kit beyond its expiration date. The expiry date is indicated on the product label.

9- The cassette (test device), the plastic pipette and the lancet are disposable and cannot be reused.

10- The kit materials that came into contact with the patient's sample must have a special destination. Lancets must be disposed of in sharps waste. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biowaste.

11- Secondary packaging (box), Instructions for Use, aluminum envelope and silica sachet can be disposed of in regular trash.

12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Service Customer Advice) by Quibasa.

14- Do not use the product if the packaging is damaged.

15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Serum or Plasma: Serum and plasma samples are stable for 3 days when kept between 2 and 8 °C and for 1 month at -20 °C, and must be thawed for one-time use only. Hemolyzed and lipemic samples should be avoided as they may cause a false positive result. Sera and plasmas containing precipitates may give inconsistent results. Prolonged storage or elevated temperature can deteriorate the cTnI protein. To obtain the plasma, perform the collection using Heparin or EDTA as an anticoagulant.

Whole Blood: Collect whole blood by venipuncture or with the aid of a lancet, using the appropriate anticoagulant (Heparin or EDTA). Whole blood can be used for immediate testing or can be stored at 2 to 8°C for up to 3 days. **Do not freeze whole blood.**



*At the time of use, samples must be stabilized at room temperature.

PROCESS DESCRIPTION

The BIOCLIN FAST TROPONIN I cassette is stable (after opening the aluminized sachet) for 2 hours. This stability may vary depending on the environmental conditions after opening the aluminized envelope.

TECHNIQUE - VENOPUNCTION









- 1- The sample must be at room temperature before starting the test.
- 2- Remove the cassette (test device) from its protective packaging, place it on a clean and level surface and identify it properly.
- 3- **For Serum, Plasma or Whole Blood:** Transfer 50 µL of serum, plasma or whole blood into the sample well.
- 4- Hold the vial vertically and apply 1 drop (approximately 50 µL) of Diluent (Reagent N° 2) in the sample well.

5- Wait for the lines to form after resting for 15 to 20 minutes. **Do not interpret after 20 minutes.**

TECHNIQUE - DIGITAL PUNCTURE / LANCET

- 1- Remove the cassette from its protective packaging, place it on a clean and level surface and identify it properly.
- 2- Perform asepsis on the tip of the finger to be perforated with 70% alcohol (w/w).
- 3- Press the tip of the finger that will be pierced by the lancet so that there is accumulation of blood in the region.
- 4- Remove the protective cap from the lancet.
- 5- Position and press the lancet firmly under the area to be lanced. Then the blood will come out through the pierced area.
- 6- Collect the blood using the disposable plastic pipette that comes with the kit. Place the pipette over the drop of blood and let the blood flow by capillarity, without pressing the bulb, up to the line marked on the disposable plastic pipette (equivalent to 50 µL of blood).
- 7- In the sample well, dispense the blood by pressing the disposable plastic pipette bulb.
- 8- Hold the flask vertically and apply 1 drop (approximately 50 µL) of Diluent (Reagent N° 2) in the sample well.
- 9- Wait for lines to form after resting for 15 to 20 minutes. **Do not interpret after 20 minutes.**


RESULT INTERPRETATION


Reagent Test for Troponin I: Formation of a red line in the test region (T) and another line in the control region (C) in the first 15 to 20 minutes. **Do not interpret after 20 minutes.**

Test Non-Reactive: Formation of a red line in the control region (C) and complete absence of a red line in the test region (T). **Do not interpret after 20 minutes.**

Invalid Result: Complete absence of line in the control region (C) with or without red line in the test region (T). Test the sample again.

**Notes:**

- 1- Non-reagent results should only be read within 20 minutes. Do not interpret after 20 minutes.
- 2- The line in the control region may appear before 15 minutes of test incubation, and this does not mean that non-reagent results can be interpreted ahead of time.
- 3- The results must not be interpreted after 20 minutes to avoid incorrect results.
- 4- Any intensity of color in the test line, even if very weak, characterizes the reagent result.

PROCESS LIMITATIONS

- 1- It is important to use the correct sample volume, as lower or higher volumes may determine erroneous results.
- 2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretation of the results.
- 3- The BIOCLIN FAST TROPONIN I kit, for cTnI, is a specific test to aid in the confirmation of a myocardial infarction. Reactive results must be confirmed with another test and available clinical information. Isolated test results cannot be used for definitive diagnosis.
- 4- The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERING

The BIOCLIN FAST TROPONIN I kit has been evaluated for possible interference. No significant interferences were found in samples with up to 0.6 mg/mL of Bilirubin, 13 mg/mL of Triglycerides, 40 mg/mL of Hemoglobin and 160 mg/mL of Albumin. Interferences with concentrations above those mentioned may cause false-positive or false-negative results. There was no evidence of interference with Rheumatoid Factor up to 980 IU/mL. No significant interferences were observed in 15 samples with reactive or non-reactive results for Troponin I.

CROSS REACTIVITY

A cross-reactivity study was evaluated using homologous proteins of Troponin I. In serum and plasma samples containing concentrations of 0, 0.5, 1.0, 2.0, 5.0, and 65 ng/ mL of cTnI, the following homologous proteins: Human Cardiac Troponin T (cTnT) 5000 ug/L, Human Skeletal Troponin (sTnI) 350 ug/L, and Human Cardiac Troponin C (cTnC) 1000 ug/L. There was no evidence of cross-reactivity with cTnT, sTnI and cTnC.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE**PRECISION****Repeatability**

For the evaluation of repeatability, 10 dosages were performed during 01 day with 04 different samples, obtaining the following results:

SAMPLES	Nº REPETITIONS	EXPECTED RESULTS	FOUND RESULTS
Serum	10	100% Reagent	100% Reagent
Serum	10	100% Non-Reagent	100% Non-Reagent
Plasma	10	100% Reagent	100% Reagent
Whole Blood	10	100% Non-Reagent	100% Non-Reagent

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 6 different samples, obtaining the following results:

SAMPLES	Nº DAYS	Nº REPETITIONS	EXPECTED RESULTS	FOUND RESULTS
Serum	03	10	100% Reagent	100% Reagent
Serum	03	10	100% Non-Reagent	100% Non-Reagent
Whole Blood	03	10	100% Non-Reagent	100% Non-Reagent
Whole Blood	03	10	100% Reagent	100% Reagent
Plasma	03	10	100% Reagent	100% Reagent
Plasma	03	10	100% Non-Reagent	100% Non-Reagent

SENSITIVITY AND CLINICAL SPECIFICITY

For the BIOCLIN FAST TROPONIN I kit, a clinical sensitivity and specificity study was carried out and 350 clinically known samples were analyzed and confirmed by the reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOCLIN FAST TROPONIN I kit is 99.00% and the clinical specificity is 99.19%.

		Expected Result	Found Result
BIOCLIN FAST TROPONIN I	Reagent	101	100
	Non-Reagent	249	247
	Total	350	347

Clinical Sensitivity: 99.00% (100/101) – 95% CI: 97.09% - 100%
Clinical Specificity: 99.19% (247/249) – 95% CI: 98.09% - 100%

ANALYTICAL SENSITIVITY

The Bioclin Fast Troponin I test was able to detect a cTnI concentration equal to or greater than 0.5 ng/mL.

DIAGNOSTIC SIGNIFICANCE

Serum levels of cardiac enzymes and isoenzymes are essential for the diagnosis or exclusion of myocardial damage. Cardiac Troponin I is specific for cardiac tissue and is detected in serum only in cases of myocardial damage. It allows the early identification of patients with symptoms suggestive of ischemia, allowing the identification of patients from 4 to 6 hours to 6 to 10 days after the infarction. It is also possible to identify false elevations of CK-MB as in rhabdomyolysis. Troponin I increases its value up to forty times in the first peak, while CK-MB only six to nine times.

BIBLIOGRAPHIC REFERENCES

1. Mair J., Genser N., Morandell D., et al. Cardiac troponin-I in the diagnosis of myocardial injury and infarction. Clin. Chem. Acta 1996;245:19-38;
2. Vallins W.J., Brand N.J., Dabhade N., et al. Molecular cloning of human cardiac troponin-I using polymerase chain reaction. FEBS Lett 1990;270:57-61;
3. Bodor G.S., Porter S., Landt Y., et al. Development of monoclonal antibodies for an assay of cardiac troponin-I and preliminary results in suspected cases of myocardial infarction. Clin. Chem. 1992;38:2203-14;
4. Zaninotto M., Altinier S., Lachin M., et al. Fluoroenzymometric method to measure cardiac troponin-I in the sera of patients with myocardial infarction. Clin. Chem. 1996;42:1460-6;
5. OH S., Foster K., Datta P., et al. Use of a dual monoclonal solid phase and polyclonal detector to create an immunoassay for the detection of human cardiac troponin-I. Clinical Biochemistry 2000;33(4):255-262.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.
7. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOCLIN FAST TROPONIN I kit: 10269360436

Review: October/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER