Bioclin

BIOCLIN POCT FIA 25-OH VITAMINA D

REF K333

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FINALIDADE

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- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- Byung-Chul, S. et al. Vitamin D side chain triazole analogs via cycloaddition 'click' chemistry. J. of Tetrahedron. Letters, v. 48, n. 45, p. 4623-4625, 2004.
- Holick, M. F. Vitamin D Deficiency. N. Engl. J. Med., v. 357, p. 266-281, 2007.
- Holick, M. F. MrOs is D-ficient. J. Clin. Endocrinol. Metab., v. 94, n. 4, p. 1092–1093, 2009.
- Holick, M. F. Vitamin D status: measurement, interpretation, and clinical application. Ann. Epidemiol., v. 19, n. 2, p. 73-79, 2009.
- Kwak, H. S. et al. Efficacy of the measurement of 25-hydroxyvitamin D₂ and D₃ levels by using PerkinElmer liquid chromatography-tandem mass spectrometry vitamin D kit compared with DiaSorin radioimmunoassay kit and Elecsys vitamin D total assay. Ann. Lab. Med., v. 35, n. 2, p. 263-265, 2015.
- Moreira, C. A. Reference values of 25-hydroxyvitamin D revisited: a position statement from the Brazilian Society of Endocrinology and Metabolism (SBEM) and the Brazilian Society of Clinical Pathology/Laboratory Medicine (SBPC). Arch. Endocrinol. Metab., v. 64, n. 4, p. 462-636, 2020.
- Bioclin – Dados de arquivos.

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- 1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 2- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- 3- Byung-Chul, S. et al. Vitamin D side chain triazole analogs via cycloaddition 'click' chemistry. J. of Tetrahedron. Letters, v. 48, n. 45, p. 4623-4625, 2004.
- 4- Holick, M. F. Vitamin D Deficiency. N. Engl. J. Med., v. 357, p. 266-281, 2007.
- 5- Holick, M. F. MrOs is D-ficient. J. Clin. Endocrinol. Metab., v. 94, n. 4, p. 1092–1093, 2009.
- 6- Holick, M. F. Vitamin D status: measurement, interpretation, and clinical application. Ann. Epidemiol., v. 19, n. 2, p. 73-79, 2009.
- 7- Kwak, H. S. et al. Efficacy of the measurement of 25-hydroxyvitamin D2 and D3 levels by using PerkinElmer liquid chromatography-tandem mass spectrometry vitamin D kit compared with DiaSorin radioimmunoassay kit and Elecsys vitamin D total assay. Ann. Lab. Med., v. 35, n. 2, p. 263-265, 2015.
- 8- Moreira, C. A. Reference values of 25-hydroxyvitamin D revisited: a position statement from the Brazilian Society of Endocrinology and Metabolism (SBEM) and the Brazilian Society of Clinical Pathology/Laboratory Medicine (SBPC). Arch. Endocrinol. Metab., v. 64, n. 4, p. 462-636, 2020.
- 9- Bioclin – Datos de archivos.

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INSTRUCTIONS FOR USE

FUNCTION

Test for the quantitative determination of 25-OH vitamin D in human biological samples of serum, plasma, and whole blood, using fluorescent immunoassay. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay.

The BIOCLIN POCT FIA 25-OH VITAMIN D kit is a fluorescent immunoassay based on the sandwich detection principle, designed for the quantitative determination of 25-OH vitamin D in human serum, plasma, or whole blood samples (collected in EDTA, citrate, or heparin). When the sample is added to its application site, the 25-OH vitamin D (25-(OH)VD) molecules present in the sample interact with the anti-25-(OH)VD antibodies conjugated to fluorescent microspheres, forming the 25-OH vitamin D-antibody conjugate complex. As the sample migrates by capillary action across the nitrocellulose membrane, this complex is captured by the test line, which also consists of anti-25-(OH)VD antibodies. The sample will continue to migrate through the membrane until it reaches the control line, which consists of anti-rabbit IgG antibodies (where rabbit IgG antibodies conjugated to fluorescent microspheres will be captured). At the end of the incubation period, the cassette is read on the immunofluorescence reader. The higher the concentration of 25-OH Vitamin D present in the sample, the more complexes are immobilized on the test line and, consequently, the higher the fluorescence signal generated and captured by the reader. The concentration of 25-OH Vitamin D in the sample, in ng/mL, is calculated according to the assay calibration curve (contained on the Identification Card).

REAGENTS

1- Cassette - Store between 2 and 30 °C. Each sachet contains:

- A)** 01 plastic device with 01 test strip containing:
 - Test line region: mouse anti-25-OH vitamin D antibody.
 - Control line region: goat anti-rabbit IgG antibody.
 - Conjugate: fluorescent microspheres conjugated to mouse anti-25-OH vitamin D antibodies and fluorescent microspheres conjugated to rabbit IgG antibodies.
- B)** 01 silica sachet.

2- Sample Diluent - Store between 2 and 30 °C. Contains: Buffer solution, stabilizers and preservatives.

3- Identification Card

PRESENTATION

REAGENTS	NUMBER OF REAGENTS PER PACKAGE
1- Cassette	25 units
2- Sample Diluent	25 units (25 x 0.160 mL)
3- Identification Card	1 unit

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation lower than 1.5%.
- 2- BIOCLIN POCT FIA Immunofluorescence Reader.
- 3- Stopwatch or clock.
- 4- Quality Control Tools.

STORAGE AND TRANSPORTATION CONDITIONS

The storage and transportation temperature should be between 2 °C and 30 °C. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- BIOCLIN POCT FIA 25-OH VITAMIN D should only be used in conjunction with the BIOCLIN POCT FIA Immunofluorescence Reader.

3- Strictly follow the proposed methodology to obtain accurate results.

4- The cassettes (test devices) do not contain contaminating substances, but since the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.

5- Open the aluminum envelope to remove the cassette only at the time of use.

6- The kit should not be frozen or exposed to high temperatures (above 45 °C), as this will cause deterioration.

7- The reagents, as well as the samples, must reach room temperature before starting the test.

8- Do not mix reagents from kits with the same batch or with different batch numbers.

9- Do not use the kit after its expiration date. The expiration date is indicated on the product label.

10- The test device (R1- Cassette) and sample dilution tubes (R2- Sample Diluent) are disposable and cannot be reused.

11- The kit materials that have come into contact with the patient's sample must be disposed of in a special way. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biological waste.

12- Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular waste.

13- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.

14- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or by requesting it through Quibasa's SAC (Customer Advisory

Service).

15- Do not use the product if the packaging is damaged.

16- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

SAMPLES

Serum or Plasma (EDTA, Citrate or Heparin)

Samples should be collected by venipuncture, in tubes with or without anticoagulant (EDTA, Citrate or Heparin). To avoid hemolysis, centrifuge the samples to obtain serum/plasma and test them as soon as possible. Do not keep the samples at room temperature for long periods. Serum and plasma samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 3 days. If the samples cannot be analyzed within 3 days, they can be stored for up to 30 days at a temperature of -20 °C. Serum and plasma samples stored for periods longer than recommended and samples collected with anticoagulants other than those mentioned should not be used. Highly hemolyzed, lipemic samples and samples containing precipitates should be avoided, as they may provide inconsistent results.

Whole Blood (EDTA, Citrate or Heparin)

Whole blood samples should be collected by venipuncture in tubes containing Heparin, Citrate or EDTA as anticoagulants. Samples should be tested immediately or as soon as possible after collection. Whole blood collected in tubes containing anticoagulants can be stored under refrigeration, between 2 and 8°C, for a maximum period of 2 days. Do not use samples collected with anticoagulants other than those mentioned and do not freeze whole blood samples.

PROCESS DESCRIPTION

Stability After Opening

The stability test results show that the BIOCLIN POCT FIA 25-OH VITAMIN D kit is stable after opening for up to 70 minutes. This stability may vary depending on the conditions of use of the product and the environment. However, for safety reasons, we recommend that the product should be used immediately after removing the cassette from the sachet.

TECHNIQUE

Check the BIOCLIN POCT FIA Immunofluorescence Reader manual for complete instructions on handling the equipment.

Before starting the test, allow all reagents (Cassette and Sample Diluent) and samples to stabilize at room temperature (15-30 °C) for at least 40 minutes.

1- Check if the information contained in the Identification Card has been properly entered into the equipment (with each new batch of the kit, the calibration and validity data contained in the Identification Card must be transferred to the equipment).

2- Remove the cassette (test device) from its protective packaging, place it on a clean, level surface and identify it appropriately.

3- Make sure that the Sample Diluent solution (Reagent No. 2) is at the bottom of the dilution tube.

4- Add 10 µL of sample to a dilution tube containing the Sample Diluent.

5- Close the tube cap and mix by inverting the sample 6 to 8 times until the sample and diluent are completely

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

homogenized.

6- Pipette 50 µL of diluted and properly mixed sample into the sample application hole of the cassette. **Apply the sample to the cassette within 15 minutes of diluting it.**

7- Option 1: Insert the cassette into the equipment immediately after applying the sample and select the “**Timing Test**” option. After 15 minutes, the equipment will take the reading and the result will appear on the display.

Option 2: After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the “**Quick Test**” option. The equipment will take the reading immediately and the result will appear on the display.

RESULTS

The BIOCLIN POCT FIA Immunofluorescence Reader calculates the 25-OH VITAMIN D concentration automatically, using the calibration curve contained in the R3 (Identification Card). The information on the Identification Card must be properly entered into the equipment for each new batch of reagents.

Quantitative results are measured in ng/mL. The measuring range of the kit is between 10 and 100 ng/mL (limit of quantification and highest point of curve, respectively). Results below the lower limit of detection are expressed as <10 ng/mL. Results above the upper limit of detection are expressed as >100 ng/mL.

PROCESS LIMITATIONS

1- It is important to use the correct sample volume, as lower or higher volumes may produce erroneous results.

2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretations of the results.

3- Reactive results must be confirmed with another test and available clinical information. The isolated test results cannot be used for definitive diagnosis.

4- It is recommended that established clinical guidelines for optimal vitamin D values be followed.

5- The results provided by this kit must be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

REFERENCE VALUES

Currently, there is no standard definition for optimal 25-OH vitamin D levels. Many experts consider commonly used population-based reference values to be too low. Health-based reference values are recommended to replace population-based reference values. Therefore, the 25-OH vitamin D level currently most widely accepted by experts is described below.

Sufficient	≥ 30 ng/mL
-------------------	------------

These values should be used as a guide, and each laboratory should create its own range of reference values according to the population served. The results provided by this kit should be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed by bilirubin up to 350 µmol/L, triglycerides up to 40 mmol/L, hemoglobin up to 2 g/L, rheumatoid factor up to 1,500 IU/mL, heterophile antibodies up to 1,000 ng/mL, and antinuclear antibodies up to 1:640.

CROSS-REACTIVITY

No false-positive results were observed for 1,25-(OH)₂ D₃ up to 100 ng/mL, 1,25-(OH)₂ D₂ (ng/mL) up to 100 ng/mL, Vitamin D₃ up to 1,000 ng/mL, and Vitamin D₂ up to 1,000 ng/mL. Despite the results found, the possibility of cross-reactivity with other analytes or with higher concentrations cannot be completely ruled out. The final diagnosis should consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. For this purpose, it is recommended to use controls, which allow the assessment of the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The BIOCLIN POCT FIA 25-OH VITAMIN D kit was compared with another commercially available method for 25-OH vitamin D measurement. 83 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 1.0259x + 0.3529$ and the correlation coefficient was 0.9702. Based on these results, it can be concluded that the kit presents good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 2 samples with different values, obtaining the following concentration results:

REPEATABILITY	SAMPLE	
	1	2
Mean	29.64	79.40
Standard Deviation	2.45	7.38
Coefficient of Variation (%)	8.3%	9.3%

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 samples with different values, obtaining the following concentration results:

REPRODUCIBILITY	SAMPLE	
	1	2
Average	30.04	78.68
Standar Deviation	2.5	6.8
Coefficient of variation (%)	8.3%	8.6%

ANALYTICAL SENSITIVITY

The Limit of Blank (LoB) of the product BIOCLIN POCT FIA 25-OH VITAMIN D is 2.0 ng/mL, the Limit of Detection (LoD) is 5.0 ng/mL and the Limit of Quantification (LoQ) is 10.0 ng/mL.

LINEARITY

The BIOCLIN POCT FIA 25-OH VITAMIN D kit is linear up to 100 ng/mL (highest point of the calibration curve).

HIGH DOSE PROZONE EFFECT

No high-dose prozone effect was observed at up to 500 ng/mL of 25-OH Vitamin D.

DIAGNOSTIC SIGNIFICANCE

The term vitamin D is used to designate the compounds vitamin D₃ (also known as cholecalciferol) and vitamin D₂ (ergocalciferol). Vitamin D is a steroid derivative belonging to the fat-soluble class of vitamins, synthesized primarily by human skin after stimulation by ultraviolet radiation. The human body is capable of synthesizing only the D₃ form of the vitamin. Therefore, a small portion of the vitamin D intake also comes from food and supplements, in the D₂ form. Vitamin D is converted to 25-OH vitamin D in the liver, which in turn is converted to the active form 1,25-OH vitamin D in the kidneys. The 25-OH vitamin D form is the main form of vitamin D found in vivo (where its concentration is approximately 1,000 times greater than the active form) and reflects the level of vitamin D storage in the human body. Consequently, 25-OH vitamin D levels are also related to the clinical symptoms of vitamin D deficiency. One of vitamin D's most important roles is to maintain skeletal calcium balance, promoting calcium absorption in the intestines, bone resorption by increasing the number of osteoclasts, maintaining calcium and phosphate levels for bone formation, and enabling the proper functioning of parathyroid hormone (PTH), also responsible for maintaining blood calcium levels. In children, severe vitamin D deficiency leads to a bone malformation known as rickets, which causes symptoms such as growth retardation, bowed legs, weakness, and pain in the spine, pelvis, and legs. In adults, especially the elderly, vitamin D deficiency causes muscle weakness, increasing the risk of falls. Furthermore, vitamin D deficiency is a common cause of secondary hyperparathyroidism, which can result in osteomalacia, increased bone turnover, and reduced bone density (osteoporosis), with a consequent increased risk of fractures. Vitamin D deficiency is also linked to diabetes, different forms of cancer, Parkinson's disease, cardiovascular disease, hypertension, chronic kidney disease, autoimmune diseases, respiratory diseases, and problems with innate immunity.

BIBLIOGRAPHIC REFERENCES

- 1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 2- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- 3- Byung-Chul, S. et al. Vitamin D side chain triazole analogs via cycloaddition 'click' chemistry. J. of Tetrahedron. Letters, v. 48, n. 45, p. 4623-4625, 2004.
- 4- Holick, M. F. Vitamin D Deficiency. N. Engl. J. Med., v. 357, p. 266-281, 2007.
- 5- Holick, M. F. MrOs is D-ficient. J. Clin. Endocrinol. Metab., v. 94, n. 4, p. 1092-1093, 2009.
- 6- Holick, M. F. Vitamin D status: measurement, interpretation, and clinical application. Ann. Epidemiol., v. 19, n. 2, p. 73-79, 2009.
- 7- Kwak, H. S. et al. Efficacy of the measurement of 25-hydroxyvitamin D₂ and D₃ levels by using PerkinElmer liquid chromatography-tandem mass spectrometry vitamin D kit compared with DiaSorin radioimmunoassay kit and Elecsys

vitamin D total assay. Ann. Lab. Med., v. 35, n. 2, p. 263-265, 2015.

8- Moreira, C. A. Reference values of 25-hydroxyvitamin D revisited: a position statement from the Brazilian Society of Endocrinology and Metabolism (SBEM) and the Brazilian Society of Clinical Pathology/Laboratory Medicine (SBPC). Arch. Endocrinol. Metab., v. 64, n. 4, p. 462-636, 2020.

9- Bioclin – Archive Data.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration number for BIOCLIN POCT FIA 25-OH VITAMIN D kit: 10269360473

Review: August/2025

UNIVERSAL SYMBOLOGY

 REF	CATALOG NUMBER	 MADE BY
 LOT	LOT NUMBER	 CONTROL
 MANUFACTURING DATE		 POSITIVE CONTROL
 VALIDITY DATE (last day of the month)		 NEGATIVE CONTROL
 TEMPERATURE LIMIT (store)		 BIOLOGICAL RISK
 CONTENT IS SUFFICIENT FOR <N> TEST		 FLAMMABLE
 SEE INSTRUCTIONS FOR USE		 CORROSIVE
 IVD IN VITRO DIAGNOSTIC PRODUCT		 TOXIC
 KEEP AWAY FROM SUNLIGHT		 DO NOT USE IF PACKAGE IS DAMAGED
 DO NOT REUSE		 PRODUCT STERILIZED
 CAUTION		 DANGER