BIOCLIN POCT FIA 25 OH VITAMINA D

REF **K279**

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FINALIDADE

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1. Borai, Anwar et al. "The pre-analytical stability of 25-hydroxyvitamin D: Storage and mixing effects."Journal of clinical laboratory analysis, vol. 34,2 (2020)
2. Byung-Chul Suh,HeungBae Jeon, Gary H. Posner and Steven M. Silverman.Vitamin D side chain triazole analogs via cycloaddition `click` chemistry.[J].J of Tetrahedron Letters, ,2004,18(45): 4623-4625.
3. Holick MF. Vitamin D Deficiency. N Engl J Med.2007;357:266–81.
4. Holick MF. MrOs is D-ficient. J Clin Endocrinol Metab.2009;94(4):1092–3.
5. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.
6. Levinson SS.The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114.
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1. Borai, Anwar et al. "The pre-analytical stability of 25-hydroxyvitamin D: Storage and mixing effects."Journal of clinical laboratory analysis, vol. 34,2 (2020)
2. Byung-Chul Suh,HeungBae Jeon, Gary H. Posner and Steven M. Silverman.Vitamin D side chain triazole analogs via cycloaddition 'click' chemistry.[J].J of Tetrahedron Letters, ,2004,18(45): 4623-4625.
3. Holick MF. Vitamin D Deficiency. N Engl J Med.2007;357:266–81.
4. Holick MF. MrOs is D-ficient. J Clin Endocrinol Metab.2009;94(4):1092–3.
5. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.
6. Levinson SS.The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114.
7. QUIBASA: Datos do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Test for quantitative determination of 25-OH D (Vitamin D) in biological samples (plasma) through fluorescent immunoassay. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay.
The BIOCLIN POCT FIA 25 OH VITAMIN D kit is a competitive test for the detection of Vitamin D in plasma samples. When sample and diluent are added, the Vitamin D in the sample dissociates from its binding protein and interacts with Anti-25-OH-D Monoclonal Antibodies conjugated to fluorescent microbead particles. Then, the sample migrates through the membrane by capillarity, up to the detection line (with 25-OH D molecules present in the "T" line). The sample will continue migrating through the membrane until it reaches the control region "C", where the conjugate containing IgG Antibodies interacts with immobilized Anti-IgG Antibodies. At the end of the incubation period, the cassette must be inserted into the reader for quantification. The Vitamin D concentration of the sample correlates inversely with the intensity of the fluorescence signal emitted by the test line of the assay.

The BIOCLIN POCT FIA reader calculates the Vitamin D concentration in the sample in ng/mL according to the assay calibration curve (contained in the Identification Chip).

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 8 °C. **Do not freeze.** Each sachet contains:

- A)** 01 plastic device with 01 test strip containing:
 - Control region (C): Anti-IgG rabbit Antibodies.
 - Test region (T): 25-OH-D molecules.
- Conjugate: Anti-25(OH)-D monoclonal antibody conjugated with fluorescent microspheres and rabbit IgG antibodies conjugated with fluorescent microspheres.
- B)** 01 sachet of silica.

Number 2 (R2) Solution A - Store between 2 and 8 °C. **Do not freeze.** Contains: Buffer, surfactant and preservative.
Number 3 (R3) Solution B - Store between 2 and 8 °C. **Do not freeze.** Contains: Buffer, surfactant and preservative.
Number 4 (R4) Diluent - Store between 2 and 8 °C. **Do not freeze.** Contains: Buffer, surfactant and preservative.
Number 5 (R5) Identification Chip.

PRESENTATION

Reagents	Number of Reagents per Package
1 - Cassette	25 units
2 - Solution A	1 unit
3 - Solution B	1 unit
4 - Diluent	25 units
5 - Identification Chip	1 unit

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table
- Instructions for use (manual)

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation less than 1.5%.
- 2- BIOCLIN POCT FIA incubator (optional).
- 3- BIOCLIN POCT FIA Reader.
- 4- Stopwatch or clock.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- **For *in vitro* diagnostic use only.**
- 2- **The BIOCLIN POCT FIA 25 OH VITAMIN D must only be used together with the BIOCLIN POCT FIA reader.**
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples used are biological, the safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Open the aluminum envelope to remove the cassette, only when using it.
- 5- The kit must not be frozen or exposed to high temperatures (45°C), as this will cause deterioration.
- 6- The reagents, as well as the samples, must reach room temperature before starting the test.
- 7- Do not mix reagents from kits with the same batch or with different batch numbers.
- 8- Do not use the kit beyond its expiration date. The expiry date is indicated on the product label.
- 9- The cassette (test device) and sample dilution tubes are disposable and cannot be reused.
- 10- The kit materials that came into contact with the patient's sample must have a special destination. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biowaste.
- 11- Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular trash.
- 12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the FISPQ (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service Customer Advice) by Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLE

Plasma: Samples must be collected by venipuncture and can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 1 day. If samples cannot be analyzed within 1 day, they can be stored for up to 3 months at -20°C (freezer)*. To obtain plasma, collect using EDTA as an anticoagulant.

*At the time of use, samples must be stabilized at room temperature.

PROCESS DESCRIPTION

Stability After Opening

The stability of the BIOCLIN POCT FIA 25 OH VITAMIN D cassette (after opening the aluminized sachet) is 1 hour. This stability may vary according to the conditions of use of the product and the environment. However, to be on the safe side, we recommend that the BIOCLIN POCT FIA 25 OH VITAMIN D kit be tested immediately after removing the cassette from the sachet.

TECHNIQUE

- 1- **The sample, Cassette and Sample Dilution Tube must be between 15-30 °C to run the test.**
- 2- Check if the identification chip is inserted in the BIOCLIN POCT FIA reader.
- 3- Pipette 40 µL of solution A (Reagent N° 2) into the sample dilution tube (Reagent N° 4). Then pipette 40 µL of Solution B (Reagent N° 3) into the Sample Dilution Tube (Reagent N° 4). Homogenize the sample dilution tube.
- 4- Transfer 20 µL of **Plasma** to the tube containing the diluent (Reagent N° 4).
- 5- Homogenize the tube with the sample.
- 6- Incubate the sample in the dilution tube at room temperature (approximately 25°C) for 10 minutes.
- 7- Remove the cassette (test device) from its protective packaging, place it on a clean and level surface and identify it properly.
- 8- Pipette **80 µL of the diluted sample** into the sample well.
- 9- Wait 10 minutes to read the test. **Do not perform the reading before 10 minutes.**

Cassette Reading

- 1- Select the test mode on the BIOCLIN POCT FIA reader. For more information on how to select the reading mode, refer to the equipment manual.
- 2- Insert the **Cassette** into the BIOCLIN POCT FIA reader and wait for the result to be issued.

NOTES

- 1- The cassette must be used immediately after opening its packaging.
- 2- The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

Calibration

The calibration of the method is carried out from the **Identification Chip**, which must be changed with each change in the lot of reagents.

PROCESS LIMITATIONS

- 1- It is important to use the correct sample volume, as lower or higher volumes may determine erroneous results.
- 2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretation of the results.
- 3- The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERING

Samples with concentrations of up to 400µmol/L of Bilirubin, 20mmol/L of Triglycerides, 15mg/mL of Hemoglobin and 3250 IU/mL of Rheumatoid Factor do not interfere with the test result.

CROSS REACTIVITY

The BIOCLIN POCT FIA 25 OH VITAMIN D kit was evaluated for cross-reactivity with different analytes. This product has 100% specificity for 25-(OH)₂ Vitamin D3 and 99% specificity for 25-(OH)₂ Vitamin D₂. No falsely elevated results were observed for concentrations up to 30 ng/mL 24,25-(OH)₂ Vitamin D3, 50 ng/mL 3-epi-25-OH Vitamin D3, 50 ng/mL 3-epi-25-OH Vitamin D₂, 100 ng/mL 1,25-(OH)₂ Vitamin D3, 100 ng/mL 1,25-(OH)₂ Vitamin D₂, 1000 ng/mL Vitamin D3 and 1000ng/mL Vitamin D₂.

Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis must consider the clinical data of the patient with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

REFERENCE VALUES

Vitamin D Level	Reference Values (ng/mL)
Deficient	< 20,0
Insufficient	20.0 to 29.0
Sufficient	30.0 to 100.0
Potentially Toxic	> 100.0

These values should be used as guidelines, and each laboratory should create its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods and Methodological Specificity

The BIOCLIN POCT FIA 25 OH VITAMIN D kit was compared with another commercially available method for measuring 25 OH D. 100 analyzes were performed and the results were evaluated. The linear equation obtained was Y=1.01X+0.69 and the correlation coefficient 0.974. With these results, it can be concluded that the kit has good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 4 samples with different concentrations, obtaining the following results:

Repeatability	Sample			
	1	2	3	4
Concentration (ng/mL)	15.9	24.0	41.2	89.9
Standard Deviation (ng/mL)	1.0	1.2	3.5	2.7
Coefficient of Variation (%)	6	5	8	3

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 4 samples with different concentrations, obtaining the following results:

Reproducibility	Sample			
	1	2	3	4
Concentration (ng/mL)	15.6	26.3	41.6	91.6
Standard Deviation (ng/mL)	0.9	1.3	2.5	2.6
Coefficient of Variation (%)	5	5	6	3

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

SENSITIVITY

The analytical sensitivity of the BIOCLIN POCT FIA 25 OH VITAMIN D kit is 4.9 ng/mL.

PROZONE EFFECT

There was no prozone effect up to a concentration of 250 ng/mL.

LINEARITY

The operating range for the BIOCLIN POCT FIA 25 OH VITAMIN D reagent is 5 ng/mL to 120 ng/mL.

DIAGNOSTIC SIGNIFICANCE

The term Vitamin D is used to designate the compounds vitamin D3 (also known as cholecalciferol) and vitamin D2 (ergocalciferol). Exposure to ultraviolet radiation (sunlight) causes the production of vitamin D3. Vitamin D is transported in the bloodstream to the liver, where it is converted into the prohormone calcifediol. Circulating calcifediol can then be converted to calcitriol, the biologically active form of vitamin D, in the kidneys. One of the most important roles of vitamin D is to maintain skeletal calcium balance by promoting calcium absorption in the intestines, promoting bone resorption by increasing the number of osteoclasts, maintaining calcium and phosphate levels for bone formation, allowing bone to function properly. adequate levels of parathyroid hormone (PTH) to maintain blood calcium levels. Vitamin D deficiency can result in lower bone mineral density and an increased risk of reduced bone density (osteoporosis) or bone fracture, as vitamin D alters mineral metabolism in the body. Thus, vitamin D is also critical for bone turnover through its role as a potent stimulator of bone resorption.

BIBLIOGRAPHIC REFERENCES

1. Borai, Anwar et al. "The pre-analytical stability of 25-hydroxyvitamin D: Storage and mixing effects."Journal of clinical laboratory analysis, vol. 34,2 (2020)

2. Byung-Chul Suh,HeungBae Jeon, Gary H. Posner and Steven M. Silverman.Vitamin D side chain triazole analogs via cycloaddition 'click' chemistry.[J].J of Tetrahedron Letters, ,2004,18(45): 4623-4625.

3. Holick MF. Vitamin D Deficiency. N Engl J Med.2007;357:266–81.

4. Holick MF. MrOs is D-ficient. J Clin Endocrinol Metab.2009;94(4):1092–3.

5. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.

6. Levinson SS.The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 – Santa Branca
CEP 31565-130 – Belo Horizonte – MG – Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOCLIN POCT FIA 25 OH VITAMIN D kit:
10269360430

Review: June/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER