Bioclin

BIOCLIN POCT FIA BETA-HCG

REF K334

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Test for quantitative determination of human chorionic gonadotropin (HCG) in human biological samples of serum, plasma and whole blood, through fluorescent immunoassay. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay.

The BIOCLIN POCT FIA BETA-HCG kit is a fluorescent immunoassay based on the "sandwich" detection principle, intended for the quantitative determination of human chorionic gonadotropin (HCG) in human serum, plasma or whole blood samples (collected in EDTA, citrate or heparin). When the sample is added to its application site, interaction occurs between the HCG molecules present in the sample and the anti-beta-HCG antibodies conjugated to fluorescent microspheres, forming the HCG-conjugated antibody complex. As the sample migrates by capillarity through the nitrocellulose membrane, this complex is captured by the test line, which also consists of anti-beta-HCG antibodies. The sample will continue to migrate through the membrane until it reaches the control line, which consists of anti-rabbit IgG antibodies (where rabbit IgG antibodies conjugated to fluorescent microspheres will be captured). At the end of the incubation time, the cassette is read in the immunofluorescence reader. The higher the concentration of beta-HCG present in the sample, the more complexes are immobilized on the test line and, consequently, the higher the fluorescence signal generated and captured by the reader. The concentration of beta-HCG in the sample, in mIU/mL, is calculated according to the assay calibration curve (contained on the Identification Card).

REAGENTS

1- Cassette - Store between 2 and 30 °C. Each sachet contains:

- A)** 01 plastic device with 01 test strip containing:
 - Test line region: mouse anti-beta-HCG antibody.
 - Control line region: goat anti-rabbit IgG antibody.
 - Conjugate: fluorescent microspheres conjugated to mouse anti-beta-HCG antibodies and fluorescent microspheres conjugated to rabbit IgG antibodies.

B) 01 silica sachet.

2- Sample Diluent - Store between 2 and 30 °C. Contains: Buffer solution, stabilizers and preservatives.

3- Identification Card

PRESENTATION

REAGENTS	NUMBER OF REAGENTS PER PACKAGE
1- Cassette	25 units
2- Sample Diluent	25 units (25 x 0.360 mL)
3- Identification Card	1 unit

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation lower than 1.5%.
- 2- BIOCLIN POCT FIA Immunofluorescence Reader.
- 3- Stopwatch or clock.
- 4- Quality Control Tools.

STORAGE AND TRANSPORTATION CONDITIONS

The storage and transportation temperature should be between 2 °C and 30 °C. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- BIOCLIN POCT FIA BETA-HCG should only be used in conjunction with the BIOCLIN POCT FIA Immunofluorescence Reader.

3- Strictly follow the proposed methodology to obtain accurate results.

4- The cassettes (test devices) do not contain contaminating substances, but since the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.

5- Open the aluminum envelope to remove the cassette only at the time of use.

6- The kit should not be frozen or exposed to high temperatures (above 45 °C), as this will cause deterioration.

7- The reagents, as well as the samples, must reach room temperature before starting the test.

8- Do not mix reagents from kits with the same batch or with different batch numbers.

9- Do not use the kit after its expiration date. The expiration date is indicated on the product label.

10- The test device (R1- Cassette) and sample dilution tubes (R2- Sample Diluent) are disposable and cannot be reused.

11- The kit materials that have come into contact with the patient's sample must be disposed of in a special way. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biological waste.

12- Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular waste.

13- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.

14- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or by requesting it through Quibasa's SAC (Customer Advisory

Service).

15- Do not use the product if the packaging is damaged.

16- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

SAMPLES

Serum or Plasma (EDTA, Citrate or Heparin)

Samples should be collected by venipuncture, in tubes with or without anticoagulant (EDTA, Citrate or Heparin). To avoid hemolysis, centrifuge the samples to obtain serum/plasma and test them as soon as possible. Do not keep the samples at room temperature for long periods. Serum and plasma samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 3 days. If the samples cannot be analyzed within 3 days, they can be stored for up to 30 days at a temperature of -20 °C. Serum and plasma samples stored for periods longer than recommended and samples collected with anticoagulants other than those mentioned should not be used. Highly hemolyzed, lipemic samples and samples containing precipitates should be avoided, as they may provide inconsistent results.

Whole Blood (EDTA, Citrate or Heparin)

Whole blood samples should be collected by venipuncture in tubes containing Heparin, Citrate or EDTA as anticoagulants. Samples should be tested immediately or as soon as possible after collection. Whole blood collected in tubes containing anticoagulants can be stored under refrigeration, between 2 and 8°C, for a maximum period of 2 days. Do not use samples collected with anticoagulants other than those mentioned and do not freeze whole blood samples.

PROCESS DESCRIPTION

Stability After Opening

The stability test results show that the BIOCLIN POCT FIA BETA-HCG kit is stable after opening for up to 70 minutes. This stability may vary depending on the conditions of use of the product and the environment. However, for safety reasons, we recommend that the product should be used immediately after removing the cassette from the sachet.

TECHNIQUE

Check the BIOCLIN POCT FIA Immunofluorescence Reader manual for complete instructions on handling the equipment.

Before starting the test, allow all reagents (Cassette and Sample Diluent) and samples to stabilize at room temperature (15-30 °C) for at least 40 minutes.

1- Check if the information contained in the Identification Card has been properly entered into the equipment (with each new batch of the kit, the calibration and validity data contained in the Identification Card must be transferred to the equipment).

2- Remove the cassette (test device) from its protective packaging, place it on a clean, level surface and identify it appropriately.

3- Make sure that the Sample Diluent solution (Reagent No. 2) is at the bottom of the dilution tube.

4- Add 50 µL of sample to a dilution tube containing the Sample Diluent.

5- Close the tube cap and mix by inverting the sample 6 to 8 times until the sample and diluent are completely homogenized.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

6- Pipette 50 µL of diluted and properly mixed sample into the sample application hole of the cassette.

7- Option 1: Insert the cassette into the equipment immediately after applying the sample and select the "Timing Test" option. After 15 minutes, the equipment will take the reading and the result will appear on the display.

Option 2: After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the "Quick Test" option. The equipment will take the reading immediately and the result will appear on the display.

RESULTS

The BIOCLIN POCT FIA Immunofluorescence Reader calculates the beta-HCG concentration automatically, using the calibration curve contained in the R3 (Identification Card). The information on the Identification Card must be properly entered into the equipment for each new batch of reagents.

Quantitative results are measured in percentage (%). The measurement range of the kit is between 2 and 200,000 mUI/mL (quantification limit and maximum point of the curve, respectively). Results below the lower detection limit are expressed as <2.0 mUI/mL. Results above the upper detection limit are expressed as >200,000 mUI/mL.

PROCESS LIMITATIONS

1- It is important to use the correct sample volume, as lower or higher volumes may produce erroneous results.

2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretations of the results.

3- Reactive results must be confirmed with another test and available clinical information. The isolated test results cannot be used for definitive diagnosis.

4- The results provided by this kit must be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

REFERENCE VALUES

Negative	<5 mUI/mL
----------	-----------

3 weeks of pregnancy	6 to 72 mUI/mL
4 weeks of pregnancy	11 to 790 mUI/mL
5 weeks of pregnancy	237 to 7,817 mUI/mL
6 weeks of pregnancy	468 to 35,035 mUI/mL
7 weeks of pregnancy	3,107 to 173,174 mUI/mL
8 weeks of pregnancy	32,451 to 156,069 mUI/mL
9 weeks of pregnancy	64,091 to 152,690 mUI/mL
10 weeks of pregnancy	47,402 to 198,989 mUI/mL
12 weeks of pregnancy	28,621 to 223,361 mUI/mL
14 weeks of pregnancy	13,977 to 65,333 mUI/mL
15 weeks of pregnancy	12,247 to 74,451 mUI/mL
16 weeks of pregnancy	9,136 to 61,802 mUI/mL
17 weeks of pregnancy	8,224 to 60,665 mUI/mL
18 weeks of pregnancy	8,156 to 63,531 mUI/mL

These values should be used as a guide, and each laboratory should create its own range of reference values according to the population served. The results provided by this kit should be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed by bilirubin up to 350 µmol/L, triglycerides up to 40 mmol/L, hemoglobin up to 2 g/L, rheumatoid factor up to 1,500 IU/mL, heterophile antibodies up to 1,000 ng/mL, and antinuclear antibodies up to 1:640.

CROSS-REACTIVITY

No false positive results were observed for FSH up to 200 mIU/mL, LH up to 200 mIU/mL and TSH up to 200 mIU/mL. Despite the results found, the possibility of cross-reactivity with other analytes or with higher concentrations cannot be completely ruled out. The final diagnosis should consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. For this purpose, it is recommended to use controls, which allow the assessment of the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The BIOCLIN POCT FIA BETA-HCG kit was compared with another commercially available method for HbA1c measurement. 95 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 0.9637x + 145.48$ and the correlation coefficient was 0.9911. Based on these results, it can be concluded that the kit presents good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following concentration results:

REPEATABILITY	SAMPLE		
	1	2	3
Mean	10.20	489.06	9,667.72
Standard Deviation	0.84	46.01	762.77
Coefficient of Variation (%)	8.2	9.4	7.9

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following concentration results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Mean	10.25	502.53	10,007.32
Standard Deviation	0.9	45.0	806.4
Coefficient of Variation (%)	8.6	9.0	8.1

ANALYTICAL SENSITIVITY

The Limit of Blank (LoB) of the BIOCLIN POCT FIA BETA-HCG product is 1.0 mUI/mL. the Limit of Detection (LoD) is 2.0 mUI/mL and the Limit of Quantification (LoQ) is 2.0 mUI/mL.

LINEARITY

The BIOCLIN POCT FIA BETA-HCG kit is linear up to 200,000 mUI/mL (highest point of the calibration curve).

HIGH DOSE PROZONE EFFECT

No high-dose prozone effect was observed at up to 240,000 mUI/mL of beta-HCG.

DIAGNOSTIC SIGNIFICANCE

Human chorionic gonadotropin (hCG) is a glycoprotein hormone that is normally produced by the placenta during pregnancy. hCG molecules are composed of an alpha and a beta subunit, with a molecular weight of approximately 38 kDa. The alpha subunit is nearly identical to the alpha subunits of other glycoprotein hormones, such as thyroid-stimulating hormone (TSH), luteinizing hormone (LH), and follicle-stimulating hormone (FSH). Differences in the beta subunit of the respective hormones account for their biological specificity and immunochemical distinction. Monoclonal antibodies that recognize unique sites on the beta chain of the hCG molecule are essential for the differentiation of hCG from LH, FSH, and TSH. After conception, approximately 5 days after implantation of the blastocyst, hCG concentrations rise rapidly until they reach a peak near the end of the first trimester. High concentrations are observed throughout pregnancy. After delivery, hCG levels fall rapidly and become undetectable after a few days.

BIBLIOGRAPHIC REFERENCES

- 1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 2- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- 3- Ross, G. T. Clinical relevance of research on the structure of human chorionic gonadotropin. Am. J. Obstet. Gynecol., v. 129, n. 7, p. 795-808, 1977.
- 4- Kadar, N.; DeVore, G.; Romero, R. Discriminatory hCG Zone: Its Use in the Sonographic Evaluation for Ectopic Pregnancy. J. Obstet. Gynecol., v. 58, p. 156-161, 1981.
- 5- Schumacher, A.; Zencussen, A.C. Human Chorionic Gonadotropin-Mediated Immune Responses That Facilitate Embryo Implantation and Placentation. Front Immunol. 2019.
- 6- Montagnana, M. Human chorionic gonadotropin in pregnancy diagnostics. Clinica Chimica Acta, 2011.
- 7- Bioclin – Dados de arquivos.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration number for BIOCLIN POCT FIA BETA-HCG kit: 10269360469

Review: July/2025

UNIVERSAL SYMBOLOGY

 REF	CATALOG NUMBER	 MADE BY
 LOT	LOT NUMBER	 CONTROL
 MANUFACTURING DATE		 POSITIVE CONTROL
 VALIDITY DATE (last day of the month)		 NEGATIVE CONTROL
 TEMPERATURE LIMIT (store)		 BIOLOGICAL RISK
 CONTENT IS SUFFICIENT FOR <N> TEST		 FLAMMABLE
 SEE INSTRUCTIONS FOR USE		 CORROSIVE
 IN VITRO DIAGNOSTIC PRODUCT		 TOXIC
 KEEP AWAY FROM SUNLIGHT		 DO NOT USE IF PACKAGE IS DAMAGED
 DO NOT REUSE		 PRODUCT STERILIZED
 CAUTION		 DANGER