Bioclin

BIOCLIN POCT FIA BNP

REF **K287**

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FINALIDADE

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- Gobinet-Georges, A et al. "Stability of brain natriuretic peptide (BNP) in human whole blood and plasma."Clinical chemistry and laboratory medicine, vol. 38,6 (2000): 519-23.
- García de Guadiana-Romualdo, Luis et al. "In vitro stability of B-type natriuretic peptide (BNP) in plasma stored under different conditions when measured with the Lumipulse assay."Scandinavian journal of clinical and laboratory investigation, vol. 79,6 (2019): 455-458.
- Alpert JS, Thygesen K, Antman E, Bassand JP. Myocardial infarction redefined--a consensus document of The Joint European Society of Cardiology/American College of Cardiology Committee for the redefinition of myocardial infarction[J]. J Am Coll Cardiol, 2000, 36(3): 959-969.
- Cummins B, et al. Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction[J]. Am.Heart J, 1987, 113: 1333.
- M.R Cowie, P Jourdain, A Maisel, U Dahlstrom, F Follath, R Isnard, A Luchner, T McDonagh, J Mair, M Nieminen, G Francis, Clinical applications of B-type natriuretic peptide (BNP) testing, European Heart Journal, Volume 24, Issue 19, 1 October 2003, Pages 1710–1718
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Gobinet-Georges, A et al. "Stability of brain natriuretic peptide (BNP) in human whole blood and plasma."Clinical chemistry and laboratory medicine, vol. 38,6 (2000): 519-23.
- García de Guadiana-Romualdo, Luis et al. "In vitro stability of B-type natriuretic peptide (BNP) in plasma stored under different conditions when measured with the Lumipulse assay."Scandinavian journal of clinical and laboratory investigation, vol. 79,6 (2019): 455-458.
- Alpert JS, Thygesen K, Antman E, Bassand JP. Myocardial infarction redefined—a consensus document of The Joint European Society of Cardiology/American College of Cardiology Committee for the redefinition of myocardial infarction[J]. J Am Coll Cardiol, 2000, 36(3): 959-969.
- Cummins B, et al. Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction[J]. Am.Heart J, 1987, 113: 1333.
- M.R Cowie, P Jourdain, A Maisel, U Dahlstrom, F Follath, R Isnard, A Luchner, T McDonagh, J Mair, M Nieminen, G Francis, Clinical applications of B-type natriuretic peptide (BNP) testing, European Heart Journal, Volume 24, Issue 19, 1 October 2003, Pages 1710–1718
- QUIBASA: Datos do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Test for quantitative determination of the analyte BNP (type B natriuretic peptide or brain natriuretic peptide), in biological samples (plasma or whole blood) using a fluorescent immunoassay test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay.
The BIOCLIN POCT FIA BNP kit is a test designed to detect the analyte BNP in plasma or whole blood samples. When a diluted sample is added, it interacts with anti-BNP monoclonal antibodies conjugated to fluorescent microsphere particles. Then, the sample migrates through the membrane by capillarity, and interacts in the detection line (with the immobilized anti-BNP monoclonal antibodies present in the "T" line). The sample will continue to migrate through the membrane until it reaches the control region "C", where the conjugate containing IgY antibodies interacts with immobilized Anti-IgY antibodies. At the end of the incubation period, the cassette must be inserted into the reader for quantification. The BNP concentration of the sample directly correlates with the intensity of the fluorescence signal emitted by the test line of the assay.

The BIOCLIN POCT FIA Reader calculates the BNP concentration in the sample in pg/mL according to the assay calibration curve (contained on the identification chip).

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 8 °C. **Do not freeze.** Each sachet contains:
A) 01 plastic device with 01 test strip containing:
- Control region (C): Anti-IgY antibody.
- Test region (T): Anti-BNP monoclonal antibody.
- Conjugate: Anti-BNP monoclonal antibody conjugated to fluorescent microspheres and IgY antibody conjugated to fluorescent microspheres.
B) 01 sachet of silica.

Number 2 (R2) Diluent - Store between 2 and 8 °C. **Do not freeze.** Contains: Buffer, surfactant and preservative.

Number 3 (R3) - Identification Chip.

PRESENTATION

Reagents	Number of Reagents per Package
1 - Cassette	25 units
2 - Diluent	25 units (25 x 0,21 mL)
3 - Identification Chip	1 unit

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table
- Instructions for use (manual)

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 5 to 100 µL with a coefficient of variation of less than 1.5%.
- 2- BIOCLIN POCT FIA incubator (optional).
- 3- BIOCLIN POCT FIA reader.
- 4- Stopwatch or clock.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- **For *in vitro* diagnostic use only.**
- 2- **The BIOCLIN POCT FIA BNP must only be used in conjunction with the BIOCLIN POCT FIA Reader.**
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Open the aluminum envelope to remove the cassette, only when using it.
- 5- The kit must not be frozen or exposed to high temperatures (above 45°C), as this will cause it to deteriorate.
- 6- The reagents, as well as the samples, must reach room temperature before starting the test.
- 7- Do not mix reagents from kits with the same batch or with different batch numbers.
- 8- Do not use a kit after its expiration date. The expiration date is indicated on the product label.
- 9- The cassette (test device) and sample dilution tubes are disposable and cannot be reused.
- 10- The kit materials that had contact with the patient's sample must have a special destination. Pipettes, cassettes (test devices) and diluent bottles must be discarded in biological waste.
- 11- Secondary packaging (box), aluminum envelope and silica sachet can be discarded in general trash.
- 12- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 13- To obtain information related to biosafety or in the event of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through a request through the SAC (Service Customer Advice) from Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Plasma: The sample must be collected by venipuncture, using collection tubes with the anticoagulant EDTA. Plasma samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 4 hours. If samples cannot be analyzed within 4 hours, they can be stored for up to 1 month at -20°C.*

Whole blood: Whole blood samples collected by venipuncture with the appropriate anticoagulant (EDTA) can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 4 hours*. **Do not freeze.**

***At the time of use, samples must be stabilized at room temperature.**

PROCESS DESCRIPTION

Stability After Opening

The stability of the BIOCLIN POCT FIA BNP cassette (after opening the aluminized sachet) is 1 hour. This stability may vary according to the conditions of use of the product and the environment. However, to be safe, we recommend that the BIOCLIN POCT FIA BNP kit should be tested immediately after removing the cassette from the sachet.

TECHNIQUE

- 1- **The sample, Cassette and Sample Dilution Tube must be between 15-30 °C to perform the test.**
- 2- Check if the identification chip is inserted in the equipment.
- 3- Remove the cassette (test device) from the protective packaging, place it on a clean and level surface and identify it appropriately.
- 4- Transfer 50 µL of **Plasma** or 80 µL of **Whole Blood** to the tube containing the diluent (Reagent No. 2).
- 5- Homogenize the tube with the sample.
- 6- Pipette 80 µL of the diluted sample into the sample well.
- 7- Wait 10 minutes to read the essay. **Do not read before 10 minutes.**

Cassette Reading

- 1- Select the test mode on the BIOCLIN POCT FIA Reader. For more information on how to select the test mode, see the equipment manual.
- 2- Insert the Cassette into the BIOCLIN POCT FIA Reader and wait for the result to be issued.

NOTES

- 1- The cassette must be used immediately after opening its packaging.
- 2- The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

Calibration

The calibration of the method is carried out using the Identification Chip, which must be changed at each batch change of reagents.

PROCESS LIMITATIONS

- 1- It is important to use the correct sample volume, as lower or higher volumes may cause erroneous results.
- 2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretations of the results.
- 3- The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

INTERFERENCES

Samples with concentrations of up to 400 µmol/L Bilirubin, 20 mmol/L Triglycerides, 15 mg/mL Hemoglobin and 3250 IU/mL Rheumatoid Factor do not affect the assay results.

CROSS REACTIVITY

The BIOCLIN POCT FIA BNP kit was evaluated for cross-reactivity with different analytes. No falsely elevated results were observed for concentrations of up to 1000 ng/L of ANP (Atrial Natriuretic Peptide) and 40000 pg/mL of NT-proBNP. Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis must consider the patient's clinical data with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and accuracy of the dosages.

REFERENCE VALUES

Patients classified in the normal group had BNP concentrations below 100 pg/mL.
These values should be used as a guideline, and each laboratory should create its range of reference values, according to the population served. The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods and Methodological Specificity

The BIOCLIN POCT FIA BNP kit was compared with another commercially available method for measuring BNP. 100 analyzes were performed and the results were evaluated. The linear equation obtained was Y=0.97x +16.91 and the correlation coefficient was 0.99. With these results, it can be concluded that the kit presents good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 4 samples with different concentrations, obtaining the following results:

Repeatability	Sample			
	1	2	3	4
Concentration (pg/mL)	31.2	56.1	2090.3	3174.3
Standard Deviation (pg/mL)	2.5	5.3	65.9	157.9
Coefficient of Variation (%)	8	10	3	5

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 4 samples with different concentrations, obtaining the following results:

Reproducibility	Sample			
	1	2	3	4
Concentration (pg/mL)	20.0	60.6	1553.5	2926.8
Standard Deviation (pg/mL)	2.0	5.0	97.9	146.8
Coefficient of Variation (%)	10	8	6	5

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

SENSITIVITY

The analytical sensitivity of the BIOCLIN POCT BNP kit is 9.61pg/mL.

PROZONE EFFECT

No prozone effect was observed up to a concentration of 15000 pg/mL.

LINEARITY

The operating range for the BIOCLIN POCT BNP reagent is 10 pg/mL to 5000 pg/mL.

DIAGNOSTIC SIGNIFICANCE

Brain Natriuretic Peptide (BNP), also known as type B natriuretic peptide, is a hormone secreted by cardiomyocytes in the heart ventricles in response to stretch caused by increased ventricular blood volume. The 134 amino acid pro-BNP precursor protein is cleaved generating BNP and NT-proBNP. BNP concentration can be used as an aid in the diagnosis and assessment of the severity of patients with heart failure and acute coronary syndromes, since an increase in the analyte is observed in these conditions. Furthermore, an increase in BNP is observed in renal failure and primary hyperaldosteronism. Medications such as diuretics and ACE inhibitors can reduce BNP concentration.

BIBLIOGRAPHIC REFERENCES

1. Gobinet-Georges, A et al. "Stability of brain natriuretic peptide (BNP) in human whole blood and plasma."Clinical chemistry and laboratory medicine, vol. 38,6 (2000): 519-23.

2. García de Guadiana-Romualdo, Luis et al. "In vitro stability of B-type natriuretic peptide (BNP) in plasma stored under different conditions when measured with the Lumpulse assay."Scandinavian journal of clinical and laboratory investigation, vol. 79,6 (2019): 455-458.

3. Alpert JS, Thygesen K, Antman E, Bassand JP. Myocardial infarction redefined--a consensus document of The Joint European Society of Cardiology/American College of Cardiology Committee for the redefinition of myocardial infarction[J]. J Am Coll Cardiol, 2000, 36(3): 959-969.

4. Cummins B, et al. Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction[J]. Am.Heart J, 1987, 113: 1333.

5. M.R Cowie, P Jourdain, A Maisel, U Dahlstrom, F Follath, R Isnard, A Luchner, T McDonagh, J Mair, M Nieminen, G Francis, Clinical applications of B-type natriuretic peptide (BNP) testing, European Heart Journal, Volume 24, Issue 19, 1 October 2003, Pages 1710–1718

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 – Santa Branca
CEP 31565-130 – Belo Horizonte – MG – Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for kit BIOCLIN POCT FIA BNP: 10269360439

Review: October/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER