# Bioclin

## BIOCLIN POCT FIA CARDÍACO

REF K341

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#### FINALIDADE

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- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- Neumeuer, D. Tissue Specific and Subcellular Distribution of Creatine Kinase Isoenzymes. In: Lang, H., editor. Creatine Kinase Isoenzymes. New York: Springer-Verlag, p. 85-109, 1981.
- Ogunro, E. A.; Hearse, D.J.; Shillingford, J.P. Creatine Kinase Isoenzymes: Their Separation and Quantitation. Cardiovasc. Res., v. 11, p. 94-102, 1977.
- Tsung, S.H. Creatine Kinase Isoenzyme Patterns in Human Tissue Obtained at Surgery. Clin. Chem., v. 22, p. 173-175, 1976.
- Tate, J. et al. Standardization of cardiac troponin I measurement: past and present. Pathology, b. 42, n. 5, p. 402-408, 2010.
- Reichlin, T. et al. Early Diagnosis of Myocardial Infarction with Sensitive Cardiac Troponin Assays. The New England Journal of Medicine, v. 361, n. 9, p. 858-867, 2009.
- Filatov, V. L. et al. Troponin: Structure, Properties, and Mechanism of Functioning. Biochemistry, v. 64, n. 9, p. 969-985, 1999.
- Wittenberg, J. B.; Wittenberg, B. a. Myoglobin function reassessed. The Journal of Experimental Biology, v. 206, n. 12, p. 2011-2020, 2003.
- Gibler, W. B. et al. Myoglobin as an early indicator of acute myocardial infarction. Annals of Emergency Medicine, v. 16, n. 8, p. 851-856, 1987.
- Varki, A. P. et al. Serum myoglobin in acute myocardial infarction: a clinical study and review of the literature. American Heart Journal, v. 96, n. 5, p. 680-688, 1978.
- European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction Redefined—A Consensus Document of The Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Journal of the American

College of Cardiology, v. 36, n. 3, p. 959-969, 2000.  
**13-** Thygesen, K. et al. Fourth Universal Definition of Myocardial Infarction (2018). Circulation, v. 138, p. e618-e651, 2018.  
**14-** Bioclin – Dados de arquivos.

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- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev 2, 2002:31.
- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- Neumeuer, D., Tissue Specific and Subcellular Distribution of Creatine Kinase Isoenzymes. In: Lang, H., editor. Creatine Kinase Isoenzymes. New York: Springer-Verlag, p. 85-109, 1981.
- Ogunro, E. A.; Hearse, D.J.; Shillingford, J.P. Creatine Kinase Isoenzymes: Their Separation and Quantitation. Cardiovasc. Res., v. 11, p. 94-102, 1977.
- Tsung, S.H. Creatine Kinase Isoenzyme Patterns in Human Tissue Obtained at Surgery. Clin. Chem., v. 22, p. 173-175, 1976.
- Tate, J. et al. Standardization of cardiac troponin I measurement: past and present. Pathology, b. 42, n. 5, p. 402-408, 2010.
- Reichlin, T. et al. Early Diagnosis of Myocardial Infarction with Sensitive Cardiac Troponin Assays. The New England Journal of Medicine, v. 361, n. 9, p. 858-867, 2009.
- Filatov, V. L. et al. Troponin: Structure, Properties, and Mechanism of Functioning. Biochemistry, v. 64, n. 9, p. 969-985, 1999.
- Wittenberg, J. B.; Wittenberg, B. a. Myoglobin function reassessed. The Journal of Experimental Biology, v. 206, n. 12, p. 2011-2020, 2003.
- Gibler, W. B. et al. Myoglobin as an early indicator of acute myocardial infarction. Annals of Emergency Medicine, v. 16, n. 8, p. 851-856, 1987.
- Varki, A. P. et al. Serum myoglobin in acute myocardial infarction: a clinical study and review of the literature. American Heart Journal, v. 96, n. 5, p. 680-688, 1978.
- European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction Redefined—A Consensus Document of The Joint European Society of

Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Journal of the American College of Cardiology, v. 36, n. 3, p. 959-969, 2000.

13- Thygesen, K. et al. Fourth Universal Definition of Myocardial Infarction (2018). Circulation, v. 138, p. e618-e651, 2018.

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### INSTRUCTIONS FOR USE

#### FUNCTION

Test for the quantitative determination of cardiac troponin I (cTnI), creatine kinase MB (CK-MB), and myoglobin in human serum, plasma, and whole blood samples using fluorescent immunoassay. For *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Fluorescent immunoassay. The BIOCLIN POCT FIA CARDIAC kit is a fluorescent immunoassay based on the sandwich detection principle, designed for the quantitative determination of cardiac troponin I (cTnI), creatine kinase (CK-MB), and myoglobin (Myo) in human serum, plasma, or whole blood samples (collected in EDTA or heparin). When the sample is added to its application site, the interaction between the cTnI, CK-MB, and Myo molecules present in the sample and the anti-cTnI, anti-CK-MB, and anti-Myo antibodies conjugated to fluorescent microspheres occurs, respectively, forming the complexes: cTnI-antibody conjugate, CK-MB-antibody conjugate, and Myo-antibody conjugate. As the sample migrates by capillary action through the nitrocellulose membrane, these complexes are captured by their respective test lines, which consist of anti-cTnI, anti-CK-MB, and anti-Myo antibodies, separately. The sample continues to migrate through the membrane until it reaches the control line, which consists of anti-rabbit IgG antibodies (where rabbit IgG antibodies conjugated to fluorescent microspheres will be captured). At the end of the incubation period, the cassette is read on the immunofluorescence reader. The higher the concentration of analytes present in the sample, the more complexes are immobilized on the test lines and, consequently, the higher the fluorescence signals generated and captured by the reader. The concentrations of cardiac troponin I, creatine kinase MB fraction, and myoglobin present in the sample, all in ng/mL, are calculated according to the assay calibration curve (contained on the Identification Card).

#### REAGENTS

- 1- Cassette** - Store between 2 and 30 °C. Each sachet contains:
- A)** 01 plastic device with 01 test strip containing:
- Test line region 1: mouse anti-CK-MB antibody <2.0 mg/mL.
  - Test line region 2: goat anti-cardiac troponin I antibody <2.0 mg/mL.
  - Test line region 3: mouse anti-myoglobin antibody <2.0 mg/mL.
  - Control line region: goat anti-rabbit IgG antibody <2.0 mg/mL.
  - Conjugate: fluorescent microspheres conjugated to mouse anti-CK-MB antibodies <0.5 mg/mL, fluorescent microspheres conjugated to mouse anti-cardiac troponin I antibodies <0.5 mg/mL, fluorescent microspheres conjugated to mouse anti-myoglobin antibodies <0.5 mg/mL, and fluorescent microspheres conjugated to rabbit IgG antibodies <0.5 mg/mL.
- B)** 01 silica sachet.
- 2- Diluent for Whole Blood Samples** - Store between 2 and 30 °C. Contains: 10 mM phosphate-buffered saline (PBS), stabilizers (casein <1%) and preservatives (ProClin™ 300 <0.5%).
- 3- Identification Card**

#### PRESENTATION

| REAGENTS                                  | NUMBER OF REAGENTS PER PACKAGE |
|---|--------------------------------|
| <b>1- Cassette</b>                        | 25 units                       |
| <b>2- Diluent for Whole Blood Samples</b> | 25 units (25 x 0.05 mL)        |
| <b>3- Identification Card</b>             | 1 unit                         |

#### EQUIPMENTS AND OPERATIONAL INPUTS

##### Materials in the kit:

- Reagents described in the table above.

##### Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation lower than 1.5%.
- 2- BIOCLIN POCT FIA Immunofluorescence Reader.
- 3- Stopwatch or clock.
- 4- Quality Control Tools.

#### STORAGE AND TRANSPORTATION CONDITIONS

The storage and transportation temperature should be between 2 °C and 30 °C. Keep away from light and avoid humidity. **Do not freeze.**

#### SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- BIOCLIN POCT FIA CARDIAC should only be used in conjunction with the BIOCLIN POCT FIA Immunofluorescence Reader.
- 3- Strictly follow the proposed methodology to obtain accurate results.
- 4- The cassettes (test devices) do not contain contaminating substances, but since the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.
- 5- Open the aluminum envelope to remove the cassette only at the time of use.
- 6- The kit should not be frozen or exposed to high temperatures (above 45 °C), as this will cause deterioration.
- 7- The reagents, as well as the samples, must reach room temperature before starting the test.
- 8- Do not mix reagents from kits with the same batch or with different batch numbers.
- 9- Do not use the kit after its expiration date. The expiration date is indicated on the product label.
- 10- The test device (R1- Cassette) and sample dilution tubes (R2- Diluent for Whole Blood Samples) are disposable and cannot be reused.
- 11- The kit materials that have come into contact with the patient's sample must be disposed of in a special way. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biological waste.
- 12- Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular waste.
- 13- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.
- 14- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or by requesting it through Quibasa's SAC (Customer Advisory Service).
- 15- Do not use the product if the packaging is damaged.
- 16- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

#### SAMPLES

##### Serum or Plasma (EDTA or Heparin)

Samples should be collected by venipuncture, in tubes with or without anticoagulant (EDTA or Heparin). To avoid hemolysis,

centrifuge the samples to obtain serum/plasma and test them as quickly as possible. Do not keep samples at room temperature for long periods. Serum and plasma samples can be stored refrigerated, between 2 and 8 °C, for a maximum of 3 days. If samples cannot be analyzed within 3 days, they can be stored for up to 30 days at -20 °C. Serum and plasma samples stored for longer than the recommended periods and samples collected with anticoagulants other than those mentioned should not be used. Highly hemolyzed, lipemic, and precipitate-containing samples should be avoided, as they may provide inconsistent results.

##### Whole Blood (EDTA or Heparin)

Whole blood samples should be collected by venipuncture in tubes containing heparin or EDTA as an anticoagulant. Samples should be tested immediately or as soon as possible after collection. Whole blood collected in tubes containing anticoagulants can be stored refrigerated at 2–8 °C for a maximum of 2 days. Do not use samples collected with anticoagulants other than those mentioned, and do not freeze whole blood samples.

#### PROCESS DESCRIPTION

##### Stability After Opening

The stability test results show that the BIOCLIN POCT FIA CARDIAC kit is stable after opening for up to 70 minutes. This stability may vary depending on the conditions of use of the product and the environment. However, for safety reasons, we recommend that the product should be used immediately after removing the cassette from the sachet

##### TECHNIQUE

**Check the BIOCLIN POCT FIA Immunofluorescence Reader manual for complete instructions on handling the equipment.**

**ATTENTION!** The techniques for serum/plasma and whole blood are different. Carefully read the instructions for use described below, according to the type of sample to be used.

##### Serum and Plasma Samples

Before starting the assay, allow all reagents (Cassette) and samples to stabilize at room temperature (15-30°C) for at least 40 minutes.

- 1- Verify that the information contained in the Identification Card has been correctly entered into the equipment (with each new batch of the kit, the calibration and validity data contained in the Identification Card must be transferred to the equipment).
- 2- Remove the cassette (test device) from its protective packaging, place it on a clean, level surface, and label it appropriately.
- 3- Pipette 50 µL of sample (serum or plasma in EDTA or Heparin) directly into the sample application port of the cassette.
- 4- **Option 1:** Insert the cassette into the equipment immediately after applying the sample and select the “**Timing Test**” option. After 15 minutes, the equipment will take the reading, and the result will appear on the display.
- Option 2:** After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the “**Quick Test**” option. The equipment will take the reading immediately, and the result will appear on the display.

##### Whole Blood Samples

Before starting the test, allow all reagents (Cassette and Whole Blood Diluent) and samples to stabilize at room temperature (15-30°C) for at least 40 minutes.

- 1- Verify that the information contained in the Identification Card has been correctly entered into the equipment (with each new batch of the kit, the calibration and expiration data contained in the Identification Card must be transferred to the equipment).
- 2- Remove the cassette (test device) from its protective packaging, place it on a clean, level surface, and label it appropriately.
- 3- Ensure that the solution of the Diluent for Whole Blood Samples (Reagent No. 2) is at the bottom of the dilution tube.
- 4- Add 50 µL of sample (whole blood collected in EDTA or

Heparin) to a dilution tube containing the Diluent for Whole Blood Samples.

- 5- Close the tube cap and mix by inverting 6 to 8 times until the sample and diluent are completely mixed.
- 6- Pipette 50 µL of the diluted and properly mixed sample into the cassette's sample application port.
- 7- **Option 1:** Insert the cassette into the equipment immediately after applying the sample and select the “**Timing Test**” option. After 15 minutes, the equipment will take the reading and the result will appear on the display.
- Option 2:** After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the “**Quick Test**” option. The equipment will take the reading immediately, and the result will appear on the display.

#### RESULTS

The BIOCLIN POCT FIA Immunofluorescence Reader automatically calculates cardiac troponin I, CK-MB and myoglobin concentrations using the calibration curve contained in the R3 (Identification Card). The information on the Identification Card must be properly entered into the equipment for each new batch of reagents.

Quantitative results are measured in ng/mL. The kit's measurement range is between 0.1 and 50.0 ng/mL for cardiac troponin, between 2.0 and 100.0 ng/mL for CK-MB, and between 10.0 and 400.0 ng/mL for myoglobin (limits of quantification and maximum curve points, respectively). Results below the lower limit of detection are expressed as <0.1 ng/mL for cardiac troponin, <2.0 ng/mL for CK-MB, and <10.0 ng/mL for myoglobin. Results above the upper limit of detection are expressed as >50 ng/mL for cardiac troponin, >100.0 ng/mL for CK-MB, and >400.0 ng/mL for myoglobin.

#### PROCESS LIMITATIONS

- 1- It is important to use the correct sample volume, as lower or higher volumes may produce erroneous results.
- 2- The reaction reading time must be followed according to the established technique to avoid misinterpretation of the results.
- 3- Reactive results must be confirmed with another test and available clinical information. Individual test results cannot be used for a definitive diagnosis.
- 4- It is recommended that established guidelines for the diagnosis of acute myocardial infarction be followed.
- 5- The results provided by this kit must be interpreted by the responsible medical professional and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

#### REFERENCE VALUES

|                           |              |
|---------------------------|--------------|
| <b>Cardiac Troponin I</b> | < 0.30 ng/mL |
| <b>CK-MB</b>              | < 5.0 ng/mL  |
| <b>Myoglobin</b>          | < 60.0 ng/mL |

These values should be used as a guide, and each laboratory should create its own range of reference values according to the population served. The results provided by this kit should be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

#### INTERFERENTS

For all analytes (cardiac troponin I, CK-MB and myoglobin), no interference was observed by Bilirubin up to 350 µmol/L, Triglycerides up to 40 mmol/L, Hemoglobin up to 2 g/L, Rheumatoid Factor up to 1,500 IU/mL, Heterophile Antibodies up to 1,000 ng/mL and Antinuclear Antibodies up to 1:640.

#### CROSS-REACTIVITY

For the analyte cardiac troponin I, no false-positive results were

observed for cardiac troponin C up to 1,000 ng/mL, cardiac troponin T up to 500 ng/mL, and skeletal troponin T up to 1,000 ng/mL. For the analyte CK-MB, no false-positive results were observed for CK-MM up to 1,000 ng/mL and CK-BB up to 1,000 ng/mL. For the analyte myoglobin, no false-positive results were observed for hemoglobin up to 2 mg/mL, cardiac troponin I up to 1,000 ng/mL, and CK-MB up to 1,000 ng/mL. Despite the results found, the possibility of cross-reactivity with other analytes or with higher concentrations cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. For this purpose, it is recommended to use controls, which allow the assessment of the precision and accuracy of the dosages.

#### PRODUCT PERFORMANCE

#### ACCURACY

##### Comparison of Methods

The BIOCLIN POCT FIA CARDIAC kit was compared with another commercially available method for measuring each of the analytes. One hundred and twenty analyses were performed, and the results were evaluated. For cardiac troponin I, the linear equation obtained was  $Y = 0.9694x + 0.32$ , with a correlation coefficient of 0.9701. For CK-MB, the linear equation obtained was  $Y = 0.9696x + 0.3608$ , with a correlation coefficient of 0.9687. For myoglobin, the linear equation obtained was  $Y = 0.9867x + 0.0056$ , with a correlation coefficient of 0.9647. Based on these results, it can be concluded that the kit presents good methodological specificity.

#### PRECISION

##### Repeatability

Repeatability was calculated from 10 successive determinations, using 2 samples with different values, obtaining the following concentration results:

| Cardiac Troponin I           |        |        |
|------------------------------|--------|--------|
| REPEATABILITY                | Sample |        |
|                              | 1      | 2      |
| Mean                         | 0.298  | 10.031 |
| Standard Deviation           | 0.02   | 0.93   |
| Coefficient of Variation (%) | 8.2    | 9.2    |

| CK-MB                        |        |       |
|------------------------------|--------|-------|
| REPEATABILITY                | Sample |       |
|                              | 1      | 2     |
| Mean                         | 4.96   | 20.02 |
| Standard Deviation           | 0.31   | 1.70  |
| Coefficient of Variation (%) | 6.2    | 8.5   |

| Myoglobin                    |        |        |
|------------------------------|--------|--------|
| REPEATABILITY                | Sample |        |
|                              | 1      | 2      |
| Mean                         | 50.88  | 201.25 |
| Standard Deviation           | 3.70   | 11.35  |
| Coefficient of Variation (%) | 7.3    | 5.6    |

#### Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 samples with different values, obtaining the following concentration results:

| Cardiac Troponin I           |        |      |
|------------------------------|--------|------|
| REPRODUCIBILITY              | Sample |      |
|                              | 1      | 2    |
| Mean                         | 0.30   | 9.98 |
| Standard Deviation           | 0.0    | 0.9  |
| Coefficient of Variation (%) | 7.7    | 9.3  |

| CK-MB                        |        |       |
|------------------------------|--------|-------|
| REPRODUCIBILITY              | Sample |       |
|                              | 1      | 2     |
| Mean                         | 4.91   | 20.48 |
| Standard Deviation           | 0.3    | 1.4   |
| Coefficient of Variation (%) | 6.5    | 7.0   |

| Myoglobin                    |        |        |
|------------------------------|--------|--------|
| REPRODUCIBILITY              | Sample |        |
|                              | 1      | 2      |
| Mean                         | 49.85  | 200.09 |
| Standard Deviation           | 3.9    | 13.1   |
| Coefficient of Variation (%) | 7.8    | 6.5    |

#### ANALYTICAL SENSITIVITY

For cardiac troponin I, the Limit of Blank (LoB) is 0.025 ng/mL, the Limit of Detection (LoD) is 0.07 ng/mL, and the Limit of Quantification (LoQ) is 0.1 ng/mL. For CK-MB, the Limit of Blank (LoB) is 0.5 ng/mL, the Limit of Detection (LoD) is 1.0 ng/mL, and the Limit of Quantification (LoQ) is 2.0 ng/mL. For myoglobin, the Limit of Blank (LoB) is 2.4 ng/mL, the Limit of Detection (LoD) is 5.0 ng/mL, and the Limit of Quantification (LoQ) is 10.0 ng/mL.

#### LINEARITY

The BIOCLIN POCT FIA CARDIAC kit is linear up to 50 ng/mL for cardiac troponin I, up to 100 ng/mL for CK-MB and up to 400 ng/mL for myoglobin (highest points of the calibration curve).

#### HIGH DOSE PROZONE EFFECT

No high-dose prozone effect was observed at up to 500 ng/mL cardiac troponin I, 800 ng/mL CK-MB, and 10,000 ng/mL myoglobin.

#### DIAGNOSTIC SIGNIFICANCE

AMI is a major cause of death and disability, and rapid diagnosis is critical for initiating effective treatment and improving prognosis. It is characterized by characteristic symptoms (such as chest pain and dyspnea), elevation of specific blood markers, and the development of a Q-wave pattern on the electrocardiogram. Acute myocardial infarction occurs due to necrosis of cardiac tissue associated with prolonged ischemia. This myocardial necrosis results in the release of several markers into the bloodstream, such as myoglobin, cardiac troponins T and I, creatine kinase, and lactate dehydrogenase. Troponins are proteins that form protein complexes present in the muscle fibers of skeletal and cardiac muscles and play important roles in the muscle contraction of these tissues. Troponins T and I have isoforms specific to cardiac muscle fibers, which are important for use as specific markers of cardiac injury, especially in the diagnosis of acute myocardial infarction (AMI). Cardiac troponin I is an important marker because it presents an elevation in blood concentration 3 to 6 hours after the infarction, peaking between 12 and 24 hours and maintaining altered levels for 5 to 14 days after the AMI episode. Although useful in the diagnosis of AMI, troponin levels increase in cases of myocardial injury without ischemia, such as myocarditis, or secondary to other diseases, such as renal failure. Creatine kinase (CK) is an enzyme present in various tissues and cell types of the human body. This

enzyme is responsible for catalyzing the conversion of creatine to phosphocreatine in different tissues. Its molecular structure consists of two subunits (M and B) that combine to produce three distinct isoenzymes: CK-MB (found in cardiac muscle), CK-BB (present in the brain and lung), and CK-MM (found in skeletal muscle). When cellular damage occurs, CK-MB concentrations increase in the blood, and as a cardiac-specific isoenzyme, CK-MB can be used to aid in the diagnosis of myocardial diseases such as acute myocardial infarction (AMI) and myocarditis. After an infarction, CK-MB concentrations increase between 4 and 6 hours after the AMI, peaking around 18 hours and returning to normal after 48 hours. In addition to these two main conditions, CK-MB levels can also be altered in other clinical conditions such as rhabdomyolysis and stroke. Myoglobin is a cytoplasmic protein of approximately 17 kDa, responsible for storing and transporting oxygen from the sarcolemma to the mitochondria in muscle cells. Like troponin and CK-MB, serum myoglobin levels increase after muscle injury. Myoglobin is the cardiac marker that changes most rapidly after an MI, exceeding the normal range one hour after the infarction, peaking 4 to 8 hours after the event, and returning to normal after approximately 24 hours. In addition to myocardial infarction, serum myoglobin concentrations may also be altered after muscle trauma, strenuous physical activity, rhabdomyolysis, kidney injury, and various myopathies.

#### BIBLIOGRAPHIC REFERENCES

- 1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 2- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- 3- Neumeuer, D. Tissue Specific and Subcellular Distribution of Creatine Kinase Isoenzymes. In: Lang, H., editor. Creatine Kinase Isoenzymes. New York: Springer-Verlag, p. 85-109, 1981.
- 4- Ogunro, E. A.; Hearse, D.J.; Shillingford, J.P. Creatine Kinase Isoenzymes: Their Separation and Quantitation. Cardiovasc. Res., v. 11, p. 94-102, 1977.
- 5- Tsung, S.H. Creatine Kinase Isoenzyme Patterns in Human Tissue Obtained at Surgery. Clin. Chem., v. 22, p. 173-175, 1976.
- 6- Tate, J. et al. Standardization of cardiac troponin I measurement: past and present. Pathology, b. 42, n. 5, p. 402-408, 2010.
- 7- Reichlin, T. et al. Early Diagnosis of Myocardial Infarction with Sensitive Cardiac Troponin Assays. The New England Journal of Medicine, v. 361, n. 9, p. 858-867, 2009.
- 8- Filatov, V. L. et al. Troponin: Structure, Properties, and Mechanism of Functioning. Biochemistry, v. 64, n. 9, p. 969-985, 1999.
- 9- Wittenberg, J. B.; Wittenberg, B. a. Myoglobin function reassessed. The Journal of Experimental Biology, v. 206, n. 12, p. 2011-2020, 2003.
- 10- Gibler, W. B. et al. Myoglobin as an early indicator of acute myocardial infarction. Annals of Emergency Medicine, v. 16, n. 8, p. 851-856, 1987.
- 11- Varki, A. P. et al. Serum myoglobin in acute myocardial infarction: a clinical study and review of the literature. American Heart Journal, v. 96, n. 5, p. 680-688, 1978.
- 12- European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction Redefined—A Consensus Document of The Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Journal of the American College of Cardiology, v. 36, n. 3, p. 959-969, 2000.
- 13- Thygesen, K. et al. Fourth Universal Definition of Myocardial Infarction (2018). Circulation, v. 138, p. e618-e651, 2018.
- 14- Bioclin – Archive data

#### QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration number for BIOCLIN POCT FIA  
CARDIAC kit: 10269360480

Review: October/2025

### UNIVERSAL SYMBOLOGY

|   |  |   |  |
|---|--|---|--|
|    | CATALOG NUMBER                           |    | MADE BY                                |
|    | LOT NUMBER                               |    | CONTROL                                |
|    | MANUFACTURING DATE                       |    | POSITIVE CONTROL                       |
|  | VALIDITY DATE<br>(last day of the month) |  | NEGATIVE CONTROL                       |
|  | TEMPERATURE LIMIT<br>(store)             |  | BIOLOGICAL RISK                        |
|  | CONTENT IS SUFFICIENT<br>FOR <N> TEST    |  | FLAMMABLE                              |
|  | SEE INSTRUCTIONS<br>FOR USE              |  | CORROSIVE                              |
|  | IN VITRO DIAGNOSTIC<br>PRODUCT           |  | TOXIC                                  |
|  | KEEP AWAY<br>FROM SUNLIGHT               |  | DO NOT USE IF<br>PACKAGE IS<br>DAMAGED |
|  | DO NOT REUSE                             |  | PRODUCT<br>STERILIZED                  |
|  | CAUTION                                  |  | DANGER                                 |