BIOCLIN POCT FIA CARDÍACO

REF K289

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FINALIDADE

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1. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
2. Albert JS, Thygesen K, Antman E, Bassand JP. Myocarddial infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. J Am Coll Cardiol. 2000 Sep;36(3):959-69
3. Cummins B, Auckland ML, Cummins P .Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction[J]. Journal of American Heart. 1987, 113:1333-1344.
4. Adams JE,et al.Diagnosis of Perioperative Myocardial Infarction with Measurement of Cardiac Troponin I[J].N.Eng.J.Med,1994,330:670-674.
5. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.
6. Levinson SS. The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Test for quantitative determination of the analytes Troponin I, CK-MB and Myoglobin in biological samples (plasma or whole blood) using a fluorescent immunoassay test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay.
The BIOCLIN POCT FIA CARDIAC kit is a test designed to detect the analyte Troponin I, CK-MB and Myoglobin in plasma or whole blood samples. This assay contains 3 distinct detection lines ("T" line).
When the sample is added, it interacts with the conjugate containing anti-cTnI, anti-CK-MB and anti-Myoglobin monoclonal antibodies conjugated to particles with fluorescent microspheres. The interaction between sample and conjugated antibodies forms 3 complexes that migrate through the membrane by capillarity to the detection lines. The first line (T1 line) is coated with the anti-cTnI monoclonal antibody. The second detection line (T2 line) is coated with the anti-CK-MB monoclonal antibody. While the third detection line (T3 line) is coated with the anti-Myoglobin monoclonal antibody. The complex formed between sample and anti-cTnI antibodies will be captured by the T1 line. The complex formed between sample and anti-CK-MB antibodies will be captured by the T2 line. Meanwhile, the complex formed by the sample and anti-Myoglobin antibodies will be captured by the T3 line. The sample will continue to migrate through the membrane until it reaches the control region "C", where the conjugate containing IgG antibodies interacts with immobilized Anti-IgG antibodies. At the end of the incubation period, the cassette must be inserted into the reader for quantification. The concentration of Troponin I, CK-MB and Myoglobin in the sample directly correlates with the intensity of the fluorescence signal emitted by the test lines of the assay (Line T1, T2 and T3).

The BIOCLIN POCT FIA Reader calculates the concentration of Troponin I, CK-MB and Myoglobin in the sample in ng/mL according to the assay calibration curve (contained on the Identification Chip).

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 8°C. **Do not freeze.** Each sachet contains:
A) 01 plastic device with 01 test strip containing:
- Control region (C): Anti-IgG antibody.
- Test region (T): Anti-cTnI monoclonal antibody, anti-CK-MB monoclonal antibody and anti-Myoglobin monoclonal antibody.
-Conjugate: Anti-cTnI monoclonal antibody conjugated to fluorescent microspheres, Anti-CK-MB monoclonal antibody conjugated to fluorescent microspheres, Anti-Myoglobin monoclonal antibody conjugated to fluorescent microspheres and IgG antibody conjugated to fluorescent microspheres.
B) 01 sachet of silica.

Number 2 (R2) Diluent - Store between 2 and 8°C. **Do not freeze.** Contains: Buffer, surfactant and preservative.

Number 3 (R3) Identification Chip.

PRESENTATION

Reagents	Number of Reagents per Package
1 - Cassette	25 units
2 - Diluent	25 units (25 x 0.21 mL)
3 - Identification Chip	1 unit

EQUIPMENT AND OPERATIONAL INPUTS
Materials contained in the kit:

- Reagents described in the previous table
- Instructions for use (manual)

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 5 to 100 µL with a coefficient of variation of less than 1.5%.
- 2- BIOCLIN POCT FIA incubator (optional).
- 3- BIOCLIN POCT FIA reader.
- 4- Stopwatch or clock.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- **For *in vitro* diagnostic use only.**
- 2- **The BIOCLIN POCT FIA CARDIAC must only be used in conjunction with the BIOCLIN POCT FIA Reader.**
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Open the aluminum envelope to remove the cassette, only when using it.
- 5- The kit must not be frozen or exposed to high temperatures (above 45°C), as this will cause it to deteriorate.
- 6- The reagents, as well as the samples, must reach room temperature before starting the test.
- 7- Do not mix reagents from kits with the same batch or with different batch numbers.
- 8- Do not use a kit after its expiration date. The expiration date is indicated on the product label.
- 9- The cassette (test device) and sample dilution tubes are disposable and cannot be reused.
- 10- The kit materials that had contact with the patient's sample must have a special destination. Pipettes, cassettes (test devices) and diluent bottles must be discarded in biological waste.
- 11- Secondary packaging (box), aluminum envelope and silica sachet can be discarded in general trash.
- 12- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 13- To obtain information related to biosafety or in the event of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through a request through the SAC (Service Customer Advice) from Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Plasma: Samples must be collected by venipuncture, using the anticoagulant EDTA or Heparin. Plasma samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 3 days. If samples cannot be analyzed within 3 days, they can be stored for up to 1 month at -20°C.*

Whole blood: Whole blood samples must be collected by venipuncture with the appropriate anticoagulant (EDTA or Heparin) and can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 3 days*. **Do not freeze.**

*** When used, samples must be stabilized at room temperature.**

PROCESS DESCRIPTION

Stability After Opening

The stability of the BIOCLIN POCT FIA CARDIAC cassette (after opening the aluminized sachet) is 1 hour. This stability may vary according to the conditions of use of the product and the environment. However, for safety, we recommend that the BIOCLIN POCT FIA CARDIAC kit should be tested immediately after removing the cassette from the sachet.

TECHNIQUE

- 1- **The Sample, Cassette and Sample Dilution Tube must be between 15-30 °C to perform the test.**
 - 2- Check if the Identification Chip is inserted in the equipment.
 - 3- Remove the cassette (test device) from the protective packaging, place it on a clean and level surface and identify it appropriately.
 - 4- **For plasma:** pipette 50µL of the sample into the sample dilution tube (Reagent No. 2).
 - 5- **For whole blood:** pipette 80µL of the sample into the sample dilution tube (Reagent No. 2).
 - 6- Homogenize the tube with the sample.
 - 7- Pipette 80 µL of the diluted sample into the sample well.
- Note:** Steps 4 to 7 must be performed within 1 minute to ensure assay accuracy.
- 8- Wait 15 minutes to read the essay. **Do not read before 15 minutes.**

Cassette Reading

- 1- Select the test mode on the BIOCLIN POCT FIA Reader. For more information on how to select the test mode, see the equipment manual.
- 2- Insert the Cassette into the BIOCLIN POCT FIA Reader and wait for the result to be issued.

NOTES

- 1- The cassette must be used immediately after opening its packaging.
- 2- The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

Calibration

The calibration of the method is carried out using the Identification Chip, which must be changed at each batch change of reagents.

PROCESS LIMITATIONS

- 1- It is important to use the correct sample volume, as lower or higher volumes may cause erroneous results.
- 2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretations of the results.
- 3- The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

INTERFERENCES

Samples with concentrations of up to 400µmol/L Bilirubin, 20mmol/L Triglycerides, 15mg/mL Hemoglobin and 3250 IU/mL Rheumatoid Factor do not affect the assay results.

CROSS REACTIVITY

The BIOCLIN POCT FIA CARDIAC kit was evaluated for cross-reactivity with different analytes.

For Troponin I, no falsely elevated results were observed for concentrations of up to 1000µg/L of skeletal human Troponin (sTnI), 1000µg/L of cardiac human Troponin C (cTnC), 1000µg/L of cardiac human Troponin T (cTnT), 1000µg /L of Myoglobin and 1000µg/L of CK-MB.

For CK-MB, no falsely elevated results were observed for concentrations of up to 10000ng/mL of Creatine kinase fraction MM(CK-MM) and 10000ng/mL of Creatine kinase fraction BB(CK-BB).

For Myoglobin, no falsely elevated results were observed for concentrations of up to 30 ng/mL of cardiac human Troponin I (cTnI), 30 ng/mL of cardiac human Troponin T (cTnT), 60 ng/mL of Creatine kinase MB fraction (CK -MB), 0.3 ng/mL of Aspirin and 200 ng/mL of Digoxin.

Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis must consider the patient's clinical data with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and accuracy of the dosages.

REFERENCE VALUES

Analyte	Reference value (ng/mL)
Troponin I	< 0.3
CK-MB	< 5.0
Myoglobin	< 55.0

These values should be used as a guideline, and each laboratory should create its range of reference values, according to the population served. The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods and Methodological Specificity

The BIOCLIN POCT FIA CARDIAC kit was compared with another commercially available method for measuring Troponin I, CK-MB and Myoglobin. 100 analyzes were performed and the results were evaluated. The linear equation obtained for Troponin I was Y=0.969x+0.008 with a correlation coefficient of 0.99, for CK-MB it was Y=1.026x+0.08 with a correlation coefficient of 0.99 and for Myoglobin Y=0.979x +2.014 with a correlation coefficient of 0.97. With these results, it can be concluded that the kit presents good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 6 samples with different concentrations, obtaining the following results:

Repeatability	Troponin I		CK-MB		Myoglobin	
	Sample		Sample		Sample	
	1	2	3	4	5	6
Concentration (ng/mL)	0.50	39.71	2.51	60.18	24.91	83.60
Standard Deviation (ng/mL)	0.04	2.89	0.19	3.17	1.97	2.63
Coefficient of Variation (%)	9%	7%	8%	5%	8%	3%

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 6 samples with different concentrations, obtaining the following results:

Reproducibility	Troponin I		CK-MB		Myoglobin	
	Sample		Sample		Sample	
	1	2	3	4	5	6
Concentration (ng/mL)	0.49	40.01	2.53	59.74	24.71	78.95
Standard Deviation (ng/mL)	0.04	3.04	0.20	2.88	1.57	3.18
Coefficient of Variation (%)	9%	8%	8%	5%	6%	4%

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

SENSITIVITY

The analytical sensitivity of the BIOCLIN POCT FIA CARDIAC kit is 0.073 ng/mL for Troponin I, 0.83 ng/mL for CK-MB and 2.18 ng/mL for myoglobin.

PROZONE EFFECT

For Troponin I, no effect was observed after the high dose zone up to a concentration of 200.00 ng/mL.

For CK-MB, no post-high dose zone effect was observed up to a concentration of 150.00 ng/mL.

For Myoglobin, no effect was observed after the high dose zone up to a concentration of 800.00 ng/mL.

LINEARITY

The operating range for the BIOCLIN POCT FIA CARDIAC reagent is 0.10 ng/mL to 50.00 ng/mL for Troponin I, 2.50 ng/mL to 80.00 ng/mL for CK-MB, and 20.00 ng/mL up to 400.00 ng/mL for Myoglobin.

DIAGNOSTIC SIGNIFICANCE

This product is used to determine the concentration of cardiac troponin I, CK-MB and myoglobin in plasma or whole blood samples. Cardiac troponin I (cTnI) has a molecular weight of 22.5 KD and combines with troponin T (TnT) and troponin C (TnC) to form the troponin complex that regulates the interaction between actin and myosin. Normally, the level of cTnI is low in the blood. However, in acute myocardial infarction (AMI), cTnI is rapidly released into the blood, increasing its plasma levels approximately 4-6 h after the onset of AMI and remaining elevated for 6-10 days. Creatine kinase (CK) is an enzyme with two M and/or B subunits present in different tissues of the human body and cell types. This enzyme is responsible for catalyzing the conversion of creatine into phosphocreatine in different tissues. Its molecular structure consists of two subunits that combine to give rise to three different isoenzymes CK-MB (found in cardiac muscle), CK-BB (present in the brain and lung) and CK-MM (found in skeletal muscle). Cellular injuries increase the concentration of CK in the blood. As a heart-specific isoenzyme, CK-MB can be used to aid in the diagnosis of myocardial diseases such as heart attack and myocarditis. Following the infarction (3 to 5 hours), the serum concentration of this enzyme rises, reaches a maximum peak after 12 - 24 hours and remains elevated for a short period (2 to 3 days). As it is a specific fraction for the heart muscle, this isoenzyme is of great importance for monitoring myocardial infarction. Myoglobin is a low molecular weight cytoplasmic protein, myoglobin levels in the blood increase in case of damage to muscle cells. Myoglobin concentration in human serum exceeds the normal range one hour after myocardial infarction reaching its peak 4-8 hours after the onset of myocardial infarction.

BIBLIOGRAPHIC REFERENCES

1. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

2. Albert JS, Thygesen K, Antman E, Bassand JP. Myocarddial infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. J Am Coll Cardiol. 2000 Sep;36(3):959-69

3. Cummins B, Auckland ML, Cummins P .Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction[J]. Journal of American Heart. 1987, 113:1333-1344.

4. Adams JE,et al.Diagnosis of Perioperative Myocardial Infarction with Measurement of Cardiac Troponin I[J].N.Eng.J.Med,1994,330:670-674.

5. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.

6. Levinson SS. The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 – Santa Branca
CEP 31565-130 – Belo Horizonte – MG – Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for kit BIOCLIN POCT FIA CARDIAC: 10269360437

Review: October/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER