BIOCLIN POCT FIA CK-MB

REF K281

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Test for quantitative determination of Creatine Kinase fraction MB (CK-MB), in biological samples (serum or plasma) using a fluorescent immunoassay test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay.
The BIOCLIN POCT FIA CK-MB kit is a test designed to detect CK-MB in human plasma and serum samples using a fluorescent immunoassay. When the sample and diluent are added, the Creatine Kinase MB fraction (CK-MB) of the sample interacts with specific anti-CK-MB monoclonal antibodies conjugated to fluorescent microspheres. The sample then migrates through the membrane by capillarity, up to the detection line (with anti-CK-MB monoclonal antibodies present in the "T" line). The sample will continue to migrate through the membrane until it reaches the control region "C", where the IgG conjugate conjugated with fluorescent microspheres interacts with Anti-IgG antibodies present in the control line. At the end of the incubation period, the cassette must be inserted into the reader for quantification. The CK-MB concentration of the sample correlates with the intensity of the fluorescence signal emitted by the test line of the assay.

The BIOCLIN POCT FIA Reader calculates the CK-MB concentration in the sample in ng/mL according to the assay calibration curve (contained on the identification chip).

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 8 °C. **Do not freeze.** Each sachet contains:
-Control region (C): Anti-IgG antibody.
-Test region (T): Anti-CK-MB specific monoclonal antibody.
-Conjugate: Anti-CK-MB monoclonal antibody conjugated to fluorescent microspheres and IgG antibodies conjugated to fluorescent microspheres.
B) 01 sachet of silica.

Number 2 (R2) Diluent - Store between 2 and 8 °C. **Do not freeze.** Contains: Buffer, surfactant and preservative.

Number 3 (R3) Identification Chip.

PRESENTATION

Reagents	Number of Reagents per Package
1 - Cassette	25 units
2 - Diluent	25 units (25 x 0.16 mL)
3 - Identification Chip	1 unit

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table
- Instructions for use (manual)

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 5 to 100 µL with a coefficient of variation of less than 1.5%.
- 2- BIOCLIN POCT FIA incubator (optional).
- 3- BIOCLIN POCT FIA reader.
- 4- Stopwatch or clock.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- *For *in vitro* diagnostic use only.*
- 2- **The BIOCLIN POCT FIA CK-MB must only be used in conjunction with the BIOCLIN POCT FIA Reader.**
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Open the aluminum envelope to remove the cassette, only when using it.
- 5- The kit must not be frozen or exposed to high temperatures (above 45°C), as this will cause it to deteriorate.
- 6- The reagents, as well as the samples, must reach room temperature before starting the test.
- 7- Do not mix reagents from kits with the same batch or with different batch numbers.
- 8- Do not use a kit after its expiration date. The expiration date is indicated on the product label.
- 9- The cassette (test device) and sample dilution tubes are disposable and cannot be reused.
- 10- The kit materials that had contact with the patient's sample must have a special destination. Pipettes, cassettes (test devices) and diluent bottles must be discarded in biological waste.
- 11- Secondary packaging (box), aluminum envelope and silica sachet can be discarded in general trash.
- 12- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 13- To obtain information related to biosafety or in the event of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through a request through the SAC (Service Customer Advice) from Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum or Plasma: Samples must be collected by venipuncture and can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 7 days. If samples cannot be analyzed within this period, they can be stored for up to 30 days at -20 °C (freezer)*. To obtain plasma, collect it using the anticoagulant heparin and separate it as soon as possible to avoid hemolysis.

***When used, samples must be stabilized at room temperature.**

PROCESS DESCRIPTION

Stability After Opening

The stability of the BIOCLIN POCT FIA CK-MB cassette (after opening the aluminized sachet) is 1 hour. This stability may vary according to the conditions of use of the product and the environment. However, to be safe, we recommend that the BIOCLIN POCT FIA CK-MB kit should be tested immediately after removing the cassette from the sachet.

TECHNIQUE

- 1- **The sample, Cassette and Sample Dilution Tube must be between 15-30 °C to perform the test.**
 - 2- Check that the identification chip is inserted in the BIOCLIN POCT FIA Reader.
 - 3- Pipette 40µL of the sample (serum or plasma) into the sample dilution tube (Reagent No. 2).
 - 4- Homogenize the tube with the sample.
 - 5- Remove the cassette (test device) from the protective packaging, place it on a clean and level surface and identify it appropriately.
 - 6- Pipette **80 µL of the diluted sample** into the sample well.
- Note:** Steps 3 to 6 must be performed within 1 minute to ensure assay accuracy.
- 7- Make sure there are no air bubbles.
 - 8- Wait 15 minutes to read the essay. **Do not read before 15 minutes.**

Cassette Reading

- 1- Select the test mode on the BIOCLIN POCT FIA Reader. For more information on how to select the test mode, see the equipment manual.
- 2- Insert the Cassette into the BIOCLIN POCT FIA Reader and wait for the result to be issued.

NOTES

- 1- The cassette must be used immediately after opening its packaging.
- 2- The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

Calibration

The calibration of the method is carried out using the Identification Chip, which must be changed at each batch change of reagents.

PROCESS LIMITATIONS

- 1- It is important to use the correct sample volume, as lower or higher volumes may cause erroneous results.
- 2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretations of the results.
- 3- The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

INTERFERENCES

Samples with concentrations of up to 400µmol/L Bilirubin, 20mmol/L Triglycerides, 15mg/mL Hemoglobin and 3250 IU/mL Rheumatoid Factor do not affect the assay results.

CROSS REACTIVITY

The BIOCLIN POCT FIA CK-MB kit was evaluated for cross-reactivity with different analytes. No falsely elevated results were observed for concentrations of up to 10000ng/mL of CK-MM (creatine kinase isoenzyme fraction MM) and 10000ng/mL of CK-BB (creatine kinase isoenzyme fraction BB).

Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis must consider the patient's clinical data with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and accuracy of the dosages.

REFERENCE VALUES

Patients classified in the normal group had CK-MB concentrations below 5 ng/mL.

These values should be used as a guideline, and each laboratory should create its range of reference values, according to the population served. The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods and Methodological Specificity

The BIOCLIN POCT FIA CK-MB kit was compared with another commercially available method for measuring CK-MB. 90 analyzes were performed and the results were evaluated. The linear equation obtained was Y=0.96x+0.10 and the correlation coefficient was 0.993. With these results, it can be concluded that the kit presents good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 4 samples with different concentrations, obtaining the following results:

Repeatability	Sample			
	1	2	3	4
Concentration (ng/mL)	1.70	5.17	50.66	67.37
Standard Deviation (ng/mL)	0.2	0.46	3.94	3.71
Coefficient of Variation (%)	12	9	8	6

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 4 samples with different concentrations, obtaining the following results:

Reproducibility	Sample			
	1	2	3	4
Concentration (ng/mL)	1.53	4.99	47.25	58.66
Standard Deviation (ng/mL)	0.17	0.43	2.48	3.39
Coefficient of Variation (%)	11	9	5	6

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

SENSITIVITY

The analytical sensitivity of the BIOCLIN POCT FIA CK-MB kit is 0.86 ng/mL.

PROZONE EFFECT

No prozone effect was observed up to a concentration of 150 ng/mL.

LINEARITY

The operating range for BIOCLIN POCT FIA CK-MB reagent is 1 ng/mL to 80 ng/mL.

DIAGNOSTIC SIGNIFICANCE

This product is used to detect the concentration of Creatine Kinase MB (CK-MB) in human serum or plasma. Creatine kinase (CK) is an enzyme with two M and/or B subunits present in different tissues of the human body and cell types. This enzyme is responsible for catalyzing the conversion of creatine into phosphocreatine in different tissues. Its molecular structure consists of two subunits that combine to give rise to three different isoenzymes CK-MB (found in cardiac muscle), CK-BB (present in the brain and lung) and CK-MM (found in skeletal muscle). Cellular injuries increase the concentration of CK in the blood. As a heart-specific isoenzyme, CK-MB can be used to aid in the diagnosis of myocardial diseases such as heart attack and myocarditis. Following the infarction (3 to 5 hours), the serum concentration of this enzyme rises, reaches a maximum peak after 12 - 24 hours and remains elevated for a short period (2 to 3 days). As it is a specific fraction for the heart muscle, this isoenzyme is of great importance for monitoring myocardial infarction.

BIBLIOGRAPHIC REFERENCES

1. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
2. Adams JE, et al. Diagnosis of Perioperative Myocardial Infarction with Measurement of Cardiac Troponin I [J]. N. Eng. J.Med, 1994, 330:670-674
3. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.
4. Levinson SS. The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 – Santa Branca
CEP 31565-130 – Belo Horizonte – MG – Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for kit BIOCLIN POCT FIA CK-MB: 10269360438

Review: October/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER