Bioclin

BIOCLIN POCT FIA D-DÍMERO

REF K332

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FINALIDADE

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1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

2- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.

3- RAMZI, D.W.; Leeper, K. V. DVT and pulmonary embolism: part I. Diagnosis. American Family Physician, v. 69, n. 12, p. 2829-2836, 2004.

4- Perrier, A. et al. D-dimer testing for suspected pulmonary embolism in outpatients. American Journal of Respiratory and Critical Care Medicine, v. 156, n. 2, p. 492-496, 1997.

5- Sakamoto, K. et al. D-Dimer is helpful for differentiating acute aortic dissection and acute pulmonary embolism from acute myocardial infarction. Hellenic J. Cardiol., v. 52, n. 2, p. 123-127, 2011.

6- Yasuka, S. Kubota, S. The value of blood D-Dimer test in the diagnosis of walk-in patients with venous thromboembolism. Vas. Health Risk Manag., v. 7, p. 125-127, 2011.

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- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- RAMZI, D.W.; Leeper, K. V. DVT and pulmonary embolism: part I. Diagnosis. American Family Physician, v. 69, n. 12, p. 2829-2836, 2004.
- Perrier, A. et al. D-dimer testing for suspected pulmonary embolism in outpatients. American Journal of Respiratory and Critical Care Medicine, v. 156, n. 2, p. 492-496, 1997.
- Sakamoto, K. et al. D-Dimer is helpful for differentiating acute aortic dissection and acute pulmonary embolism from acute myocardial infarction. Hellenic J. Cardiol., v. 52, n. 2, p. 123-127, 2011.
- Yasuka, S. Kubota, S. The value of blood D-Dimer test in the diagnosis of walk-in patients with venous thromboembolism. Vas. Health Rish Manag., v. 7, p. 125-127, 2011.
- Bioclin – Datos de archivos.

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INSTRUCTIONS FOR USE

FUNCTION

Test for the quantitative determination of D-dimer in human biological samples of plasma and whole blood through a fluorescence immunoassay. It is for *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay.

The BIOCLIN POCT FIA D-Dimer kit is a fluorescence immunoassay based on the "sandwich" detection principle, designed for the quantitative determination of D-dimer in human plasma or whole blood samples (collected with heparin, EDTA, or citrate). When the sample is added to the application site, the D-dimer molecules present interact with anti-D-dimer antibodies conjugated to fluorescent microspheres, forming a D-dimer-antibody conjugate complex. As the sample migrates by capillarity through the nitrocellulose membrane, this complex is captured by the test line, which also contains anti-D-dimer antibodies. The sample continues migrating until reaching the control line, consisting of anti-rabbit IgG antibodies, where rabbit IgG antibodies conjugated to fluorescent microspheres are captured. At the end of the incubation period, the cassette is read by the immunofluorescence reader. The higher the D-dimer concentration in the sample, the more complexes are immobilized on the test line, and consequently, the stronger the fluorescence signal generated and detected by the reader. The concentration of D-dimer in the sample, expressed in mg/L FEU, is calculated according to the calibration curve provided in the Identification Card.

REAGENTS

1- Cassette - Store between 2 and 30 °C. Each sachet contains:

- A) 01 plastic device with 01 test strip containing:
 - Test line region: mouse anti-D-dimer antibody.
 - Control line region: goat anti-rabbit IgG antibody.
 - Conjugate: fluorescent microspheres conjugated to mouse anti-D-dimer antibodies and fluorescent microspheres conjugated to rabbit IgG antibodies.

B) 01 silica sachet.

2- Sample Diluent - Store between 2 and 30 °C. Contains: 10 mM phosphate-buffered saline (PBS), stabilizers (casein <1%) and preservatives (ProClin™ 300 <0.5%).

3- Identification Card

PRESENTATION

REAGENTS	NUMBER OF REAGENTS PER PACKAGE
1- Cassette	25 units
2- Sample Diluent	25 units (25 x 0.160 mL)
3- Identification Card	1 unit

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation lower than 1.5%.
- 2- BIOCLIN POCT FIA Immunofluorescence Reader.
- 3- Stopwatch or clock.
- 4- Quality Control Tools.

STORAGE AND TRANSPORTATION CONDITIONS

The storage and transportation temperature should be between 2 °C and 30° C. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- BIOCLIN POCT FIA D-DIMER should only be used in conjunction with the BIOCLIN POCT FIA Immunofluorescence Reader.

3- Strictly follow the proposed methodology to obtain accurate results.

4- The cassettes (test devices) do not contain contaminating substances, but since the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.

5- Open the aluminum envelope to remove the cassette only at the time of use.

6- The kit should not be frozen or exposed to high temperatures (above 45°C), as this will cause deterioration.

7- The reagents, as well as the samples, must reach room temperature before starting the test.

8- Do not mix reagents from kits with the same batch or with different batch numbers.

9- Do not use the kit after its expiration date. The expiration date is indicated on the product label.

10- The test device (R1- Cassette) and sample dilution tubes (R2- Sample Diluent) are disposable and cannot be reused.

11- The kit materials that have come into contact with the patient's sample must be disposed of in a special way. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biological waste.

12- Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular waste.

13- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.

14- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or by requesting it through Quibasa's SAC (Customer Advisory

Service).

15- Do not use the product if the packaging is damaged.

16- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

SAMPLES

Plasma (Heparin, EDTA, or Citrate)

Samples should be collected by venipuncture in tubes containing anticoagulant (Heparin, EDTA, or Citrate). To avoid hemolysis, centrifuge the samples to obtain the plasma and test them as quickly as possible. Do not keep samples at room temperature for long periods. Plasma samples can be stored refrigerated at 2–8°C for a maximum of 3 days. If samples cannot be analyzed within 3 days, they can be stored for up to 30 days at -20°C. Plasma samples stored for longer than the recommended periods and samples collected with anticoagulants other than those mentioned should not be used. Highly hemolyzed, lipemic, and precipitate-containing samples should be avoided, as they may provide inconsistent results.

Whole Blood (Heparin, EDTA, or Citrate)

Whole blood samples should be collected by venipuncture in tubes containing heparin, EDTA, or citrate as an anticoagulant. Samples should be tested immediately or as soon as possible after collection. Whole blood collected in tubes containing anticoagulants can be stored refrigerated at 2–8°C for a maximum of 2 days. Do not use samples collected with anticoagulants other than those mentioned, and do not freeze whole blood samples

PROCESS DESCRIPTION

Stability After Opening

The stability test results show that the BIOCLIN POCT FIA D-DIMER kit is stable after opening for up to 70 minutes. This stability may vary depending on the conditions of use of the product and the environment. However, for safety reasons, we recommend that the product should be used immediately after removing the cassette from the sachet.

TECHNIQUE

Check the BIOCLIN POCT FIA Immunofluorescence Reader manual for complete instructions on handling the equipment.

Before starting the test, allow all reagents (Cassette and Sample Diluent) and samples to stabilize at room temperature (15-30 °C) for at least 40 minutes.

1- Check if the information contained in the Identification Card has been properly entered into the equipment (with each new batch of the kit, the calibration and validity data contained in the Identification Card must be transferred to the equipment).

2- Remove the cassette (test device) from its protective packaging, place it on a clean, level surface and identify it appropriately.

3- Make sure that the Sample Diluent solution (Reagent No. 2) is at the bottom of the dilution tube.

4- Add 50 µL of sample to a dilution tube containing the Sample Diluent.

5- Close the tube cap and mix by inverting the sample 6 to 8 times until the sample and diluent are completely homogenized.

6- Pipette 50 µL of diluted and properly mixed sample into the

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST,

CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

sample application hole of the cassette.

7- Option 1: Insert the cassette into the equipment immediately after applying the sample and select the "**Timing Test**" option. After 15 minutes, the equipment will take the reading and the result will appear on the display.

Option 2: After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the "**Quick Test**" option. The equipment will take the reading immediately and the result will appear on the display.

RESULTS

The BIOCLIN POCT FIA Immunofluorescence Reader calculates the cardiac troponin I concentration automatically, using the calibration curve contained in the R3 (Identification Card). The information on the Identification Card must be properly entered into the equipment for each new batch of reagents.

Quantitative results are measured in mg/L FEU. The kit's measurement range is between 0.1 and 10 mg/L FEU (limit of quantification and maximum point of the curve, respectively). Results below the lower limit of detection are expressed as <0.1 mg/L FEU. Results above the upper limit of detection are expressed as >10 mg/L FEU.

PROCESS LIMITATIONS

1- It is important to use the correct sample volume, as lower or higher volumes may produce erroneous results.

2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretations of the results.

3- Reactive results must be confirmed with another test and available clinical information. The isolated test results cannot be used for definitive diagnosis.

4- The results provided by this kit must be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

REFERENCE VALUES

Normal	< 0.5 mg/L FEU
--------	----------------

Conversion: ng/mL = mg/L x 1.000

These values should be used as a guide, and each laboratory should create its own range of reference values according to the population served. The results provided by this kit should be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed by bilirubin up to 350 µmol/L, triglycerides up to 40 mmol/L, hemoglobin up to 2 g/L, rheumatoid factor up to 1,500 IU/mL, heterophile antibodies up to 1,000 ng/mL, and antinuclear antibodies up to 1:640.

CROSS-REACTIVITY

No false-positive results were observed for fibrin degradation products up to 1,000 mg/L. Despite the results, the possibility of cross-reactivity with other analytes or higher concentrations cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along

with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. For this purpose, it is recommended to use controls, which allow the assessment of the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The BIOCLIN POCT FIA D-DÍMER kit was compared with another commercially available D-Dimer assay method. 107 analyses were performed, and the results were evaluated. The linear equation obtained was $Y = 0.9904x + 0.0114$, and the correlation coefficient was 0.9739. Based on these results, it can be concluded that the kit presents good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 2 samples with different values, obtaining the following concentration results:

REPEATABILITY	SAMPLE	
	1	2
Mean	0.53	4.84
Standard Deviation	0.04	0.47
Coefficient of Variation (%)	6.8	9.7

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 samples with different values, obtaining the following concentration results:

REPRODUCIBILITY	SAMPLE	
	1	2
Mean	0.51	4.94
Standard Deviation	0.0	0.4
Coefficient of Variation (%)	6.7	8.9

ANALYTICAL SENSITIVITY

The Limit of Blank (LoB) of the BIOCLIN POCT FIA D-DIMER product is 0.03 mg/L FEU, the Limit of Detection (LoD) is 0.05 mg/L FEU and the Limit of Quantification (LoQ) is 0.1 mg/L FEU.

LINEARITY

The BIOCLIN POCT FIA D-DIMER kit is linear up to 10 mg/L FEU (highest point of the calibration curve).

HIGH DOSE PROZONE EFFECT

No high-dose prozone effect was observed at up to 60 mg/L FEU of D-Dimer.

DIAGNOSTIC SIGNIFICANCE

The blood clotting process produces fibrin as a final product, which, together with platelets, forms a network responsible

for stopping bleeding, protecting the wound from infection, and providing a matrix for new cell growth. After the damaged tissue has healed, fibrin is degraded in a process called fibrinolysis, triggered by the enzyme plasmin. Plasmin degradation generates fibrin degradation products (FDPs), which include D-dimer. D-dimer is not normally present in human blood plasma except when the coagulation system, and consequently the fibrinolytic process, has been severely activated. Therefore, it is considered a specific marker for thromboembolic events, such as deep vein thrombosis, pulmonary embolism, and disseminated intravascular coagulation. Deep vein thrombosis is a condition in which a clot forms in a deep vein, most commonly in the lower limbs. The main symptoms are pain, redness, and swelling at the affected site. Pulmonary embolism occurs when a clot blocks an artery in the lung, primarily causing shortness of breath and chest pain. Finally, disseminated intravascular coagulation is characterized by severe bleeding due to the depletion of clotting factors and platelets due to excessive clot formation in various vessels throughout the body. It occurs most frequently after obstetric complications, gram-negative infections, cancer, and shock. It is a complex condition that can be fatal if not properly identified and treated. While it is important to rule out the conditions listed above, there are other clinical conditions where increased D-dimer concentrations can be observed, such as postoperatively, during pregnancy, peripheral vascular disease, cancer, kidney failure, sepsis, and other inflammatory diseases, such as COVID-19.

BIBLIOGRAPHIC REFERENCES

- 1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 2- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- 3- RAMZI, D.W.; Leeper, K. V. DVT and pulmonary embolism: part I. Diagnosis. American Family Physician, v. 69, n. 12, p. 2829-2836, 2004.
- 4- Perrier, A. et al. D-dimer testing for suspected pulmonary embolism in outpatients. American Journal of Respiratory and Critical Care Medicine, v. 156, n. 2, p. 492-496, 1997.
- 5- Sakamoto, K. et al. D-Dimer is helpful for differentiating acute aortic dissection and acute pulmonary embolism from acute myocardial infarction. Hellenic J. Cardiol., v. 52, n. 2, p. 123-127, 2011.
- 6- Yasuka, S. Kubota, S. The value of blood D-Dimer test in the diagnosis of walk-in patients with venous thromboembolism. Vas. Health Rish Manag., v. 7, p. 125-127, 2011.
- 7- Bioclin – Dados de arquivos.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration number for BIOCLIN POCT FIA D-DIMER kit: 10269360476

Review: September/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER