# Bioclin

## BIACLIN POCT FIA TROPONINA I

REF K339

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#### FINALIDADE

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**1-** WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

**2-** Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.

**3-** Cummins, B.; Auckland, M. L.; Cummins, P. Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction. American Heart Journal, v. 113, n. 6, p. 1333-1344, 1987.

**4-** Tate, J. et al. Standardization of cardiac troponin I measurement: past and present. Pathology, b. 42, n. 5, p. 402-408, 2010.

**5-** Reichlin, T. et al. Early Diagnosis of Myocardial Infarction with Sensitive Cardiac Troponin Assays. The New England Journal of Medicine, v. 361, n. 9, p. 858-867, 2009.

**6-** Paneghini, M. et al. Standardization of troponin I measurements: an update. Clin. Chem. Lab. Med., v. 46, n. 11, p. 1501-1506, 2008.

**7-** Filatov, V. L. et al. Troponin: Structure, Properties, and Mechanism of Functioning. Biochemistry, v. 64, n. 9, p. 969-985, 1999.

**8-** Keller, T. et al. Sensitive Troponin I Assay in Early Diagnosis of Acute Myocardial Infarction. The New England Journal of Medicine, v. 361, n. 9, p. 868-877, 2009.

**9-** European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction Redefined—A Consensus Document of The Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Journal of the American College of Cardiology, v. 36, n. 3, p. 959-969, 2000.

**10-** Thygesen, K. et al. Fourth Universal Definition of Myocardial Infarction (2018). Circulation, v. 138, p. e618-e651, 2018.

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### INSTRUCTIONS FOR USE

#### FUNCTION

Test for the quantitative determination of cardiac troponin I (cTnI) in human biological samples of serum, plasma, and whole blood, using fluorescent immunoassay. For *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Fluorescent immunoassay.

The BIOCLIN POCT FIA TROPONIN I kit is a fluorescent immunoassay based on the sandwich detection principle, designed for the quantitative determination of cardiac troponin I in human serum, plasma, or whole blood samples (collected in EDTA or heparin). When the sample is added to its application site, the interaction between the troponin I (cTnI) molecules present in the sample and the anti-troponin I antibodies conjugated to fluorescent microspheres occurs, forming the cTnI-antibody conjugate complex. As the sample migrates by capillary action through the nitrocellulose membrane, this complex is captured by the test line, which consists of anti-troponin I antibodies. The sample will continue to migrate through the membrane until it reaches the control line, which consists of anti-rabbit IgG antibodies (where rabbit IgG antibodies conjugated to fluorescent microspheres will be captured). At the end of the incubation period, the cassette is read in the immunofluorescence reader. The higher the troponin I concentration in the sample, the more complexes are immobilized on the test line and, consequently, the higher the fluorescence signal generated and captured by the reader. The troponin I concentration in the sample, in ng/mL, is calculated according to the assay calibration curve (contained on the Identification Card).

#### REAGENTS

**1- Cassette** - Store between 2 and 30 °C. Each sachet contains:

- A) 01 plastic device with 01 test strip containing:
  - Test line region: mouse anti-troponin I antibody.
  - Control line region: goat anti-rabbit IgG antibody.
  - Conjugate: fluorescent microspheres conjugated to mouse anti-troponin I antibodies and fluorescent microspheres conjugated to rabbit IgG antibodies.
- B) 01 silica sachet.

**2- Diluent for Whole Blood Samples** - Store between 2 and 30 °C. Contains: 10 mM phosphate-buffered saline (PBS), stabilizers (casein <1%) and preservatives (ProClin™ 300 <0.5%).

#### 3- Identification Card

#### PRESENTATION

REAGENTS	NUMBER OF REAGENTS PER PACKAGE
1- Cassette	25 units
2-Diluent for Whole Blood Samples	25 units (25 x 0.05 mL)
3- Identification Card	1 unit

#### EQUIPMENTS AND OPERATIONAL INPUTS

**Materials in the kit:**

- Reagents described in the table above.

#### Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation lower than 1.5%.
- 2- BIOCLIN POCT FIA Immunofluorescence Reader.
- 3- Stopwatch or clock.
- 4- Quality Control Tools.

#### STORAGE AND TRANSPORTATION CONDITIONS

The storage and transportation temperature should be between 2 °C and 30 °C. Keep away from light and avoid humidity. **Do not freeze.**

#### SPECIAL CARE

**1- For *in vitro* diagnostic use only.**

**2- BIOCLIN POCT FIA TROPONIN I should only be used in conjunction with the BIOCLIN POCT FIA Immunofluorescence Reader.**

**3-** Strictly follow the proposed methodology to obtain accurate results.

**4-** The cassettes (test devices) do not contain contaminating substances, but since the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.

**5-** Open the aluminum envelope to remove the cassette only at the time of use.

**6-** The kit should not be frozen or exposed to high temperatures (above 45 °C), as this will cause deterioration.

**7-** The reagents, as well as the samples, must reach room temperature before starting the test.

**8-** Do not mix reagents from kits with the same batch or with different batch numbers.

**9-** Do not use the kit after its expiration date. The expiration date is indicated on the product label.

**10-** The test device (R1- Cassette) and sample dilution tubes (R2- Diluent for Whole Blood Samples) are disposable and cannot be reused.

**11-** The kit materials that have come into contact with the patient's sample must be disposed of in a special way. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biological waste.

**12-** Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular waste.

**13-** We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.

**14-** To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or by

requesting it through Quibasa's SAC (Customer Advisory Service).

**15-** Do not use the product if the packaging is damaged.

**16-** It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

#### SAMPLES

**Serum or Plasma (EDTA or Heparin)**

Samples should be collected by venipuncture, in tubes with or without anticoagulant (EDTA or Heparin). To avoid hemolysis, centrifuge the samples to obtain serum/plasma and test them as quickly as possible. Do not keep samples at room temperature for long periods. Serum and plasma samples can be stored refrigerated, between 2 and 8°C, for a maximum of 3 days. If samples cannot be analyzed within 3 days, they can be stored for up to 30 days at -20°C. Serum and plasma samples stored for longer than the recommended periods and samples collected with anticoagulants other than those mentioned should not be used. Highly hemolyzed, lipemic, and precipitate-containing samples should be avoided, as they may provide inconsistent results.

#### Whole Blood (EDTA or Heparin)

Whole blood samples should be collected by venipuncture in tubes containing heparin or EDTA as an anticoagulant. Samples should be tested immediately or as soon as possible after collection. Whole blood collected in tubes containing anticoagulants can be stored refrigerated at 2–8°C for a maximum of 2 days. Do not use samples collected with anticoagulants other than those mentioned, and do not freeze whole blood samples.

#### PROCESS DESCRIPTION

##### Stability After Opening

The stability test results show that the BIOCLIN POCT FIA TROPONIN I kit is stable after opening for up to 70 minutes. This stability may vary depending on the conditions of use of the product and the environment. However, for safety reasons, we recommend that the product should be used immediately after removing the cassette from the sachet.

#### TECHNIQUE

**Check the BIOCLIN POCT FIA Immunofluorescence Reader manual for complete instructions on handling the equipment.**

**ATTENTION!** The techniques for serum/plasma and whole blood are different. Carefully read the instructions for use described below, according to the type of sample to be used.

#### Serum and Plasma Samples

Before starting the assay, allow all reagents (Cassette) and samples to stabilize at room temperature (15-30°C) for at least 40 minutes.

**1-** Verify that the information contained in the Identification Card has been correctly entered into the equipment (with each new batch of the kit, the calibration and validity data contained in the Identification Card must be transferred to the equipment).

**2-** Remove the cassette (test device) from its protective packaging, place it on a clean, level surface, and label it appropriately.

**3-** Pipette 50 µL of sample (serum or plasma in EDTA or Heparin) directly into the sample application port of the cassette.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

**4- Option 1:** Insert the cassette into the equipment immediately after applying the sample and select the “**Timing Test**” option. After 15 minutes, the equipment will take the reading, and the result will appear on the display.

**Option 2:** After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the “**Quick Test**” option. The equipment will take the reading immediately, and the result will appear on the display.

#### Whole Blood Samples

Before starting the test, allow all reagents (Cassette and Whole Blood Diluent) and samples to stabilize at room temperature (15-30°C) for at least 40 minutes.

**1-** Verify that the information contained in the Identification Card has been correctly entered into the equipment (with each new batch of the kit, the calibration and expiration data contained in the Identification Card must be transferred to the equipment).

**2-** Remove the cassette (test device) from its protective packaging, place it on a clean, level surface, and label it appropriately.

**3-** Ensure that the solution of the Diluent for Whole Blood Samples (Reagent No. 2) is at the bottom of the dilution tube.

**4-** Add 50 µL of sample (whole blood collected in EDTA or Heparin) to a dilution tube containing the Diluent for Whole Blood Samples.

**5-** Close the tube cap and mix by inverting 6 to 8 times until the sample and diluent are completely mixed.

**6-** Pipette 50 µL of the diluted and properly mixed sample into the cassette's sample application port.

**7- Option 1:** Insert the cassette into the equipment immediately after applying the sample and select the “**Timing Test**” option. After 15 minutes, the equipment will take the reading and the result will appear on the display.

**Option 2:** After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the “**Quick Test**” option. The equipment will take the reading immediately, and the result will appear on the display.

#### RESULTS

The BIOCLIN POCT FIA Immunofluorescence Reader automatically calculates cardiac troponin I concentrations using the calibration curve contained in the R3 (Identification Card). The information on the Identification Card must be properly entered into the equipment for each new batch of reagents.

Quantitative results are measured in ng/mL. The kit's measurement range is between 0.1 and 50 ng/mL (limit of quantification and maximum point of the curve, respectively). Results below the lower limit of detection are expressed as <0.1 ng/mL. Results above the upper limit of detection are expressed as >50 ng/mL.

#### PROCESS LIMITATIONS

**1-** It is important to use the correct sample volume, as lower or higher volumes may produce erroneous results.

**2-** The reaction reading time must be followed according to the established technique to avoid misinterpretation of the results.

**3-** The BIOCLIN POCT Troponin I kit for cTnI is a specific test for aiding in the diagnosis of acute myocardial infarction.

Reactive results must be confirmed with another test and available clinical information. Test results alone cannot be used for a definitive diagnosis.

**4-** It is recommended that established guidelines for the diagnosis of acute myocardial infarction be followed.

**5-** The results provided by this kit must be interpreted by the responsible medical professional and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

#### REFERENCE VALUES

Normal	<0.3 ng/mL
--------	------------

These values should be used as a guide, and each laboratory should create its own range of reference values according to the population served. The results provided by this kit should be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

#### INTERFERENTS

No interference was observed by bilirubin up to 350 µmol/L, triglycerides up to 40 mmol/L, hemoglobin up to 2 g/L, rheumatoid factor up to 1,500 IU/mL, heterophile antibodies up to 1,000 ng/mL, and antinuclear antibodies up to 1:640.

#### CROSS-REACTIVITY

No false-positive results were observed for cardiac troponin C up to 1,000 ng/mL, cardiac troponin T up to 500 ng/mL, and skeletal troponin T up to 1,000 ng/mL. Despite these results, the possibility of cross-reactivity with other analytes or higher concentrations cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. For this purpose, it is recommended to use controls, which allow the assessment of the precision and accuracy of the dosages.

#### PRODUCT PERFORMANCE

##### ACCURACY

##### Comparison of Methods

The BIOCLIN POCT FIA TROPONIN I kit was compared with another commercially available troponin I assay method. 104 analyses were performed, and the results were evaluated. The linear equation obtained was  $Y = 0.9848x + 0.1186$  and the correlation coefficient was 0.974. Based on these results, it can be concluded that the kit presents good methodological specificity.

##### PRECISION

##### Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following concentration results:

REPEATABILITY	SAMPLE		
	1	2	3
Mean	0.74	1.00	19.45
Standart Deviation	0.02	0.1	1.74
Coefficient of Variation	3.3	7.7	8.9

##### Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following concentration results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Mean	0.76	0.98	19.93
Standart Deviation	0.1	0.1	1.8
Coefficient of Variation	7.4	8.3	8.8

##### ANALYTICAL SENSITIVITY

The Limit of Blank (LoB) of the BIOCLIN POCT FIA TROPONIN I product is 0.02 ng/mL, the Limit of Detection (LoD) is 0.07 ng/mL and the Limit of Quantification (LoQ) is 0.1 ng/mL.

##### LINEARITY

The BIOCLIN POCT FIA TROPONIN I kit is linear up to 50 ng/mL (highest point of the calibration curve).

##### HIGH DOSE PROZONE EFFECT

No high-dose prozone effect was observed at up to 500 ng/mL of cardiac troponin I.

##### DIAGNOSTIC SIGNIFICANCE

Troponins, classified as troponin C (TnC), troponin T (TnT), and troponin I (TnI), are proteins that form protein complexes present in the muscle fibers of skeletal and cardiac muscles and play important roles in muscle contraction in these tissues. While TnC is similar in both types of muscle tissue, troponins T and I have isoforms specific to cardiac muscle fibers. This characteristic is important in the use of TnT and TnI as specific markers of cardiac injury, especially in the diagnosis of acute myocardial infarction (AMI). AMI is a major cause of death and disability, and rapid diagnosis is critical for initiating effective treatment and improving prognosis. It is characterized by characteristic symptoms (such as chest pain and dyspnea), elevation of specific markers in the blood, and the development of a Q-wave pattern on the electrocardiogram. Acute myocardial infarction occurs due to necrosis of cardiac tissue associated with prolonged ischemia. This myocardial necrosis results in the release of several markers into the bloodstream, such as myoglobin, cardiac troponins T and I, creatine kinase, and lactate dehydrogenase. Cardiac troponin I is an important marker because its blood concentration rises 3–6 hours after the infarction, peaking between 12–24 hours, and maintaining altered levels for 5–14 days after the AMI episode. Although troponin T is also cardiospecific, troponin I measurement is preferred because it often presents higher levels and rises more rapidly after an AMI. Although useful in the diagnosis of AMI, troponin levels rise in cases of myocardial injury without ischemia, such as myocarditis, or secondary to other diseases, such as renal failure.

#### BIBLIOGRAPHIC REFERENCES

- 1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 2- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- 3- Cummins, B.; Auckland, M. L.; Cummins, P. Cardiac-specific troponin-I radioimmunoassay in the diagnosis of acute myocardial infarction. American Heart Journal, v. 113, n. 6, p. 1333-1344, 1987.
- 4- Tate, J. et al. Standardization of cardiac troponin I measurement: past and present. Pathology, b. 42, n. 5, p. 402-408, 2010.
- 5- Reichlin, T. et al. Early Diagnosis of Myocardial Infarction with Sensitive Cardiac Troponin Assays. The New England Journal of Medicine, v. 361, n. 9, p. 858-867, 2009.
- 6- Paneghini, M. et al. Standardization of troponin I measurements: an update. Clin. Chem. Lab. Med., v. 46, n. 11, p. 1501-1506, 2008.
- 7- Filatov, V. L. et al. Troponin: Structure, Properties, and Mechanism of Functioning. Biochemistry, v. 64, n. 9, p. 969-985, 1999.
- 8- Keller, T. et al. Sensitive Troponin I Assay in Early Diagnosis of Acute Myocardial Infarction. The New England Journal of Medicine, v. 361, n. 9, p. 868-877, 2009.
- 9- European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction Redefined—A Consensus Document of The Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Journal of the American College of Cardiology, v. 36, n. 3, p. 959-969, 2000.
- 10- Thygesen, K. et al. Fourth Universal Definition of Myocardial Infarction (2018). Circulation, v. 138, p. e618-e651, 2018.
- 11- BIOCLIN – Dados de arquivos.

#### QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

**ANVISA registration number for BIOCLIN POCT FIA TROPONIN I kit: 10269360479**

Review: October/2025

#### UNIVERSAL SYMBOLOGY

 REF	CATALOG NUMBER	 MADE BY
 LOT	LOT NUMBER	 CONTROL
 MANUFACTURING DATE		 POSITIVE CONTROL
 VALIDITY DATE (last day of the month)		 NEGATIVE CONTROL
 TEMPERATURE LIMIT (store)		 BIOLOGICAL RISK
 CONTENT IS SUFFICIENT FOR <N> TEST		 FLAMMABLE
 SEE INSTRUCTIONS FOR USE		 CORROSIVE
 IN VITRO DIAGNOSTIC PRODUCT		 TOXIC
 KEEP AWAY FROM SUNLIGHT		 DO NOT USE IF PACKAGE IS DAMAGED
 DO NOT REUSE		 PRODUCT STERILIZED
 CAUTION		 DANGER