BIOCLIN POCT FIA PCR/PCR-US

REF **K282**

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FINALIDADE

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1. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
2. Milano SS; et al. C-reactive Protein is a Predictor of Mortality in ST-segment Elevation Acute Myocardial Infarction. Int J Cardiovasc Sci, 2019, 32(2)):118–24.
3. Kao, Pai C et al. "Serum C-reactive protein as a marker for wellness assessment." Annals of clinical and laboratory science vol. 36,2 (2006): 163-9.
4. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.
5. Levinson SS.The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Test for the quantitative determination of C-Reactive Protein (CRP) and Ultrasensitive C-Reactive Protein (US-CRP) in biological samples (serum, plasma or whole blood) through fluorescent immunoassay. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay
The BIOCLIN POCT FIA CRP/US-CRP kit is a test designed to detect CRP in serum, plasma and whole blood samples. When sample and diluent are added, C-Reactive Protein (CRP) in the sample interacts with Anti-CRP monoclonal Antibodies conjugated to fluorescent microbead particles. Then, the sample migrates through the membrane by capillarity, up to the detection line (with specific Anti-CRP monoclonal Antibodies present in the "T" line). The sample will continue migrating through the membrane until it reaches the "C" control region, where the Conjugate interacts with Anti-IgG Antibodies present in the control line. At the end of the incubation period, the cassette must be inserted into the reader for quantification. The CRP concentration of the sample correlates with the intensity of the fluorescence signal emitted by the test line of the assay.

The BIOCLIN POCT FIA Reader calculates the CRP concentration in the sample in mg/L according to the assay calibration curve (contained on the Identification Chip).

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 8 °C. **Do not freeze.** Each sachet contains:
A) 01 plastic device with 01 test strip containing:
- Control region (C): Anti-IgG Antibodies.
- Test region (T): Anti-CRP monoclonal Antibodies.
- Conjugate: Anti-CRP monoclonal Antibody conjugated with fluorescent microspheres.
B) 01 sachet of silica.

Number 2 (R2) Diluent - Store between 2 and 8 °C. **Do not freeze.** Contains: Buffer, surfactant and preservative.

Number 3 (R3) Identification Chip.

PRESENTATION

Reagents	Number of Reagents per Packaging
1 - Cassette	25 units
2 - Diluent	25 units
3 - Chip Identification	1 unit

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table
- Instructions for use (manual)

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation less than 1.5%.
- 2- BIOCLIN POCT FIA incubator (optional).
- 3- BIOCLIN POCT FIA Reader.
- 4- Stopwatch or clock.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- **For *in vitro* diagnostic use only.**
- 2- **The BIOCLIN POCT FIA CRP/US-CRP must only be used together with the BIOCLIN POCT FIA reader.**
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples used are biological, the safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Open the aluminum envelope to remove the cassette, only when using it.
- 5- The kit must not be frozen or exposed to high temperatures (45°C), as this will cause deterioration.
- 6- The reagents, as well as the samples, must reach room temperature before starting the test.
- 7- Do not mix reagents from kits with the same batch or with different batch numbers.
- 8- Do not use the kit beyond its expiration date. The expiry date is indicated on the product label.
- 9- The cassette (test device) and sample dilution tubes are disposable and cannot be reused.
- 10- The kit materials that came into contact with the patient's sample must have a special destination. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biowaste.
- 11- Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular trash.
- 12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the FISPQ (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service Customer Advice) by Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Serum or Plasma: sample collection must be performed by venipuncture using collection tubes with EDTA anticoagulant (to obtain plasma samples) or tube without anticoagulant (to collect serum), avoiding prolonged tourniquet. The samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 7 days. If the samples cannot be analyzed within this period, they can be stored for up to 6 months at -20 °C (freezer)*. Separate plasma as soon as possible to avoid hemolysis.

Whole Blood: Whole blood collected by venipuncture with EDTA anticoagulant, the sample can be kept under refrigeration, between 2 and 8 °C, for a maximum period of 7 days*. **Do not freeze.**

*At the time of use, samples must be stabilized at room temperature.

PROCESS DESCRIPTION

Stability After Opening

The stability of the BIOCLIN POCT FIA CRP/US-CRP cassette (after opening the aluminized sachet) is 1 hour. This stability may vary according to the conditions of use of the product and the environment. However, to be on the safe side, we recommend that the BIOCLIN POCT FIA CRP/US-CRP kit be tested immediately after removing the cassette from the sachet.

TECHNIQUE

- 1- **The Sample, Cassette and Sample Dilution Tube must be between 15-30 °C to run the test.**
 - 2- Check if the identification chip is inserted in the BIOCLIN POCT FIA Reader.
 - 3- **For Serum or Plasma:** pipette 3 µL of the sample into the sample dilution tube (Reagent N° 2).
 - 4- **For Whole Blood:** pipette 5 µL of the sample into the sample dilution tube (Reagent N° 2).
 - 5- Homogenize the tube with the sample.
 - 6- Remove the cassette (test device) from its protective packaging, place it on a clean and level surface and identify it properly.
 - 7- Pipette **80 µL of the diluted sample** into the sample well.
- Note:** Steps 2 through 7 must be performed within 1 minute to ensure assay accuracy.
- 8- Make sure there are no air bubbles.
 - 9- **Wait 3 minutes to read the test. Do not perform the reading before 3 minutes.**

Cassette Reading

- 1- Select the test mode on the BIOCLIN POCT FIA reader. For more information on how to select the reading mode, refer to the equipment manual.
- 2- Insert the Cassette into the BIOCLIN POCT FIA reader and wait for the result to be issued.

NOTES

- 1- The cassette must be used immediately after opening its packaging.
- 2- The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

Calibration

The calibration of the method is carried out from the Identification Chip, which must be changed with each change in the lot of reagents.

PROCESS LIMITATIONS

- 1- It is important to use the correct sample volume, as lower or higher volumes may determine erroneous results.
- 2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretation of the results.
- 3- The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERING

Samples with concentrations of up to 400 µmol/L of Bilirubin, 20 mmol/L of Triglycerides, 15 mg/mL of Hemoglobin and 3250 IU/mL of Rheumatoid Factor do not interfere with the test result.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

REFERENCE VALUES

Analyte	Value in Reference (mg/L)
Ultrasensitive CRP	< 1.0
CRP	< 10.0

These values should be used as guidelines, and each laboratory should create its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods and Methodological Specificity

The BIOCLIN POCT FIA CRP/US-CRP kit was compared with another commercially available CRP assay method. 115 analyzes were performed and the results were evaluated. The linear equation obtained was Y= 0.90x + 0.26 and the correlation coefficient 0.997. With these results, it can be concluded that the kit has good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 4 samples with different concentrations, obtaining the following results:

Repeatability	Sample			
	1	2	3	4
Concentration (mg/L)	2.55	15.16	47.56	82.19
Standard Deviation (mg/L)	0.19	1.00	1.46	1.94
Coefficient of Variation (%)	7	7	3	2

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 4 samples with different concentrations, obtaining the following results:

Reproducibility	Sample			
	1	2	3	4
Concentration (mg/L)	2.30	16.02	47.52	82.83
Standard Deviation (mg/L)	0.29	1.22	2.60	3.50
Coefficient of Variation (%)	13	8	5	4

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

SENSITIVITY

The analytical sensitivity of the BIOCLIN POCT FIA CRP/US-CRP kit is 0.42 mg/L.

PROZONE EFFECT

There was no prozone effect up to a concentration of 500 mg/L.

LINEARITY

The operating range for the BIOCLIN POCT FIA CRP/US-CRP reagent is 0.5 mg/L to 200 mg/L.

DIAGNOSTIC SIGNIFICANCE

C-Reactive Protein (CRP) is an annular (ring-shaped) pentameric protein found in the blood whose levels increase in response to inflammation. It is an acute phase protein of hepatic origin that increases after secretion of interleukin-6 by macrophages and T cells. In healthy adults, normal CRP concentrations are less than 10 mg/L. When there is a stimulus, the CRP level can increase rapidly. CRP levels are applied to monitor infectious and non-infectious inflammatory processes. In addition, the CRP concentration may be increased in conditions such as neoplastic diseases, intense trauma and burns. While the evaluation of hs-CRP (ultrasensitive CRP) can help in assessing cardiac risk, since CRP levels above the reference value, in apparently normal patients, can signal coronary artery disease (indicator for AMI and CVA). In this scenario, Ultrasensitive CRP can be used to assess cardiovascular diseases if the patient has no inflammatory conditions, since, in these situations, its values are significantly increased.

REFERENCES BIBLIOGRAPHICS

1. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

2. Milano SS; et al. C-reactive Protein is a Predictor of Mortality in ST-segment Elevation Acute Myocardial Infarction. Int J Cardiovasc Sci, 2019, 32(2)):118–24.

3. Kao, Pai C et al. "Serum C-reactive protein as a marker for wellness assessment." Annals of clinical and laboratory science vol. 36,2 (2006): 163-9.

4. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.

5. Levinson SS.The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 – Santa Branca
CEP 31565-130 – Belo Horizonte – MG – Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOCLIN POCT FIA CRP/US-CRP kit:
10269360428

Review: June/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER