BIOCLIN POCT FIA PROCALCITONINA

REF K286

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FINALIDADE

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BIOCLIN POCT FIA PROCALCITONIN

REF K286
INSTRUCTIONS FOR USE

FUNCTION
Test for quantitative determination of the analyte Procalcitonin (PCT) in samples biological samples (serum, whole blood and plasma) through fluorescent immunoassay. Only for use *in vitro* diagnosis use.

PRINCIPLE OF ACTION
Methodology: Immunoassay fluorescent.
The BIOCLIN POCT FIA PROCALCITONIN kit is a test intended for the detection of the analyte Procalcitonin in plasma, whole blood or serum samples. When a sample is added, it interacts with Anti-Procalcitonin monoclonal Antibodies conjugated to fluorescent beads. Then, the sample migrates through the membrane by capillarity, and interacts up to the detection line (with the specific Anti-Procalcitonin monoclonal Antibodies present in the “T” test line). The sample will continue migrating through the membrane until it reaches the control region “C”, where the Conjugate containing IgG Antibodies with fluorescent microspheres interacts with immobilized Anti-IgG Antibodies. At the end of the incubation period, the cassette must be inserted into the reader for quantification. The Procalcitonin concentration of the sample directly correlates with the intensity of the fluorescence signal emitted by the test line of the assay.

The BIOCLIN POCT FIA Reader calculates the concentration of Procalcitonin in the sample in ng/mL according to the assay calibration curve (contained on the Identification Chip).

The BIOCLIN POCT FIA Reader calculates the concentration of Procalcitonin in the sample in ng/mL according to the assay calibration curve (contained on the Identification Chip).

REAGENTS
Number 1 (R1) Cassette - Store between 2 and 8 °C. **Do not freeze.** Each sachet contains:
A) 01 plastic device with 01 test strip containing:
-Control region (C): Goat Anti-IgG Antibody.
-Test region (T): Anti-Procalcitonin specific monoclonal Antibodies.
-Conjugate: Anti-Procalcitonin monoclonal Antibody conjugated with fluorescent microspheres and goat IgG antibody conjugated with fluorescent microspheres.
B) 01 sachet of silica.

Number 2 (R2) Diluent - Store between 2 and 8 °C. **Do not freeze.** Contains: Buffer, surfactant and preservative.

Number 3 (R3) Identification Chip.

Reagents	Number of Reagents per Package
1 - Cassette	25 units
2 - Diluent	25 units
3 - Identification Chip	1 unit

EQUIPMENT AND OPERATIONAL INSTRUCTIONS
Materials contained in the kit:
- Reagents described in the previous table
- Instructions for use (manual)

Required materials not contained in the kit:
1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation less than 1.5%.
2- BIOCLIN POCT FIA incubator (optional).
3- BIOCLIN POCT FIA Reader.
4- Stopwatch or clock.

STORAGE AND TRANSPORT CONDITIONS
The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture.
Do not freeze.

SPECIAL CARES
1- For *in vitro* diagnostic use only.
2- The BIOCLIN POCT FIA PROCALCITONIN must only be used together with the BIOCLIN POCT FIA reader.
3- The cassettes (test devices) do not contain contaminating substances, but as the samples used are biological, the safety routines must be strictly followed and the use of disposable gloves is essential.
4- Open the aluminum envelope to remove the cassette, only when using it.
5- The kit must not be frozen or exposed to high temperatures (45°C), as this will cause deterioration.
6- The reagents, as well as the samples, must reach room temperature before starting the test.
7- Do not mix reagents from kits with the same batch or with different batch numbers.
8- Do not use the kit beyond its expiration date. The expiry date is indicated on the product label.
9- The cassette (test device) and sample dilution tubes are disposable and cannot be reused.
10- The kit materials that came into contact with the patient's sample must have a special destination. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biowaste.
11- Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular trash.
12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
13- To obtain information related to biosafety or in case of accidents with the product, consult the FISPQ (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service Customer Advice) by Quibasa.
14- Do not use the product if the packaging is damaged.
15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLE
Serum or plasma: Sample collection must be performed by venipuncture using collection tubes with heparin anticoagulant (to obtain plasma samples) or tube without anticoagulant (to collect serum), avoiding prolonged tourniquet. Samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, store for up to 3 months at -20 °C.*

Whole blood: Sample collection must be performed by venipuncture using collection tubes with heparin anticoagulant (to obtain whole blood samples), avoiding prolonged tourniquet. The samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 3 days*. **Do not freeze.**

*At the time of use, samples must be stabilized at room temperature.

PROCESS DESCRIPTION
Stability After Opening
The stability of the BIOCLIN POCT FIA PROCALCITONIN cassette (after opening the aluminized sachet) is 1 hour. This stability may vary according to the conditions of use of the product and the environment. However, to be on the safe side, we recommend that the BIOCLIN POCT FIA PROCALCITONIN kit be tested immediately after removing the cassette from the sachet.

TECHNIQUE
1- The Sample, Cassette and Sample Dilution Tube must be between 15-30 °C to run the test.
2- Check if the identification chip is inserted in the equipment.
3- Remove the cassette (test device) from its protective packaging, place it on a clean and level surface and identify it properly.
4- **For serum or plasma:** Transfer 50 µL of the sample to the tube containing the diluent (Reagent N° 2).
5- **For whole blood:** Transfer 80 µL of the sample to the tube containing the diluent (Reagent N° 2).
6- Homogenize the tube with the sample.
7- Pipette **80 µL of the diluted sample** into the sample well.
Note: Steps 4-7 must be completed within 1 minute.
8- Make sure there are no air bubbles.
9- **Wait 15 minutes to read the test. Do not perform the reading before 15 minutes.**

Cassette Reading
1- Select the test mode on the BIOCLIN POCT FIA reader. For more information on how to select the reading mode, refer to the equipment manual.
2- Insert the **Cassette** into the BIOCLIN POCT FIA reader and wait for the result to be issued.

NOTES
1- The cassette must be used immediately after opening its package.
2- The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

Calibration
The calibration of the method is carried out from the **Identification Chip**, which must be changed with each change in the lot of reagents.

PROCESS LIMITATIONS
1- It is important to use the correct sample volume, as lower or higher volumes may determine erroneous results.
2- The reaction reading time must be followed according to the established technique, in order to avoid false interpretation of the results.
3- The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERING
Samples with concentrations of up to 400µmol/L of Bilirubin, 20mmol/L of Triglycerides, 15mg/mL of Hemoglobin and 3250 IU/mL of Rheumatoid Factor do not interfere with the test result.

CROSS REACTIVITY
The BIOCLIN POCT FIA PROCALCITONIN kit was evaluated for cross-reactivity with different analytes. No falsely elevated results were observed for concentrations up to 100 g/L albumin, 100 ng/mL human calcitonin, 10 µg/mL alpha calcitonin gene-related peptide (α-CGRP) and 10 µg/mL beta peptide related to the calcitonin gene (β-CGRP). Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis must consider the clinical data of the patient with other laboratory data.

INTERNAL QUALITY CONTROL
The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

REFERENCE VALUES
Patients classified in the normal group had Procalcitonin concentrations below 0.5 ng/mL.

These values should be used as guidelines, and each laboratory should create its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY
Comparison of Methods and Methodological Specificity
The BIOCLIN POCT FIA PROCALCITONIN kit was compared with another commercially available method for measuring procalcitonin. 138 analyzes were performed and the results were evaluated. The linear equation obtained was Y= 0.97x - 0.02 and the correlation coefficient 0.995. With these results, it can be concluded that the kit has good methodological specificity.

PRECISION
Repeatability
Repeatability was calculated from 10 successive determinations, using 4 samples with different concentrations, obtaining the following results:

Repeatability	Sample			
	1	2	3	4
Concentration (ng/mL)	0.8	25.1	50.3	73.2
Standard Deviation (ng/mL)	0.1	1.9	4.1	2.3
Coefficient of Variation (%)	9	8	8	3

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

Reproducibility
The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 4 samples with different concentrations, obtaining the following results:

Reproducibility	Sample			
	1	2	3	4
Concentration (ng/mL)	0.8	26.3	51.2	68.9
Standard Deviation (ng/mL)	0.1	2.0	3.1	3.2
Coefficient of Variation (%)	7	8	6	5

The product presents results within the expected variation range (coefficient of variation ≤ 15%).

SENSITIVITY
The analytical sensitivity of the BIOCLIN POCT FIA PROCALCITONIN Kit is 0.8 ng/mL.

PROZONE EFFECT
There was no prozone effect up to a concentration of 500 ng/mL.

LINEARITY
The operating range for the BIOCLIN POCT FIA PROCALCITONIN reagent is 0.1 ng/mL to 100 ng/mL.

DIAGNOSTIC SIGNIFICANCE
Procalcitonin (PCT) is a precursor peptide of calcitonin synthesized by thyroid C cells, liver, kidneys, lungs, intestine and leukocytes. Procalcitonin has a molecular weight of 13000 Dalton, consisting of 116 amino acids. PCT is stable in vivo, with a half-life of 25 to 30 hours. Normally, PCT levels in the blood are low, however, in bacterial infections, an increase in this protein is observed. PCT levels can be used to monitor patients with severe bacterial infections with suspected sepsis, allowing assessment of treatment efficacy. In addition, there is an increase in PCT levels related to non-infectious causes such as: trauma, burns, carcinomas, cardiogenic shock, cirrhosis, in patients undergoing peritoneal dialysis and immunomodulatory therapy.

BIBLIOGRAPHIC REFERENCES

1. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
2. Oberhoffer M, Stonans I, Russwurm S, et al. Procalcitonin expression in Human PeriPheral blood mononuclear cells and its mudalution by lipopolysaccharides and sepsis related cytokines in vitro[J]. Lab ClinMed. 1999, 134(1): 49-55.
3. Tranca, SD;et al;. Biomarkers in polytrauma induced systemic inflammatory response syndrome and sepsis - a narrative review. Rom J Anaesth Intensive Care. 2014 ;21(2):118- 122.
4. Castro,CM; Silva,RM. Procalcitonina no diagnóstico de doenças infecciosas.Rev Bras Clin Med, São Paulo, 2012,10(6):521-8..
5. Tate J, Ward G. Interferences in immunoassay. Clin Biochem Rev. 2004;25(2):105-120.
6. Levinson SS.The Nature of Heterophilic Antibodies and the Role in Immunoassay Interference[J].J of Clin Immunoassay,1992,15:108-114.
7. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 – Santa Branca
CEP 31565-130 – Belo Horizonte – MG – Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 – Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for kit BIOCLIN POCT FIA PROCALCITONIN:
10269360429

Review: June/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER