Bioclin

BIOCLIN POCT FIA TROPONINA I-US

REF K331

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FINALIDADE

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- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- Apple, F. S. IFCC educational materials on selected analytical and clinical applications of high sensitivity cardiac troponin assays. Clinical Biochemistry, v. 48, p. 201-203, 2015.
- Thygesen, K. et al. Fourth Universal Definition of Myocardial Infarction (2018). Circulation, v. 138, p. e618-e651, 2018.
- Neumann, J. T. Application of High-Sensitivity Troponin in Suspected Myocardial Infarction. The New England Journal of Medicine, v. 380, n. 26, p. 2529-2540, 2009.
- Apple, F. S. Implementation of High-Sensitivity and Point-of-Care Cardiac Troponin Assays into Practice: Some Different Thoughts. Clinical Chemistry, v. 67, n. 1, p. 70-78, 2021.
- Lazar, D. R. High-Sensitivity Troponin: A Review on Characteristics, Assessment, and Clinical Implications. Disease Markers, v. 2022.
- Chapman, A. R. High-Sensitivity Cardiac Troponin and the Universal Definition of Myocardial Infarction. Circulation, v. 141, p. 161-171, 2020.
- Twerenbold, R. Clinical Use of High-Sensitivity Cardiac Troponin in Patients with Suspected Myocardial Infarction. Journal of the American College of Cardiology, v. 70, n. 8, 2017.
- Bioclin – Dados de arquivos.

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- 1- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 2- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- 3- Apple, F. S. IFCC educational materials on selected analytical and clinical applications of high sensitivity cardiac troponin assays. Clinical Biochemistry, v. 48, p. 201-203, 2015.
- 4- Thygesen, K. et al. Fourth Universal Definition of Myocardial Infarction (2018). Circulation, v. 138, p. e618-e651, 2018.
- 5- Neumann, J. T. Application of High-Sensitivity Troponin in Suspected Myocardial Infarction. The New England Journal of Medicine, v. 380, n. 26, p. 2529-2540, 2009.
- 6- Apple, F. S. Implementation of High-Sensitivity and Point-of-Care Cardiac Troponin Assays into Practice: Some Different Thoughts. Clinical Chemistry, v. 67, n. 1, p. 70-78, 2021.
- 7- Lazar, D. R. High-Sensitivity Troponin: A Review on Characteristics, Assessment, and Clinical Implications. Disease Markers, v. 2022.
- 8- Chapman, A. R. High-Sensitivity Cardiac Troponin and the Universal Definition of Myocardial Infarction. Circulation, v. 141, p. 161-171, 2020.
- 9- Twerenbold, R. Clinical Use of High-Sensitivity Cardiac Troponin in Patients with Suspected Myocardial Infarction. Journal of the American College of Cardiology, v. 70, n. 8, 2017.
- 10- Bioclin – Datos de archivos.

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INSTRUCTIONS FOR USE

FUNCTION

Test for ultrasensitive quantitative determination of cardiac troponin I (cTnI) in human biological samples of serum, plasma and whole blood, through fluorescent immunoassay. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Fluorescent immunoassay.

The BIOCLIN POCT FIA TROPONIN I-US kit is a fluorescent immunoassay based on the "sandwich" type detection principle, intended for the ultra-sensitive quantitative determination of cardiac troponin I in human samples of serum, plasma or whole blood (collected in EDTA or heparin). When the sample is added at its application site, an interaction occurs between the troponin I (cTnI) molecules present in the sample and the anti-troponin I antibodies conjugated to fluorescent microspheres, forming the cTnI-antibody conjugated complex. As the sample migrates by capillarity through the nitrocellulose membrane, this complex is captured by the test line, which consists of anti-troponin I antibodies. The sample will continue to migrate through the membrane until reaching the control line, which consists of anti-rabbit IgG antibodies (where rabbit IgG antibodies conjugated to fluorescent microspheres will be captured). At the end of the incubation time, the reading of the cassette is carried out with the immunofluorescence reader. The higher the concentration of troponin I present in the sample, the more complexes are immobilized in the test line and, consequently, the higher the fluorescein signal is generated and captured by the reader. The troponin concentration in the sample, in ng/mL, is calculated according to the assay calibration curve (contained on the Identification Card).

REAGENTS

1- Cassette - Store between 2 and 30 °C. Each sachet contains:

A) 01 plastic device with 01 test strip containing:

- Test line region: mouse anti-cardiac troponin I antibody.

- Control line region: goat anti-rabbit IgG antibody.

- Conjugate: fluorescent microspheres conjugated to mouse anti-cardiac troponin I antibodies and fluorescent microspheres conjugated to rabbit IgG antibodies.

B) 01 silica sachet.

2- Whole Blood Diluent - Store between 2 and 30 °C. Contains: 10 mM phosphate-buffered saline (PBS), stabilizers (casein <1%) and preservatives (ProClin™ 300 <0.5%).

3- Identification Card

PRESENTATION

REAGENTS	NUMBER OF REAGENTS PER PACKAGE
1- Cassette	25 units
2- Whole Blood Diluent	25 units (25 x 0.05 mL)
3- Identification Card	1 unit

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

1- Pipettes capable of dispensing volumes from 5 to 100 µL with a coefficient of variation lower than 1.5%.

2- BIOCLIN POCT FIA Immunofluorescence Reader.

3- Stopwatch or clock.

4- Quality Control Tools.

STORAGE AND TRANSPORTATION CONDITIONS

The storage and transportation temperature should be between 2 °C and 30° C. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- BIOCLIN POCT FIA TROPONIN I-US should only be used in conjunction with the BIOCLIN POCT FIA Immunofluorescence Reader.

3- Strictly follow the proposed methodology to obtain accurate results.

4- The cassettes (test devices) do not contain contaminating substances, but since the samples used are biological, safety routines must be strictly followed and the use of disposable gloves is essential.

5- Open the aluminum envelope to remove the cassette only at the time of use.

6- The kit should not be frozen or exposed to high temperatures (above 45°C), as this will cause deterioration.

7- The reagents, as well as the samples, must reach room temperature before starting the test.

8- Do not mix reagents from kits with the same batch or with different batch numbers.

9- Do not use the kit after its expiration date. The expiration date is indicated on the product label.

10- The test device (R1- Cassette) and sample dilution tubes (R2- Whole Blood Diluent) are disposable and cannot be reused.

11- The kit materials that have come into contact with the patient's sample must be disposed of in a special way. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biological waste.

12- Secondary packaging (box), aluminum envelope and silica sachet can be disposed of in regular waste.

13- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.

14- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data

Sheet) available on the website www.bioclin.com.br or by requesting it through Quibasa's SAC (Customer Advisory Service).

15- Do not use the product if the packaging is damaged.

16- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Samples should be collected by venipuncture, in tubes with or without anticoagulant (EDTA or Heparin). To avoid hemolysis, centrifuge the samples to obtain serum/plasma and test them as soon as possible. Do not keep samples at room temperature for long periods. Serum and plasma samples can be stored refrigerated at 2–8°C for a maximum of 3 days. If samples cannot be analyzed within 3 days, they can be stored for up to 30 days at -20°C. Serum and plasma samples stored for longer than the recommended periods and samples collected with anticoagulants other than those mentioned should not be used. Highly hemolyzed, lipemic, and precipitate-containing samples should be avoided, as they may provide inconsistent results.

Whole Blood (EDTA or Heparin)

Whole blood samples should be collected by venipuncture in tubes containing heparin or EDTA as an anticoagulant. Samples should be tested immediately or as soon as possible after collection. Whole blood collected in tubes containing anticoagulants can be stored refrigerated at 2–8°C for a maximum of 2 days. Do not use samples collected with anticoagulants other than those mentioned, and do not freeze whole blood samples.

PROCESS DESCRIPTION

Stability After Opening

The stability test results show that the BIOCLIN POCT FIA TROPONIN I-US kit is stable after opening for up to 70 minutes. This stability may vary depending on the conditions of use of the product and the environment. However, for safety reasons, we recommend that the product should be used immediately after removing the cassette from the sachet.

TECHNIQUE

Check the BIOCLIN POCT FIA Immunofluorescence Reader manual for complete instructions on handling the equipment.

ATTENTION! The techniques for serum/plasma and whole blood are different. Carefully read the instructions for use described below, according to the type of sample to be used.

Serum and Plasma Samples

Before starting the assay, allow the Cassette and samples to stabilize at room temperature (15-30°C) for at least 40 minutes.

1- Check if the information contained in the Identification Card has been properly entered into the equipment (with each new batch of the kit, the calibration and validity data contained in the Identification Card must be transferred to the equipment).

2- Remove the cassette (test device) from its protective packaging, place it on a clean, level surface and identify it appropriately.

3- Pipette 50 µL of sample (serum or plasma in EDTA or Heparin) directly into the cassette's sample application port.

4- Option 1: Insert the cassette into the equipment immediately after applying the sample and select the "Timing Test" option. After 15 minutes, the equipment will take the reading, and the result will appear on the display.

Option 2: After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the "Quick Test" option. The equipment will take the reading immediately, and the result will appear on the display.

Whole Blood Samples

Before starting the test, allow all reagents (Cassette and Whole Blood Diluent) and samples to stabilize at room temperature (15-30 °C) for at least 40 minutes.

1- Check if the information contained in the Identification Card has been properly entered into the equipment (with each new batch of the kit, the calibration and validity data contained in the Identification Card must be transferred to the equipment).

2- Remove the cassette (test device) from its protective packaging, place it on a clean, level surface and identify it appropriately.

3- Make sure that the Whole Blood Diluent solution (Reagent No. 2) is at the bottom of the dilution tube.

4- Add 50 µL of sample to a dilution tube containing the Whole Blood Diluent.

5- Close the tube cap and mix by inverting the sample 6 to 8 times until the sample and diluent are completely homogenized.

6- Pipette 50 µL of diluted and properly mixed sample into the sample application hole of the cassette.

7- Option 1: Insert the cassette into the equipment immediately after applying the sample and select the "Timing Test" option. After 15 minutes, the equipment will take the reading and the result will appear on the display.

Option 2: After adding the sample, wait 15 minutes. After the time has elapsed, insert the cassette into the equipment and select the "Quick Test" option. The equipment will take the reading immediately and the result will appear on the display.

RESULTS

The BIOCLIN POCT FIA Immunofluorescence Reader calculates the cardiac troponin I concentration automatically, using the calibration curve contained in the R3 (Identification Card). The information on the Identification Card must be properly entered into the equipment for each new batch of reagents.

Quantitative results are measured in ng/mL. The measuring range of the kit is between 0.03 and 20 ng/mL (limit of quantification and highest point of curve, respectively). Results below the lower limit of detection are expressed as <0.03 ng/mL. Results above the upper limit of detection are expressed as >20 ng/mL.

PROCESS LIMITATIONS

1- It is important to use the correct sample volume, as lower or higher volumes may produce erroneous results.

2- The reaction reading time must be followed according to the established technique to avoid misinterpretation of the results.

3- The BIOCLIN POCT TROPONIN I-US kit for cTnI is a

specific test for aiding in the diagnosis of acute myocardial infarction. Reactive results must be confirmed with another test and available clinical information. Test results alone cannot be used for a definitive diagnosis.

4- It is recommended that established guidelines for the diagnosis of acute myocardial infarction be followed.

5- The results provided by this kit should be interpreted by the responsible medical professional and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

REFERENCE VALUES

Normal	< 0.04 ng/mL
--------	--------------

These values should be used as a guide, and each laboratory should create its own range of reference values according to the population served. The results provided by this kit should be interpreted by the responsible medical professional, and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

Cardiac troponin I assays are defined as highly sensitive when they detect cTnI levels in >50% of the apparently healthy population (i.e., detect troponin levels above the method's Limit of Quantification - 0.03 ng/mL) and have a coefficient of variation <10% at the 99th percentile of samples from healthy individuals. Of 160 samples from healthy individuals analyzed, 80 had values ≥ 0.03 ng/mL, representing 50% of the sample. Furthermore, at the 99th percentile (0.04 ng/mL), a CV of 8.4% was observed. Therefore, the BIOCLIN POCT FIA TROPONIN I-US kit is considered a highly sensitive assay

INTERFERENTS

No interference was observed for bilirubin up to 350 $\mu\text{mol/L}$, triglycerides up to 40 mmol/L, hemoglobin up to 2 g/L, rheumatoid factor up to 1,500 IU/mL, heterophile antibodies up to 1,000 ng/mL, and antinuclear antibodies up to 1:640.

CROSS-REACTIVITY

No false-positive results were observed for cardiac troponin C up to 1,000 ng/mL, cardiac troponin T up to 500 ng/mL, and skeletal troponin T up to 1,000 ng/mL. Despite these results, the possibility of cross-reactivity with other analytes or higher concentrations cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. For this purpose, it is recommended to use controls, which allow the assessment of the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The BIOCLIN POCT FIA TROPONIN I-US kit was compared with another commercially available method for troponin I measurement. 86 analyses were performed and the results were evaluated. The linear equation obtained was $Y = 0.9663x + 0.1184$ and the correlation coefficient was 0.9742. Based on these results, it can be concluded that the kit presents good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following concentration results:

REPEATABILITY	SAMPLE		
	1	2	3
Mean	0.10	0.97	5.05
Standard Deviation	0.01	0.05	0.38
Coefficient of Variation (%)	6.4	5.5	7.6

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following concentration results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Mean	0.10	0.99	5.02
Standard Deviation	0.0	0.1	0.4
Coefficient of Variation (%)	7.9	8.3	8.3

ANALYTICAL SENSITIVITY

The Limit of Blank (LoB) of the BIOCLIN POCT FIA TROPONIN I-US product is 0.004 ng/mL, the Limit of Detection (LoD) is 0.01 ng/mL and the Limit of Quantification (LoQ) is 0.03 ng/mL.

LINEARITY

The BIOCLIN POCT FIA TROPONIN I-US kit is linear up to 20 ng/mL (highest point of the calibration curve).

HIGH DOSE PROZONE EFFECT

No high-dose prozone effect was observed at up to 800 ng/mL of troponin I.

DIAGNOSTIC SIGNIFICANCE

Troponins, classified as troponin C (TnC), troponin T (TnT), and troponin I (TnI), are proteins that form protein complexes present in the muscle fibers of skeletal and cardiac muscles and play important roles in muscle contraction in these tissues. While TnC is similar in both types of muscle tissue, troponins T and I have isoforms specific to cardiac muscle fibers. This characteristic is important in the use of TnT and TnI as specific markers of cardiac injury, especially in the diagnosis of acute myocardial infarction (AMI). AMI is a major cause of death and disability, and rapid diagnosis is critical for initiating effective treatment and improving prognosis. It is characterized by characteristic symptoms (such as chest pain and dyspnea), elevation of specific markers in the blood, and the development of a Q-wave pattern on the electrocardiogram. Acute myocardial infarction occurs due to necrosis of cardiac tissue associated with prolonged ischemia. This myocardial necrosis results in the release of several markers into the bloodstream, such as myoglobin, cardiac troponins T and I, creatine kinase, and lactate dehydrogenase. Cardiac troponin I is an important marker because its blood concentration rises 3–6 hours after the infarction, peaking between 12–24 hours, and maintaining altered levels for 5–14 days after the AMI episode. Although troponin T is also cardiospecific, troponin I measurement

is preferred because it often presents higher levels and rises more rapidly after AMI. The latest troponin detection assays, called ultrasensitive or high-sensitivity, are capable of detecting lower levels than previous-generation assays, providing greater sensitivity for detecting AMI upon patient admission, speeding diagnosis, and improving efficiency in emergency departments. Although useful in diagnosing AMI, troponin levels rise in cases of non-ischemic myocardial injury, such as myocarditis, or secondary to other diseases, such as renal failure.

BIBLIOGRAPHIC REFERENCES

- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Tate, J.; Ward, G. Interferences in immunoassay. Clin. Biochem. Rev., v. 25, n. 2, p. 105-120, 2004.
- Apple, F. S. IFCC educational materials on selected analytical and clinical applications of high sensitivity cardiac troponin assays. Clinical Biochemistry, v. 48, p. 201-203, 2015.
- Thygesen, K. et al. Fourth Universal Definition of Myocardial Infarction (2018). Circulation, v. 138, p. e618-e651, 2018.
- Neumann, J. T. Application of High-Sensitivity Troponin in Suspected Myocardial Infarction. The New England Journal of Medicine, v. 380, n. 26, p. 2529-2540, 2009.
- Apple, F. S. Implementation of High-Sensitivity and Point-of-Care Cardiac Troponin Assays into Practice: Some Different Thoughts. Clinical Chemistry, v. 67, n. 1, p. 70-78, 2021.
- Lazar, D. R. High-Sensitivity Troponin: A Review on Characteristics, Assessment, and Clinical Implications. Disease Markers, v. 2022.
- Chapman, A. R. High-Sensitivity Cardiac Troponin and the Universal Definition of Myocardial Infarction. Circulation, v. 141, p. 161-171, 2020.
- Twerenbold, R. Clinical Use of High-Sensitivity Cardiac Troponin in Patients with Suspected Myocardial Infarction. Journal of the American College of Cardiology, v. 70, n. 8, 2017.
- Bioclin – Archive Files.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration number for BIOCLIN POCT FIA TROPONIN I-US kit: 10269360475

Review: August/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER