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REF K137

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FINALIDADE

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1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. National Committee For Clinical Laboratory Standards: Colleccion, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
4. International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. National Committee For Clinical Laboratory Standards: Coleccion, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
4. International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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USAGE INSTRUCTIONS

FUNCTION

Reagent used as control in coagulation tests. For *in vitro* diagnostic use only.

REAGENT

Number 1 – Plasma Control Normal - store between 2 and 8°C.

Contains: Plasma control containing 3 analytes incorporated in a lyophilized protein matrix and preservative. **Potentially Infectious.**

Attention: the analyte concentration vary on each lot refer to table sent along with the product..

PRESENTATION

Presentation	Volume
1	1 x 1 mL
2	3 x 1 mL
3	5 x 1 mL
4	10 x 1 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Automated and semi-automated biochemistry analyzers, watches or stopwatches, pipettes, test tubes, water bath at 37°C and Bioclin's Fibrinogen, TP and TTPA kits. Can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transportation, at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Follow strictly the proposed methodology to obtain results accurate.

3- The water used in material cleansing must be recent and free of contaminants.

4- The reagent must be handle cautiously because it is susceptible to biologic contamination.

5- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the product is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).

6- Handle with caution, the reagent contain sodium azide, irritating to mucous and skin.

7- The reagent was tested for anti-HIV, anti-HCV antibodies and antigen HBs using methods of recent

latest generation and present negative results. The risk of infection can not be excluded and the reagent must be handle with the same caution taken for the patient sera.

Potentially infectious.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

PROCESS DESCRIPTION

Reconstitution of the Plasma Control

Carefully open the bottle and add 1 mL of deionized water.

Mix. Leave the bottle resting for 30 minutes. On every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid foam formation.

Attention: The reconstituted control is stable for 6 hours between 2 and 8°C. After reconstitution, freezing is not recommended.

REFERENCE VALUES

The values in the table enclosed were obtained through repeated measurements using the indicated method.

The values provided should be used for guidelines, each laboratory should establish its own limits of precision.

The laboratory's average should, however, be within the acceptable average for the given parameters

See Reference Values table enclosed.

PRODUCT PERFORMANCE

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Analyte	Method	Average	SD	CV	Un.
Fibrinogen	Coagulation	252.05	0.60	0.24	mg/dL
TP	Coagulation	12.15	0.49	4.03	seconds
TTPA	Coagulation	32.15	0.59	1.83	seconds

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Analyte	Methodology	Average	SD	CV	Un.
Fibrinogen	Coagulation	264.06	0.56	0.21	mg/dL
TP	Coagulation	11.35	0.49	4.31	seconds
TTPA	Coagulation	33.70	0.47	1.40	seconds

BIBLIOGRAPHIC REFERENCES

1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. National Committee For Clinical Laboratory Standards: Collection, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
4. International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

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CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA resgistration for the Biocontrol N Coagulation Kit: 10269360223.

Review: August/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER