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REF K138

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FINALIDADE

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1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. National Committee For Clinical Laboratory Standards: Colleccion, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
4. International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. National Committee For Clinical Laboratory Standards: Coleccion, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
4. International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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INSTRUCTIONS FOR USE

FUNCTION

Reagent used as control in coagulation tests. For *in vitro* diagnostic use only.

REAGENT

Number 1 - Plasma Control Pathological - store between 2 and 8°C.

Contains: Plasma control containing 3 analytes incorporated in a lyophilized protein matrix and preservative. **Potentially Infectious.**

Attention: the analyte concentration vary on each lot. Refer to table sent along with the product.

PRESENTATION

Presentation	Volume
1	1 x 1 mL
2	3 x 1 mL
3	5 x 1 mL
4	10 x 1 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Automated and semi-automated biochemistry instruments, watches or stopwatches, pipettes, test tubes, water bath at 37°C and Bioclin's Fibrinogen, TP and TTPA kits. Can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transportation at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The water used to clean the material must be recent and free of contaminating agents.
- 4- The reagent must be handled with caution, as it is liable to biological contamination.
- 5- To avoid possible analytical errors, it is advisable that the distilled water used in the reconstitution of the product has adequate quality. For this, it is recommended that the water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- 6- Handle with care, as the reagent contains Sodium Azide, which is irritating to the skin and mucous membranes.

7- The reagent was tested for anti-HIV, anti-HCV and HBs antigen antibodies using state-of-the-art methods and showed negative results. The risk of infection cannot be excluded and the reagent must be handled with the same care as observed for patient serum. **Potentially infectious.**

8- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product if the packaging is damaged.

PROCESS DESCRIPTION

Reconstitution of the Plasma Control

Carefully open the bottle and add 1mL of deionized water. Mix. Leave the bottle resting for 30 minutes. On every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid foam formation.

Caution: The reconstituted Control is stable for 6 hours at 2 to 8°C. After freezing at -20°C, the Control is stable for up to 72 hours (3 days).

REFERENCE VALUES

The values in the table enclosed were obtained through repeated measurements using the indicated method. The values provided should be used for guidelines, each laboratory should establish its own limits of precision. The laboratory's average should, however, be within the acceptable average for the given parameters.

See Reference Values table enclosed.

PRODUCT PERFORMANCE

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Analyte	Method	Average	SD	CV	Un.
Fibrinogen	Coagulation	482.50	0.61	0.13	mg/dL
TP	Coagulation	34.20	0.52	1.53	seconds
TTPA	Coagulation	57.55	0.60	1.05	seconds

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Analyte	Methodology	Average	SD	CV	Un.
Fibrinogen	Coagulation	484.50	0.76	0.16	mg/dL
TP	Coagulation	32.64	0.67	2.05	seconds
TTPA	Coagulation	57.10	0.55	0.97	seconds

BIBLIOGRAPHIC REFERENCES

1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. National Committee For Clinical Laboratory Standards: Collection, Transport And Processing Of Blood Specimens For Coagulation Test And Performance Of Coagulation Assays. Approved Guideline Fourth Edition. NCCLS document H12-A4, 2003.
4. International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
5. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for the Biocontrol Coagulation Family kit: 10269360223

Review: May/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER