BIOCONTROL HEMATO 3P

REF K142

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FINALIDADE

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1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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USAGE INSTRUCTIONS

FUNCTION

Reference hematology control used to monitor the performance of automated hematology analyzers. Only for *in vitro* diagnostic use.

REAGENT

Number 1 - Level 1 - Store between 2 and 8°C. Contains: suspension with white cells, platelets, stabilized erythrocyte, dextrose and preservatives. **Potentially infective.**

Number 2 - Level 2 - Store between 2 and 8°C. Contains: suspension with white cells, platelets, stabilized erythrocyte, dextrose and preservatives. **Potentially infective.**

Number 3 - Level 3 - Store between 2 and 8°C. Contains: suspension with white cells, platelets, stabilized erythrocyte, dextrose and preservatives. **Potentially infective.**

Attention: The concentration of parameters varies with each batch.

PRESENTATION

Presentation	Nivel 1	Nivel 2	Nivel 3
1	1 x 3.0 mL	1 x 3.0 mL	1 x 3.0 mL
27	1 x 3.0 mL	---	---
28	---	1 x 3.0 mL	---
29	---	---	1 x 3.0 mL

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. Do not freeze.

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the proposed methodology to obtain exact results.

3- The water used to clean the material must be recent and free from contaminating agents.

4- The control must be handled carefully, as it is susceptible to biological contamination.

5- The control was tested for anti-HIV antibodies, anti-HCV and HBs antigen using state-of-the-art methods and showed negative results.

The risk of infection cannot be excluded and the reagent should be handled with the same care as for patient serum. **Potentially Infectious.**

6- If the test results are outside the upper and lower limits shown in the quality control charts, this indicates reagent deterioration or a problem with the analyzer. In this case, perform the dosage again. If the test results are still outside the limits, use a new control bottle. If the problem persists, check the analyzer for system errors.

7- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the FISPQ (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Safety Service Customer Advisory).

9- Do not use the product in case of damage to the packaging.

10- It is essential that the instruments and equipment used are properly calibrated and subject to periodic maintenance.

PROCESS DESCRIPTION

Control ready to use.

Before using the control, remove it from the refrigerator and let it stabilize at room temperature (15 - 30°C) for 10 to 30 minutes. Homogenize gently for 20 to 30 seconds, avoiding the formation of bubbles.

Repeat the homogenization 20 times.

Carry out the test according to the equipment instruction manual.

After the test, clean the vial, replace the cap and leave it stored at a temperature of 2 to 8°C.

Once opened, the control is stable for 14 days.

REFERENCE VALUES

All the values assigned to the parameters are obtained through repeated dosages of the Control in the hematology analyzer calibrated by reference methods. The values provided should be used as a guide and each laboratory should establish its own limits of accuracy. The laboratory average must, however, be within the acceptable average for the parameters provided.

BIBLIOGRAPHIC REFERENCES

- Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsgesellschaft; 1998. p. 1393-1401
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiration date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of the Biocontrol Hemato 3P kit at ANVISA:
10269360230.

Review: September/2022

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER