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## VETERINARY USE

**REF** K074

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

## USAGE INSTRUCTIONS LEVEL P - PATHOLOGICAL

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

### FUNCTION

Reagent used to monitor the accuracy and precision of analytical testing in manual and automated methodologies. Has the concentration of analytes representing the pathological parameters. For use in *In vitro* Diagnostics only.

### REAGENTS

**Number 1 - Control Serum** - Store between 2 and 8 °C. Contains: Control serum containing approximately 30 analytes embedded in a lyophilized matrix and preservative. **Potentially infectious.**

**Warning:** The concentration of analytes varies each lot. See table that accompanies the product.

### PRESENTATION

Presentation	Volume
K074-3-VET	1 x 5 mL

### EQUIPMENTS AND OPERATIONAL INPUTS

Automated or semi-automated biochemistry analyzers, watches or stopwatches, pipettes, test tubes, kits and water bath at 37 °C. These items are found at markets specialized on Clinical Analysis Laboratories.

### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8 °C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture.

### SPECIAL CARE

- 1- **For *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- The reagent should be handled cautiously since it is subject to biological contamination.
- 5- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the product is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- 6- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- 7- The reagent has a contaminating biological matrix. Although it was negative for the main human and veterinary contaminating pathogens, the risk of infection cannot be excluded and the reagent should be handled with the same care observed for patient serum. **Potentially Infectious.**
- 8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 9- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

- 10- Do not use the product in case of damaged packaging.
- 11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

### PROCESS DESCRIPTION

#### RECONSTITUTION OF CONTROL SERUM

**Before reconstitute the Control Serum should be noted the volume indicated on the bottle and proceed as described below.**

**Presentation 1 and 2:** Open carefully the bottle and add exactly **3 mL of distilled water**, as volume described in the label of the bottle.

**Presentation 3 and 4:** Open carefully the bottle and add exactly **5 mL of distilled water**, as volume described in the label of the bottle.

Mix. Leave the bottle resting for 30 minutes. Every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid formation of foam.

After mixing, observe the stability according to the table below.

	-20 °C	+4 °C	+25 °C
<b>Bilirubin (protected from the light)</b>	14 days	6 hours	2 hours
<b>CK and CK-MB</b>	30 days	6 hours	2 hours
<b>Other parameters</b>	30 days	5 days	12 hours

#### Warning:

- The reconstituted Control Serum should not be thawed more than once. Therefore, for minor routines, fractionate the material into aliquots before freezing.
- The Control Serum should be used in an identical way as patient's serum.
- After pipetting, if there is leftover serum control, do not return to the original bottle.

#### REFERENCE VALUES

The values enclosed were obtained through repeated measurements using the method indicated. The values provided should be used for orientation, each laboratory should establish its own limits of precision. The laboratory's average should, however, be within the acceptable average for the given parameters.

**See table that accompanies the product.**

## BIBLIOGRAPHIC REFERENCES

1. Röhle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt:TH-Books Verlagsge-sellschaft; 1998. p. 1393-1401.
2. Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services, Washington 1993 (HHS Publication No. [CDC] 93-8395).
3. Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien. Deutsches Ärzteblatt 1988;85: B519-B532.
4. Richtlinie der Bundesärztekammer zur Qualitätssicherung quantitativer laboratoriums-medizinischer Untersuchungen. Deutsches Ärzteblatt 2002;98:A 2747-59.
5. QUIBASA: Dados do departamento de Pesquisa e Desenvolvimento.

## QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

## CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

**Registration Number:** Product exempt from registration on the Ministry of Agriculture, Livestock and Supply.

**Technical Responsible:** Camila Eckstein (CRMV/MG 20.611)

**Review:** June/2025

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER