BIOLÁTEX PCR

REF K044

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INSTRUCTIONS FOR USE

FUNCTION

Method for determination qualitative and semiquantitative of C Reactive Protein (CRP) in biological samples or serum, by agglutination of latex particles, without previous dilution of the sample. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Latex

This method is based on a particle agglutination reaction of latex particles covered with Gamma-Globulin anti-CRP, especially treated to prevent nonspecific agglutination. Agglutination is visible in samples with concentrations of CRP equal to or greater than 6 mg/L, according to the references established by the WHO International Standards.

REAGENTS

- Number 1 - Latex PCR** - Store between 2 and 8°C. **Do not freeze.** Contains: Latex particles in suspension sensitized with Antibody Anti-PCR. Homogenize before using it.
- Number 2 - Positive Control** - Store between 2 and 8°C. Contains: Latex Binder Solution, 0.9% Saline Solution, stabilizer and preservative. **Potentially infectious.**
- Number 3 - Negative Control** - Store between 2 and 8°C. Contains: 0.9% Saline solution and preservative.

PRESENTATION

Reagent	1 (Economic Packing)	2 (Normal Packing)
Nº 1	2 mL	2 mL
Nº 2	--	1 mL
Nº 3	--	1 mL

The normal package is accompanied by test cards.

EQUIPMENTS AND OPERATIONAL INPUTS

Cards, slides or dark bottom plates, spatulas, pipettes and watch or stopwatch. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- 2-** Strictly follow the methodology proposed to obtain exact results.
- 3-** Water used in material cleaning must be recent and free of contaminants.
- 4- Do not freeze reagents.**
- 5-** Always use reagents from same lot.
- 6-** Do not use lipemic serum. Do not use plasma.

- 7-** For cards reuse, wash immediately after use in distilled or deionized water until all residue is removed. If the washing can not be immediate, wash the cards with neutral detergent and rinse thoroughly in distilled or deionized water until all the residue is removed.
- 8-** Handle with care Reagents Nº 2 and 3, they contain Sodium Azide.
- 9-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 10-** To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 11-** Do not use the product in case of damaged packaging.
- 12-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum, with no previous dilution. The analyte is stable for 8 days between 2 and 8°C or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing.

PROCESS DESCRIPTION

TECHNIQUE

QUALITATIVE TEST

In each circle of the card, place the following:

	Circle Nº 1	Circle Nº 2	Circle Nº 3
Negative Control	20 µL	--	--
Positive Control	--	20 µL	--
Sample	--	--	20 µL
Reagent Nº 1*	20 µL	20 µL	20 µL

*Previously homogenized

Homogenize with the aid of a spatula using the entire length of each circle of the card. Then, agitate the card with circular motions for 2 minutes. Make immediately reading with artificial light.

A clear agglutination indicates the presence of C Reactive Protein at a concentration equal to or greater than 6 mg/L. In this case, perform the semiquantitative test.

SEMIQUANTITATIVE TEST

- 1-** Perform sample dilutions with saline, starting from initial sample (1:2, 1:4, 1:8, 1:16, 1:32, etc.).
- 2-** Follow the process as described in qualitative test for each of the dilutions.
- It will be considered as title the biggest serum dilution that presents agglutination.

RESULTS

- Positive:** Clear presence of agglutination.
- Negative:** Absence of agglutination (homogenized suspension).

CALCULATIONS

Sample	Concentration (mg/L)
No dilution	6
1/2	12
1/4	24
1/8	48
1/16	96
1/32	192

Results can be expressed in titles or in mg/L.

mg/L = 6 x title from last dilution (dilution nº)

Negative Tests: express the results as negative or minor than 6 mg/L.

PROCEDURE LIMITATIONS

Do not use plasma, hemolysate serum or limpemic, as they may produce nonspecific agglutination.

By correlating methods for determining C Reactive Protein, to verify the sensitivity of the reagents. The obtained results should only be compared when expressed in mg/L.

The detergent residue on the card can interfere with the test and generate a false result.

INTERFERENT

No interference was observed for Bilirubin up to 20 mg/dL, Hemoglobin up to 10 g/L and Lipids up to 10 g/L. Rheumatoid Factor up to 100 IU/mL does not interfere with the results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's traceability was determined through the reference material NIBSC 85/506 (Human C-Reactive Protein 1st International Standard - WHO International Standard).

REFERENCE VALUES

Up to 6 mg/L

These values should be used as guidance, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Biolatex CRP kit was compared with another method for confirmation of inflammatory process. According to results from 100 clinical samples, all methods showed excellent correlation. With these results we can conclude the kit has good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 3 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Obtained Results
01	20	Positive 1/4	Positive 1/4
02	20	Positive 1/64	Positive 1/64
03	20	Negative	Negative

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Obtained Results
01	20	Positive 1/4	Positive 1/4
02	20	Positive 1/64	Positive 1/64
03	20	Negative	Negative

ANALYTICAL SENSITIVITY

The analytical sensitivity study of the Kit Biolatex CRP was performed by dilution of a positive sample with a known concentration of 30 IU/mL for Anti-Streptolysin O. The analytical sensitivity was 6 IU/mL.

HIGH DOSE HOOK EFFECT

Hook effect wasn't observed with a high concentration of PCR up to 1600 mg/L.

DIAGNOSTIC SIGNIFICANCE

C Reactive Protein is a useful indicator of inflammatory process in activity, whether it is infectious (pneumonia, tuberculosis) or non-infectious (active rheumatic fever, rheumatoid arthritis, lupus erythematosus). Is also present in several other pathological conditions such as in acute myocardial infarction, neoplasms, intense trauma, viruses and burns. The determination of plasma concentrations constitutes an effective test to monitor the therapy used and prognosis of inflammation.

BIBLIOGRAPHIC REFERENCES

- 1 - WARWORTH, E. Wadsworth, Ch. Clin. Chim. Acta, 138, 1984.
- 2 - PEPYS, M. B.; DASH, A. C.; ASHLEY M. J., Clin. Exp. Immunol, 30, 32-37, 1977.
- 3 - DEYO, R. A.; POPE, R. M., PERSELLIN, R. H.; J. Rheumatol, 279, 1980.
- 4 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Biolatex CRP kit: 10269360093

Review: August/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER