BIOLISA ANTI HBs

^{REF} K121

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FINALIDADE

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1- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44. 2017 AACR.

2- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

3- Frank Fenner and David O. White, Medical Virology, 4th Edition, Academic Press, 1994.

4- Richman, D., R. Whitley, F. Hayden. Clinical Virology. New York: Churchill Livingstone Inc., 1997.

5- WHO. Guidelines on Hepatitis B and C Testing. February 2017.

6- WHO. Guidelines for the Prevention, Care and Tretatment of Persons with Chronic Hepatitis B Infection. March 2015.

7- WHO. WHO Expert Committee on Biological Standardization: Sixty-first report. 2013.

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INSTRUCTIONS FOR USE

FUNCTION

Test for the quantitative determination of antibodies against the surface antigen of the hepatitis B virus (anti-HBs) in biological samples of serum or plasma, through the enzyme immunoassay test. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The Biolisa Anti HBs kit is a solid-phase enzyme immunoassay based on the “sandwich” principle for the quantitative detection of antibodies against hepatitis B virus surface antigen (anti-HBs) in human serum or plasma samples. The microplate is coated with hepatitis B virus specific surface antigens (HBs). At first, the sample is incubated with the conjugate, made up of specific HBs antigens conjugated to peroxidase. Anti-HBs antibodies present in the sample bind to the specific HBs antigens immobilized on the plate and to the specific HBs antigens present in the conjugate, forming immune complexes. After the initial incubation, the microplate is washed to remove unbound materials. Then the substrate is added and incubated, producing a blue color that indicates the detection of anti-HBs antibodies in the sample. Stop solution is added to stop the reaction, with a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1- Reference Standard (A - E)** - Store between 2 and 8 °C. Contains: Five (5) vials of Reference Standards (A - E), containing anti-HBs antibodies in different concentrations, buffer solution, stabilizers, preservatives, surfactant and dye. **Concentrations vary from batch to batch (from 0 to 500 mIU/mL), see vial label. Potentially infectious.**
- 2- Conjugate** - Store between 2 and 8 °C. Contains: Buffer Solution containing Hepatitis B virus-specific surface antigens linked to Peroxidase, stabilizers, preservatives, surfactant and dye.
- 3- Sensitized Plate** - Store between 2 and 8 °C. Contains: Plate sensitized with specific surface antigens for Hepatitis B virus.
- 4- Concentrated Washing** - Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 5- Substrate** - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8 °C. Contains: 1N Sulfuric Acid.
- 7- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Reference Standards (A – E)	5 Vials x 1 mL	10 Vials x 1 mL	25 Vials x 1 mL
2- Conjugate	1 Vial x 6 mL	2 Vials x 6 mL	5 Vials x 6 mL
3- Sensitized Plate	1 Unit	2 Units	5 Units
4- Concentrated Washing	1 Vial x 50 mL	2 Vials x 50 mL	5 Vials x 50 mL
5- Substrate	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
6- Stop Solution	1 Vial x 12 mL	2 Vials x 12 mL	2 Vials x 12 mL
7- Plate Sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the previous item.

Materials needed but not contained in the Kit:

- 1- Pipettes capable of dispensing volumes of 5 to 500 µL with a variation coefficient of less than 1.5%.
- 2- Repipettor for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer.
- 4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5- Absorbent paper to dry the wells.
- 6- Stopwatch or watch.
- 7- Vial to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator at 37 °C ± 2 °C.

TRANSPORTATION AND STORAGE CONDITIONS

Storage temperature should be between 2 and 8 °C. Transport at temperatures up to 30°C should not exceed 5 days. Keep out of the light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The sachet containing the microplate must be opened only after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be fresh and free from contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Sulfuric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw material of the product is tested and must be non-reactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not offer complete security in the absence of infectious agents. The handling of any product that contains serum is potentially capable of transmitting disease. Therefore, it is necessary to take due care of biosafety in the handling of these products.
- 8- Pipette the reagents always in the same order to minimize the difference in reaction time between the wells.
- 9- As a protection measure, the plate must be covered during the reaction.
- 10- Make sure that the bottom of the cavity is clean and dry and that there are no bubbles on the liquid surface before reading the plate. Do not allow the wells to dry out during the test.
- 11- Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.
- 12- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service of Chemicals) Customer Service) of Quibasa.
- 14- Do not use the product in case of damage to the packaging.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolyzed and highly lipemic samples should not be used. The samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C. Plasma samples stored for periods longer than recommended and samples collected with anticoagulants other than those mentioned should not be used.

Decreased anti-HBs concentration values can be obtained in plasma samples collected in EDTA and Heparin compared to serum samples.

DESCRIPTION OF PROCESS

Stability After Open

The results of the stability test show that the Biolisa Anti HBs kit is stable after being opened up to 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using internal controls of the kit and the criteria for validating the technique.

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of Reagent N° 4 (Concentrated Washing) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

Substrate

The Substrate is ready for use.

TECHNIQUE

For use in automatic equipment, consult the Customer Service Department (SAC).

Before starting the assay, allow all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

- 1- Separate the cavities to be used considering: Reference Standards (from A to E) and Samples. It is recommended to test in duplicate. Return unused microplate strips to the original sealed package.
- 2 - Separate the first cavity for the Blank (OPTIONAL).
- 3- Pipette 50 µL of each of the Reference Standards (A-E) and Samples into the previously determined wells.
- 4- Pipette 50 µL of Conjugate into all wells, **except** for the Blank well.
- 5- Homogenize gently for ± 10 seconds. Cover the wells with plate sealer.
- 6- Incubate for 60 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 7- Remove the sealer from the cavities.
- 8- Discard the contents of the wells by aspiration in a Microplate Washer. Use approximately 500 µL of previously prepared Wash Solution and perform a total of six (6) wash cycles. To ensure plate drying, at the end of the wash, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause poor results.
- 9- Pipette 100 µL of Substrate into all wells, including the Blank well.
- 10- Homogenize gently for ± 10 seconds. Cover the wells with the plate sealer.
- 11- Incubate for 15 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 12- Remove the sealer from the cavities.
- 13- Pipette 100 µL of Stop Solution into all wells, including the Blank well.
- 14- Homogenize gently for ± 30 seconds.
- 15- Read using double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and Reference Standards are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	< 0.100
Reference Standard A	< 0.100
Reference Standard B	Entre 0.050 y 0.200
Reference Standard C	> B y < D
Reference Standard D	> C y < E
Reference Standard E	> 1.500

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS

A calibration curve is used to determine the concentration of Anti-HBs in unknown samples.

Preparation of Calibration Curve

If they were performed in duplicate, calculate the averages of the absorbances obtained from the microplate reader of each of the Reference Standards (A-E), as in the example:

REFERENCE STANDARD	CONCENTRATION (SEE LABEL)	ABSORBANCE	AVERAGE ABSORBANCE
A	0.0 mUI/mL	0.020	0.024 mUI/mL
		0.028	
B	10.0 mUI/mL	0.083	0.082 mUI/mL
		0.081	
C	100.0 mUI/mL	0.512	0.525 mUI/mL
		0.538	
D	300.0 mUI/mL	1.331	1.337 mUI/mL
		1.343	
E	500.0 mUI/mL	2.248	2.271 mUI/mL
		2.294	

Plot the average absorbance of each Reference Standard (A-E) versus the corresponding concentration in mIU/mL on graph paper or through computer programs. Trace the curve point by point.

Calculation of Anti-HBs Concentration of Samples

Calculate the concentration of each analyzed sample from the Calibration Curve made on millimeter paper, interpolating the absorbance value obtained from the sample reading with the corresponding concentration value (mUI/mL). The calculation of the concentration of the tested samples can also be performed using suitable computer programs, through linear regression.

The measurement range of the kit is between 2 and 500 mIU/mL (analytical sensitivity and maximum point of the curve, respectively). Results below the lower limit of detection should be expressed as < 2 mIU/mL. Results above the upper detection limit should be expressed as > 500 mIU/mL. Samples that present values greater than the last point on the curve can be diluted and reanalyzed, as described below.

Sample Dilution

Samples with a concentration > 500 mIU/mL can be diluted with 0.9% sodium chloride at a ratio of 1:10 and reanalyzed. To do so, add 10 µL of sample to 90 µL of 0.9% sodium chloride. The concentration obtained from the diluted sample must be multiplied by the dilution factor (if diluted 1:10, multiply by 10). If the concentration of the diluted sample is less than the cut-off point (10 mIU/mL), a smaller dilution should be made. If the sample concentration is still greater than 500 mIU/mL, a higher dilution should be performed.

Notes:

- 1-** The data presented are for illustration purposes only and cannot be used as a substitute for the calibration curve, which must be constructed in the laboratory. Reference Standards should always be tested with each new assay.
- 2-** Due to the diversity of available antigens used in the manufacture of each test, the concentrations of anti-HBs antibodies may vary between different immunoassay manufacturers.

INTERPRETATION OF RESULTS

According to the recommendation of the World Health Organization (WHO) for the quantification of anti-HBs in samples from individuals vaccinated or infected with the hepatitis B virus, a concentration of 10 IU/L (or 10 mIU/L) is recommended as a cut-off point. Thus, after calculating the concentration of anti-HBs in the samples, consider the values below to determine qualitative results.

Result	QUALITATIVE
	Anti-HBs Concentration
Negative (No Reactive)	< 10.0 mIU/mL
Positive (Reactive)	> 10.0 mIU/mL

Note: In the case of a result between 9 and 11 mIU/mL, the sample must be reanalyzed. Samples that repeatedly get results between 9 and 11 mIU/mL should be retested using an alternative method. If the results are confirmed, a new sample should be collected in two weeks. The result of the last sample collected must prevail. The results provided by this kit must be interpreted by the responsible medical professional and are not the only criterion for determining the patient's diagnosis and/or treatment.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before the descriptive diagnosis of the disease.

INTERFERENTS

No interference was observed for Triglycerides 1500 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Acid Oxalic 60 mg/dL, C-Reactive Protein 41.2 mg/dL and anti-Streptolysin O 1023 IU/mL. High concentrations of Rheumatoid Factor (above 500 IU/mL) may interfere with the result obtained in samples where the concentration of anti-HBs is close to the cut-off point (10 mIU/mL). Plasma samples stored for periods longer than recommended (30 days) may show fibrin and fibronectin precipitation that may interfere with the test. Plasma samples collected in EDTA and Heparin can be used in the product, but decreased values of anti-HBs concentration compared to serum can be obtained.

CROSS REACTIVITY

A cross-reactivity study was performed, evaluating 75 serum and plasma samples negative for anti-HBs but positive for other infections. Among them, 5 samples positive for HIV, 5 samples positive for HCV, 5 samples positive for HSV, 5 samples positive for COVID-19, 5 samples positive for Dengue, 10 samples positive for Rubella, 15 samples positive for CMV, 5 samples positive for Chagas Disease, 5 samples positive for Syphilis and 15 samples positive for Toxoplasmosis. No false positive results were observed for the listed pathogens. Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPEATABILITY	SAMPLE		
	1	2	3
Average	1.750	0.375	0.034
Standard Deviation	0.075	0.025	0.005
Coefficient of Variation (%)	4.31	6.69	14.42

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	1.715	0.368	0.033
Standard Deviation	0.144	0.039	0.005
Coefficient of Variation (%)	8.41	10.58	15.13

ANALYTICAL SENSITIVITY

The analytical sensitivity of the Biolisa Anti HBs kit is 2 mIU/mL.

CLINICAL SENSITIVITY AND SPECIFICITY

The Biolisa Anti HBs kit was used for the analysis of clinical samples previously characterized and confirmed by another reference method. The results show that the Biolisa Anti HBs kit has a clinical sensitivity of 95.0% and a clinical specificity of 99.2%.

		REFERENCE RESULT		
		Positive	Negative	Total
BIOLISA ANTI HBS	Positive	115	1	116
	Negative	6	121	127
	Total	121	122	243

Clinical Sensitivity: 95.0% (115/121) - IC 95%: 89.5 – 98.2%

Clinical Specificity: 99.2% (121/122) - IC 95%: 95.5 – 99.9%

LINEARITY

The reaction is linear up to the concentration at the highest point of the calibration curve. For samples with higher values, the samples must be diluted as described in "Sample Dilution", repeat the dosage and multiply the result obtained by the dilution factor.

HIGH DOSE HOOK EFFECT

No high-dose hook effect was observed up to a concentration of 6,300 mIU/mL of anti-HBs. However, a reduction in absorbance can be observed from a concentration of 12,500 mIU/mL. Also, due to a high-dose hook effect, samples with an anti-HBs concentration ≥ 100,000 mIU/mL may result in results lower than the upper measurement limit (500 mIU/mL).

DIAGNOSTIC SIGNIFICANCE

Hepatitis B is an infection caused by the hepatitis B virus (HBV), considered a public health problem, responsible for a variety of liver diseases, from acute to chronic forms, including cirrhosis and hepatocellular carcinoma. In 2019, the World Health Organization (WHO) estimated that 296 million people were living with the virus worldwide, with 1.5 million new infections and more than 800,000 deaths each year. Hepatitis B is spread through sexual contact, blood exposure, mother-to-child transmission during childbirth, or sharing sharps. In some cases, the virus can persist for a lifetime. HBV belongs to the Hepadnaviridae family and has a protein nucleocapsid (core protein), involving a double-stranded DNA genetic material. The virus is enveloped in a lipoprotein envelope containing surface proteins (HBsAg), one of the first markers to appear in the blood after infection with the Hepatitis B virus. The immune response to HBV infection includes the development of antibodies specific to HBsAg, called anti-HBs. These antibodies appear a few weeks after HBsAg is detected in the blood, and are associated with recovery from infection. These antibodies also appear after the immunization process through hepatitis B vaccines. Therefore, the detection and monitoring of antibodies against HBsAg has become an important tool in the screening and follow-up of infected individuals, as well as an indicator of success of the vaccinated against HBV.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of Biolisa Anti HBs kit at ANVISA: 10269360289

Review: October/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER