BIOLISA CHAGAS RECOMBINANTE DBS

REF **K272**

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FINALIDADE

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- 1- Burns Jr. JM, ET al. Identification and synthesis of a major conserved antigenic epitope of Trypanosoma cruzi. Proc Natl Acad Sci USA 89: 1239-1243, 1992.
- 2- Gruber, A. anda Zingales, B. Trypanosoma cruzi: Characterization of two recombinant antigens with potencial application in the diagnosis of Chagas Disease. Exp. Parasitol. 76:1-12, 1993.
- 3- Ibañez, C.F., Affranchino, J.L., Medina, R.A., Reyes, M.B., Leguizamon, S., Camargo, M.E., Aslund, L., Petteron, U. and Frasch, A.C.C. Multiple Trypanossoma Cruzi antigens containing tabdemly repeated amino acid sequence motifs. Mol. Biochem. Parasitol 30:27-34, 1998.
- 4- Peralta JM, et al. Serodiagnosis of Chagas' disease by enzyme-linked immunosorbent assay using two synthetic peptides as antigens. J Clin Microbiology 32(4): 971-974, 1994.
- 5- Spencer, H.C., Allain, D.S., Sulzer, A.J., Collins, W.E. Evaluation of the Micro Enzyme Lynked Immunosorbent Assay for Antibodies to Trypanosoma cruzi. Am. J. Trop. Med. Hyg., 29 (2): 179.182, 1980.
- 6- Umezawa ES, ET al. Evaluation of recombinant antigen for serodiagnosis of Chagas' disease in South and Central America, J Clin Microbiol 37(5): 1554-1560, 1999.
- 7- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- 8- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 9- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- 10- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1- Burns Jr. JM, ET al. Identification and synthesis of a major conserved antigenic epitope of Trypanosoma cruzi. Proc Natl Acad Sci USA 89: 1239-1243, 1992.
- 2- Gruber, A. anda Zingales, B. Trypanosoma cruzi: Characterization of two recombinant antigens with potencial application in the diagnosis of Chagas Disease. Exp. Parasitol. 76:1-12, 1993.
- 3- Ibañez, C.F., Affranchino, J.I., Medina, R.A., Reyes, M.B., Leguizamon, S., Camargo, M.E., Aslund, L., Petteron, U. and Frasch, A.C.C. Multiple Trypanossoma Cruzi antigens containing tabdemly repeated amino acid sequence motifs. Mol. Biochem. Parasitol 30:27-34, 1998.
- 4- Peralta JM, et al. Serodiagnosis of Chagas' disease by enzyme-linked immunosorbent assay using two synthetic peptides as antigens. J Clin Microbiology 32(4): 971-974, 1994.
- 5- Spencer, H.C., Allain, D.S., Sulzer, A.J., Collins, W.E. Evaluation of the Micro Enzyme Lynked Inmunosorbent Assay for Antibodies to Trypanosoma cruzi. Am. J. Trop. Med. Hyg., 29 (2): 179.182, 1980.
- 6- Umezawa ES, ET al. Evaluation of recombinant antigen for serodiagnosis of Chagas' disease in South and Central America, J Clin Microbiol 37(5): 1554-1560, 1999.
- 7- Cancer Epidemiol Biomarkers Prev; 26(8); 1337-44.2017 AACR
- 8- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 9- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- 10- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Indirect immunoenzymatic test (ELISA) for the determination of IgG Antibodies against *Trypanosoma cruzi*, the protozoan that causes Chagas disease, in **dried blood spot collected on filter paper (DBS)**. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

**Methodology:** Enzyme immunoassay or enzyme immunoassay.  
The BIOLISA CHAGAS RECOMBINANT DBS kit is a solid phase immunoenzymatic assay based on the principle of qualitative detection of IgG Antibodies against *Trypanosoma cruzi* in dried human blood spot collected on filter paper. Specific IgG Antibodies for *T. cruzi*, present in the dried blood spot are eluted and bind to recombinant *T. cruzi* Antigens coated on the microplate forming Antigen-Antibody immune complexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated Anti-IgG Antibodies are added to the microplate, which is then incubated. Peroxidase-conjugated Anti-IgG Antibodies bind to plate-bound Antigen-Antibody complexes. A new wash is carried out to remove the excesses. After this step, the Substrate is added and incubated producing a blue color, which indicates the amount of Anti-*T. cruzi* present in the samples. Stop Solution is added to stop the reaction with a color change from blue to yellow. Color intensity is measured at 450/620 nm.

REAGENTS

- 1- Sensitized Plate** – Store between 2 and 8 °C. Contains: Plate sensitized with *Trypanosoma cruzi* recombinant Antigens and preservative.
- 2- Conjugate** – Store between 2 and 8°C. Contains: Buffer Solution containing Anti-IgG linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** – Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4- Substrate** – Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 5- Stop Solution** – Store between 2 and 8°C. Contains: Hydrochloric Acid 1 M.
- 6- Elution Buffer** – Store between 2 and 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 7- DBS Negative Control** – Store between 2 and 8 °C. Contains: Non-reactive sample for Anti *T. cruzi* IgG Antibodies impregnated on filter paper. **Potentially infectious.**
- 8- DBS Positive Control** – Store between 2 and 8 °C. Contains: Sample reactive for IgG Anti-*T. cruzi*, impregnated on filter paper. **Potentially infectious.**

PRESENTATION

| REAGENTS                | 1              | 2               | 3               |
|-------------------------|----------------|-----------------|-----------------|
|                         | 96 Cavities    | 192 Cavities    | 480 Cavities    |
| 1- Sensitized Plate     | 1 Unit         | 2 Units         | 5 Units         |
| 2- Conjugate            | 1 Vial x 12 mL | 1 Vial x 24 mL  | 1 Vial x 60 mL  |
| 3- Concentrated Washing | 1 Vial x 50 mL | 1 Vial x 100 mL | 1 Vial x 250 mL |
| 4- Substrate            | 1 Vial x 12 mL | 1 Vial x 24 mL  | 1 Vial x 60 mL  |
| 5- Stop Solution        | 1 Vial x 12 mL | 1 Vial x 24 mL  | 1 Vial x 60 mL  |
| 6-Elution Buffer        | 1 Vial x 20 mL | 1 Vial x 40 mL  | 1 Vial x 100 mL |
| 7- DBS Negative Control | 1 Unit         | 2 Units         | 3 Units         |
| 8- DBS Positive Control | 1 Unit         | 2 Units         | 3 Units         |

EQUIPMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 100 to 500 µL with a coefficient of variation less than 1.5%.
- 2- Repeater for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer for technique on filter paper.
- 4- ELISA reader with absorbance capacity at 450 and 620 nm of wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or clock.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator 37 ± 2 °C.
- 11- Paper shredder (diameter 3 ± 0.2 mm).

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be between 2 and 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For in vitro diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate should only be opened after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be recent and free of contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw material of the product is tested and non-reactive for HBsAg and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. Handling any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take the proper biosafety precautions when handling these products.
- 8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microwells.

- 9- As a protection measure, the plate must be covered during the reaction.
- 10- It must be ensured that the bottom of the well is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the assay.
- 11- Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.
- 12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) by Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

**Dried Blood collected on filter paper (DBS) - EDTA or Puncture**  
The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use.  
Blood samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8 °C. Storage for a period longer than 2 years must be carried out at a temperature of -20°C.

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test show that the BIOLISA CHAGAS RECOMBINANTE DBS kit is stable after opening for up to 30 days. This stability may vary according to test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of vial N° 3 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Washing solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30°C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolution.

All other reagents are ready to use.

TECHNIQUE

**For use in automatic equipment, consult the SAC (Customer Advisory Service).**

- Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
- 1- Separate the wells to be used considering: DBS Positive Control, DBS Negative Control and dried blood spot collected on filter paper (it is recommended to test in duplicate). Return unused microplate strips to the original sealed package.
- 2- Separate the first well for the Blank (OPTIONAL).
- 3- Add a 3mm disk of DBS Positive Control, DBS Negative Control and dried blood spot collected on filter paper, **previously perforated**, in the determined wells.
- 4- Pipette 100 µL of Elution Buffer into all wells, **including the Blank well**.
- 5- Gently homogenize for ± 30 seconds, cover the wells with plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 7- After incubation, discard the contents of the wells by aspiration (Washer for filter paper technique). Use approximately 300 µL of **previously prepared** Washing Solution and perform a total of five (5) wash cycles with shake of 5 seconds. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause poor results.
- 8- Pipette 100 µL of Conjugate into all wells, **including the Blank well**.

- 9- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- 10- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 11- Remove the plate sealer from the wells.
- 12- Repeat item 7.
- 13- Pipette 100 µL of Substrate into all wells. **Including the Blank well**.
- 14- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- 15- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 16- Remove the plate sealer from the wells.
- 17- Pipette 100 µL of Stop Solution into all wells.
- 18- Gently homogenize for ± 30 seconds.
- 19- Read using a double filter: 450 nm / 620 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

| ITEM                 | ABSORBANCE |
|----------------------|------------|
| Blank                | < 0.150    |
| DBS Negative Control | < 0.150    |
| DBS Positive Control | > 0.500    |

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula:

Cut Off = Mean Absorbance of DBS Negative Control + 0.400.

Example:

| ITEM                                                        | ABSORBANCE                            |
|-------------------------------------------------------------|---------------------------------------|
| DBS Negative Control                                        | 0.020                                 |
|                                                             | 0.022                                 |
| Cut Off = (Mean Absorbance of DBS Negative Control) + 0.400 | ((0.020 + 0.022) / 2) + 0.400 = 0.421 |

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

| ITEM                           | ABSORBANCE           |
|--------------------------------|----------------------|
| Sample                         | 1.074                |
| Cut Off Value                  | 0.421                |
| Index = Sample / Cut Off Value | 1.074 / 0.421 = 2.55 |

**Note:** Data shown in examples is for illustration only and cannot be used for calculation of results.

Each laboratory must validate the Cut-off according to the instrumentation used and the researched population.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

| RESULTS                | QUALITATIVE         |
|------------------------|---------------------|
|                        | INDEX               |
| Negative (Non-Reagent) | ≤ 0.8               |
| Undetermined           | Between 0.8 and 1.2 |
| Positive (Reagent)     | ≥ 1.2               |

Samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative, and samples with an index between 0.8 and 1.2 are considered indeterminate. Thus, we have that the sample mentioned in the example, whose absorbance was 1.074 and index of 2.55, presents a positive result (index ≥ 1.2).

**Notes:** In case of indeterminate results, the sample must be reanalyzed in duplicate. Samples that give repeatedly indeterminate results should be retested using an alternate method. If results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected shall prevail. The interpretation of a diagnostic test should not be based on a single test. Additional confirmatory tests must be included before a specimen is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, prior to making a descriptive diagnosis of the disease. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS REACTIVITY

A study was carried out with 72 dried blood samples collected on filter paper negative for Anti-Trypanosoma cruzi IgG antibodies, but positive for other infections, in order to evaluate the possibility of cross-reactivity of these interferents in the result of BIOLISA CHAGAS RECOMBINANT DBS. Among them 10 positive samples for Syphilis, 10 positive samples for HCV, 9 positive samples for HBsAg, 8 positive samples for Zika, 10 positive samples for CMV, 10 positive samples for Rubella, 5 positive samples for HIV and 10 positive samples for Toxoplasmosis. No cross-reactivity was observed with positive samples for Syphilis, HCV, HBsAg, Zika, CMV, Rubella, HIV and Toxoplasmosis. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

| Repeatability                | Sample |        |        |
|------------------------------|--------|--------|--------|
|                              | 1      | 2      | 3      |
| Average                      | 0.780  | 0.103  | 0.116  |
| Standard Deviation           | 0.092  | 0.013  | 0.019  |
| Coefficient of Variation (%) | 11.855 | 12.887 | 16.394 |

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

| Reproducibility              | Sample |        |        |
|------------------------------|--------|--------|--------|
|                              | 1      | 2      | 3      |
| Average                      | 0.740  | 0.141  | 0.092  |
| Standard Deviation           | 0.074  | 0.014  | 0.015  |
| Coefficient of Variation (%) | 9.910  | 10.703 | 16.268 |

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA CHAGAS RECOMBINANT DBS kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOLISA CHAGAS RECOMBINANT DBS kit is > 99.9% and the clinical specificity is 99.01%.

|                                |          | Result Reference |          |       |
|--------------------------------|----------|------------------|----------|-------|
|                                |          | Positive         | Negative | Total |
| BIOLISA CHAGAS RECOMBINANT DBS | Positive | 104              | 0        | 104   |
|                                | Negative | 1                | 100      | 101   |
|                                | Total    | 105              | 100      | 205   |

Clinical Sensitivity: >99.9% (104/104) CI 95% = 96.50 to 100%  
Clinical Specificity: 99.01% (104/105) CI 95% = 94.60 to 100%

DIAGNOSTIC SIGNIFICANCE

Chagas disease is a chronic disease caused by infection with the protozoan *Trypanosoma cruzi* endemic to Latin America. The parasite is transmitted to humans by a group of insects of the *Reduviidae* family, with the kissing bug (*Triatoma infestans*) being the main vector. The infection can also be transmitted congenitally, through blood transfusion or organ transplants. This infection affects various organs to different degrees and systems, especially the heart and the gastrointestinal tract.

BIBLIOGRAPHIC REFERENCE

- 1- Burns Jr. JM, ET al. Identification and synthesis of a major conserved antigenic epitope of Trypanosoma cruzi. Proc Natl Acad Sci USA 89: 1239-1243, 1992.
- 2- Gruber, A. anda Zingales, B. Trypanosoma cruzi: Characterization of two recombinant antigens with potencial application in the diagnosis of Chagas Disease. Exp. Parasitol. 76:1-12, 1993.
- 3- Ibañez, C.F., Affranchino, J.I., Medina, R.A., Reyes, M.B., Leguizamon, S., Camargo, M.E., Aslund, L., Petteron, U. and Frasch, A.C.C. Multiple Trypanossoma Cruzi antigens containing tabdemyly repeated amino acid sequence motifs. Mol. Biochem. Parasitol 30:27-34, 1998.
- 4- Peralta JM, et al. Serodiagnosis of Chagas' disease by enzyme-linked immunosorbent assay using two synthetic peptides as antigens. J Clin Microbiology 32(4): 971-974, 1994.
- 5- Spencer, H.C., Allain, D.S., Sulzer, A.J., Collins, W.E. Evaluation of the Micro Enzyme Lynked Immunosorbent Assay for Antibodies to Trypanosoma cruzi. Am. J. Trop. Med. Hyg., 29 (2): 179.182, 1980.
- 6- Urmezawa ES, ET al. Evaluation of recombinant antigen for serodiagnosis of Chagas' disease in South and Central America, J Clin Microbiol 37(5): 1554-1560, 1999.
- 7- Cancer Epidemiol Biomarkers Prev; 26(8): 1337-44.2017 AACR
- 8- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 9- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- 10- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

 QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service  
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of the BIOLISA CHAGAS RECOMBINANT DBS kit at ANVISA: 10269360426

Review: July/2023

UNIVERSAL SYMBOLOGY

|                                                                                       |                                          |                                                                                       |                                        |
|---------------------------------------------------------------------------------------|------------------------------------------|---------------------------------------------------------------------------------------|----------------------------------------|
|    | CATALOG NUMBER                           |    | MADE BY                                |
|    | LOT NUMBER                               |    | CONTROL                                |
|    | MANUFACTURING DATE                       |    | POSITIVE CONTROL                       |
|    | VALIDITY DATE<br>(last day of the month) |    | NEGATIVE CONTROL                       |
|   | TEMPERATURE LIMIT<br>(store)             |   | BIOLOGICAL RISK                        |
|  | CONTENT IS SUFFICIENT<br>FOR <N> TEST    |  | FLAMMABLE                              |
|  | SEE INSTRUCTIONS<br>FOR USE              |  | CORROSIVE                              |
|  | IN VITRO DIAGNOSTIC<br>PRODUCT           |  | TOXIC                                  |
|  | KEEP AWAY<br>FROM SUNLIGHT               |  | DO NOT USE IF<br>PACKAGE IS<br>DAMAGED |
|  | DO NOT REUSE                             |  | PRODUCT<br>STERILIZED                  |
|  | CAUTION                                  |  | DANGER                                 |