Bioclin

BIOLISA CHAGAS RECOMBINANTE

REF **K365**

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FINALIDADE

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1. Cancer Epidemiol Biomarkers Prev; 26(8): 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. Burns Jr. JM, ET al. Identification and synthesis of a major conserved antigenic epitope of Trypanosoma cruzi. Proc Natl Acad Sci USA 89: 1239-1243, 1992.
4. Gruber, A. anda Zingales, B. Trypanosoma cruzi: Characterization of two recombinant antigens with potencial application in the diagnosis of Chagas Disease. Exp. Parasitol. 76:1-12, 1993.
5. Ibañez, C.F., Afranchino, J.I., Medina, R.A., Reyes, M.B., Leguizamon, S., Camargo, M.E., Aslund, L., Petteron, U. and Frascch, A.C.C. Multiple Trypanosoma Cruzii antigens containing tabdemly repeated amino acid sequence motifs. Mol. Biochem. Parasitol 30:27-34, 1998.
6. Peralta JM, et al. Serodiagnosis of Chagas' disease by enzyme-linked immunosorbent assay using two synthetic peptides as antigens. J Clin Microbiology 32(4): 971-974, 1994.
7. Spencer, H.C., Allain, D.S., Sulzer, A.J., Collins, W.E. Evaluation of the Micro Enzyme Lynked Immunosorbent Assay for Antibodies to Trypanosoma cruzi. Am. J. Trop. Med. Hyg., 29 (2): 179.182, 1980.
8. Umezawa ES, ET al. Evaluation of recombinant antigen for serodiagnosis of Chagas' disease in South and Central America, J Clin Microbiol 37(5): 1554-1560, 1999.
9. Bioclin - Dados de Arquivo.

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1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. Burns Jr. JM, ET al. Identification and synthesis of a major conserved antigenic epitope of Trypanosoma cruzi. Proc Natl Acad Sci USA 89: 1239-1243, 1992.
4. Gruber, A. anda Zingales, B. Trypanosoma cruzi: Characterization of two recombinant antigens with potencial application in the diagnosis of Chagas Disease. Exp. Parasitol. 76:1-12, 1993.
5. Ibañez, C.F., Afranchino, J.I., Medina, R.A., Reyes, M.B., Leguizamón, S., Camargo, M.E., Aslund, L., Petteron, U. and Frascch, A.C.C. Multiple Trypanosoma Cruzi antigens containing tabdemly repeated amino acid sequence motifs. Mol. Biochem. Parasitol 30:27-34, 1998.
6. Peralta JM, et al. Serodiagnosis of Chagas' disease by enzyme-linked immunosorbent assay using two synthetic peptides as antigens. J Clin Microbiology 32(4): 971-974, 1994.
7. Spencer, H.C., Allain, D.S., Sulzer, A.J., Collins, W.E. Evaluation of the Micro Enzyme Lynked Immunosorbent Assay for Antibodies to Trypanosoma cruzi. Am. J. Trop. Med. Hyg., 29 (2): 179.182, 1980.
8. Umezawa ES, ET al. Evaluation of recombinant antigen for serodiagnosis of Chagas' disease in South and Central America, J Clin Microbiol 37(5): 1554-1560, 1999.
9. Bioclin - Datos de Archivo.

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INSTRUCTIONS FOR USE

FUNCTION

Test for the qualitative determination of IgG antibodies against *Trypanosoma cruzi*, the protozoan that causes Chagas disease, in biological samples (human serum and plasma) using an enzyme immunoassay. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic

The BIOLISA CHAGAS RECOMBINANT kit is a solid-phase enzyme immunoassay based on the principle of qualitative detection of IgG antibodies against *Trypanosoma cruzi* in human serum or plasma samples. Specific IgG antibodies for *T. cruzi*, present in the sample, bind to recombinant *T. cruzi* antigens coated on the microplate, forming antigen-antibody immune complexes. After initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated anti-IgG antibodies are added to the microplate, which is then incubated. The peroxidase-conjugated anti-IgG antibodies bind to the antigen-antibody complexes on the plate. A further wash is performed to remove excess material. After this step, the substrate is added and incubated, producing a blue color, which indicates the detection of anti-*T. cruzi* IgG antibodies present in the samples. The Stopping Solution is added to interrupt the reaction, resulting in a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** – Store between 2 and 8°C. Contains: Plate sensitized with recombinant *Trypanosoma cruzi* antigen and preservative.
- 2- Conjugate** – Store between 2 and 8°C. Contains: Buffer solution containing anti-IgG bound to Peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Wash** – Store between 2 and 8°C. Contains: Buffer solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4- Sample Diluent** – Store between 2 and 8°C. Contains: Buffer solution, stabilizers, surfactant, dye and preservative.
- 5- Substrate** – Store between 2 and 8°C. Contains: Buffer solution containing urea peroxide, tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** – Store between 2 and 8°C. Contains: 1 M hydrochloric acid.
- 7- Negative Control** – Store between 2 and 8°C. Contains: Buffer solution, stabilizer, surfactant and preservative. Potentially infectious.
- 8- Positive Control** – Store between 2 and 8°C. Contains: Buffer solution, anti-*T. cruzi* IgG antibodies, dye, stabilizers, surfactant and preservative. Potentially infectious.

PRESENTATION

REAGENTS	1	2	3
1- Sensitized Plate	1 units (96 cavities)	2 units (192 cavities)	5 units (480 cavities)
2- Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3- Concentrated Wash	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
4- Sample Diluent	1 Vial x 25 mL	1 Vial x 50 mL	1 Vial x 110 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
7- Negative Control	1 Vial x 300 µL	1 Vial x 300 µL	1 Vial x 1 mL
8- Positive Control	1 Vial x 300 µL	1 Vial x 300 µL	1 Vial x 1 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 300 µL with a coefficient of variation lower than 1.5%.
- 2- Repipette for repetitive pipetting of volumes of 300 µL with a coefficient of variation lower than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5- Absorbent paper to dry the microwells.
- 6- Stopwatch or clock.
- 7- Vial to store the Wash Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator 37 °C ± 2 °C.

STORAGE AND TRANSPORTATION CONDITIONS

The storage temperature should be between 2 and 8 °C. Transportation at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For professional *in vitro* diagnostic use only.

- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate should only be opened after it has reached room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be fresh and free of contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid, which is a strong acid. Therefore, handle it with due care.
- 7- All raw materials of the product are tested and found to be non-reactive for HBsAg, Anti-HIV and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. Handling any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take the appropriate biosafety precautions when handling these products.
- 8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microwells.
- 9- As a protective measure, the plate must be covered during the reaction.
- 10- Ensure that the bottom of the well is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.
- 11- Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during the storage or incubation stages.
- 12- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or by requesting it through Quibasa's SAC (Customer Advisory Service).
- 14- Do not use the product if the packaging is damaged.
- 15- It is essential that the instruments and equipment.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolyzed or highly lipemic samples should not be used. Samples can be stored under refrigeration, between 2 and 8° C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at a temperature of -20 °C. Plasma samples stored for periods longer than the recommended period (30 days) and samples collected with anticoagulants other than those mentioned should not be used.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

ITEM	ABSORBANCE
Blank	< 0.100
Negative Control	< 0.100
Positive Control	> 1.000

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

QUALITATIVE

Calculate the Cut-Off using the following formula:

Cut-Off = (Average Absorbance Obtained from the Positive Control x 0.100) + 0.200.

Example:

ITEM	ABSORBANCE
Positive Control	A1 = 2.051
	A2 = 2.034
Cut Off: (Average Absorbance Obtained with Positive Control x 0.100) + 0.200	((2.051 + 2.034) / 2) x 0.100 + 0.200 = 0.404

Calculate the Index by dividing the absorbance of the sample by the Cut-Off value.

Example:

ITEM	ABSORBANCE
Sample	1.525
Cut-Off Value	0.404
Index = Sample / Cut-Off Value	1.525 / 0.404 = 3.78

Note:The data presented in the examples are for illustrative purposes only and cannot be used to calculate the results.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indices below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (Non-Reactive)	< 0.9
Indeterminate	Between 0.9 and 1.1
Positive (Reactive)	> 1.1

Samples with an index greater than 1.1 are considered positive, samples with an index less than 0.9 are considered negative, and samples with an index between 0.9 and 1.1 are considered indeterminate. Therefore, the sample cited in the example, with an absorbance of 1.525 and an index of 3.78, shows a positive result (index > 1.1).

Note: In the case of an indeterminate result, the sample should be reanalyzed. Samples that repeatedly obtain indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected should prevail. The interpretation of a diagnostic test should not be established based on a single assay. Other confirmatory tests should be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test show that the BIOLISA CHAGAS RECOMBINANTE kit is stable after opening for up to 30 days. This stability may vary depending on the test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of Vial No. 3 (Concentrated Wash) in a ratio of 1:20 with distilled or deionized water; for example: 50 mL of Concentrated Washing Solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

All other reagents are ready for use.

TECHNIQUE

For use in automatic equipment, consult SAC (Customer Advisory Service).

Before starting the assay, allow all reagents, samples, and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

- 1- Separate the wells to be used, considering: Controls and Samples (it is recommended to test in duplicate). Return unused strips from the microplate to the original sealed packaging.
- 2- Separate the first well for the Blank (OPTIONAL).
- 3- Pipette 100 µL of Sample Diluent. Including into the well for the Blank.
- 4- Pipette 5 µL of Controls and Samples into the previously determined wells. Observe the color change of the diluent at the moment of sample addition. The color change indicates that the sample was added to the well correctly.
- 5- Gently homogenize for ± 30 seconds, cover the wells with plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 7- Remove the sealer from the wells.
- 8- Discard the contents of the wells by aspiration (Washer) or by decantation (Manual). Use approximately 300 µL of Washing Solution, previously prepared, and perform a total of five (5) washing cycles with 5-second shaking. To ensure the plate is dry, at the end of the washing, tap the plate for a few seconds on absorbent paper.
- Note: Inadequate washing/drying may cause inadequate results.
- 9- Pipette 100 µL of Conjugate into all wells, including the Blank well.
- 10- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 12- Remove the plate sealer from the wells.
- 13- Repeat item 8.
- 14- Pipette 100 µL of Substrate into all wells, including the Blank well.
- 15- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- 16- Incubate for 15 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 17- Remove the plate sealer from the wells.
- 18- Pipette 100 µL of Stop Solution into all wells, including the Blank well.
- 19- Gently homogenize for ± 30 seconds.
- 20- Read using a dual filter: 450 nm / 630 nm for up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check whether the results obtained for reading the Blank and Controls are compatible with the values presented below:

other available clinical information before a descriptive diagnosis of the disease.

The results provided by this kit should be interpreted by the responsible medical professional and should not be the sole criterion for determining the diagnosis and/or treatment of the patient.

LIMITATIONS OF THE PROCESS

The interpretation of a diagnostic test should not be established based on a single assay. Other confirmatory tests should be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other available clinical information before a descriptive diagnosis of the disease.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS-REACTIVITY

A cross-reactivity study was conducted, evaluating 100 serum and plasma samples negative for Chagas disease but positive for other infections. Among them, 10 samples were positive for HIV, 10 samples positive for Zika, 10 samples positive for CMV, 10 samples positive for Syphilis, 10 samples positive for Toxoplasmosis, 10 samples positive for HBsAg, 10 samples positive for HCV, 10 samples positive for Chikungunya, 10 samples positive for Rubella, and 10 samples positive for Human Visceral Leishmaniasis. No cross-reactivity was observed with samples positive for HIV, Zika, CMV, Syphilis, Toxoplasmosis, HBsAg, HCV, Chikungunya, Rubella, and Human Visceral Leishmaniasis. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits, and tolerances for variations are clearly established. It is important to emphasize that all measurement systems exhibit a characteristic analytical variability, which must be monitored by the laboratories themselves. Therefore, the use of controls is recommended, allowing for the evaluation of the precision and accuracy of dosages.

PRODUCT PERFORMANCE

Precision

REPEATABILITY

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	3
Average	0.780	1.810	0.116
Standard Deviation	0.092	0.146	0.019
Coefficient of Variation (%)	11.855	8.053	16.394

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Average	0.740	1.734	0.141
Standard Deviation	0.074	0.142	0.014
Coefficient of Variation (%)	9.910	8.140	10.703

Clinical Sensitivity and Specificity

The BIOLISA CHAGAS RECOMBINANT kit was used to analyze clinical samples previously characterized and confirmed by another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA CHAGAS RECOMBINANT kit is >99.9% and the clinical specificity is >99.9%.

METHOD	Biolisa Chagas Recombinant		TOTAL	
	Positive	Negative		
Kit Reference	Positive	157	0	157
	Negative	0	172	172
Total Result		157	172	329

Clinical Sensitivity: >99.9% (157/157) - 95% CI (97.7 – 100%)

Clinical Specificity: >99.9% (172/172) - 95% CI (97.9 – 100%)

DIAGNOSTIC SIGNIFICANCE

Chagas disease is a chronic disease caused by infection with the protozoan *T. cruzi*. The parasite is transmitted to humans by a group of insects of the family Roduvidae, with the kissing bug (*Triatoma infestans*) being the main vector. The infection can also be transmitted congenitally, through blood transfusion or organ transplants. This infection affects various organs and systems to varying degrees, especially the heart and gastrointestinal tract.

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. Burns Jr. JM, ET al. Identification and synthesis of a major conserved antigenic epitope of Trypanosoma cruzi. Proc Natl Acad Sci USA 89: 1239-1243, 1992.
4. Gruber, A. anda Zingales, B. Trypanosoma cruzi: Characterization of two recombinant antigens with potencial application in the diagnosis of Chagas Disease. Exp. Parasitol. 76:1-12, 1993.
5. Ibañez, C.F., Afranchino, J.I., Medina, R.A., Reyes, M.B., Leguizamón, S., Camargo, M.E., Aslund, L., Petteron, U. and Frasc, A.C.C. Multiple Trypanosoma Cruzi antigens containing tabdelym repeated amino acid sequence motifs. Mol. Biochem. Parasitol 30:27-34, 1998.
6. Peralta JM, et al. Serodiagnosis of Chagas' disease by enzyme-linked immunosorbent assay using two synthetic peptides as antigens. J Clin Microbiology 32(4): 971-974, 1994.
7. Spencer, H.C., Allain, D.S., Sulzer, A.J., Collins, W.E. Evaluation of the Micro Enzyme Lynked Immunosorbent Assay for Antibodies to Trypanosoma cruzi. Am. J. Trop. Med. Hyg., 29 (2): 179.182, 1980
8. Umezawa ES, ET al. Evaluation of recombinant antigen for serodiagnosis of Chagas' disease in South and Central America, J Clin Microbiol 37(5): 1554-1560, 1999.
9. Bioclin - Dados de Arquivo.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA BIOLISA CHAGAS RECOMBINANTE:
10269360489

Review: April/2026

UNIVERSAL SYMBOLOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER