BIOLISA CHAGAS RECOMBINANTE

REF K180

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Test for qualitative determination of anti-*Trypanosoma cruzi* IgG antibodies in serum or human plasma (EDTA or Heparin) by microplate enzyme immunoassay. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The Biolisa Chagas Recombinant kit is a solid phase immunoenzymatic assay based on the principle of qualitative detection of antibodies against *T. cruzi* in human serum or plasma. Specific *T. cruzi* antibodies, present in the sample, bind to recombinant *T. cruzi* antigens coated on the microplate forming antigen-antibody complexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated anti-IgG antibodies are added to the microplate, which is then incubated. Anti-Peroxidase conjugated anti-IgG antibodies bind to antigen-antibody complexes. A new wash is performed to remove surplus. After this step, the Substrate is added and incubated producing a blue color, which indicates the amount of anti-T IgG antibodies. *T. Cruzi* present in the samples. The Stop Solution is added to stop the reaction with a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1 - Sensitized plate** - Store between 2 and 8°C. Microplate impregnated with recombinant *T. cruzi* antigen.
- 2 - Conjugate** - Store between 2 and 8°C. Peroxidase bound anti-IgG conjugate and stabilizer.
- 3 - Concentrated Washing** - Store between 2 and 8°C. Solution Buffer, surfactant and preservative.
- 4 - Sample Diluents** - Store between 2 and 8°C. Solution Buffer, stabilizer and surfactant.
- 5 - Substrate** - Store between 2 and 8°C. Solution containing Hydrogen Peroxide and Tetramethylbenzidine (TMB).
- 6 - Stop Solution** - Store between 2 and 8°C. Acid Solution 1N.
- 7 - Negative Control** - Store between 2 and 8°C. Serum inactivated, not reactive for *T. cruzi* and preservative. **Potentially infective.**
- 8 - Positive Control** - Store between 2 and 8°C. Inactivated serum containing anti-*T. cruzi* antibodies and preservative. **Potentially infective.**
- 9 - Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
1- Sensitized Plate	1 units (96 cavities)	2 units (192 cavities)	5 units (480 cavities)
2- Conjugate	15 mL	2 x 15 mL	5 x 15 mL
3- Concentrated Washing	50 mL	2 x 50 mL	5 x 50 mL
4- Sample Diluents	30 mL	2 x 30 mL	5 x 30 mL
5- Substrate	15 mL	2 x 15 mL	5 x 15 mL
6- Stop Solution	10 mL	2 x 10 mL	5 x 10 mL
7- Negative Control	0.5 mL	2 x 0.5 mL	5 x 0.5 mL
8- Positive Control	0.5 mL	2 x 0.5 mL	5 x 0.5 mL
9- Plate Sealers	3 units	5 units	10 units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Instructions for use (manual)

Required materials not contained in the kit:

- 1-** Pipettes capable of dispensing volumes of 10, 50, 100 and 200 µL with coefficient of variation smaller than 1.5%.
- 2-** Pipettes for repetitive pipetting of volumes of 100 µL and 300 µL, with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3-** Microplate washer (optional).
- 4-** ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5-** Absorbent paper to dry the wells.
- 6-** Stopwatch or clock.
- 7-** Bottle to stock the wash solution after dilution.
- 8-** Distilled or deionized water.
- 9-** Quality Control Tools.
- 10-** Incubator of 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- 2-** Strictly follow the methodology proposed to obtain accurate results.
- 3-** The envelope containing the strips should only be opened after it reaches room temperature. Place the strips with unused cavities in the aluminum bag, seal and store between 2 and 8°C.
- 4-** The water used in material cleaning must to be recent and free of contaminants.
- 5-** Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6-** The Stop Solution contains a strong acid. Therefore, handle it with care.
- 7-** All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1 & 2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manual manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.
- 8-** Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 9-** As a safety measure, you should cover the plate during the reaction.
- 10-** Yu must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fluid before reading the plate. Do not let the cavities run dry during the test.
- 11-** Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13-** In order to obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or at the request of the SAC (Customer Assistance) of Quibasa.
- 14-** Do not use the product in case of damage to the packaging.
- 15-** It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin).
Hemolyzed or highly lipemic samples should not be used. Samples may be stored under refrigeration at 2-8°C for a maximum period of 7 days. If samples can not be analyzed within 7 days, they can be stored for up to 30 days at -20 ° C (freezer).

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of vial N°3 (Concentrated Wash) into 1000 mL of distilled or deionized water. After preparation, the solution may be stored at 2-8°C until the expiration date printed on the original vial. Can be stored at room temperature. If crystallization occurs, heat to 37°C until dissolution.

STABILITY AFTER OPEN

The stability test results show that the Biolisa Chagas Recombinante kit is stable after opening for up to 30 days. This stability may vary according to test conditions and environment. Therefore, it is suggested to monitor product performance using internal kit controls and technique validation criteria.

TECHNIQUE

- Before starting the assay, place all Reagents, Samples and Controls to stabilize at room temperature (15-30°C) for at least 40 minutes.
- 1-** Select the cavities to be used considering: Controls and Samples (They can be tested in duplicates). Return the strips of the plate will not be used for the original sealed packaging.
- 2-** Select the first cavity for Blank (OPTIONAL).
- 3-** Pipette 200 µL of Sample Diluent into all the cavities.
- 4-** Pipette 10 µL of Negative Control, Positive Control and Sample into the cavities previously determined.
- 5-** Homogenize gently for ± 30 seconds. Cover the cavities with sealer.
- 6-** Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 7-** Remove the sealing from the cavities.
- 8-** Discard the contents of the cavities by aspiration (Washer) or by decanting (manual). Use approximately 300 µL of Wash Solution, **previously diluted***, to perform a total of five (5) washing cycles. To ensure the drying of the plate at the end of washing, beat it for a few seconds on absorbent paper.
- Note:** Poor washing and drying can cause inadequate results.
- 9-** Pipette 100 µL of Conjugate into each cavities. **Including in the Blank (if you made this choice).**
- 10-** Homogenize gently for ± 30 seconds. Cover the cavities with sealer.
- 11-** Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 12-** Remove the sealing from the cavities.
- 13-** Repeat item 8.
- 14-** Pipette 100 µL of Substrate in all cavities.
- 15-** Homogenize gently for ± 30 seconds. Cover the cavities with sealer.
- 16-** Incubate for exactly 30 minutes at **room temperature**.
- 17-** Remove the sealing from the cavities.
- 18-** Pipette 50 µL Stop Solution into all cavities.
- 19-** Homogenize gently for ± 30 seconds.
- 20-** Read: 450 nm (primary filter)/630 nm (secondary filter) up to 30 minutes maximum.
- Note:** The reaction changes color at each step of the procedure in all phases. Reagents are ready to use, except Wash Solution.

TECHNIQUE VERIFICATION

Verify if the results obtained for reading the Blank and Controls are compatible with the values shown below:

ITEM	ABSORBANCE
Blank	< 0.15
Negative Control	< 0.20
Positive Control	≥ 1.00

If the values are outside the expected values, the technique must be repeated.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

CALCULATIONS QUALITATIVE

Sample serum and plasma

Calculate the Cut Off accoding to the formula below:

Cut Off = (Average Abs. Positive Control x 0,1) + 0,200

Example:

ITEM	ABSORBANCE
Positive Control	A1 = 2.051
	A2 = 2.034
Cut Off = (Average Absorbance Positive Control x 0.1) + 0.200	Cut Off = ((2.051+2.034)/2 x 0.1) + 0.200 = 0.404

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1.525
Cut Off Value	0.404
Index: Sample/ Cut Off Value	1.525 / 0.404 = 3.774

INTERPRETATION OF RESULTS

RESULTS	QUALITATIVE
	INDEX
Negative	< 0.9
Positive	≥ 1.1
Undetermined	≥ 0.9 e < 1.1

Note: In case of undetermined result, the sample should be retested. Samples that obtain repeatedly indeterminate results should be retested using an alternative method. If the results remain undetermined, a new sample should be collected within two weeks. The result of the last sample collected should prevail.

The results provided by this kit should be interpreted by the responsible medical professional and not the only criterion for determining the diagnosis and / or treatment of the patient.

Note: The data presented in the examples are for ilustration only and can not be used for calculation of the results.

PROCEDURE LIMITATIONS

All indeterminate or positive samples should be duplicated. Repeatedly positive samples should be confirmed by techniques such as IFI, Western Blot, Xenodiagnosis and PCR. Samples with microbial contamination or people with other parasites or autoimmune disease may produce false results.

A negative result does not exclude the possibility of *T. cruzi* infection. A humoral immune response is not detected during the first few weeks after infection by any method on the market. The diagnosis of Chagas' disease should be based on the combination of results, including clinical history, diagnosis of infection and subsequent serological tests. The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests should be included before a sample is considered positive.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, norms, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For this purpose, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

PRODUCT PERFORMANCE
ACCURACY

Repeatability

Repeatability was calculated from 10 successive determinations using 3 different samples, yielding the following results:

Repeatability	SAMPLE		
	1	2	3
Average	5.28	4.36	2.63
Standard Deviation	0.14	0.27	0.20
Coefficient of Variation (%)	2.63	6.10	7.70

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 different samples, obtaining the following results:

Reproducibility	SAMPLE		
	1	2	3
Average	5.30	4.38	2.67
Standard Deviation	0.02	0.02	0.04
Coefficient of Variation (%)	0.38	0.46	1.50

CLINICAL SENSITIVITY AND SPECIFICITY

The Biolisa Chagas Recombinant kit analyzed clinical specimens compared to the reference kit. The results show that the clinical sensitivity of the Biolisa Chagas Recombinant Kit is > 99.9%, and the clinical specificity is 99.3%.

METHOD		Biolisa Chagas Recombinant		TOTAL
		Positive	Negative	
Kit Reference	Positive	209	0	209
	Negative	7	940	947
Total Result		216	940	1156

Clinical Sensitivity: > 99.9% (209 / 209)

Clinical Specificity: 99.3% (940 / 947)

DIAGNOSTIC SIGNIFICANCE

Chagas disease is a chronic disease caused by infection with the protozoan *T. cruzi*. The parasite is transmitted to humans by a group of insects of the family Roduvidae, being the barber (*Triatoma infestans*), the main vector. The infection can also be congenitally transmitted by blood transfusion or organ transplants. This infection affects several organs in varying degrees and systems, especially the heart and gastrointestinal tract.

BIBLIOGRAPHIC REFERENCES

- 1- Burns Jr. JM, ET al. Identification and synthesis of a major conserved antigenic epitope of Trypanosoma cruzi. Proc Natl Acad Sci USA 89: 1239-1243, 1922.
- 2- Gruber, A. anda Zingales, B. Trypanosoma cruzi: Characterization of two recombinant antigens with potencial application in the diagnosis of Chagas Disease. Exp. Parasitol. 76:1-12, 1993.
- 3- Ibañez, C.F., Affranchino, J.I., Medina, R.A., Reyes, M.B., Leguizamón, S., Camargo, M.E., Aslund, L., Petteron, U. and Frascch, A.C.C. Multiple Trypanossoma Cruzi antigens containing tabdemyly repeated amino acid sequence motifs. Mol. Biochem. Parasitol 30:27-34, 1998.
- 4- Peralta JM, et al. Serodiagnosis of Chagas' disease by enzyme-linked immunosorbent assay using two synthetic peptides as antigens. J Clin Microbiology 32(4): 971-974, 1994.
- 5- Spencer, H.C., Allain, D.S., Sulzer, A.J., Collins, W.E. Evaluation of the Micro Enzyme Lynked Immunosorbent Assay for Antibodies to Trypanosoma cruzi. Am. J. Trop. Med. Hyg., 29 (2): 179.182, 1980.
- 6- Umezawa ES, ET al. Evaluation of recombinant antigen for serodiagnosis of Chagas' disease in South and Central America, J Clin Microbiol 37(5): 1554-1560, 1999.
- 7- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Biolisa Chagas Recombinant kit: 10269360306

Review: July/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER