Bioclin

BIOLISA CMV IgG DBS

REF K274

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FINALIDADE

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- Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
- Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after "Silent" Congenital Cytomegalovirus Infection. N. Engl. J. Med. (1976) 295:468-470.
- Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovirus. N. Engl. J. Med. (1974) 790:291-296.
- Stern, H. Cytomegalovirus Vaccine: Justification and Problems. In: Waterson AP (ed.) Recent Advances in Clinical Virology (1977) 117-134.

5. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

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INSTRUCTIONS FOR USE

FUNCTION

Indirect enzyme-linked immunosorbent assay (ELISA) for the quantitative and qualitative determination of cytomegalovirus (CMV) IgG antibodies in **dried blood samples collected on filter paper (DBS)**. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA CMV IgG DBS kit is a solid-phase immunoassay based on the principle of indirect qualitative and quantitative detection of IgG antibodies to CMV in dried human blood samples collected on filter paper. Anti-CMV IgG antibodies present in the dried blood sample are eluted and bind to the Antigens coated on the microplate forming Antigen-IgG Antibody immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated Human Anti-IgG Antibodies are added to the microplate, which is then incubated. The enzyme-conjugated Human Anti-IgG Antibodies bind to the IgG Antibodies present on the Antigen-coated plate. Another wash is carried out to remove any excess. After this step, the Substrate is added and incubated, producing a blue color, which indicates the amount of Anti-CMV IgG Antibodies present in the sample. The Stop Solution is added to stop the reaction and the color changes from blue to yellow. The intensity of the color is measured at 450/620 nm.

REAGENTS

1- Reference Standards (A - E) - Store between 2 and 8°C. Contains: Five (5) vials (A - E) of Reference Standards containing Anti-CMV IgG Antibodies in different concentrations in a buffer solution containing surfactant, stabilizers, dye and preservative. The concentrations of the Reference Standards (A-E) vary with each batch. See bottle labels. **Potentially infectious.**

2- Conjugate - Store between 2 and 8°C. Contains: Buffer solution, Peroxidase-linked Human Anti-IgG Antibody, surfactant, stabilizers, dye and preservative.

3- Sensitized plate - Store between 2 and 8°C. Contains: Sensitized plate with purified CMV antigens.

4- Concentrated Wash - Store between 2 and 8°C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.

5- Substrate - Store between 2 and 8°C. Contains: Buffer solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.

6- Stop Solution - Store between 2 and 8°C. Contains: 1 M Hydrochloric Acid.

7- Elution Buffer - Store between 2 and 8°C. Contains: Buffer solution, surfactant, stabilizer and preservative.

8- DBS Negative Control - Store between 2 and 8°C. Contains: Non-reactive sample for IgG Anti-CMV antibodies impregnated on filter paper. **Potentially infectious.**

9- DBS Positive Control - Store between 2 and 8°C. Contains: Reactive sample for IgG Anti-CMV antibodies impregnated on filter paper. **Potentially infectious.**

PRESENTATION

REAGENTS	1	2	3
	96 Cavidades	192 Cavidades	480 Cavidades
1- Reference Standards (A-E)	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 2 mL
2- Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3- Sensitized Plate	1 Unit	2 Units	5 Units
4- Concentrated Wash	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
7- Elution Buffer	1 Vial x 20 mL	1 Vial x 40 mL	1 Vial x 100 mL
8- DBS Negative Control	1 Unit	2 Units	3 Units
9- DBS Positive Control	1 Unit	2 Units	3 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table

Required materials not contained in the kit:

- Pipettes capable of dispensing volumes from 100 to 500 µL with a coefficient of variation of less than 1.5%.
- Repipettor for repetitive pipetting of volumes of 300 µL with a coefficient of variation of less than 1.5% or multichannel pipette (optional).
- Microplate washer for filter paper technique.
- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- Absorbent paper to dry the microcavities.
- Stopwatch or clock.
- Bottle to store the washing solution after dilution.
- Distilled or deionized water.
- Quality control tools.
- Incubator 37 °C ± 2 °C.
- Paper mincer (diameter 3 ± 0.2 mm).

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8 °C. Transportation at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the proposed methodology in order to obtain accurate results.

3- The sachet containing the microplate should only be opened after it has reached room temperature. Replace the unused microcavity strips in the sachet, seal and store between 2 and 8 °C.

4- The water used to clean the material must be fresh and free of contaminants.

5- Saturated deionizing columns release alkaline water, various ions, oxidizing and reducing agents which can significantly alter the results.

6- The Stop Solution contains hydrochloric acid, which is a strong acid. Therefore, handle it with due care.

7- All the product's raw materials are tested as non-reactive for HBsAg, anti-HIV and Anti-HCV. However, these tests do not provide complete assurance

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

of the absence of infectious agents. The handling of any product containing a biological sample is potentially capable of transmitting diseases. Biosafety precautions must therefore be taken when handling these products.

8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microcavities.

9- As a protective measure, cover the plate during the reaction.

10- Ensure that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.

11- Do not expose reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.

12- We recommend applying local, state and federal environmental protection regulations so that reagents and biological material are disposed of in accordance with current legislation.

13- For information related to biosafety or in the event of accidents with the product, consult the SDS (Safety Data Sheet) available at www.bioclin.com.br or by contacting Quibasa's Customer Service Center.

14- Do not use the product if the packaging is damaged.

15- It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

SAMPLES

Dried blood collected on filter paper (DBS) - EDTA or Puncture

The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must check the quality of the samples before use. Stability studies carried out by Bioclin's Research and Development department show that dried blood samples on filter paper can be stored for up to 35 days at room temperature, 2 years at 2-8 °C or 2 years at -20°C. The Manual "Coleta de sangue: Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB)"¹⁰, reports that dried blood samples on filter paper can be stored at room temperature, as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, samples should be kept refrigerated between 2 and 8 °C. Storage for longer than 2 years should be at -20 °C.

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test show that the BIOLISA CMV IgG DBS kit is stable after opening for up to 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the product's performance using the kit's internal controls and the technique's validation criteria.

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of vial No. 4 (Concentrated Wash) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Wash solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

All other reagents are ready for use.

TECHNIQUE

For use in automatic equipment, consult the Customer Service Department (SAC).

Before starting the test, leave all reagents, samples and controls to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.

1- Separate the wells to be used considering: Reference Standards (A - E), DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper (it is recommended to test in duplicate). Return the unused microplate strips to their original sealed packaging.

2- Separate the first well for the blank (OPTIONAL).

3- Add a 3 mm disk of DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper, previously minced, to the determined wells.

4- Pipette 100 µL of Reference Standards (A - E) into the previously determined wells.

5- Pipette 100 µL of Elution Buffer, **including** into the Blank well. **Do not**

pipette into the Reference Standards (A - E) wells.

6- Homogenize gently for ± 30 seconds, cover the wells with a plate sealer.

7- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

8- Remove the sealer from the wells.

9- After incubation, discard the contents of the wells by aspiration (washer for filter paper technique). Use approximately 300 µL of **previously prepared** washing solution and carry out a total of five (5) washing cycles with a 5-second shake. To ensure that the plate is dry, at the end of the wash, tap the plate for a few seconds on absorbent paper.

Note: Poor washing/drying can cause inadequate results.

10- Pipette 100 µL of Conjugate into all the wells, **including** the blank well.

11- Homogenize gently for ± 30 seconds. Cover the wells with the plate sealer.

12- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

13- Remove the plate sealer from the wells.

14- Repeat item 9.

15- Pipette 100 µL of Substrate into all wells, **including** the blank well.

16- Homogenize gently for ± 30 seconds. Cover the wells with the plate sealer.

17- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

18- Remove the plate sealer from the wells.

19- Pipette 100 µL of Stop Solution into all wells.

20- Homogenize gently for ± 30 seconds.

21- Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check that the results obtained for the Blank, Controls and Reference Standards are compatible with the values shown below:

ITEM	ABSORBANCE
Blank	< 0.100
Reference Standard A	< 0.100
Reference Standard B	> 0.200 and < 0.600
Reference Standard C	> B and < D
Reference Standard D	> C and < E
Reference Standard E	> 1.400
DBS Negative Control	< 0.100
DBS Positive Control	> 0.500

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS

QUALITATIVE

Consider the average absorbance obtained with Reference Standard B as the Cut Off.

Example:

ITEM	ABSORBANCE
Reference Standard B	0.412
	0.411
Cut Off = Average Absorbance of Reference Standard B	$(0.412 + 0.411) / 2 + 0.200 = 0.412$

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCIA
Sample	1.254
Cut-Off Value	0.412
Index = Sample / Cut-Off Value	$1.254 / 0.412 = 3.04$

Note: The data presented in the examples is for illustrative purposes only and cannot be used to calculate results. Each laboratory must validate the cut-off according to the instrumentation used and the population surveyed.

QUANTITATIVE

A calibration curve is used to determine the concentration of CMV IgG antibodies in unknown samples.

Preparation of the Calibration Curve

Record the absorbances obtained on the microplate reader, as shown in example below. Calculate the averages of the duplicates (if duplicates are made). Plot the average absorbances of each Reference Standard (A-E) versus the corresponding Concentration in IU/mL (described on the label) on graph paper and draw the curve.

REFERENCE STANDARD	CONCENTRATION (SEE LABEL)	ABSORBANCE	DUPLICATES AVERAGE
A	0.0 IU/mL	0.005	0.006
		0.006	
B	1.2 IU/mL	0.412	0.411
		0.409	
C	4.5 IU/mL	0.902	0.906
		0.910	
D	9.0 IU/mL	1.301	1.305
		1.308	
E	18.0 IU/mL	2.121	2.122
		2.123	

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	QUALITATIVE	QUANTITATIVE
	INDEX	CONCENTRATION
Negative (Non Reactive)	$\leq 0,8$	$\leq 1,2$
Indeterminate	Between 0,8 and 1,2	Between 1,2 and 1,4
Positive (Reactive)	$\geq 1,2$	$\geq 1,4$

Samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative and samples with an index between 0.8 and 1.2 are considered indeterminate.

Thus, the sample cited in the example, whose absorbance was 1.254 and index 3.04, shows a positive result (index ≥ 1.2).

The concentration of the samples tested can also be calculated using suitable computer programs, using linear regression.

SAMPLE	ABSORBANCE	CALCULATED CONCENTRATION (UI/mL)	QUANTITATIVE RESULT
1	0.125	Below the lower limit of detection	< 1,2 UI/mL
2	2.942	Above the upper limit of detection	> 18,0 UI/mL
3	0.685	1,66	1,66 UI/mL

Observations: In the event of an indeterminate result, the sample must be retested in duplicate. Samples with repeatedly indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be taken within two weeks. The result of the last sample taken should prevail. All results must be interpreted in conjunction with other available clinical information before a descriptive diagnosis of the disease can be made. A negative result does not exclude the possibility of exposure. The results provided by this kit must be interpreted by the medical professional responsible and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

LIMITATIONS OF THE PROCESS

The interpretation of a diagnostic test should not be based on a single test. Other confirmatory tests should be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before a descriptive diagnosis of the disease can be made.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS-REACTIVITY

A study was carried out with 73 dried blood samples collected on filter paper that were negative for CMV IgG but positive for other infections, in order to assess the possibility of cross-reactivity of these interferents in the BIOLISA CMV IgG DBS result. These included 7 samples positive for HTLV, 10 samples positive for Syphilis, 6 samples positive for HBsAg, 10 samples positive for HCV, 7 samples positive for Chagas Disease, 6 samples positive for Zika, 10 samples positive for Toxoplasmosis, 10 samples positive for Rubella and 7 samples positive for HIV. No cross-reactivity was observed with samples positive for HTLV, Syphilis, HBsAg, HCV, Chagas Disease, Zika, Toxoplasmosis, Rubella and HIV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must take into account the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is advisable to use controls to assess the precision and accuracy of dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	3
Average	2.422	1.151	0.119
Standard Deviation	0.100	0.059	0.011
Coefficient of Variation (%)	4.123	5.101	8.931

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Average	2.441	1.159	0.119
Standard Deviation	0.075	0.040	0.010
Coefficient of Variation (%)	3.051	3.441	8.556

Considering the variation of all samples tested in the precision study, it was observed that the kit presented a coefficient of variation of less than 20%.

Analytical Sensitivity

The analytical sensitivity of the BIOLISA CMV IgG DBS kit is 0.08 IU/mL.

Clinical Sensitivity and Specificity

The BIOLISA CMV IgG DBS kit was used for the analysis of clinical samples previously characterized and confirmed using another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA CMV IgG DBS kit is 98.5% and the clinical specificity is 98.0%.

BIOLISA CMV IgG DBS	REFERENCE		
	Positive	Negative	Total
Positive	65	1	66
Negative	1	49	50
Total	66	50	116

Clinical Sensitivity: 98.5% (65/66) - IC 95% (91.8 a 99.9%)

Clinical Specificity: 98.0% (49/50) - IC 95% (89.4 a 99.9%)

DIAGNOSTIC SIGNIFICANCE

Cytomegalovirus (CMV) is a member of the herpesvirus family that includes herpes simplex virus (HSV) 1 and 2, varicella zoster virus (VZV) and Epstein-Barr virus (EBV). It is a human pathogen that is transmitted through saliva, sexual contact, perinatal care, organ transplants or blood transfusions. In most cases, the infection remains asymptomatic. However, CMV infection can cause severe disease in newborns and immunosuppressed individuals, such as patients with AIDS, cancer or who have received organ transplants. During immunosuppressive therapy, reactivation of latent virus or primary infection often occurs. CMV infections can be acquired before birth, during delivery and later in life. They can cause severe congenital anomalies such as microcephaly, motor disability and mental retardation. Therefore, primary determination of maternal infection and its distinction from latent infection are of great importance. The presence of IgM antibodies indicates the presence of primary infection, while the presence of IgG antibodies indicates the immune status of the patients.

BIBLIOGRAPHIC REFERENCES

- Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
- Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after "Silent" Congenital Cytomegalovirus Infection. N. Engl. J. Med. (1976) 295:468-470.
- Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovirus. N. Engl. J. Med. (1974) 790:291-296.
- Stern, H. Cytomegalovirus Vaccine: Justification and Problems. In: Waterson AP (ed.) Recent Advances in Clinical Virology (1977) 117-134.
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2,

2002:31.

- Cancer Epidemiol Biomarkers Prev; 26(8); 1337-44.2017 AACR
- Coleta de sangue: Diagnóstico e monitoramento das DST, Aids e Hepatites Virais: Brasília, Ministério da Saúde, Departamento de DST, Aids e Hepatites Virais. 2010, 98 p. (Série TELELAB).
- Bioclin – Archive Data.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for CMV IgG DBS kit: 10269360460

Review: March/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER