BIOLISA CMV IgM DBS

REF **K273**

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FINALIDADE

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- Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
- Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after “Silent” Congenital Cytomegalovirus Infection. N. Engl. J. Med. (1976) 295:468-470.
- Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovuruss. N. Engl. J. Med. (1974) 790:291-296.

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- Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
- Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after "Silent" Congenital Cytomegalovirus Infection. N. Engl. J. Med. (1976) 295:468-470.
- Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovirus. N. Engl. J. Med. (1974) 790:291-296.

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REF K273

INSTRUCTIONS FOR USE

FUNCTION
 Enzymatic immunoassay capture test (ELISA) for determining anti-Cytomegalovirus (CMV) class IgM antibodies in biological samples of dried blood collected on filter paper (DBS), for Pre-Natal and Neonatal screening, using an enzyme immunoassay test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION
Methodology: Enzymeimmunoassay or enzyme immunoassay.
 The BIOLISA CMV IgM DBS kit is a solid-phase enzyme immunoassay based on the principle of qualitative detection by capture of anti-Cytomegalovirus (CMV) IgM class antibodies in human dried blood samples collected on filter paper. IgM antibodies present in the dried blood sample are eluted and bind to the anti-IgM antibodies coated on the microplate forming immune complexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated CMV antigens are added to the microplate, which is then incubated. The enzyme-conjugated Antigens bind to the anti-CMV IgM Antibodies present, attached to the plate coated with anti-IgM antibodies. New washing is carried out to remove excess. After this step, the Substrate is added, producing a blue color that indicates the amount of anti-CMV IgM Antibodies present in the sample. The Stop Solution is added to stop the reaction and there is a color change from blue to yellow. Color intensity is measured at 450/620 nm.

REAGENTS
1- Sensitized Plate – Store between 2 and 8°C. Contains: Plate sensitized with anti-IgM antibodies.
2- Conjugate – Store between 2 and 8°C. Contains: Buffer solution, CMV antigens conjugated to peroxidase, surfactant, stabilizers, dye and preservative.
3- Concentrated Washing – Store between 2 and 8°C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
4- Substrate – Store between 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
5- Stop Solution – Store between 2 and 8°C. Contains: 1 M Hydrochloric Acid.
6- Elution Buffer – Store between 2 and 8°C. Contains: Buffer solution, surfactant, stabilizer and preservative.
7- DBS Negative Control – Store between 2 and 8°C. Contains: Non-reactive sample for IgM Anti-CMV antibodies impregnated in filter paper. **Potentially infectious.**
8- DBS Positive Control – Store between 2 and 8°C. Contains: Reactive sample for IgM Anti-CMV antibodies impregnated in filter paper. **Potentially infectious.**

REAGENTS	1	2	3
	96 cavities	192 cavities	480 cavities
1- Sensitized Plate	1 Unit	2 Units	5 Units
2- Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3- Concentrated Washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
4- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
5- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6- Elution Buffer	1 Vial x 20 mL	1 Vial x 40 mL	1 Vial x 100 mL
7- DBS Negative Control	1 Unit	2 Units	3 Units
8- DBS Positive Control	1 Unit	2 Units	3 Units

EQUIPMENTS AND OPERATIONAL INPUTS
 - Reagents described in the previous table

Required materials not contained in the kit:
1- Pipettes capable of dispensing volumes of 100 to 500 µL with a coefficient of variation of less than 1.5%.
2- Repipettor for repetitive pipetting of 300 µL volumes with a coefficient of variation less than 1.5% or multichannel pipette (optional).
3- Microplate washer for filter paper technique.
4- ELISA reader with absorbance capacity at 450 and 620 nm wavelength.
5- Absorbent paper to dry the microcavities.
6- Stopwatch or clock.
7- Bottle to store the Washing Solution after dilution.
8- Distilled or deionized water.
9- Quality Control Tools.
10- Incubator 37 °C ± 2 °C.
11- Paper shredder (diameter 3 + 0.2 mm).

STORAGE AND TRANSPORT CONDITIONS
 The storage temperature should be 2 to 8 °C. Transport at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES
1- For *in vitro* diagnostic use only.
2- Strictly follow the proposed methodology to obtain accurate results.
3- The sachet containing the microplate must only be opened after reaching room temperature. Return the unused microcavity strips to the sachet, seal and store at 2 to 8 °C.
4- The water used to clean the material must be recent and free from contaminants.
5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
7- All raw materials for the product are tested and are non-reactive for HBsAg and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. The handling of any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take appropriate biosafety precautions when handling these products.

8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microcavities.
9- As a protective measure, the plate must be covered during the reaction.
10- It must be ensured that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.
11- Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation stages.
12- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or through a request through SAC (Customer Advisory Service) from Quibasa.
14- Do not use the product if the packaging is damaged.
15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES
Dried blood collected on filter paper (DBS) – EDTA, Puncture or Heel (in neonates)
 The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use.
 Stability studies, carried out by Bioclin's Research and Development department, show that dried blood samples on filter paper can be stored for up to 35 days at room temperature, 2 years at 2-8 °C or 2 years at -20° W. The Manual "Blood Collection: Diagnosis and monitoring of STDs, AIDS and Viral Hepatitis: Brasília, Ministry of Health, Department of STD, AIDS and Viral Hepatitis. 2010, 98 p. (TELELAB Series)" ⁵, reports that dried blood samples on filter paper can be stored at room temperature, as long as they are protected from direct sunlight and in low humidity, in hermetically sealed plastic bags with or without desiccant. For storage of up to 2 years, samples must remain refrigerated, between 2 and 8 °C. Storage for more than 2 years must be carried out at a temperature of -20 °C.

PROCESS DESCRIPTION
Stability After Opening
 The stability test results prove that the BIOLISA CMV IgM DBS kit is stable after opening for up to 30 days. This stability may vary depending on test and environmental conditions. Therefore, it is suggested to monitor the product's performance using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS
Washing Solution
 Dilute the contents of bottle No. 3 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Wash solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

All other reagents are ready to use.

TECHNIQUE
For use in automatic equipment, consult SAC (Customer Advisory Service).
 Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
1- Separate the cavities to be used considering: DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper (it is recommended to test in duplicate). Return unused microplate strips to the original sealed packaging.
2- Separate the first cavity for the Blank (OPTIONAL).
3- Add a 3mm disc of DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper, previously perforated, into the determined cavities.

4- Pipette 100 µL of Elution Buffer into all wells, **including the Blank well.**
5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
7- After incubation, discard the contents of the cavities by aspiration (Washer for filter paper technique). Use approximately 300 µL of Washing Solution, **previously prepared**, and perform a total of five (5) washing cycles with 5-second shaking. To ensure the plate dries, at the end of washing, tap the plate for a few seconds on absorbent paper.
Note: Poor washing/drying may cause inadequate results.
8- Pipette 100 µL of Conjugate into all wells, **including the Blank well.**
9- Gently homogenize for ± 30 seconds. Cover the cavities with plate sealer.
10- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
11- Remove the plate sealer from the cavities.
12- Repeat item 7.
13- Pipette 100 µL of Substrate into all wells. **Including in the Blank cavity.**
14- Gently homogenize for ± 30 seconds. Cover the cavities with plate sealer.
15- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
16- Remove the plate sealer from the cavities.
17- Pipette 100 µL of Stop Solution into all wells.
18- Gently homogenize for ± 30 seconds.
19- Read using a double filter: 450 nm / 620 nm in up to 15 minutes (maximum).

TECHNIQUE CHECK
 Check whether the results obtained for reading the Blank and Reference Standards are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	< 0.150
DBS Negative Control	< 0.150
DBS Positive Control	> 0.500

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS
QUALITATIVE
 Calculate Cut Off according to the following formula:

Cut Off = (DBS Negative Control Mean Absorbance x 0.05) + 0.220.

Example:

ITEM	ABSORBANCES
DBS Negative Control	0.040
	0.041
Cut Off = (DBS Negative Control Mean Absorbance x 0.05) + 0.220.	((0.040 + 0.041) / 2 x 0.05) + 0.220 = 0.222

Calculate the Index by dividing the sample absorbance by the Cut Off value. Example:

ITEM	ABSORBANCES
Sample	0.532
Cut Off Value	0.222
Index = Sample / Cut Off Value	0.532 / 0.222 = 2.40

Note: The data presented in the examples are for illustration only and cannot be used to calculate results. Each laboratory must validate the cut-off according to the instrumentation used and the population researched.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (Non-Reactive)	≤ 0.8
Undetermined	0.8 and 1.2
Positive (Reagent)	≥ 1.2

Samples with an index greater than or equal to 1.2 are considered positive, samples with an index less than or equal to 0.8 are considered negative and samples with an index between 0.8 and 1.2 are considered indeterminate. Therefore, the sample mentioned in the example, whose absorbance was 0.532 and index of 2.40, presents a positive result (index ≥ 1.2).

Observations: In case of an indeterminate result, the sample must be re-analyzed in duplicate. Samples that yield repeatedly indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected within two weeks. The result of the last sample collected must prevail. The interpretation of a diagnostic test should not be established based on a single assay. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before descriptive diagnosis of the disease. The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis. The use of dried blood samples on filter paper that have been collected, dried and stored unsatisfactorily may generate inadequate results. The results provided by this kit must be interpreted by the responsible medical professional and are not the only criteria for determining the diagnosis and for patient treatment.

INTERFERENCES

No interference was observed for concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS REACTIVITY

A study was carried out with 48 whole blood samples on filter paper that were negative for CMV IgM, but positive for other infections, in order to evaluate the possibility of cross-reactivity of these interferences in the BIOLISA CMV IgM DBS result. Among them 5 positive samples for HTLV, 6 positive samples for Syphilis, 6 positive samples for HBsAg, 5 positive samples for HCV, 7 positive samples for Chagas Disease, 6 positive samples for Zika, 4 positive samples for Toxoplasmosis, 4 positive samples for Rubella and 5 positive samples for HIV. No cross-reactivity was observed with samples positive for HTLV, Syphilis, HBsAg, HCV, Chagas Disease, Zika, Toxoplasmosis, Rubella and HIV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	3
Average	2.624	0.075	0.026
Standard Deviation	0.053	0.005	0.004
Coefficient of Variation (%)	2.04	6.45	14.38

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Average	2.296	0.074	0.038
Standard Deviation	0.059	0.005	0.005
Coefficient of Variation (%)	2.58	6.68	13.46

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA CMV IgM DBS kit was used for the analysis of clinical samples previously characterized and confirmed using another reference enzyme immunoassay method. The samples analyzed were from neonatal and adult patients. The results show that the clinical sensitivity of the BIOLISA CMV IgM DBS kit is 98.0% and the clinical specificity is 96.4%.

		Reference Result		
		Positive	Negative	Total
BIOLISA CMV IgM DBS	Positive	98	2	100
	Negative	5	137	142
	Total	103	137	242

Clinical Sensitivity: 98.0% (98/100) IC 95% = 95.3 to 100.0 %

Clinical Specificity: 96.5% (137/142) IC 95% = 93.5 to 99.5%

DIAGNOSTIC MEANING

Cytomegalovirus (CMV) is a member of the Herpes Virus family that includes Herpes Simplex Virus (HSV) 1 and 2, Varicella Zoster Virus (VZV), and Epstein-Barr Virus (EBV). It is a human pathogen transmitted through saliva, sexual contact, perinatal contact, organ transplantation or blood transfusion. In most cases, the infection remains asymptomatic. However, CMV infection can cause serious illness in newborns and immunocompromised individuals, such as patients with AIDS, cancer, or patients who have received organ transplants. During immunosuppressive therapy, reactivation of the latent virus or primary infection often occurs. CMV infections can be acquired before birth, during childbirth, and later in life. It can cause serious congenital anomalies, such as microcephaly, motor disability and mental retardation. Therefore, the primary determination of maternal infection and its distinction from latent infection are of great importance. The presence of IgM antibodies indicates the presence of primary infection, while the presence of IgG antibodies indicates the immunological status of the patients. Therefore, the Biolisa CMV IgM DBS kit is a valuable tool in monitoring Pre-Natal and Neonatal screening.

BIBLIOGRAPHIC REFERENCE

1. Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
2. Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after "Silent" Congenital Cytomegalovirus Infection. N. Engl. J. Med. (1976) 295:468-470.
3. Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovovirus. N. Engl. J. Med. (1974) 790:291-296.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA CMV IgM DBS kit: 10269360447

Review: February/2024

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER