BIOLISA CMV IgM

REF K123

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FINALIDADE

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REFERÊNCIA BIBLIOGRÁFICA
1. Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
2. Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after "Silent" Congenital Cytomegalovirus Infection. N. Engl. J. Med.(1976) 295:468-470.
3. Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovirus. N. Engl. J. Med. (1974) 790:291-296.
4. Stern, H. Cytomegalovirus Vaccine: Justification and Problems. In: Waterson AP (ed.) Recent Advances in Clinical Virology (1977) 117-134.
5. TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
- Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after "Silent" Congenital Cytomegalovirus Infection. N. Engl. J. Med.(1976) 295:468-470.
- Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovurus. N. Engl. J. Med. (1974) 790:291-296.
- Stern, H. Cytomegalovirus Vaccine: Justification and Problems. In: Waterson AP (ed.) Recent Advances in Clinical Virology (1977) 117-134.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Test for qualitative determination of IgM antibodies to Cytomegalovirus (CMV) in human samples of serum, plasma and whole blood on microplate enzyme immunoassay filter paper. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA CMV IgM kit is a solid phase immunoenzymatic assay based on the principle of qualitative detection by capture of IgM antibodies to CMV in human serum, plasma and whole blood samples on filter paper. IgM antibodies present in the sample bind to the anti IgM antibodies coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Peroxidase-conjugated CMV antigens are added to the microplate which is then incubated. The enzyme-conjugated antigens bind to the IgM Anti-CMV antibodies present, attached to the plate coated with Anti IgM antibodies. A new wash is performed to remove the surplus.
After this stage, the Substrate is added and incubated, producing a blue color that indicates the amount of IgM Anti-CMV Antibodies present in the sample. The Stop Solution is added to stop the reaction with a color change from blue to yellow, measured on a microplate reader.

REAGENTS

- 1 - Sensitized Plate** - Store at 2 and 8°C. Contains: Plate sensitized with Anti-IgM antibodies.
2- Conjugate - Store between 2 and 8°C. Contains: Buffer Solution, Peroxidase-conjugated CMV antigen, surfactant, stabilizers, dye and preservative.
3- Concentrated Washing - Store at 2 and 8°C. Contains: Buffer Solution, surfactant and preservative.
4- Sample Diluent - Store at 2 and 8°C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
5- Substrate - Store at 2 and 8°C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
6- Stop Solution - Store at 2 and 8°C. Contains: Hydrochloric Acid 1 M.
7- Negative Control - Store at 2 and 8°C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infectious.**
8- Positive Control - Store between 2 and 8°C. Contains: Antibodies IgM Anti-CMV, Buffer Solution, dye, stabilizers, surfactant and preservative. **Potentially infectious.**
9- Plate Sealers.
10- Filter Paper Elution Buffer - Store at 2 and 8°C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
11- Negative Extraction Control - Store at 2 and 8°C. Contains: Non-reactive sample for Anti-CMV IgM antibodies impregnated in Filter Paper. **Potentially infectious.**
12- Positive Extraction Control - Store at 2 and 8°C. Contains: Reactive sample for Anti-CMV IgM antibodies impregnated in filter paper. **Potentially infectious.**

PRESENTATION

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1 – Sensitized Plate	1 Unit (96 Cavities)	2 Units (192 Cavities)	5 Units (480 Cavities)
2 - Conjugate	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
3 –Concentrate Washing	1 Flask x 50mL	2 Flasks x 50mL	5 Flasks x 50mL
4 – Sample Diluent	1 Flask x 42mL	2 Flasks x 42mL	5 Flasks x 42mL
5 – Substrate	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
6 – Stop Solution	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
7 –Negative Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
8 –Positive Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
9 – Plate Sealers	3 Units	6 Units	15 Units
10 – Filter Paper Elution Buffer	1 Flask x 20 mL	2 Flasks x 20 mL	5 Flasks x 20 mL
11 – Negative Extraction Control	1 unit	2 units	5 units
12 – Positive Extraction Control	1 unit	2 units	5 units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table

Required materials not contained in the Kit:

- 1-** Pipette capable of dispensing volumes from 5 to 300 µL with lower coefficient of variation than 1.5%.
2- Re-pipettor for repetitive pipetting volumes of 300 µL with lower coefficient of variation than 1.5% or multicontrol pipette (Optional).
3- Microplate washer (optional).
4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
5- Paper towel to dry cavities
6- Stopwatch or watch.
7- Flask to store the Washing Solution after diluted.
8- Distilled or deionized water.
9- Tools of Quality Control.
10- Incubator 37°C ± 2°C.
11- Paper shredder (3 mm diameter) for filter paper technique.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
2- Strictly follow the methodology proposed to obtain accurate results.
3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store at 2 to 8°C.
4- The water used in material cleaning must be recent and free of contaminants.

- 5-** Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can signifi cantly alter the results.
6- Stop Solution contains Chloridric Acid, which is a strong acid. Handle it with care.
7- All the raw material of product is tested and should be nonreactive for HBsAg, HIV 1 & 2 and Anti HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the biosafety of these products.
8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
9- As a safety measure, you should cover the plate during the reaction.
10- You must ensure that the bottom of the cavity is clean and dry and there are no bubbles on the surface fl uid before reading the plate. Do not let the cavities run dry during the test.
11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
13- To obtain information related to biosafety or in case of accidents with he product, consult the SDS (Safety Data Sheet) available on the website www. bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
14- Do not use the product in case of damaged packaging.
15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

Whole Blood on Filter Paper

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Samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8°C. For storage longer than 2 years, samples must be stored at a temperature of -20°C.

DESCRIPTION OF PROCESS

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of the Reagent No3 (Concentrated Washing Solution) in 1000 mL of distilled water or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

STABILITY AFTER OPEN

The results of the stability test show that the Biolisa kit CMV IgM is stable after opening for up to 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

TECHNIQUE

Samples of Serum or Plasma

Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.

- 1-** Select the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used, for the original sealed packaging.
2- Select the first cavity for Blank (OPTIONAL).
3- Pipette 100 µL of Sample Diluent, included in the Blank cavity.
4- Pipette 5 µL of Controls and Samples into previously determined wells. Observe the color change of the diluent at the time of sample addition. The color change indicates that the sample was added to the well correctly.
5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
7- Remove the sealing of wells.
8- Discard the contents of the wells by aspiration (Washer) or by decanting (manual).
Use approximately 300 µL of Washing Solution, previously prepared to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
Note: Poor washing and drying can cause inadequate results.
9- Pipette 100 µL of Conjugate into each well except in the Blank cavity (if you made this option).
10- Mix gently for ± 30 seconds. Cover cavities with plate sealer.
11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
12- Remove the sealer from plate cavities.
13- Repeat item 8.
14- Pipette 100 µL of Substrate into all wells, included in the Blank cavity.
15- Shake gently for ± 30 seconds. Cover wells with plate sealer.
16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
17- Remove the sealer from plate cavities.
18- Pipette 100 µL of Stop Solution into all wells, included in the Blank cavity.
19- Mix gently for ± 30 seconds.
20- Read using double filter: 450nm/630nm within up to 15 minutes (maximum).

Whole Blood Sample on Filter Paper

Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.

- 1-** Select the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.
2- Select the first cavity for Blank (OPTIONAL).
3- Add a 3 mm disk of Positive Extraction Control, Negative Extraction Control and Whole Blood Sample on Filter Paper, previously perforated, in the determined wells.
4- Pipette 100 µL of Elution Buff er in Filter Paper, included in the Blank cavity.
5- Pipette 5 µL of Controls into previously determined wells.
6- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
7- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
8- Remove the sealing of wells.
9- Discard the contents of the wells by decanting (manual). Remove the paper discs, if necessary, with the aid of a needle or tip.
Use approximately 300 µL of Washing Solution, previously prepared to perform a total of five (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
Note: Poor washing and drying can cause inadequate results.
10- Pipette 100 µL of Conjugate into each well, included in the Blank cavity.
11- Mix gently for ± 30 seconds. Cover cavities with plate sealer.
12- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
13- Remove the sealer from plate cavities.
14- Repeat item 9.

- 15- Pipette 100 µL of Substrate into all wells, included in the Blank cavity.
16- Shake gently for ± 30 seconds. Cover wells with plate sealer.
17- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
18- Remove the sealer from plate cavities.
19- Pipette 100 µL of Stop Solution into all wells, included in the Blank cavity.
20- Mix gently for ± 30 seconds.
21- Read using double filter: 450nm/630nm within up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Samples of Serum, Plasma and Filter Paper Whole Blood
Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	< 0.150
Negative Control	< 0.150
Positive Control	> 1.000
Negative Extraction Control	< 0.300
Positive Extraction Control	> 1.000

If the values are out of the expected values, you must repeat the technique.

DESCRIPTION OF CALCULATIONS

QUALITATIVE

Sample of Serum or Plasma

Calculate Cut Off according to the formula below:

Cut Off = (Positive Control Mean Abs. x 0.1) + 0.300

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N°8)	A1 = 1.964
	A2= 1.960
Cut Off = (Positive Control Mean Absorbance x 0.1) + 0.300	Cut Off = ((1.964 + 1.960)/2 x 0.1) + 0.300 Cut Off = 0.496

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1.279
Cut Off Value	0.496
Index = Sample/ Cut Off Value	1.279 / 0.496 = 2.579

Filter Paper Whole Blood Sample

Calculate Cut Off according to the formula below:

Cut Off = (Positive Control Mean Absorbance x 0.05) + 0.120

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N°8)	A1 = 1.964
	A2 = 1.960
Cut Off = (Positive Control Mean Absorbance x 0.05) + 0.120	Cut Off = ((1.964 + 1.960)/2 x 0.05) + 0.120 Cut Off = 0.218

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1.279
Cut Off Value	0.218
Index = Sample/ Cut Off Value	1.279 / 0.218 = 5.867

INTERPRETATION OF RESULTS

Samples of Serum, Plasma and Filter Paper Whole Blood

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	ABSORBANCE
	INDEX
Negative (No Reagent)	< 0.9
Undetermined	0.9 – 1.1
Positive (Regent)	> 1.1

The results provided by this kit should be interpreted by the responsible medical professional and not the only criterion for determining the diagnosis and/or treatment of the patient.

PROCEDURE LIMITATIONS

IgM antibodies are present in the serum, plasma and whole blood for a short period of time, disappearing up to six months after infection, while IgG antibodies remain present for a long period of time. However, the high sensitivity of the new diagnostic methodologies, in some cases, allows to find very low levels of IgM antibodies, called residual IgM, for a longer period of time. Detection of these residual IgM antibodies hinders clinical interpretation over the period of infection. In these cases, to confirm the result, an IgG avidity test is recommended. The interpretation of a diagnostic test, should not be based on a single run. This should include confirmation of other tests before that a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other clinical information available before the descriptive diagnosis of the disease.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE

ACCURACY

Repeatability

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	3
Average	4.067	1.267	0.083
Standard Deviation	0.015	0.025	0.011
Coefficient of Variation (%)	0.360	2.013	13.426

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Average	3.943	1.259	0.073
Standard Deviation	0.160	0.005	0.001
Coefficient of Variation (%)	4.069	0.377	1.821

CLINICAL SENSITIVITY AND SPECIFICITY

Serum and Plasma Sample

BIOLISA CMV IgM Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA CMV IgM kit is > 99.99% and clinical specificity is > 99.99%.

	Expected Result	Biolisa CMV IgM
Positive Sample	53	53
Negative Sample	94	94
Total Tested Sample	147	

Clinical sensitivity: > 99.99% (53/53)

Clinical specificity: > 99.99% (94/94)

Whole Blood Samples on Filter Paper

BIOLISA CMV IgM Kit analyzed clinical samples in comparison with other methods of EIA. The results show that the clinical sensitivity of the BIOLISA CMV IgM kit is > 99.99% and clinical specificity is > 99.99%.

	Expected Result	Biolisa CMV IgM
Positive Sample	50	50
Negative Sample	52	52
Total Tested Samples	102	

Clinical sensitivity: > 99.99% (50/50)

Clinical specificity: > 99.99% (52/52)

DIAGNOSTIC SIGNIFICANCE

Cytomegalovirus (CMV) is a member of the Herpes Virus family that includes Herpes Simplex viruses (HSV) 1 and 2, Varicella Zoster virus (VZV) and Epstein-Barr virus (EBV). It is a human pathogen transmitted through saliva, sexual contact, perinatal, organ transplantation or blood transfusion. In most cases, the infection remains asymptomatic.

However, CMV infection can cause serious illness in newborns and immunocompromised individuals, such as patients with AIDS, cancer or patients who have received organ transplants. During immunosuppressive therapy, latent virus reactivation or primary infection often occurs. CMV infections can be acquired before birth, during childbirth and later in life. It can cause serious congenital abnormalities, such as microcephaly, motor disability and mental retardation. Therefore, the primary determination of maternal infection and its distinction from latent infection are of great importance.

The presence of IgM antibodies indicates the presence of primary infection, while the presence of IgG antibodies indicates the patients' immune status.

BIBLIOGRAPHIC REFERENCES

- Hodinka, RL, and Friedman, HM. Human Cytomegalovirus. In: Manual of Clinical Microbiology 6th Edition (1995) 884-894.
- Hanshaw, JB, Scheiner, AP, Moxley, AW, Gaev, L, Abel, V, and Scheiner, B. School Failure and Deafness after "Silent" Congenital Cytomegalovirus Infection. N. Engl. J. Med. (1976) 295:468-470.
- Reynolds, DW, Stagno, S, Stubbs, KG, Dabte, AJ, Livingston, NM, Saxon, SS, Alford, CA. Inapparent Congenital Cytomegalovirus. N. Engl. J. Med. (1974) 790:291-296.
- Stern, H. Cytomegalovirus Vaccine: Justification and Problems. In: Waterson AP (ed.) Recent Advances in Clinical Virology (1977) 117-134.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA CMV IgM kit: 10269360189

Review: August/2024

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER