BIOLISA COVID-19 Ac TOTAL

REF K236-1

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FINALIDADE

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Departamento de Pesquisa e Desenvolvimento.

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods;approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Departamento de Pesquisa e Desenvolvimento.

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REF **K236-1**

INSTRUCTIONS FOR USE

FUNCTION
Test for qualitative determination of total antibodies (IgG, IgM and IgA) for SARS-CoV-2 (virus that causes COVID-19 disease) in biological samples (serum, plasma or whole blood on filter paper) through enzyme immunoassay test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION
Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA COVID-19 TOTAL Ac kit is an immunoenzymatic assay in solid phase based on the "sandwich" principle for the detection of total antibodies (IgG, IgM and IgA) for SARS-CoV-2 in serum, plasma and human whole blood on paper filter. IgG, IgM and IgA antibodies present in the sample bind to the recombinant SARS-COV-2 antigens coated on the microplate forming immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Recombinant antigens conjugated to peroxidase are added to the microplate which is then incubated. The enzyme-conjugated antigens bind to the anti-SARS-CoV-2 IgG, IgM and IgA antibodies present, attached to the antigen-coated plate. A new wash is performed to remove the surplus. After this stage, the Substrate is added and incubated, producing a blue color that indicates the amount of IgG, IgM and IgA anti-SARS-CoV-2 antibodies present in the sample. The Stop Solution is added to stop the reaction with a color change from blue to yellow, measured on a microplate reader.

REAGENTS
1- Sensitized Plate - Store between 2 and 8 °C. Contains: Plate Sensitized with recombinant SARS-CoV-2 antigen.
2- Conjugate - Store at 2 to 8 °C. Contains: Buffer Solution containing recombinant SARS-CoV-2 antigen linked to Peroxidase, surfactant, stabilizers, dye and preservative.
3- Concentrated Washing Solution - Store between 2 and 8 °C. Contains: Buffer Solution, surfactant and preservative.
4- Sample Diluent - Store at 2 to 8 °C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
5- Substrate - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
6- Stop Solution - Store at 2 to 8 °C. Contains: 1M Hydrochloric Acid.
7- Negative Control - Store between 2 and 8 °C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infective.**
8- Positive Control - Store between 2 and 8 °C. Contains: Buffer Solution, IgG, IgM and IgA anti-SARS-CoV-2 antibodies, dye, stabilizers, surfactant and preservative. **Potentially infective.**
9- Plate Sealers
10- Filter Paper Elution Buffer - Store at 2 to 8 °C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
11- Negative Extraction Control - Store between 2 and 8 °C. Contains: Non-reactive sample for IgG, IgM and IgA anti-SARS-CoV-2 antibodies impregnated on filter paper. **Potentially infective.**
12- Positive Extraction Control - Store between 2 and 8 °C. Contains: Reactive sample for IgG, IgM and IgA anti-SARS-CoV-2 antibodies impregnated in filter paper. **Potentially infective.**

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 – Sensitized Plate	1 Unit (96 cavities)	2 Units (192 cavities)	5 Units (480 cavities)
2 – Conjugate	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
3 – Concentrated Washing	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
4 – Sample Diluent	1 Flask x 42 mL	2 Flasks x 42 mL	5 Flasks x 42 mL
5 - Substrate	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
9 – Stop Solution	1 Flask x 12mL	2 Flasks x 12mL	5 Flasks x 12mL
7 –Negative Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
8 –Positive Control	1 Flask x 300 µL	2 Flasks x 300 µL	5 Flasks x 300 µL
9 – Plate Sealers	3 Units	6 Units	15 Units
10 – Filter Paper Elution Buffer	1 Flask x 12 mL	2 Flasks x 12 mL	5 Flasks x 12 mL
11 – Negative Extraction Control	1 unit	2 units	5 units
12 – Positive Extraction Control	1 unit	2 units	5 units

EQUIPMENTS AND OPERATIONAL INPUTS
Materials in the kit:
- Reagents described in the above table
- Instructions for use (manual)

Necessary materials not contained in the kit:
1- Pipettes capable of dispensing volumes of 5 to 500 µL with a variation coefficient of less than 1.5%.
2- Repipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
3- Microplate washer (optional).
4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
5- Absorbent paper to dry the wells.
6- Stopwatch or watch.
7- Bottle to store the Washing Solution after dilution.
8- Distilled or deionized water.
9- Quality Control Tools.
10- Incubator at 37 °C ± 2 °C.
11- Paper shredder (3 mm diameter) for filter paper technique.

TRANSPORTATION AND STORAGE CONDITIONS
The storage and transport temperature should be 2 and 8 °C. Keep out of the light and avoid humidity. **Do not freeze.**

SPECIAL CARE
1- For *in vitro* diagnostic use only.
2- Strictly follow the proposed methodology to obtain exact results.
3- The sachet containing the microplate must be opened only after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
4- The water used to clean the material must be fresh and free from contaminants.
5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.

7- All raw material of the product is tested and must be non-reactive for HBsAg, Anti-HIV 1 & 2 and Anti HCV. However, these tests do not offer complete security in the absence of infectious agents. The handling of any product that contains serum is potentially capable of transmitting disease. Therefore, it is necessary to take due care of biosafety in the handling of these products.
8- Pipette the reagents always in the same order to minimize the difference in reaction time between the wells.
9- As a protection measure, the plate must be covered during the reaction.
10- Make sure that the bottom of the cavity is clean and dry and that there are no bubbles on the liquid surface before reading the plate. Do not allow the wells to dry out during the test.
11- Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.
12- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.
13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Service Customer Service) of Quibasa.
14- Do not use the product in case of damage to the packaging.
15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES
Serum or Plasma (EDTA or Heparin)
Hemolyzed or highly lipemic samples should not be used.The samples can be kept refrigerated, between 2 and 8 °C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at a temperature of -20 °C.

Whole Blood on Filter Paper (EDTA or Puncture)
Dry samples on filter paper can be stored at room temperature as long as they are protected from direct sunlight and low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 to 8 °C. For storage for more than 2 years, the samples must be stored at a temperature of -20 °C.

PROCESS DESCRIPTION
Stability After Open
The results of the stability test show that the BIOLISA COVID-19 TOTAL Ac kit is stable after being opened for at least 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using internal controls of the kit and the criteria for validating the technique.

Preparation of Work Reagents
Washing Solution
Dilute the contents of vial N° 3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for at least 30 days. It can be stored at room temperature. If crystallization occurs, heat to 37 °C until dissolved.

Substrate
The Substrate is ready for use.

TECHNIQUE
For use in automatic equipment, consult the Customer Service Department (SAC)

Serum or Plasma samples
Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.
1- Separate the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the unused strips of the microplate to the original sealed package.
2- Separate the first cavity for Blank (OPTIONAL).
3- Pipette 100 µL of Sample Diluent, including in the Blank cavity.
4- Pipette 10 µL of Sample and Controls into the previously determined wells. Observe the color change of the diluent when adding the sample. The color change indicates that the sample has been added to the well correctly.
5- Homogenize gently for ± 30 seconds and cover the wells with plate sealer.
6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
7- Remove the sealer from the cavities.
8- Discard the contents of the cavities by aspiration (Washer) or by decantation (Manual).
Use approximately 300 µL of Wash Solution, previously prepared, and perform a total of five (5) wash cycles. To guarantee the drying of the plate, at the end of the washing, beat the plate for a few seconds on absorbent paper.
Note: Poor washing / drying can cause inadequate results.
9- Pipette 100 µL of Conjugate into all wells, including the Blank cavity (OPTIONAL).
10- Homogenize gently for ± 30 seconds. Cover the cavities with the plate sealer.
11- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
12- Remove the plate sealer from the cavities.
13- Repeat item 8.
14- Pipette 100 µL of Substrate into all wells.
15- Gently homogenize for ± 30 seconds. Cover the cavities with the plate sealer.
16- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
17- Remove the plate sealer from the cavities.
18- Pipette 100 µL of Stop Solution into all wells.
19- Gently mix for ± 30 seconds.
20- Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

Samples of Whole Blood on Filter Paper
Before starting the test, place all reagents, samples and controls to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.
1- Separate the cavities to be used considering: Positive Control, Negative Control, Positive Extraction Control, Negative Extraction Control and Whole Blood Samples on Filter Paper (it is recommended to test in duplicate). Return the unused strips of the microplate to the original sealed package.
2- Separate the first cavity for Blank (OPTIONAL).
3- Add a 3 mm disc of Positive Extraction Control, Negative Extraction Control and Whole Blood Samples on Filter Paper, previously perforated, in the determined wells.
4- Pipette 100 µL of the Filter Paper Elution Buffer, including in the well for the Blank.
5- Pipette 10 µL of Controls into the previously determined wells.
6- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
7- Incubate for 60 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
8- Remove the sealer from the cavities.
9- Discard the contents of the cavities by aspiration (Washer) or by decantation (Manual). Remove the paper discs, if necessary, with the aid of a needle or tip. Use approximately 300 µL of Wash Solution, previously prepared, and perform a total of five (5) wash cycles. To guarantee the drying of the plate, at the end of the washing, beat the plate for a few seconds on absorbent paper.
Note: Poor washing / drying can cause inadequate results.
10- Pipette 100 µL of Conjugate into all wells, including the Blank cavity (OPTIONAL).
11- Homogenize gently for ± 30 seconds. Cover the cavities with the plate sealer.
12- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
13- Remove the plate sealer from the cavities.
14- Repeat item 9.

- 15- Pipette 100 µL of Substrate into all wells.
16- Gently homogenize for ± 30 seconds. Cover the cavities with the plate sealer.
17- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
18- Remove the plate sealer from the cavities.
19- Pipette 100 µL of Stop Solution into all wells.
20- Gently mix for ± 30 seconds.
21- Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Serum, Plasma and Whole Blood on Filter Paper samples

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	<0.250
Negative Control	<0.250
Positive Control	> 1.000
Negative Extraction Control	<0.250
Positive Extraction Control	> 1.000

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS

QUALITATIVE

Serum and Plasma samples

Calculate Cut Off according to the following formula:
Cut Off = (Average Positive Control Absorbance x 0.05) + 0.060

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N° 8)	A1 = 2.453 A2 = 2.470
Cut Off = (Average Positive Control Absorbance x 0.05) + 0.060	((2.453 + 2.470) / 2) x 0.05 + 0.060 = 0.183

Calculate the Index by dividing the Sample's absorbance by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	1.245
Cut Off Value	0.183
Index = Sample / Cut Off Value	1.245 / 0.183 = 6.70

Whole blood on filter paper samples

Calculate Cut Off according to the following formula:
Cut Off = (Average Positive Control Absorbance x 0.05) + 0.040

Example:

ITEM	ABSORBANCE
Positive Control (Reagent N° 8)	A1 = 2.325 A2 = 2.390
Cut Off = (Average Positive Control Absorbance x 0.05) + 0.040	((2.235 + 2.390) / 2) x 0.05 + 0.040 = 0.156

Calculate the Index by dividing the Sample's absorbance by the Cut Off value.

Example:

ITEM	ABSORBANCE
Sample	0.752
Cut Off Value	0.156
Index = Sample / Cut Off Value	0.752 / 0.156 = 4.82

Note: The data presented in the examples are for illustration only and cannot be used to calculate the results.

INTERPRETATION OF RESULTS

Serum, Plasma and Whole Blood on Filter Paper Samples

After calculating the sample index, consider the indices below to determine the results.

RESULTS	QUALITATIVE FOR SERUM SAMPLE, PLASMA AND TOTAL BLOOD ON FILTER PAPER
	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	0.9 – 1.1

Note: In the case of an undetermined result, the sample must be re-analyzed. Samples that obtain repeatedly indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected in two weeks. The result of the last collected sample must prevail. The results provided by this kit must be interpreted by the responsible medical professional, and it is not the only criterion for determining the diagnosis and / or treatment of the patient.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before the descriptive diagnosis of the disease.

INTERFERENT

No interference was observed for Triglycerides up to 1500 mg/dL, Acetylsalicylic Acid up to 20 mg/dL, Ascorbic Acid up to 2 g/dL, Creatine up to 200 mg/dL, Bilirubin up to 1 g/dL, Albumin up to 2 g/dL, Hemoglobin up to 1000 mg/dL, Oxalic Acid up to 60 mg/dL, Rheumatoid Factor up to 980 IU/mL, C Reactive Protein up to 41.2 mg/dL and Anti Streptolysin O up to 1023 IU/mL.

CROSS REACTIVITY

Serum and Plasma Cross Reactivity

A cross-reactivity study was carried out, evaluating 162 samples negative for SARS-CoV-2, but positive for other infections. Among them 12 positive samples for Influenza, 27 positive for Rhinovirus, 13 positive samples for Respiratory Syncytial Virus, 14 positive samples for Cytomegalovirus (CMV), 14 positive samples for Human Immunodeficiency Virus (HIV), 10 positive samples for HBV, 12 samples positive for HCV, 10 positive samples for Zika, 14 positive samples for Dengue, 10 positive samples for Chikungunya, 13 positive samples for Toxoplasmosis and 13 positive samples for Rubella. No crossreactivity was observed with the samples evaluated. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

Whole Blood on Filter Paper Cross Reactivity

A cross-reactivity study was carried out, evaluating 81 samples negative for SARS-CoV-2, but positive for other infections. Among them 10 positive samples for Cytomegalovirus (CMV), 5 positive samples for Human Immunodeficiency Virus (HIV), 9 positive samples for HBV, 10 positive samples for HCV, 7 positive samples for Zika, 10 positive samples for Toxoplasmosis, 10 positive samples for Rubella, 10 positive samples for Syphilis and 10 positive samples for Chagas Disease. No cross reactivity was observed with the samples evaluated. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

PRUDCT PERFORMANCE

ACCURACY

Repeitbility

Reproducibility was calculated from 10 successive determinations using 3 samples with different values, obtaining the following absorbance results:

Repeitbility	Sample		
	1	2	3
Average	0.906	0.955	0.044
Standard Deviation	0.079	0.153	0.006
Coefficient of Variation (%)	8.75	16.07	14.37

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	1.021	1.013	0.043
Standard Deviation	0.178	0.149	0.006
Coefficient of Variation (%)	17.42	14.67	13.49

CLINICAL SENSITIVITY AND SPECIFICITY

Serum and Plasma samples

The BIOLISA COVID-19 TOTAL Ac kit was used to analyze clinical samples. The results show that the BIOLISA COVID-19 TOTAL Ac kit has a clinical sensitivity of 96.89% and a clinical specificity of 97.63%.

	Expected Result	BIOLISA COVID-19 Ac TOTAL
Positive Sample	354	343
Negative Sample	295	288
Total Tested Sample	649	631

Clinical Sensitivity: 96.89% (343/354) - IC95%: 94.5% - 98.4%
Clinical Specificity: 97.63 (288/295) - IC95%: 95.2% - 99.0%

Whole blood on filter paper samples

The BIOLISA COVID-19 TOTAL Ac kit was used to analyze clinical samples. The results show that the BIOLISA COVID-19 TOTAL Ac kit has a clinical sensitivity of 93.62% and a clinical specificity of >99.9%.

	Expected Result	BIOLISA COVID-19 Ac TOTAL
Positive Sample	47	44
Negative Sample	162	162
Total Tested Sample	209	206

Clinical Sensitivity: 93.62% (44/47) - IC95%: 82.5% - 98.7%
Clinical Specificity: >99.9% (162/162) - IC95%: 97.7% - 100.0%

CLINICAL SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family, and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are: fever, cough, dyspnoea, myalgia or fatigue. It can cause diseases ranging from the common cold to more serious illnesses, such as the Middle East Respiratory Syndrome (MERS-CoV) and the Severe Acute Respiratory Syndrome (SARS-CoV). In December 2019, a series of pneumonia cases of unknown cause appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples it was shown to be Coronavirus infection, which was called the new coronavirus (2019 -nCoV). In approximately 150 countries on all continents, cases of infection have already been reported.

BIBLIOGRAPHIC REFERENCES

- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods;approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019- nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration number for BIOLISA COVID-19 Ac Total kit: 10269360339

Review: August/2024

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER