# Bioclin

## BIOLISA COVID-19 IgA

REF **K235-1**

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#### FINALIDADE

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
- Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
- Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497-506.
- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
- QUIBASA: Departamento de Pesquisa e Desenvolvimento.

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### USE INSTRUCTIONS

#### FUNCTION

Test for qualitative determination of IgA antibodies to SARS-CoV-2 (virus that causes COVID-19 disease) in biological samples (serum, plasma or whole blood on filter paper) through enzyme immunoassay test. Only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Enzyme immunoassay or immunoenzymatic  
The BIOLISA COVID-19 IgA kit is an immunoenzymatic assay in solid phase based on the principle of qualitative detection of IgA antibodies for SARSCoV-2 in serum, plasma and whole blood on filter paper human. IgA antibodies present in the sample bind to the SARS-COV-2 recombinant antigens coated on the microplate forming immunocomplexes. After incubation the microplate is washed to remove unbound materials. Peroxidase-conjugated Anti-IgA antibodies are added to the microplate which is then incubated. The enzyme-conjugated antibodies bind to the IgA Anti-SARS-CoV-2 Antibodies present, attached to the plate coated with anti-IgA antibodies. New wash is performed to remove surpluses. After this stage, the substrate is added and incubated, producing a blue color that indicates the amount of IgA Antibodies SARS-CoV-2 present in the sample. Stop Solution is added to stop the reaction with a color change from blue to yellow, measured in a microplate reader.

#### REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Sensitized Plate with recombinant SARS-CoV-2 antigen and preservative.
- 2- Conjugate** - Store between 2 and 8 °C. Contains: Buffer Solution containing human anti-IgA antibody conjugated to peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** - Store at 2 to 8 °C. Contains: Buffer Solution (Phosphate <0.5 mol/L, Potassium Chloride <100 mmol/L, Sodium Chloride <5 mol/L), surfactant and preservative.
- 4- Sample Diluent** - Store between 2 and 8 °C. Contains: Buffer Solution, stabilizers, surfactant and preservative.
- 5- Substrate** - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8 °C. Contains: Hydrochloric Acid 1 M.
- 7- Negative Control** - Store between 2 and 8 °C. Contains: Buffer Solution, stabilizer, surfactant and preservative. **Potentially infective.**
- 8- Positive Control** - Store between 2 and 8 °C. Contains: Buffer Solution, IgA Anti - SARS-CoV-2 antibodies, dye, stabilizers, surfactant and preservative. **Potentially infective.**
- 9- Plate Sealers**
- 10- Filter Paper Elution Buff er** - Store between 2 and 8 °C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 11- Negative Extraction Control** - Store between 2 and 8 °C. Contains: Non-reactive sample for IgA Anti-SARS-CoV-2 antibodies impregnated in Filter Paper. **Potentially infective.**
- 12- Positive Extraction Control** - Store between 2 and 8 °C. Contains: Reactive sample for IgA Anti - SARS-CoV-2 antibodies impregnated in Filter Paper. **Potentially infective.**

#### PRESENTATION

| REAGENTS                         | 1                       | 2                         | 3                         |
|----------------------------------|-------------------------|---------------------------|---------------------------|
|                                  | 96 Cavities             | 192 Cavities              | 480 Cavities              |
| 1 – Sensitized Plate             | 1 Unit<br>(96 cavities) | 2 Units<br>(192 cavities) | 5 Units<br>(480 cavities) |
| 2 – Conjugate                    | 1 Flask x 12 mL         | 2 Flasks x 12 mL          | 5 Flasks x 12 mL          |
| 3 – Concentrate Washing          | 1 Flask x 50 mL         | 2 Flasks x 50 mL          | 5 Flasks x 50 mL          |
| 4 – Sample Diluent               | 1 Flask x 42 mL         | 2 Flasks x 42mL           | 5 Flasks x 42 mL          |
| 5 - Substrate                    | 1 Flask x 12mL          | 2 Flasks x 12mL           | 5 Flasks x 12mL           |
| 9 – Stop Solution                | 1 Flask x 12mL          | 2 Flasks x 12mL           | 5 Flasks x 12mL           |
| 7 –Negative Control              | 1 Flask x 300 µL        | 2 Flasks x 300 µL         | 5 Flasks x 300 µL         |
| 8 –Positive Control              | 1 Flask x 300 µL        | 2 Flasks x 300 µL         | 5 Flasks x 300 µL         |
| 9 – Plate Sealers                | 3 Units                 | 6 Units                   | 15 Units                  |
| 10 – Filter Paper Elution Buffer | 1 Flask x 12 mL         | 2 Flasks x 12 mL          | 5 Flasks x 12 mL          |
| 11 – Negative Extraction Control | 1 unit                  | 2 units                   | 5 units                   |
| 12 – Positive Extraction Control | 1 unit                  | 2 units                   | 5 units                   |

#### EQUIPMENTS AND OPERATIONAL INPUTS

##### Materials in the kit:

- Reagents described in the above table
- Operating instructions (manual)

##### Required materials not contained in the kit:

- 1- Pipette capable of dispensing volumes of 5 to 500 µL with lower coefficient of variation than 1.5%.
- 2- Re-pipettor for repetitive pipetting volumes of 500 µL, with lower coefficient of variation than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Paper towel to dry cavities
- 6- Stopwatch or watch.
- 7- Flask to store the washing solution after diluted.
- 8- Distilled or deionized water.
- 9- Tools of Quality Control.
- 10- Incubator 37 °C ± 2 °C.
- 11- Paper shredder (3 mm diameter) for filter paper technique.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature should be 2 and 8 °C. Keep out of the light and avoid humidity.**Do not freeze.**

#### SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can significantly alter the results.
- 6- Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.
- 7- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.

- 8- Always add reagents in the same order to minimize the difference in reaction time between the cavities.
- 9- As a safety measure, you should cover the plate during the reaction.
- 10- You must ensure that the bottom of the cavity is clean and dry, and there are no bubbles on the surface fl uid before reading the plate. Do not let the cavities run dry during the test.
- 11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

##### Serum or Plasma (EDTA or Heparin)

Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8 °C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20 °C.

##### Whole Blood on Filter Paper (Puncture or EDTA)

Dry samples on filter paper can be stored at room temperature as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 to 8 °C. For storage for more than 2 years, the samples must be stored at a temperature of -20 °C.

#### PROCESS DESCRIPTION

##### Stability After Open

The results of the stability test show that the Biolisa kit COVID-19 IgA is stable after opening for at least 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

##### Washing Solution

Dilute the contents of the Reagent N° 3 (Concentrate Washing) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30 °C until expiration date printed on the original bottle. In case of crystallization, heat it at 37 °C until dissolution.

##### Substrate

The Substrate is ready for use.

#### TECHNIQUE

##### Samples of Serum or Plasma

- Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
- 1- Select the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the strips of the microplate will not be used, for the original sealed packaging.
- 2- Select the fi rst cavity for Blank (OPTIONAL).
- 3- Pipette 100 µL of Sample Diluent, included in the cavity for Blank.
- 4- Pipette 5 µL of Sample and Controls into the previously determined wells.. Observe the color change of the diluent when adding the sample. A Color change indicates that the sample was successfully added to the well.
- 5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 7- Remove the sealing of wells.
- 8- Discard the contents of the wells by aspiration (Washer) or by decanting (manual). Use approximately 300 µL of Washing Solution, previously prepared to perform a total of fi ve (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
- Note:** Poor washing/drying can cause inadequate results.

- 9- Pipette 100 µL of Conjugate into each well.
- 10- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 12- Remove the sealer from plate cavities.
- 13- Repeat item 8.
- 14- Pipette 100 µL of Substrate into all wells.
- 15- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 16- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 17- Remove the sealer from plate cavities.
- 18- Pipette 100 µL of Stop Solution into all wells.
- 19- Mix gently for ± 30 seconds.
- 20- Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

##### Whole Blood Sample on Filter Paper

- Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.
- 1- Select the wells to be used considering Positive Extraction Control, Negative Extraction Control and Whole Blood Samples on Filter Paper (it is recommended to test in duplicate). Return the strips of the microplate will not be used for the original sealed packaging.
- 2- Select the first cavity for Blank (OPTIONAL).
- 3- Add a 3 mm disk of Positive Extraction Control, Negative Extraction Control and Whole Blood Sample on Filter Paper, previously perforated, in the determined wells.
- 4- Pipette 100 µL of Filter Paper Elution Buff er into each well. Included in the Blank cavity.
- 5- Pipette 5 µL of Controls into previously determined wells.
- 6- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 7- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 8- Remove the sealing of wells.
- 9- Discard the contents of the cavities by aspiration (Washing machine) or by decanting (Manual). Remove the paper discs, if necessary, with the aid of a needle or tip. Use approximately 300 µL of Washing Solution, previously prepared to perform a total of fi ve (5) washing cycles. To ensure drying of the plate, at the end of the wash, hit the board for a few seconds on absorbent paper.
- Note:** Poor washing/drying can cause inadequate results.
- 10- Pipette 100 µL of Conjugate into all wells.
- 11- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 12- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 13- Remove the sealer from plate cavities.
- 14- Repeat item 9.
- 15- Pipette 100 µL of Substrate into each well.
- 16- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.
- 17- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 18- Remove the sealer from plate cavities.
- 19- Pipette 100 µL of Stop Solution into all wells.
- 20- Mix gently for ± 30 seconds.
- 21- Read using double filter: 450nm / 630nm within up to 15 minutes (maximum).

#### TECHNIQUE VERIFICATION

Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

| ITEM                        | ABSORBANCE |
|-----------------------------|------------|
| Blank                       | < 0.200    |
| Negative Control            | < 0.200    |
| Positive Control            | > 1.000    |
| Negative Extraction Control | < 0.200    |
| Positive Extraction Control | > 0.300    |

If the values are out the expected values, you must repeat the technique.

#### CALCULATIONS

##### QUALITATIVE

##### Samples of Serum or Plasma

Calculate Cut Off according to the following formula:

Cut Off = (Average Absorbance Positive Control x 0,05) + 0,140

Example:

| ITEM                                                           | ABSORBANCE                                 |
|----------------------------------------------------------------|--------------------------------------------|
| Positive Control<br>(Reagent N° 8)                             | A1 = 2.512                                 |
|                                                                | A2 = 2.454                                 |
| Cut Off = (Average Absorbance Positive Control x 0.05) + 0.140 | ((2.512 + 2.454/2) x 0.05) + 0.140 = 0.264 |

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

| ITEM                           | ABSORBANCE            |
|--------------------------------|-----------------------|
| Sample                         | 1.321                 |
| Cut Off Value                  | 0.264                 |
| Index = Sample / Cut Off Value | 1.321 / 0.264 = 5.003 |

**Note:** The data presented in the examples are for illustration only and cannot be used to calculate results.

#### Whole Blood Sample on Filter Paper

Calculate Cut Off according to the following formula:  
Cut Off = (Average Absorbance of Positive Control x 0.01) + 0.09

Example:

| ITEM                                                             | ABSORBANCE                              |
|------------------------------------------------------------------|-----------------------------------------|
| Positive Control<br>(Reagent N° 8)                               | A1 = 2.856                              |
|                                                                  | A1 = 2.858                              |
| Cut Off = (Average Absorbance of Positive Control x 0.01) + 0.09 | (2.856 + 2.858)/2 x 0.01 + 0.09 = 0,118 |

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

| ITEM                           | ABSORBANCE          |
|--------------------------------|---------------------|
| Sample                         | 0.936               |
| Cut Off Value                  | 0.118               |
| Index = Sample / Cut Off Value | 0.936/ 0.118= 7.932 |

#### INTERPRETATION OF RESULTS

##### Samples of Serum, Plasma and Filter Paper Whole Blood

After calculating the index of the samples, consider the indices below to determine the results.

| RESULTS      | QUALITATIVE FOR SERUM SAMPLE, PLASMA AND TOTAL BLOOD ON FILTER PAPER |
|--------------|----------------------------------------------------------------------|
|              | INDEX                                                                |
| Negative     | < 0.90                                                               |
| Positive     | > 1.1                                                                |
| Undetermined | 0.9 – 1.1                                                            |

**Note:** In the case of an undetermined result, the sample must be reanalyzed. Samples that obtain results repeatedly Indeterminate conditions must be retested using an alternative method. If results remain indeterminate, a new sample in two weeks. The result of the last sample must prevail collected. The results provided by this kit must be interpreted responsible medical professional, and it is not the only criterion for the determining the diagnosis and/or treatment of the patient.

#### PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be established with basis of a single trial. Other confirmatory tests should be included, before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before the definitive diagnosis of the disease.

#### INTERFERENT

No interference was observed for Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rematoid Factor 980 IU/mL, C Reactive Protein 41.2 mg/dL and Anti Streptolysin O 1023 IU/ml.

#### CROSS REACTIVITY

##### Samples of Serum or Plasma

A cross-reactivity study was carried out, evaluating 130 samples negative for SARS-CoV-2, but positive for other infections. Among them 10 positivesamples for Influenza, 27 positive for Rhinovirus, 13 positive samples for Respiratory Syncytial Virus, 10 positive samples for Cytomegalovirus, 10 positive samples for HIV, 10 positive samples for Zika, 10 positive samples for HBV, 10 positive samples for HCV, 10 positive samples for Dengue, 10 positive samples for Chikungunya and 10 samples for Toxoplasmosis. It was not observed cross-reactivity with positive samples for Influenza, Rhinovirus, Respiratory syncytial virus, Cytomegalovirus, HIV, Zika, HCV, Chikungunya and Toxoplasmosis. However, 1 cross-reactivity was observed for HBV and 1 for Dengue, with 90% specificities for HBV and Dengue. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data together with other laboratory data.

##### Whole Blood Sample on Filter Paper

A cross-reactivity study was carried out, evaluating 82 samples negative for SARS-CoV-2, but positive for other infections. Among them 9 positive samples for Zika, 11 positive samples for Toxoplasmosis, 13 positive samples for Rubella, 8 positive samples for Chagas, 5 positive samples for HIV, 10 positive samples for HCV, 10 positive samples for Syphilis, 7 positive samples for Cytomegalovirus and 9 HBV positive samples. It was not observed cross-reactivity with positive samples forZika, Toxoplasmosis, Rubella, Chagas, HIV, HCV, Syphilis, Cytomegalovirus and HBV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data together with other laboratory data.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### PRUDUCT PERFORMANCE

##### Accuracy REPETIBILITY

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

| REPETIBILITY                 | SAMPLE |       |        |
|------------------------------|--------|-------|--------|
|                              | 1      | 2     | 3      |
| Average                      | 0.312  | 0.965 | 0.138  |
| Standard Deviation           | 0.044  | 0.078 | 0.019  |
| Coefficient of Variation (%) | 14.107 | 8.056 | 13.591 |

##### REPRODUCIBILITY

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

| REPRODUCIBILITY              | SAMPLE |       |        |
|------------------------------|--------|-------|--------|
|                              | 1      | 2     | 3      |
| Average                      | 0.302  | 10.37 | 0.140  |
| Standard Deviation           | 0.038  | 0.084 | 0.018  |
| Coefficient of Variation (%) | 12.879 | 8.112 | 12.653 |

##### Clinical Sensitivity and Specificity

###### Serum and Plasma Sample

The BIOLISA COVID-19 IgA kit was used to analyze clinical samples. The results show that the clinical sensitivity of the BIOLISA COVID-19 IgA kit is 96,97% and clinical specificity is 95,83%.

|                     | Expected Result | BIOLISA COVID-19 IgA |
|---------------------|-----------------|----------------------|
| Positive Sample     | 66              | 64                   |
| Negative Sample     | 144             | 138                  |
| Total Tested Sample | 210             | 202                  |

Clinical Sensitivity: 96,97% (64/66) – IC 95%: 92,90% - 100%  
Clinical Specificity: 95,83% (138/144) – IC 95%: 91,71% - 98,57%

#### Whole Blood Samples on Filter Paper

|                     | Expected Result | BIOLISA COVID-19 IgG |
|---------------------|-----------------|----------------------|
| Positive Sample     | 50              | 48                   |
| Negative Sample     | 154             | 141                  |
| Total Tested Sample | 204             | 189                  |

Clinical Sensitivity: 96,00% (48/50) IC 95%: 86,30% - 99,50%  
Clinical Specificity: 91,56 % (141/154) IC 95%: 86,00% - 95,40%

#### DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family, and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are: fever, cough, dyspnoea, myalgia or fatigue. It can cause diseases ranging from the cold common to more serious diseases, such as the Eastern Respiratory Syndrome (MERS-CoV) and Severe Acute Respiratory Syndrome (SARSCoV). In December 2019, a series of cases of pneumonia caused disease appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples was found to be Coronavirus infection, which it was called the new coronavirus (2019-nCoV). In approximately 150 countries, on all continents, cases of infection have already been reported.

#### BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR.
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
4. Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
5. Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.
6. TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
7. QUIBASA: Departamento de Pesquisa e Desenvolvimento.

#### QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

##### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454  
E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454  
E-mail: sac@bioclin.com.br  
ANVISA registration for BIOLISA COVID-19 IgA kit: 10269360338

**Review:** June/2021

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY  
(last day of month)

TEMPERATURE LIMITATION  
(store at)

CONTAINS SUFFICIENT  
FOR <N> TESTS

CONSULT INSTRUCTIONS  
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED  
REPRESENTATIVE

KEEP AWAY  
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF  
PACKAGE IS  
DAMAGED