BIOLISA COVID-19 ANTICORPO NEUTRALIZANTE

^[REF] **K243**

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FINALIDADE

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1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
4. Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019- nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
5. Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.
6. Huo, J., Le Bas, A., Ruza, R.R. et al. Neutralizing nanobodies bind SARS-CoV-2 spike RBD and block interaction with ACE2. Nat Struct Mol Biol 27, 846–854 (2020).
7. Quibasa: Datos do Departamento de Pesquisa e Desenvolvimento.

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BIOLISA COVID-19 ANTIBODY NEUTRALIZER

REF K243

INSTRUCTIONS FOR USE

FUNCTION

Test for qualitative determination of the neutralization capacity of total antibodies anti-SARS-CoV-2 in biological samples of serum or plasma (EDTA or Heparin) by inhibiting the formation of the S antigen and ACE2 protein complex (S-ACE2) by competition test. For use only *in vitro* diagnosis.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzyme
The BIOLISA COVID-19 ANTIBODY NEUTRALIZER kit is a solid phase immunoassay based on the principle of qualitative detection by competition of total antibodies anti-SARS-CoV-2 capable of inhibiting the S-ACE2 complex in serum and human plasma samples.
In the initial incubation, Conjugate 1, composed of the protein ACE2 conjugated to Biotin, binds to coated Streptavidin on the microplate. After the initial incubation, the microplate is washed to remove unbound materials. In the second incubation, serum and plasma samples containing total antibodies anti-SARS-CoV-2 compete with the ACE2 protein for binding with the Spike Antigen conjugated to Peroxidase present in Conjugate 2, preventing the formation of the S-ACE2 complex. A New wash is performed to remove the total antibodies anti-SARS-CoV-2, linked or not linked to Spike Antigen. After this step, the substrate is added and incubated. An blue color for samples indicates the absence of Total Antibodies anti-SARS-CoV-2 and complex formation between the Antigen Spike conjugated to Peroxidase and ACE2. On the other hand, a lower intensity of the blue color indicates a greater amount of Total anti-SARS-CoV-2 antibodies present in the sample prevented the formation of the S-ACE2 complex. The Stop Solution is added to stop the reaction if there is a color change from blue to yellow, measured on a reader microplate.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8°C. Contains: Plate sensitized with Streptavidin and preservative.
- 2- Conjugate 1** - Store between 2 and 8°C. Contains: Buffer Solution containing Biotin-bound recombinant ACE2 protein, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** - Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate <0.5 mol/L, Potassium Chloride <100 mmol/L, Sodium Chloride <5 mol/L), surfactant and preservative.
- 4- Sample Diluent** - Store at 2 to 8°C. Contains: Buffer Solution, stabilizers, surfactant, dye and preservative.
- 5- Conjugate 2** - Store between 2 and 8°C. Contains: Buffer Solution containing Spike Antigen linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- 6- Substrate** - Store between 2 and 8°C. Contains: Solution Buffer containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 7- Stop Solution** - Store at 2 to 8°C. Contains: Acid Hydrochloride 1 M.
- 8- Negative Control** - Store at 2 to 8°C. Contains: Buffer, stabilizer, surfactant and preservative solution. **Potentially infective.**
- 9- Positive Control** - Store between 2 and 8°C. Contains: Buffer Solution, Total Anti-SARS-CoV-2 Antibodies, Dye, stabilizers, surfactant and preservative. **Potentially infecting.**
- 10- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1 - Sensitized Plate	1 Unit	2 Units	5 Units
2- Conjugate 1	1 Vial x 6 mL	2 Vials x 6 mL	5 Vials x 6 mL
3- Concentrated Washing	1 Vial x 50 mL	2 Vials x 50 mL	5 Vials x 50 mL
4- Sample Diluent	1 Vial x 15 mL	2 Vials x 15 mL	5 Vials x 15 mL
5- Conjugate 2	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
6- Substrate	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
7- Stop Solution	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
8- Negative Control	1 Vial x 300 µL	2 Vials x 300 µL	5 Vials x 300 µL
9- Positive Control	1 Vial x 300 µL	2 Vials x 300 µL	5 Vials x 300 µL
10- Plate Sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table
- Instructions for use (manual)

Required materials not contained in the kit:

- 1- Pipette capable of dispensing volumes of 10 to 100 µL with lower coefficiente of variation than 1.5%.
- 2- Re-pipettor for repetitive pipetting volumes of 10 to 100 µL, with lower coefficient of variation than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader capable of absorbance at 450 and 630 nm wavelength.
- 5- Paper towel to dry cavities
- 6- Stopwatch or watch.
- 7- Flask to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Tools of Quality Control.
- 10- Incubator 37°C ± 2°C.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature and transport should be 2 to 8°C. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The sachet containing the microplate should be opened only after it reaches room temperature. Place the strip with unused cavities in the sachet, seal and store be 2 to 8°C.
- 4- The water used in material cleaning must to be recent and free of contaminants.
- 5- Deionized and saturated columns release alkaline water, several ions and oxidizing and reducing agents that can signifi cantly alter the results.
- 6- Stop Solution contains Hydrochloric Acid which is a strong acid. Handle it with care.
- 7- All the raw material of product is tested and should be nonreactive for HBsAg, Anti-HIV 1&2 and Anti-HCV. However, these tests do not provide total assurance of the absence of infectious agents. The manipulation of any product containing human serum is potentially capable of transmitting diseases. Therefore, we must take due care in handling the bio safety of these products.

- 8- Always add reagents in the same order to minimize the diff erence in reaction time between the cavities.
- 9- As a safety measure, you should cover the plate during the reaction.
- 10- You must ensure that the bottom of the cavity is clean and dry, and there are no bubbles on the surface fl uid before reading the plate. Do not let the cavities run dry during the test.
- 11- Do not expose reagents, especially the Substrate, to strong light or Hypochlorite fumes during storage or incubation steps.
- 12- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Service) of Quibasa.
- 14- Do not use the product in case of damaged packaging.
- 15- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolysed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C.

PROCESS DESCRIPTION

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of the Reagent N° 3 (Concentrate Washing) in 1000 mL of distilled or deionized water. After preparation the solution may be stored at 2 to 30°C until expiration date printed on the original bottle. In case of crystallization, heat it at 37°C until dissolution.

Substrate

The Substrate is ready for use.

STABILITY AFTER OPEN

The results of the stability test show that the KIT BIOLISA COVID-19 ANTIBODY NEUTRALIZER is stable after opening for up to 30 days. This stability may vary according to the conditions of the test and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and technique validation criteria.

TECHNIQUE

For use in automatic equipment, consult the SAC (Customer Service).

- Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 - 30°C) for at least 40 minutes.
- 1- Separate the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate).
- Return unused strips from the microplate to the packaging sealed original.
- 2- Separate the first cavity for Blank (OPTIONAL).
- 3- Pipette 50 µL of Conjugate 1 in all wells, including in the Blank cavity.
- 4- Homogenize gently for ± 30 seconds, cover the cavities with plate sealer.
- 5- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 6- Discard the contents of the cavities by aspiration (Washer) or by decanting (manual). Use approximately 500 mL of Wash Solution, **previously prepared**, for a total of five (5) washing with a 30-second sauce per cycle.
- Note:** Evaluate the suction and dispensing capacity of the equipment. To guarantee the drying of the plate, at the end of the washing, tap the plate for a few seconds on absorbent paper as washing/Poor drying can cause inadequate results.

- 7- Pipette 50 µL of Sample Diluent into all wells.
- 8- Pipette 10 µL of Sample or Controls into the wells previously determined, then pipette 100 µL of the Conjugate 2 in all wells.
- 9- Homogenize gently for ± 30 seconds. Cover the cavities with the plate sealer.
- 10- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 11- Remove the plate sealer from the cavities.
- 12- Repeat item 6.
- 13- Pipette 100 µL of Substrate into all wells.
- 14- Gently homogenize for ± 30 seconds. Cover the cavities with the plate sealer.
- 15- Incubate for 15 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
- 16- Remove the plate sealer from the cavities.
- 17- Pipette 100 µL of Stop Solution into all wells.
- 18- Gently mix for ± 30 seconds.
- 19- Read using double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Verify if the results obtained by the reading of the Blank and Controls are compatible to the values below:

ITEM	ABSORBANCE
Blank	> 1.000
Negative Control	> 1.200
Positive Control	< 1,000

If the values are out the expected values, you must repeat the technique.

DESCRIPTION OF CALCULATIONS

QUALITATIVE

Reaction Inhibition Percentage Calculation:
% inhibition of reaction=[1-(Sample/NegativeControl)] x100
Example:

ITEM	ABSORBANCES
Negative Control	Negative Control = 2,181
Sample	Sample = 1,013
%=[1-(Sample/NegativeControl)] x100	[1-(1,013/2,181)]x100=53,55%

INTERPRETATION OF RESULTS

After calculating the index of the samples, consider the indices below to determine the results.

RESULTS	QUALITATIVE
Negative	> 30%
Positive	> 35%
Undetermined	Entre 30 or 35%

Note: In the case of an undetermined result, the sample must be re-analyzed. Samples that obtain results repeatedly indeterminate must be retested using an alternative method. If the results remain indeterminate, a new sample must be collected in two weeks. The result of the last sample collected must prevail.. The results provided by this kit should be interpreted by the responsible medical professional, and the sole criterion for determining diagnosis and/or treatment of the patient.

Note: The data presented in the examples are only for illustration and cannot be used to calculate the results.

PROCEDURE LIMITATIONS

The kit is a qualitative assay with the possibility of correlation of the inhibition capacity of the reaction by the analysis in percentage.

There are still no clinical data that establish a relationship between the ability to inhibit the reaction presented in a sample and the clinical status of the patient. The results obtained with the BIOLISA COVID-19 ANTIBODY NEUTRALIZING kit demonstrate only that patients with results < 30% do not have neutralizing antibodies, and samples > 35% do have these antibodies.

Results <30% do not exclude SARS-CoV-2 infection.

Results > 35% may be due to infection by other coronaviruses, although no cross-reactions were evidenced in studies performed during development.

This product simulates the immunodominant inhibition mechanism based on viral neutralization by neutralizing antibodies, blocking the formation of the S-ACE2 immune complex. Patients with positive virus may develop other neutralizing or even non-neutralizing antibodies against SARS-CoV-2, both undetectable by the proposed methodology. All results should be interpreted in conjunction with other clinical information available prior to definitive diagnosis of the disease.

INTERFERENT

No interference was observed for Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rematoid Factor 980 IU/mL, C Reactive Protein 41.2 mg/dL and Anti Streptolysin O 1023 IU/ml.

CROSS REACTIVITY

A study was carried out with 140 negative samples for COVID-19, but positive for other infections. Amongst them, 10 positive samples for Influenza, 27 positive samples for Rhinovirus, 13 positive samples for Syncytial Virus Respiratory, 10 Zika positive samples, 10 samples positive for Dengue, 10 positive samples for Chikungunya, 10 positive samples for Toxoplasmosis, 10 positive samples for Cytomegalovirus, 10 positive samples for Rubella, 10 HIV positive samples, 10 positive samples for HBV and 10 HCV positive samples. It was not observed cross-reactivity with positive samples for Zika, Dengue, Chikungunya, Toxoplasmosis, Cytomegalovirus, Rubella, HIV, HBV and HCV. Despite the results found, it is not possible to completely rule out the possibility of cross-reactivity.

The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRUDUCT PERFORMANCE**ACCURACY****Repeatability**

The repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Samples		
	1	2	3
Average	0.262	0.084	1.306
Standard Deviation	0.024	0.011	0.083
Coefficient of Variation (%)	9.141	12.586	6.357

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Samples		
	1	2	3
Average	0.262	0.079	1.206
Standard Deviation	0.025	0.007	0.078
Coefficient of Variation (%)	10.125	8.936	6.446

CLINICAL SENSITIVITY AND SPECIFICITY

To carry out the sensitivity and specificity study 156 clinically known samples were analyzed, 96 positive and 60 negative. The results can be observed in the table below:

	Expected Result	BIOLISA COVID-19 ANTIBODY NEUTRALIZER
Positive Sample	96	93
Negative Sample	60	57
Total Samples Tested	156	

Clinical Sensitivity: 96.87% (93/96) - 95% CI: 88.9% - 100%

Clinical Specificity: 95% (57/60) - 95% CI: 89.62% - 100%

A new clinical sensitivity study was carried out with 70 positive samples for the PRNT gold standard, which are considered positive by the method when they present titer $\geq 1/20$.

A new clinical specificity study was performed with 100 pre-pandemic samples.

Sensitivity (68/70) = 97.1% 95% CI: 93.3 - 100%

Specificity (92/100) = 92.0% 95%CI: 86.9 - 97.1%

The study of detection of neutralizing antibodies in different groups of samples was also carried out. Below are the new data:

Study by sample group:

Grup*	Total evaluated samples	Positive Neutralizing Ac
Group of infected individuals with symptoms > 14 days	37	28
Group of uninfected individuals vaccinated with CORONAVAC (> 30 days after 2nd dose)	73	64

*The data presented above do not represent the sensitivity of the product, but the expression of neutralizing antibodies in different situations per group evaluated in the tested panel.

DIAGNOSTIC SIGNIFICANCE

Coronaviruses (CoV) belong to the Coronaviridae family, and are widely distributed infecting humans and other mammals. In humans, the most common symptoms presented are: fever, cough, dyspnoea, myalgia or fatigue. It can cause diseases ranging from the cold common to more serious diseases, such as the Eastern Respiratory Syndrome (MERS-CoV) and Severe Acute Respiratory Syndrome SARS-CoV. In December 2019, a series of cases of pneumonia caused disease appeared in Wuhan in China, with clinical symptoms very similar to viral pneumonia, after complete sequencing of the respiratory samples was found to be Coronavirus infection, which it was called the new coronavirus (2019-nCoV). In approximately 150 countries on all continents have already been reported cases of infection.

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
4. Corman VM, Landt O, et al. Detection of 2019 novel coronavirus (2019- nCoV) by real-time RT-PCR. Euro Surveill. 2020 Jan;25(3).
5. Huang C, Wang Y, et al. Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. Lancet. 2020 Feb 15;395(10223):497- 506.
6. Huo, J., Le Bas, A., Ruza, R.R. et al. Neutralizing nanobodies bind SARS-CoV-2 spike RBD and block interaction with ACE2. Nat Struct Mol Biol 27, 846–854 (2020).
7. Quibasa: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte/MG - Brasil
Tel.: +55 31 3439 5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Indústria Brasileira

CUSTOMER SERVICE

Customer Advisory Service

Phone: +55 31 3439 5454 | E-mail: customerservice@bioclin.com.br

ANVISA registration for BIOLISA COVID-19 ANTIBODY NEUTRALIZER kit: 10269360337

Review: June/2023

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER