BIOLISA DENGUE NS1

[REF] **K253**

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FINALIDADE

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- REFERÊNCIA BIBLIOGRÁFICA**
1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
 2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
 3. Halstead SB, Selective primary health care: strategies for control of disease in the develop-ing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
 4. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
 - 5.Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detect-ed by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
 6. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
 7. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Ge-neva: World Health Organization.
 8. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neu-tralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
 9. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephali-tis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
 10. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Micro-biol Immunol. 1994; 38(12): 951-5.
 11. Dussart P, et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
 12. Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
 13. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- Halstead SB, Selective primary health care: strategies for control of disease in the develop-ing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detect-ed by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Ge-neva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neu-tralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephali-tis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Micro-biol Immunol. 1994; 38(12): 951-5.
- Dussart P., et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
- Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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FUNCTION

Test for qualitative determination of Dengue virus NS1 antigen in biological samples (serum or plasma) through enzyme-immunoassay. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or Immunoenzyme

The BIOLISA DENGUE NS1 kit is a solid-phase enzyme immunoassay based on the “sandwich” principle for the qualitative detection of Dengue Virus NS1 antigen in human serum or plasma samples. The microplate is coated with anti-NS1 monoclonal antibodies specific for Dengue subtypes I, II, III and IV. At first, the sample is incubated with the Conjugate, which consists of anti-NS1 antibodies conjugated to peroxidase. Dengue Virus NS1 antigens, present in the sample, bind to specific anti-Dengue antibodies immobilized on the plate and to specific anti-Dengue antibodies present in the conjugate, forming immune complexes. After the initial incubation, the microplate is washed to remove unbound materials. Then the Substrate is added and incubated, producing a blue color that indicates the detection of Dengue NS1 antigens in the specimen. Stopping Solution is added to stop the reaction and there is a color change from blue to yellow as measured on a microplate reader.

REAGENTS

- 1 - Sensitized Plate** - Store between 2 and 8°C. Contains: Plate sensitized with monoclonal antibodies anti-NS1 from Dengue.
- 2 - Conjugate** - Store between 2 and 8°C. Contains: Buffer solution containing monoclonal antibodies to Dengue's NS1 linked to peroxidase, surfactant, stabilizers, dye and preservative.
- 3 - Concentrated Wash** - Store between 2 and 8 °C. Contains: Contains: Buffer solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4 - Substrate** - Store between 2 and 8 °C. Contains: Buffer solution containing urea peroxide, tetramethylbenzidine (TMB) and preservative.
- 5 - Stop Solution** – Store between 2 and 8°C. Contains: 1 M Hydrochloric Acid Solution.
- 6 - Negative Control** – Store between 2 and 8°C. Contains: Buffer solution, stabilizer, surfactant and preservative. **Potentially infectious.**
- 7 - Positive Control** – Store between 2 and 8°C. Contains: Buffer solution, NS1 antigen for Dengue, dye, stabilizers, surfactant and preservative. **Potentially infectious.**
- 8 - Plate Sealers.**

PRESENTATION

REAGENTES	1	2	3
1- Sensitized Plate	1 Unit (96 cavities)	2 Units (96 cavities)	5 Units (96 cavities)
2 - Conjugate	1 Vial x 6 mL	2 Vials x 6 mL	5 Vials x 6 mL
3 - Concentrated Washing Solution	1 Vial x 50 mL	2 Vials x 50 mL	5 Vials x 50 mL
4 - Substrate	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
5 - Stop Solution	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
6 - Negative Control	1 Vial x 1 mL	2 Vials x 1 mL	5 Vials x 1 mL
7 - Positive Control	1 Vial x 1 mL	2 Vials x 1 mL	5 Vials x 1 mL
8 - Plate Sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INPUTS

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Necessary materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 5 to 300 µL with a coefficient of variation lower than 1.5%.
- 2- Repeater for repetitive pipetting of 300 µL volumes with a coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or clock.
- 7- Bottle to store the Wash Solution after diluted.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator at 37 °C ± 2 °C.

TRANSPORTATION AND STORAGE CONDITIONS

Storage temperature should be 2 to 8 °C. Transport, at room temperature (up to 30 °C), should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate must be opened only after reaching room temperature. Replace strips from unused microwells in the sachet, seal and store at 2 and 8 °C.
- 4- The water used to clean the material must be fresh and free from contaminants.
- 5- Saturated deionizer columns release alkaline water, different ions, and oxidizing and reducing agents that can significantly alter the results.

- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. So handle it with due care.
- 7- All raw material of the product is tested and must be non-reactive for HBsAg and Anti-HIV 1&2. However, these tests do not offer complete assurance of the absence of infectious agents. The handling of any product that contains serum is potentially capable of transmitting disease.
- Therefore, it is necessary to take proper biosafety precautions when handling these products.
- 8- Pipette the reagents always in the same order to minimize the difference in reaction time between the microwells.
- 9- As a protection measure, the plate must be covered during the reaction.
- 10- Ensure that the bottom of the cavity is clean and dry, and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow wells to dry during the assay.
- 11- Do not expose reagents, especially Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.
- 12- We recommend applying local, state and federal regulations for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service of Customer Service) by Quibasa.
- 14- Do not use the product in case of damage to the packaging.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolyzed and highly lipemic specimens should not be used. The samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C. Plasma samples stored for periods longer than recommended (30 days) may present fibrin and fibronectin precipitation, which could interfere with the test. Samples collected with other anticoagulants, other than those mentioned, should not be used.

PROCESSS DESCRIPTION

Stability After Open

The results of the stability test prove that the BIOLISA DENGUE NS1 kit is stable after being opened for up to 30 days. This stability may vary depending on test conditions and environment. Therefore, it is suggested to monitor the product performance using the kit's internal controls and the technique's validation criteria.

Preparation of Work Reagents

WASHING SOLUTION

Dilute the contents of vial No3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for at least 30 days. It can be stored at room temperature. If crystallization occurs, heat to 37 °C until dissolved.

SUBSTRATE

The Substrate is ready for use.

TECHNIQUE

For use in automatic equipment, consult the (Customer Advisory Service) SAC

Before starting the assay, allow all reagents, samples and controls to stabilize at room temperature (15 – 30°C) for at least 40 minutes. Return unused microplate strips to the original sealed packaging.

- 1- Separate the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return unused microplate strips to the original sealed packaging.
- 2- Separate the first cavity for White (OPTIONAL).
- 3- Pipette 50 µL of the Samples, Positive Control and Negative Control in the previously determined wells.
- 4- Pipette 50 µL of the Conjugate into all wells, except for the Blank well.
- 5- Homogenize gently for ± 10 seconds. Cover wells with plate sealer.
- 6- Incubate for 90 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 7- Remove the sealant from the cavities.
- 8- Discard the contents of the cavities by aspiration (Washer) or by decantation (Manual). Use approximately 300 µL of Wash Solution, previously prepared, and carry out a total of five (5) wash cycles, with shaking (shake) for 5 seconds. To guarantee the plate drying, at the end of washing, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying may cause inappropriate results.
- 9- Pipette 100 µL of Substrate into all wells, including the Blank well.
- 10- Mix gently for ± 30 seconds. Cover the wells with the plate sealer.
- 11- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 12- Remove the plate sealer from the wells.
- 13- Pipette 100 µL of Stop Solution into all wells, including the Blank well.
- 14- Mix gently for ± 30 seconds.
- 15- Read using a double filter: 450 nm / 630 nm within 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0.100
Negative Control	< 0.100
Positive Control	> 1.000

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS

Qualitative

Calculate Cut Off according to the following formula:
Cut Off = (Average Positive Control Absorbance x 0.02)
+ 0.180.

Example:

ITEM	ABSORBANCE
Positive Control	A1 = 2.101
	A2 = 2.104
Cut Off: (Average Positive Control Absorbance x 0.02 + 0.180	((2.101 + 2.104) / 2 x 0.02) + 0.180=0.222

Calculate the Index by dividing the Sample's absorbance by the Cut Off value.

Example:

ITEM	ABSORBÂNCIA
Sample	2.675
Cut Off Value	0.222
Index = Samples / Cut Off Value	2.675 / 0.222 = 12.05

Note: The data presented in the examples are for illustration only and cannot be used to calculate the results.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indices below to determine the results.

RESULTS	QUALITATIVE FOR SERUM AND PLASMA SAMPLES
	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	Entre 0.9 - 1.1

Note: In case of indeterminate result, the sample must be reanalyzed. Samples that repeatedly obtain indeterminate results should be retested using an alternative method. If results remain indeterminate, a new sample should be collected within 9 days of symptom onset. The result of the last sample collected must prevail. The results provided by this kit must be interpreted by the responsible medical professional and are not the sole criterion for determining the patient's diagnosis and/or treatment.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before the descriptive diagnosis of the disease.

INTERFERENTS

No interference was observed by Triglycerides 1500 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Acid Oxalic 60 mg/dL, Rheumatoid Factor 980 UI/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 UI/mL. Plasma samples stored for periods longer than recommended (30 days) may present fibrin and fibronectin precipitation, which could interfere with the test. These samples should not be used, as well as samples collected with other anticoagulants than those mentioned (EDTA and heparin).

CROSS REACTIVITY

A cross-reactivity study was performed, evaluating 147 serum and plasma samples negative for Dengue NS1 but positive for other infections. Among them 10 HIV positive samples, 10 HTLV positive samples, 10 HBsAg positive samples, 10 anti-HBs positive samples, 10 Syphilis positive samples, 10 Yellow Fever positive samples, 10 Zika positive samples, 5 positive samples for Chikungunya, 10 samples positive for Chagas Disease, 10 samples positive for CMV, 10 samples positive for Rubella, 10 samples positive for Toxoplasmosis, 7 samples positive for H1N1, 10 samples positive for Sars-CoV-2 and 15 samples positive for HCV . No cross-reactivity was observed with positive samples for HIV, HTLV, HBsAg, anti-HBs, Syphilis, Yellow Fever, Zika, Chikungunya, Chagas Disease, CMV, Rubella, Toxoplasmosis, H1N1, Sars-Cov-2 and HCV. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

Accuracy

REPEATABILITY

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

REPEATABILITY	MUESTRA		
	1	2	3
Average	1.244	3.334	0.075
Standard Deviation	0.028	0.135	0.005
Coefficient of Variation (%)	2.29	4.04	6.45

REPRODUTIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

REPRODUTIBILITY	MUESTRA		
	1	2	3
Average	1.287	3.301	0.075
Standard Deviation	0.066	0.077	0.007
Coefficient of Variation (%)	5.12	2.33	9.53

Clinical Sensitivity and Specificity

The BIOLISA DENGUE NS1 kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA DENGUE NS1 kit is 96% and the clinical specificity is 96.3%.

CLINICAL SENSITIVITY AND SPECIFICITY	EXPECTED RESULT	BIOLISA DENGUE NS1
Positive Sample	50	48
Negative Sample	54	52
Total Tested Samples	104	100

Clinical Sensitivity: 96% (48/50) - IC 95%: 83.3% - 99.5%

Clinical Specificity: 96.3% (52/54) - IC 95%: 87.3% - 99.5%

CLINICAL SIGNIFICANCE

Dengue virus is a Flavivirus transmitted by the Aedes aegypti mosquito. It is distributed throughout the tropical and subtropical areas of the world, causing more than 100 million infections annually. The classic dengue infection is characterized by high fever (usually between 38 and 40°C) with an abrupt onset, malaise, anorexia (poor appetite), headache, muscle pain and eye pain. Cases of Hemorrhagic Dengue can cause gingivorrhages and epistaxis, internal hemorrhages and disseminated intravascular coagulation, with damage to several organs, which can cause death. NS1 is a highly conserved glycoprotein that is present in high concentrations in the serum of dengue-infected patients. The NS1 antigen is found from 1st to 9th days after the onset of fever in samples from patients with primary and secondary infection with the dengue virus. However, its concentration decreases over the days of infection, considerably reducing the chances of detecting the antigen. Articles report a sensitivity of less than 50% for detection after the 4th day of infection by various diagnostic methods. Primary dengue infection causes a rise in IgM antibodies 3 to 5 days after the onset of fever.

IgM antibodies normally remain in circulation for 30 to 90 days. Patients from endemic regions may present with secondary infections, which result in elevated levels of IgG antibodies, either alone or simultaneously with an IgM response. In case of secondary infection, the NS1 antigen may not be detected due to the possibility of lower viremia.

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. Halstead SB, Selective primary health care: strategies for control of disease in the develop-ing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
4. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- 5.Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detect-ed by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
6. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
7. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Ge-neva: World Health Organization.
8. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neu-tralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
9. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephali-tis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
10. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus infections. Micro-biol Immunol. 1994; 38(12): 951-5.
11. Dussart P., et al. Evaluation of an enzyme immunoassay for detection of dengue virus NS1 antigen in human serum. Clinical and vaccine immunology. 2006; p. 1185–1189.
12. Dussart P., et al. Evaluation of two new commercial tests for the diagnosis of acute dengue virus infection using NS1 antigen detection in human serum. Plos – Neglected Tropical Diseases. 2008; Volume 2; Issue 8; e280.
13. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439 5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454
E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA DENGUE NS1 kit:10269360348

Revision: January/2023

REF

CATALOG NUMBER

LOT

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER