Bioclin

BIOLISA HBsAg DBS

REF **K252**

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FINALIDADE

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1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
3. ALLEN JUNIOR, L. V.; ANSEL, H. C. Ansel's Pharmaceutical Dosage Forms and Drug Delivery Systems. 10. ed. Baltimore: Wolters Kluwer, 2014. 809 p.
4. CLINICAL AND LABORATORY STANDARD INSTITUTE. EP-25A: Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. Wayne: Clinical And Laboratory Standard Institute, 2009. 52 p.
5. DESHPANDE, S.S. Enzyme Imunoassays: from concept to development. New York: Chapman & Hall, 1996. 472 p.
6. TSITSILONIS, Ourania E.; THRASYVOULIDES, Apollon; BALAFAS, Apostolos; VOUTSAS, John F.; PAPAMICHAIL, Michael; LYMBERI, Peggy. Serological detection of hepatitis B viral infection by a panel of solid-phase enzyme-linked immunosorbent assays (ELISA). Journal Of Pharmaceutical And Biomedical Analysis, [S.L.], v. 34, n. 4, p. 811-822, mar. 2004. Elsevier BV. http://dx.doi.org/10.1016/s0731-7085(03)00563-6.
7. TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
8. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- REFERENCIAS BIBLIOGRÁFICAS**
- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
 - ALLEN JUNIOR, L. V.; ANSEL, H. C. Ansel's Pharmaceutical Dosage Forms and Drug Delivery Systems. 10. ed. Baltimore: Wolters Kluwer, 2014. 809 p.
 - CLINICAL AND LABORATORY STANDARD INSTITUTE. EP-25A: Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. Wayne: Clinical And Laboratory Standard Institute, 2009. 52 p.
 - DESHPANDE, S.S. Enzyme Immunoassays: from concept to development. New York: Chapman & Hall, 1996. 472 p.
 - TSITSILONIS, Ourania E.; THRASYVOULIDES, Apollon; BALAFAS, Apostolos; VOUTSAS, John F.; PAPAMICHAIL, Michael; LYMBERI, Peggy. Serological detection of hepatitis B viral infection by a panel of solid-phase enzyme-linked immunosorbent assays (ELISA). Journal Of Pharmaceutical And Biomedical Analysis, [S.L.], v. 34, n. 4, p. 811-822, mar. 2004. Elsevier BV. [http://dx.doi.org/10.1016/s0731-7085\(03\)00563-6](http://dx.doi.org/10.1016/s0731-7085(03)00563-6).
 - TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Immunoenzymatic test for qualitative detection of the presence of Hepatitis B virus surface Antigen in dried blood spot collected on filter paper (DBS). Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic. The BIOLISA HBsAg DBS kit is a solid-phase enzyme immunoassay based on the "sandwich" principle for the qualitative detection of HBsAg in human samples of dried blood spot collected on filter paper. The microplate is coated with monoclonal Antibodies specific for the Ad and Ay subtypes of HBsAg. At first, the dried blood spot impregnated on filter paper is incubated with the Elution Buffer. If the sample contains HBsAg, it will bind to the Antibodies coated on the microplate. After the initial incubation, the microplate is washed to remove unbound materials and then Peroxidase-conjugated Anti-HBs Antibodies are added. The microplate is again incubated and, in the presence of HBsAg in the sample, it will form immunocomplexes Antibody-HBsAg-Conjugate. A new wash is carried out to remove the excess.

After this step, the Substrate is added and incubated, producing a blue color that indicates the presence of Hepatitis B Antigen in the sample. The Stop Solution is added to stop the reaction, there is a color change from blue to yellow as measured on a microplates.

REAGENTS

1- Sensitized Plate - Store at 2 to 8°C. Contains: Microplate coated with Antibodies (mouse monoclonal) Anti-HBs (< 100 ng/well).

2- Conjugate - Store between 2 and 8°C. Contains: Anti-HBs Antibodies linked to Peroxidase, surfactant, stabilizers, colorant and preservative.

3- Concentrated Washing - Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.

4- Substrate - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.

5- Stop Solution – Store between 2 and 8°C. Contains: 1 M Hydrochloric Acid Solution.

6 - Plate Sealers.

7 - Elution Buffer - Store between 2 and 8 °C. Contains: Buffer Solution, surfactant, stabilizer and preservative.

8 - DBS Negative Control - Store at 2 to 8 °C. Contains: Non-reactive sample for HBsAg, HCV, HIV-1 and HIV-2 impregnated in filter paper. **Potentially infectious.**

9 - DBS Positive Control - Store at 2 to 8 °C. Contains: Sample reactive for HBsAg and negative for HCV, HIV-1 and HIV-2 impregnated in filter paper. **Potentially infectious.**

PRESENTATION

REAGENT	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Sensitized Plate	1 Unit (96 Cavities)	2 Units (192 Cavities)	5 Units (480 Cavities)
2 - Conjugate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
3 - Concentrated Washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
4 - Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
5 - Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6 - Plate Sealers	3 Units	1 Unit	1 Unit
7- Elution Buffer	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
8- DBS Negative Control	1 Unit	2 Units	3 Units
9- DBS Positive Control	1 Unit	2 Units	3 Units

EQUIPAMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous item.

Necessary materials not contained in the kit:

1- Pipettes capable of dispensing volumes of 5 to 300 µL with a variation coefficient of less than 1.5%.

2- Repipettor for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).

3- Microplate washer for filter paper technique.

4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.

5- Absorbent paper to dry the wells.

6- Stopwatch or watch.

7- Vial to store the Washing Solution after dilution.

8- Distilled or deionized water.

9- Quality Control Tools.

10- Incubator at 37 °C ± 2 °C.

11- Paper shredder (3 mm diameter).

TRANSPORTATION AND STORAGE CONDITIONS

Storage temperature should be between 2 and 8 °C. Transport at temperatures up to 30°C should not exceed 5 days. Keep out of the light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the proposed methodology to obtain exact results.

3- The sachet containing the microplate must be opened only after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.

4- The water used to clean the material must be fresh and free from contaminants.

5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.

6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.

7- All raw material of the product is tested and must be non-reactive for HBsAg and Anti-HIV 1 & 2. However, these tests do not offer complete security in the absence of infectious agents. The handling of any product that contains serum is potentially capable of transmitting disease. Therefore, it is necessary to take due care of biosafety in the handling of these products.

8- Pipette the reagents always in the same order to minimize the difference in reaction time between the wells.

9- As a protection measure, the plate must be covered during the reaction.

10- Make sure that the bottom of the cavity is clean and dry and that there are no bubbles on the liquid surface before reading the plate. Do not allow the wells to dry out during the test.

11- Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.

12- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Customer Service Department) of Quibasa.

14- Do not use the product in case of damage to the packaging.

15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Dried Blood collected on filter paper (DBS) - EDTA or Puncture
Samples dried on filter paper can be stored at room temperature as long as they are protected from direct sunlight and low humidity. For storage of up to 2 years, samples must be kept refrigerated at between 2 and 8 °C. Storage for more than 2 years must be carried out at a temperature of -20 °C.

PROCESS DESCRIPTION

Stability After Open

The results of the stability test show that the BIOLISA HBsAg DBS kit is stable after being opened for 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using internal controls of the kit and the criteria for validating the technique.

PREPARATION OF WORK REAGENTS

Washing Solution

Dilute the contents of vial Nº 3 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Washing solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30°C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolution.

Substrate

The Substrate is ready for use.

TECHNIQUE

For use in automatic equipment, consult the Customer Service Department (SAC)

Before starting the assay, allow all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

1- Separate the cavities to be used considering: DBS Positive Control, DBS Negative Control and Samples (it is recommended to test in duplicate). Return unused microplate strips to the original sealed packaging.

2- Separate the first cavity for the Blank (OPTIONAL).

3- Add a 3mm disk of DBS Positive Control, DBS Negative Control and Dried Blood Spot collected on filter paper, previously perforated, in the determined cavities.

4- Pipette 100 µL of the Elution Buffer into the wells previously determined for filter paper, including the well for the White.

5- Homogenize gently for ± 30 seconds, cover the wells with plate sealer.

6- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

7- Remove the sealant from the cavities.

8- Discard the contents of the cavities by aspiration (Washer for filter paper technique). Use approximately 300 µL of Wash Solution, previously prepared, and make a total of five (5) wash cycles with shaking (*shake*) of 5 seconds. To guarantee the drying of the plate, at the end of the washing, tap the plate for a few seconds on absorbent paper. Note: Poor washing/drying may cause inappropriate results.

9- Pipette 100 µL of Conjugate into all wells, including the Blank well.

10- Mix gently for ± 30 seconds. Cover the wells with the plate sealer.

11- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

12- Remove the plate sealer from the wells.

13- Repeat item 8.

14- Pipette 100 µL of Substrate into all wells.

15- Mix gently for ± 30 seconds. Cover the wells with the plate sealer.

16- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

17- Remove the plate sealer from the wells.

18- Pipette 100 µL of Stop Solution into all wells.

19- Mix gently for ± 30 seconds.

20- Read using a double filter: 450 nm / 630 nm within 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	< 0.100
DBS Positive Control	< 0.100
DBS Positive Control	> 0.500

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula:
Cut Off = Average DBS Negative Control Absorbance + 0.100
Example:

ITEM	ABSORBANCE
DBS Negative Control	A1= 0.020
	A2= 0.021
Cut Off = Average DBS Negative Control Absorbance + 0.100	(0.020 + 0.021) / 2) + 0.100 = 0.121

Calculate the Index by dividing the Sample's absorbance by the Cut Off value.
Example:

ITEM	ABSORBANCE
Sample	2.145
Cut Off Value	0.121
Índice = Samples / Cut Off Value	2.145 / 0.121 = 17.73

Note: The data presented in the examples are for illustration only and cannot be used to calculate the results.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indices below to determine the results.

RESULTS	QUALITATIVE FOR WHOLE BLOOD ON FILTER PAPER SAMPLES
	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	Between 0.9 - 1.1

Note: In the case of an undetermined result, the sample must be re-analyzed. Samples that obtain repeatedly indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected in two weeks. The result of the last collected sample must prevail. The results provided by this kit must be interpreted by the responsible medical professional, and it is not the only criterion for determining the diagnosis and / or treatment of the patient.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before the descriptive diagnosis of the disease.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 IU/mL, Protein C Reactive 41.2 mg/dL and Anti-Streptolysin O 1023 IU/mL.

CROSS REACTIVITY

A cross-reactivity study was performed, evaluating 106 HBsAg-negative but positive whole blood samples on filter paper for other infections. Among them 5 HIV positive samples, 9 HTLV positive samples, 9 Syphilis positive samples, 33 samples positive for HCV, 10 samples positive for Chagas Disease, 10 samples positive for Zika, 10 samples positive for CMV, 10 samples positive for Rubella and 10 positive samples for Toxoplasmosis. No cross-reactivity was observed with positive samples for HIV, HTLV, Syphilis, HCV, Chagas Disease, Zika, CMV, Rubella and Toxoplasmosis. Despite the results found, the possibility of cross-eactivity. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Sample		
	1	2	3
Average	1.312	3.897	0.115
Standard Deviation	0.087	0.101	0.013
Coefficient of Variation (%)	6.65	2.59	11.58

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	1.258	3.813	0.111
Standard Deviation	0.089	0.119	0.010
Coefficient of Variation (%)	7.11	3.13	9.42

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA HBsAg DBS kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA HBsAg DBS kit is >99.9% and the clinical specificity is >99.9%.

CLINICAL SENSITIVITY AND SPECIFICITY	Resultado Esperado	BIOLISA HBSAG DBS
Positive Sample	50	50
Negative Sample	50	50
Total Tested Samples	100	100

Clinical Sensitivity: >99.9% (50/50) - IC 95%: 92.9% - 100,0%
Clinical Specificity: >99.9% (50/50) - IC 95%: 92.9% - 100,0%

CLINICAL SIGNIFICANCE

Hepatitis B is an infectious disease caused by the HBV virus that affects the liver. The HBV virus is transmitted through sexual contact, blood exposure (including sharing sharps) and mother-to-child transmission via the placenta or during childbirth. HBV infection results in an apparent number of serological markers and one of the first of these markers is the Surface Antigen of Hepatitis B (HBsAg). The main HBsAg subtypes include ad and ay, all sharing the common determinant 'a'. HBsAg appears between 1 and 10 weeks after infection and before biochemical evidence of liver disease or jaundice. However, for most people, due to its high production, HBsAg can be detected after 30 days of infection. This period is known as the diagnostic window. In the bearer state chronic, HBsAg persists for more than 6 months. Antibodies against HBV antigens appear between 30 and 60 days after infection, varying between individuals. This is the period known as the immune window. HBV infection can cause self-limiting acute infection, hepatitis fulminant, chronic hepatitis with progression to cirrhosis and liver failure, and asymptomatic chronic carrier status. The HBV virus in infected people, persists for the rest of their lives and can be transmitted to other people. Thus, hepatitis B has become a problem of public health. Therefore, screening for HBsAg is recommended for all donors, pregnant women and people in high-risk groups.

BIBLIOGRAPHIC REFERENCES

1. Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002: 31.
3. ALLEN JUNIOR, L. V.; ANSEL, H. C. Ansel's Pharmaceutical Dosage Forms and Drug Delivery Systems. 10. ed. Baltimore: Wolters Kluwer, 2014. 809 p.
4. CLINICAL AND LABORATORY STANDARD INSTITUTE. EP-25A: Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. Wayne: Clinical And Laboratory Standard Institute, 2009. 52 p.
5. DESHPANDE, S.S. Enzyme Imunoassays: from concept to development. New York: Chapman & Hall, 1996. 472 p.
6. TSITSILONIS, Ourania E.; THRASYVOULIDES, Apollon; BALAFAS, Apostolos; VOITSAS, John F.; PAPAMICHAIL, Michael; LYMBERI, Peggy. Serological detection of hepatitis B viral infection by a panel of solid-phase enzyme-linked immunosorbent assays (ELISA). Journal Of Pharmaceutical And Biomedical Analysis, [S.L.], v. 34, n. 4, p. 811-822, mar. 2004. Elsevier BV. [http://dx.doi.org/10.1016/s0731-7085\(03\)00563-6](http://dx.doi.org/10.1016/s0731-7085(03)00563-6).
7. TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.
8. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of the BIOLISA HBsAg DBS kit at ANVISA: 10269360349

Review: May/2023

UNIVERSAL SYMBOLOGY

REF	CATALOG NUMBER	MADE BY
LOT	LOT NUMBER	CONTROL
MANUFACTURING DATE		POSITIVE CONTROL
VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
SEE INSTRUCTIONS FOR USE		CORROSIVE
IN VITRO DIAGNOSTIC PRODUCT		TOXIC
KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
DO NOT REUSE		PRODUCT STERILIZED
CAUTION		DANGER