# Bioclin

## BIOLISA HBsAg

REF **K120**

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#### FINALIDADE

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- ALLEN JUNIOR, L. V.; ANSEL, H. C. Ansel's Pharmaceutical Dosage Forms and Drug Delivery Systems. 10. ed. Baltimore: Wolters Kluwer, 2014. 809 p.
- CLINICAL AND LABORATORY STANDARD INSTITUTE. EP-25A: Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. Wayne: Clinical And Laboratory Standard Institute, 2009. 52 p.
- DESPANDE, S.S. Enzyme Immunoassays: from concept to development. New York: Chapman & Hall, 1996. 472 p.
- TSITSILONIS, Ourania E.; THRASYVOULIDES, Apollon; BALAFAS, Apostolos; VOUTSAS, John F.; PAPAMICHAIL, Michael; LYMBER I, Peggy. Serological detection of hepatitis B viral infection by a panel of solid-phase enzyme-linked immunosorbent assays (ELISA). Journal Of Pharmaceutical And Biomedical Analysis, [S.L.], v. 34, n. 4, p. 811-822, mar. 2004. Elsevier BV. [http://dx.doi.org/10.1016/s0731-7085\(03\)00563-6](http://dx.doi.org/10.1016/s0731-7085(03)00563-6).
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

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### INSTRUCTIONS FOR USE

#### FUNCTION

Test for qualitative detection of the presence of Hepatitis B virus surface antigen in biological samples (serum and plasma collected with EDTA or Heparin), through enzyme immunoassay. For *in vitro* diagnostic use only.

#### PRINCIPLE OF ACTION

**Methodology:** Enzyme immunoassay or immunoenzymatic assay based on the "sandwich" principle for the detection of HBsAg in human serum or plasma (EDTA or Heparin). The microplate is coated with monoclonal antibodies specific for the Ad and Ay subtypes of HBsAg. First, the sample and conjugate are added to the antibody-coated microplate and incubated. If the sample contains HBsAg, it will bind to the antibodies coated on the microplate and, simultaneously, will bind to the conjugate to form Antibody-HBsAg-Conjugate complexes. After incubation, the microplate is washed to remove unbound materials. After this step, the Substrate is added and incubated, producing a blue color that indicates the presence of Hepatitis B antigen in the sample. Stop Solution is added to stop the reaction, resulting in a color change from blue to yellow, measured on a microplate reader.

#### REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Microplate coated with anti-HBs antibodies (mouse monoclonal).
- 2- Conjugate** - Store between 2 and 8 °C. Contains: Anti-HBs antibodies linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Wash** - Store between 2 and 8 °C. Contains: Buffer solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4- Substrate** - Store between 2 and 8 °C. Contains: Buffer solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 5- Stop Solution** - Store between 2 and 8 °C. Contains: 1 N sulfuric acid.
- 6- Negative Control** - Store between 2 and 8 °C. Contains: Non-reactive solution for HBsAg, HCV, HIV-1 and HIV-2, stabilizer and preservative. **Potentially infectious.**
- 7- Positive Control** - Store between 2 and 8 °C. Contains: Reactive solution for HBsAg and non-reactive for HCV, HIV-1 and HIV-2, dye, stabilizer and preservative. **Potentially infectious.**
- 8- Plate Sealers.**

#### PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
<b>1- Sensitized Plate</b>	1 Unit	2 Units	5 Units
<b>2- Conjugate</b>	1 Vial x 8 mL	2 Vials x 8 mL	5 Vials x 8 mL
<b>3- Concentrated Washing</b>	1 Vial x 50 mL	2 Vials x 50 mL	5 Vials x 50 mL
<b>4- Substrate</b>	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
<b>5- Stop Solution</b>	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
<b>6- Negative Control</b>	1 Vial x 1 mL	2 Vials x 1 mL	5 Vials x 1 mL
<b>7- Positive Control</b>	1 Vial x 1 mL	2 Vials x 1 mL	5 Vials x 1 mL
<b>8- Plate Sealers</b>	3 Units	5 Units	10 Units

#### EQUIPMENTS AND OPERATIONAL INPUTS

##### Materials in the kit:

- Reagents described in the above table

##### Required materials not contained in the kit:

- Pipettes capable of dispensing volumes from 5 to 300 µL with a coefficient of variation lower than 1.5%.
- Repipette for repetitive pipetting of volumes of 300 µL with a coefficient of variation lower than 1.5% or multichannel pipette (optional).
- Microplate washer (optional).
- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- Absorbent paper to dry the microwells.
- Timer or clock.
- Vial to store the Washing Solution after dilution.
- Distilled or deionized water.
- Quality Control Tools.
- Incubator at 37 °C ± 2 °C.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8 °C. Transportation at room temperature (up to 30 °C) should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

#### SPECIAL CARE

- For professional *in vitro* diagnostic use only.**
- Strictly follow the proposed methodology to obtain accurate results.
- The sachet containing the microplate should only be opened after it has reached room temperature. Replace the unused microwell strips in the sachet, seal and store between 2 and 8 °C.
- The water used to clean the material must be fresh and free of contaminants.
- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- The Stop Solution contains Sulfuric Acid, which is a strong acid. Therefore, handle it with due care.
- All raw materials of the product are tested and must be non-reactive for Anti-HIV 1&2 and Anti-HCV. However, these tests

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

do not offer total assurance of the absence of infectious agents. Handling any product containing serum is potentially capable of transmitting diseases. Therefore, it is necessary to take appropriate biosafety precautions when handling these products.

**8-** Always pipette the reagents in the same order to minimize the difference in reaction time between the microcavities.

**9-** As a protective measure, the plate should be covered during the reaction.

**10-** Ensure that the bottom of the well is clean and dry, and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.

**11-** Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation stages.

**12-** We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is done in accordance with current legislation.

**13-** To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or through a request through SAC (Customer Advisory Service) from Quibasa.

**14-** Do not use the product if the packaging is damaged.

**15-** It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

#### SAMPLES

##### Serum or plasma (EDTA or Heparin).

Hemolyzed and highly lipemic samples should not be used. Samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at a temperature of -20 °C. Plasma samples stored for periods longer than the recommended period (30 days) may present fibrin and fibronectin precipitation that may interfere with the test. Samples collected with anticoagulants other than those mentioned should not be used.

#### PROCESS DESCRIPTION

##### Stability After Opening

The results of the stability test show that the BIOLISA HBsAg kit is stable after opening for up to 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

#### PREPARATION OF WORKING REAGENT

##### Washing Solution

Dilute the contents of Vial No. 3 (Concentrated Wash) in a ratio of 1:20 with distilled or deionized water; for example: 50 mL of Concentrated Washing Solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

**Note:** Small gelatinous aggregates may form in reagent No. 3 (Concentrated Wash), which are easily dissolved after dilution. Therefore, if this characteristic is observed in the reagent, proceed normally as described above.

**All other reagents are ready for use.**

#### TECHNIQUE

##### For use in automatic equipment, consult the Customer Service Department (SAC).

Before starting the test, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

- Separate the wells to be used considering: Controls and

Samples (it is recommended to test in duplicate). Return the unused strips of the microplate to the original sealed packaging.

**2-** Separate the first well for the Blank (OPTIONAL).

**3-** Pipette 50 µL of Negative Control, Positive Control and Samples into the previously determined wells.

**4-** Pipette 50 µL of Conjugate into all wells, except the Blank well (if you have chosen to use it).

**5-** Gently mix for ± 30 seconds. Cover the wells with plate sealer.

**6-** Incubate for 80 minutes ± 2 minutes at 37 °C ± 2 °C.

**7-** Remove the plate sealer from the wells.

**8-** Discard the contents of the wells by aspiration (Washer) or by decantation (Manual). Use approximately 350 µL of Washing Solution, previously prepared, and perform a total of eight (8) washing cycles with shaking for 5 seconds.

To ensure that the plate is dry, at the end of the wash, tap the plate for a few seconds on absorbent paper.

**Note:** Inadequate washing/drying may cause inadequate results.

**9-** Pipette 100 µL of Substrate into all wells, including the Blank well.

**10-** Gently mix for ± 30 seconds. Cover the wells with the plate sealer.

**11-** Incubate for 10 minutes ± 2 minutes at 37 °C ± 2 °C.

**12-** Remove the plate sealer from the wells.

**13-** Pipette 100 µL of Stop Solution into all wells.

**14-** Gently mix for ± 30 seconds.

**15-** Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

#### TECHNIQUE VERIFICATION

Check whether the results obtained for reading the Controls are compatible with the values presented below:

ITEM	ABSORBANCE
Blank	< 0.100
Negative Control	< 0.100
Positive Control	> 1.000

**Note:** If the values are found outside the specified range, the technique must be repeated. The values of the Controls and Samples must not be subtracted from the Blank value, if the option to use it has been chosen.

#### CALCULATIONS

##### QUALITATIVE

Calculate the Cut-Off according to the following formula:

Cut-Off: Average Absorbance Obtained with the Negative Control + 0.050.

Example:

ITEM	ABSORBANCE
Negative Control	A1 = 0.040
	A2 = 0.041
Cut-Off: (Average Absorbance of Negative Control) + 0.050	(0.040 + 0.041) / 2 + 0.050 = 0.091

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

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#### INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	INDEX
Negative (No Reactive)	< 0.9
Indeterminate	Between 0.9 and 1.1
Positive (Reactive)	> 1.1

#### NO REACTIVE

A sample with absorbance lower than the cut-off or index lower than 0.9 is considered non-reactive for Hepatitis B antigen and can be considered negative.

#### REACTIVE

A sample with an absorbance higher than the cut-off or an index higher than 1.1 is initially considered reactive for Hepatitis B antigen. The sample should be reanalyzed in duplicate before the end of the interpretation. A sample that is reactive in at least one of the reanalysis is presumed to be positive and should be confirmed using a confirmatory method.

**Nota:** In the case of an indeterminate result, the sample should be reanalyzed. Samples that repeatedly yield indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected should prevail. The interpretation of a diagnostic test should not be based on a single assay. Other confirmatory tests should be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results should be interpreted in conjunction with other available clinical information before a descriptive diagnosis of the disease is made.

**Nota:** The data presented in the examples are for illustration purposes only and cannot be used to calculate results.

#### LIMITATIONS OF THE PROCESS

The interpretation of a diagnostic test should not be based on a single assay. All results should be interpreted in conjunction with other available clinical information before a definitive diagnosis is made. The results provided by this kit should be interpreted by the responsible medical professional and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

**Nota:** Be careful not to touch the inner surface of the well with the pipette tip when the conjugate is dispensed. Contamination may occur.

#### INTERFERENTS

No interference was observed for Triglycerides 1500 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatinine 200 mg/dL, Bilirubin 1 g/dL, Albumin 10 g/dL, Hemoglobin 1000mg/dL, Oxalic Acid 60 mg/dL, Rheumatoid Factor 980 UI/mL, C-Reactive Protein 41.2 mg/dL and Anti-Streptolysin O 1023 UI/mL. Plasma samples collected in EDTA and Heparin can be used. Samples collected with anticoagulants other than those mentioned cannot be used. Plasma samples stored for periods longer than recommended (30 days) may present fibrin and fibronectin precipitation that can interfere with the test.

#### CROSS-REACTIVITY

A cross-reactivity study was performed, evaluating 08 samples positive for HTLV, 08 samples positive for HIV, 08 samples

positive for CMV, 08 samples positive for Chagas disease, 08 samples positive for Zika, 08 samples positive for Rubella, 08 samples positive for Toxoplasmosis and 08 samples positive for Syphilis. No cross-reactivity was observed with samples positive for HTLV, HIV, CMV, Chagas disease, Zika, Rubella, Toxoplasmosis and Syphilis. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data together with other laboratory data.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to emphasize that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow the assessment of the precision and accuracy of the dosages.

#### PRODUCT PERFORMANCE

##### PRECISION

##### Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	3
Average	1.709	1.500	0.019
Standard Deviation	0.075	0.119	0.002
Coefficient of Variation (%)	4.40	7.94	11.66

##### Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Average	1.575	1.670	0.007
Standard Deviation	0.057	0.089	0.001
Coefficient of Variation (%)	3.6	5.5	15.1

##### Analytical Sensitivity

The analytical sensitivity of the BIOLISA HBsAg kit was determined using the International Standard (HBV genotype B4, HBsAg subtypes ayw1/adw2) WHO NIBSC 12/226. The analytical sensitivity found was 0.05 UI/mL or 0.2 ng/mL.

##### Clinical Sensitivity and Specificity

The BIOLISA HBsAg kit analyzed clinical samples in comparison with another EIA kit. The results show that the clinical sensitivity of the BIOLISA HBsAg kit is >99.9% and the clinical specificity is >99.9%.

RESULT		REFERENCE		
		Positive	Negative	Total
BIOLISA HBsAg	Positive	96	0	96
	Negative	0	154	154
	Total	96	154	250

Clinical Sensitivity: >99.9% (96/96) - IC 95% (96.2% - 100%)

Clinical Specificity: >99.9% (154/154) - IC 95% (97.6% - 100%)

#### DIAGNOSTIC SIGNIFICANCE

Hepatitis B is an infectious disease caused by the HBV virus that affects the liver. The HBV virus is transmitted through sexual contact, blood-borne exposure (including sharing sharp objects), and mother-to-child transmission via the placenta or during childbirth. HBV infection results in a number of apparent serological markers, and one of the earliest of these markers is Hepatitis B Surface Antigen (HBsAg). The major subtypes of HBsAg include ad and ay, all sharing the common determinant 'a'. HBsAg appears between 1 and 10 weeks after infection and before biochemical evidence of liver disease or jaundice. However, in most people, due to its high production, HBsAg can be detected after 30 days of infection. This period is known as the diagnostic window. In the chronic carrier state, HBsAg persists for more than 6 months. Antibodies against HBV antigens appear between 30 and 60 days after infection, varying between individuals. This is the period known as the immunological window. HBV infection can cause acute self-limiting infection, fulminant hepatitis, chronic hepatitis with progression to cirrhosis and liver failure, and asymptomatic chronic carrier states. The HBV virus in infected individuals persists for the rest of their lives and can be transmitted to others. Thus, hepatitis B has become a public health problem. Therefore, screening for HBsAg is recommended for all donors, pregnant women, and individuals in high-risk groups.

#### BIBLIOGRAPHIC REFERENCES

- ALLEN JUNIOR, L. V.; ANSEL, H. C. Ansel's Pharmaceutical Dosage Forms and Drug Delivery Systems. 10. ed. Baltimore: Wolters Kluwer, 2014. 809 p.
- CLINICAL AND LABORATORY STANDARD INSTITUTE. EP-25A: Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline. Wayne: Clinical And Laboratory Standard Institute, 2009. 52 p.
- DESHPANDE, S.S. Enzyme Immunoassays: from concept to development. New York: Chapman & Hall, 1996. 472 p.
- TSITSILONIS, Ourania E.; THRASYVOULIDES, Apollon; BALAFAS, Apostolos; VOUTSAS, John F.; PAPAMICHAIL, Michael; LYMBERI, Peggy. Serological detection of hepatitis B viral infection by a panel of solid-phase enzyme-linked immunosorbent assays (ELISA). Journal Of Pharmaceutical And Biomedical Analysis, [S.L.], v. 34, n. 4, p. 811-822, mar. 2004. Elsevier BV. [http://dx.doi.org/10.1016/s0731-7085\(03\)00563-6](http://dx.doi.org/10.1016/s0731-7085(03)00563-6).
- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento

#### QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 (31) 3439.5454 | E-mail: [bioclin@bioclin.com.br](mailto:bioclin@bioclin.com.br)  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

ANVISA registration for BIOLISA HBsAg kit: 10269360288

Review: February/2025

### UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER