BIOLISA HCV DBS

[REF] K245

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Test for qualitative determination of Total Antibodies against the Hepatite C Virus (HCV) in **dried blood spot collected on filter paper (DBS)**. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzyme.
BIOLISA HCV DBS is a solid-phase immunoenzymatic assay, based on the "sandwich" principle, for the qualitative detection of Antibodies against Hepatitis C Virus (HCV) in dried blood spot collected on filter paper. The microplate is coated with HCV Recombinant Antigens (Core, NS3, NS4 and NS5). At first, the dried blood spot is eluted by the Elution Buffer and incubated with Conjugate 1, consisting of HCV Recombinant Antigens conjugated to Biotin. Anti-HCV Antibodies present in the sample bind to HCV Antigens immobilized on the microplate and to biotinylated Antigens present in Conjugate 1, forming the Antigen-Antibody-biotinylated Antigen immunocomplex. After the initial incubation, the microplate is washed to remove unbound materials. Conjugate 2, consisting of Streptavidin conjugated to Peroxidase, is added to the microplate, which is then incubated. Streptavidin binds to biotinylated Antigens, in turn bound to Anti-HCV Antibodies. A new wash is carried out to remove the excesses. After this step, the Substrate is added and incubated, producing a blue color that indicates the detection of Anti-HCV Antibodies present in the sample. Stop Solution is added to stop the reaction, with a color change from blue to yellow measured on a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Plate sensitized with HCV Recombinant Antigens.
- 2- Conjugate 1** - Store between 2 and 8 °C. Contains: Buffer Solution containing HCV Recombinant Antigens linked to Biotin, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** - Store between 2 and 8 °C. Contains: Buffer Solution (Phosphate < 0.5 mol/L, Potassium Chloride < 100 mmol/L, Sodium Chloride < 5 mol/L), surfactant and preservative.
- 4- Conjugate 2** - Store between 2 and 8 °C. Contains: Buffer Solution containing Streptavidin linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- 5- Substrate** - Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8 °C. Contains: 1 M Hydrochloric Acid.
- 7- Elution Buffer on Filter Paper** - Store between 2 and 8 °C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 8- DBS Negative Control** - Store between 2 and 8 °C. Contains: Non-Reactive Sample for Total Anti-HCV Antibody impregnated on filter paper. **Potentially infective.**
- 9- DBS Positive Control** - Store between 2 and 8 °C. Contains: Reagent Sample for Total Anti-HCV Antibody impregnated on filter paper. **Potentially infective.**

PRESENTATION

REAGENTS	1	2	3	4
	96 Cavities	192 Cavities	480 Cavities	960 Cavities
1- Sensitized Plate	1 Unit	2 Units	5 Units	10 Units
2-Conjugated 1	1 Vial x 6 mL	1 Vial x 12 mL	1 Vial x 30 mL	1 Vial x 60 mL
3-Concentrated Washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL	1 Vial x 500 mL
4- Conjugated 2	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL	1 Vial x 120 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL	1 Vial x 120 mL
6- Stop Solution	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL	1 Vial x 120 mL
7-Elution Buffer on Filter Paper	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL	1 Vial x 120 mL
8- DBS Negative Control	1 Unit	2 Units	3 Units	5 Units
9- DBS Positive Control	1 Unit	2 Units	3 Units	5 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

- Pipettes capable of dispensing volumes from 50 to 500 µL with a coefficient of variation less than 1.5%.
- Repeater for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- Microplate washer for technique on filter paper.
- ELISA reader with absorbance capacity at 450 and 620 nm of wavelength.
- Absorbent paper to dry the microcavities.
- Stopwatch or watch.
- Bottle to store the Washing Solution after dilution.
- Distilled or deionized water.
- Quality Control Tools.
- Incubator 37 ± 2 °C.
- Paper shredder (diameter 3 ± 0.2 mm).

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be between 2 and 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- For *in vitro* diagnostic use only.**
- Strictly follow the proposed methodology to obtain accurate results.
- The sachet containing the microplate should only be opened after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- The water used to clean the material must be recent and free of contaminants.
- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.

- All raw material of the product is tested, being non-reactive for HBsAg and Anti-HIV. However, these tests do not offer complete assurance of the absence of infectious agents. Handling any product that contains a biological sample is potentially capable of transmitting diseases. Therefore, it is necessary to take the proper biosafety precautions when handling these products.
- Always pipette the reagents in the same order to minimize the difference in reaction time between the microwells.
- As a protection measure, the plate must be covered during the reaction.
- It must be ensured that the bottom of the well is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the assay.
- Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.
- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by QUIBASA's SAC (Customer Advisory Service).
- Do not use the product if the packaging is damaged.
- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Dried Blood collected on filter paper (DBS) - EDTA or Puncture

The filter paper used must be specific for collecting dried blood spot, such as the S&S 903 and Ahlstrom 226 models. The laboratory must certify the quality of the samples before use.
Blood samples dried on filter paper can be stored for up to 35 days at room temperature, provided they are protected from direct sunlight and in low humidity. For storage of up to 2 years, the samples must remain refrigerated, between 2 and 8 °C. Storage for a period longer than 2 years must be carried out at a temperature of -20 °C.⁷

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test prove that the BIOLISA HCV DBS kit is stable after opening for up to 30 days. This stability may vary according to test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of vial Nº 3 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Washing solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30°C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolution.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, consult the SAC (Customer Advisory Service).

Before starting the assay, place all Reagents, Samples and Controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

- Separate the wells to be used considering: DBS Positive Control, DBS Negative Control and dried blood spot collected on filter paper (it is recommended to test in duplicate). Return unused microplate strips to the original sealed package.
- Separate the first cavity for the Blank (OPTIONAL).
- Add a 3 mm disk of DBS Positive Control, DBS Negative Control and dried blood spot collected on filter paper, previously perforated, in the determined wells.

- Pipette 50 µL of Elution Buffer into all wells, **including the Blank well.**
- Pipette 50 µL of Conjugate 1, **including in the Blank well.**
- Gently homogenize for ± 30 seconds. Cover wells with plate sealer.
- Incubate for 60 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- After incubation, discard the contents of the wells by aspiration (Washer suitable for the filter paper technique). Use approximately 300 µL of Washing Solution, **previously prepared**, and perform a total of five (5) wash cycles. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
- Poor washing/drying may cause poor results.
- Pipette 100 µL of Conjugate 2 into all wells, **including the Blank well.**
- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- Remove the plate sealer from the wells.
- Repeat item 8.
- Pipette 100 µL of Substrate in all wells, **including the Blank well.**
- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.
- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- Remove the plate sealer from the wells.
- Pipette 100 µL of Stop Solution into all wells.
- Gently homogenize for ± 30 seconds.
- Read using double filter: 450 nm / 620 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	< 0.120
DBS Negative Control	< 0.120
DBS Positive Control	> 0.500

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula:

Cut Off = Mean Absorbance of DBS Negative Control + 0.070.

Example:

ITEM	ABSORBANCES
DBS Negative Control	0.054
	0.057
Cut Off = (DBS Negative Control Mean Absorbance) + 0.070	((0.054 + 0.057) / 2) + 0.070 = 0.126

Calculate the Index by dividing the Sample absorbance by the Cut Off value.

Example:

ITEM	ABSORBANCES
Sample	1.023
Cut Off Value	0.126
Index = Sample / Cut Off Value	1.023 / 0.126 = 8.12

Note: Data shown in examples is for illustration only and cannot be used for calculation of results. Each laboratory must validate the cut-off according to the instrumentation used and the researched population.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (Non-Reagent)	≤ 0.8
Undetermined	0.8 and 1.2
Positive (Reagent)	≥ 1.2

Notes: In case of indeterminate results, the sample must be reanalyzed in duplicate. Samples that give repeatedly indeterminate results should be retested using an alternate method. If results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected shall prevail. The interpretation of a diagnostic test should not be based on a single test. Additional confirmatory tests must be included before a specimen is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, prior to making a descriptive diagnosis of the disease. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERING

No interference was observed for Triglycerides up to 1500 mg/dL, Acetylsalicylic Acid up to 20 mg/dL, Ascorbic Acid up to 2 g/dL, Creatine up to 200 mg/dL, Bilirubin up to 1 g/dL, Albumin up to 10 g/dL, Hemoglobin up to 1000 mg/dL, Oxalic Acid up to 60 mg/dL, Rheumatoid Factor up to 980 IU/mL, C-Reactive Protein up to 41.2 mg/dL and Anti Streptolysin O up to 1023 IU/mL.

CROSS REACTIVITY

A cross-reactivity study was performed, evaluating 88 samples negative for HCV but positive for other infections. Among them 5 positive samples for HIV, 9 positive samples for HTLV, 9 positive samples for Syphilis, 9 positive samples for HBsAg, 10 positive samples for Chagas Disease, 6 positive samples for Zika, 10 positive samples for CMV, 10 positive samples for Rubella, 10 positive samples for Toxoplasmosis and 10 positive samples for SARS-CoV-2. No cross-reactivity was observed with positive samples for HIV, HTLV, Syphilis, HBsAg, Chagas Disease, Zika, CMV, Rubella, Toxoplasmosis and SARS-CoV-2. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Sample		
	1	2	3
Average	2.159	1.709	0.038
Standard Deviation	0.049	0.067	0.007
Coefficient of Variation (%)	2.26	3.94	18.51

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	2.340	1.899	0.042
Standard Deviation	0.135	0.081	0.006
Coefficient of Variation (%)	5.68	4.33	14.86

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA HCV DBS kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOLISA HCV DBS kit is >99.9% and the clinical specificity is 99.43%.

		Result Reference		
		Positive	Negative	Total
BIOLISA HCV DBS	Positive	83	1	84
	Negative	0	175	175
	Total	83	176	259

Clinical Sensitivity: >99.9% (83/83) - CI 95% = 95.7 to 100%.

Clinical Specificity: 99.4% (175/176) - CI 95% = 96.9 to 100%.

DIAGNOSTIC SIGNIFICANCE

Hepatitis C, formerly called non-A and non-B hepatitis, is a disease caused by parenteral transmission of the Hepatitis C Virus (HCV). HCV is an enveloped virus consisting of only one RNA molecule in its genome, which encodes 3 structural proteins (a capsid protein, or core, and two envelope proteins, E1 and E2) and 7 non-structural proteins (NS1, NS2, NS3, NS4A, NS4B, NS5A and NS5B). The main routes of transmission of the virus are through transfusion of blood and blood products, organ transplantation and sharing of contaminated needles and syringes. Most cases of HCV infection are asymptomatic and about 15 to 20% have nonspecific symptoms such as fatigue, loss of appetite, weight loss, fever and abdominal pain. About 50 to 80% of cases of HCV infection become chronic and may progress to liver fibrosis, cirrhosis and hepatocellular carcinoma. Serological diagnosis of Hepatitis C is the first step towards identifying infected individuals, and early diagnosis is the main means of reducing morbidity and mortality caused by HCV.

BIBLIOGRAPHIC REFERENCES

- 1- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR
- 2- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 3- Warkad, S. D. et al. Developments in the HCV Screening Technologies Based on the Detection of Antigens and Antibodies. Sensors. 2019;19:4257.
- 4- Ministério da Saúde. Manual Técnico para o Diagnóstico da Infecção pelas Hepatites Virais. 2018.
- 5- Van der Poel, C. L., H.T.M. Cuypers, H.W. Reesink, and P.N.Lelie. Confirmation of Hepatitis C Virus Infection by New Four-antigen Recombinant Immunoblot Assay. Lancet. 1991;337:317.
- 6- Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.
- 7- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST,Aids e Hepatites Virais.
- 8- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of the BIOLISA HCV DBS kit at ANVISA:
10269360435

Review: July/2023

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER