BIOLISA HCV

REF K128

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FINALIDADE

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1- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR

2- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.

3- Warkad, S. D. et al. Developments in the HCV Screening Technologies Based on the Detection of Antigens and Antibodies. Sensors. 2019;19;4257.

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6- Wilber, J.C. Development and Use of Laboratory Tests for Hepatitis C Infection: A Review. J. Clinical Immunoassay. 1993;16:204.

7- TELELAB. Manual de Coleta de Sangue - Diagnóstico e monitoramento das DST, Aids e Hepatites Virais.

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INSTRUCTIONS FOR USE

FUNCTION

Test for qualitative determination of total antibodies against Hepatitis C Virus (HCV) in biological samples (serum or plasma) through fourth-generation enzyme immunoassay test. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA HCV is a fourth-generation solid-phase enzyme immunoassay based on the "sandwich" principle for the qualitative detection of total antibodies against Hepatitis C Virus (HCV) in human serum and plasma samples. The microplate is coated with recombinant HCV antigens (Core, NS3, NS4 and NS5). At first, the sample is incubated with Conjugate 1, consisting of recombinant HCV antigens conjugated to biotin. Anti-HCV antibodies present in the sample bind to the HCV antigens immobilized on the microplate and to the biotinylated antigens present in Conjugate 1, forming the antigen-antibody-biotinylated antigen immune complex. After the initial incubation, the microplate is washed to remove unbound materials. Conjugate 2, consisting of Streptavidin conjugated to Peroxidase, is added to the microplate, which is then incubated. Streptavidin binds to biotinylated antigens bound to anti-HCV antibodies. A new wash is carried out to remove the excess. After this step, the Substrate is added and incubated, producing a blue color that indicates the detection of anti-HCV antibodies present in the sample. The Stop Solution is added to stop the reaction, with a color change from blue to yellow, measured in a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Plate sensitized with recombinant HCV antigens.
- 2- Conjugate 1** - Store between 2 and 8 °C. Contains: Buffer solution containing recombinant HCV antigens linked to biotin, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing Solution** - Store between 2 and 8 °C. Contains: Buffer solution, surfactant and preservative.
- 4- Conjugate 2** - Store between 2 and 8 °C. Contains: Buffer solution containing streptavidin linked to peroxidase, surfactant, stabilizers, colorant and preservative.
- 5- Substrate** - Store between 2 and 8 °C. Contains: Buffer solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store at 2 to 8 °C. Contains: 1 M Hydrochloric Acid.
- 7- Negative Control** - Store between 2 and 8 °C. Contains: Buffer solution, stabilizer, surfactant and preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8 °C. Contains: Buffer solution, total anti-HCV antibodies, colorant, stabilizers, surfactant and preservative. **Potentially infectious.**
- 9- Plate Sealers**

PRESENTATION

REAGENTS	1	2	3
	96 CAVITIES	192 CAVITIES	480 CAVITIES
1- Sensitized Plate	1 Unit (96 cavities)	2 Units (96 cavities)	5 Units (96 cavities)
2- Conjugate 1	1 Vial x 6 mL	2 Vials x 6 mL	5 Vials x 6 mL
3 - Concentrated Washing Solution	1 Vial x 50 mL	2 Vials x 50 mL	5 Vials x 50 mL
4- Conjugate 2	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
6- Substrate	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
7- Stop Solution	1 Vial x 12 mL	2 Vials x 12 mL	5 Vials x 12 mL
8- Negative Control	1 Vial x 1 mL	2 Vials x 1 mL	5 Vials x 1 mL
9- Positive Control	1 Vial x 1 mL	2 Vials x 1 mL	5 Vials x 1 mL
10- Plate Sealers	3 Units	6 Units	15 Units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the previous item.

Necessary materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes of 5 to 500 µL with a variation coefficient of less than 1.5%.
- 2- Repipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer (optional).
- 4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5- Absorbent paper to dry the wells.
- 6- Stopwatch or watch.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator at 37 °C ± 2 °C.

TRANSPORTATION AND STORAGE CONDITIONS

Storage temperature should be between 2 and 8 °C. Transport at temperatures up to 30°C should not exceed 5 days. Keep out of the light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- For in vitro diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The sachet containing the microplate must be opened only after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be fresh and free from contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw material of the product is tested and must be non-reactive for HBsAg and Anti-HIV 1 & 2. However, these tests do not offer complete security in the absence of infectious agents. The handling of any product that contains serum is potentially capable of transmitting disease. Therefore, it is necessary to take due care of biosafety in the handling of these products.
- 8- Pipette the reagents always in the same order to minimize the difference in reaction time between the wells.

- 9- As a protection measure, the plate must be covered during the reaction.
- 10- Make sure that the bottom of the cavity is clean and dry and that there are no bubbles on the liquid surface before reading the plate. Do not allow the wells to dry out during the test.
- 11- Do not expose the reagents, especially the Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.
- 12- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Service of Chemicals) Customer Service) of Quibasa.
- 14- Do not use the product in case of damage to the packaging.
- 15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolyzed or highly lipemic samples should not be used. The samples can be kept refrigerated, between 2 and 8 °C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at a temperature of -20 °C.

PROCESSS DESCRIPTION

Stability After Open

The results of the stability test show that the BIOLISA HCV kit is stable after being opened for up to 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the performance of the product using internal controls of the kit and the criteria for validating the technique.

PREPARATION OF WORK REAGENTS

WASHING SOLUTION

Dilute the contents of vial N°3 (Concentrated Washing Solution) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C for at least 30 days. It can be stored at room temperature. If crystallization occurs, heat to 37 °C until dissolved.

SUBSTRATE

The Substrate is ready for use.

TECHNIQUE

For use in automatic equipment, consult the Customer Service Department (SAC)

- Before starting the assay, allow all reagents, Controls and Samples to stabilize at room temperature (15 – 30°C) for at least 40 minutes.
- 1- Separate the wells to be used considering: Controls and Samples (it is recommended to test in duplicate). Return unused microplate strips to the original sealed packaging.
- 2- Separate the first cavity for Blank (OPTIONAL).
- 3- Pipette 50 µL of the Samples, Positive Control and Negative Control in the previously determined wells.
- 4- Pipette 50 µL of Conjugate 1, including in the Blank cavity.
- 5- Homogenize gently for ± 10 seconds. Cover wells with plate sealer.
- 6- Incubate for 60 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 7- Remove the sealant from the cavities.
- 8- Discard the contents of the cavities by aspiration (Washer) or by decantation (Manual). Use approximately 300 µL of previously prepared Wash Solution and perform a total of five (5) wash cycles. To guarantee the plate drying, at the end of washing, tap the plate for a few seconds on absorbent paper.

Note: Poor washing/drying may cause inappropriate results.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

- 9- Pipette 100 µL of Conjugate 2 into all wells, including the Blank well.
- 10- Mix gently for ± 10 seconds. Cover the wells with the plate sealer.
- 11- Incubate for 30 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 12- Remove the plate sealer from the wells.
- 13- Repeat item 8.
- 14- Pipette 100 µL of Substrate into all wells, including the Blank well.
- 15- Mix gently for ± 30 seconds. Cover the cavities with the plate sealer.
- 16- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.
- 17- Remove the plate sealer from the wells.
- 18- Pipette 100 µL of Stop Solution into all wells, including the Blank well.
- 19- Mix gently for ± 30 seconds.
- 20- Read using a double filter: 450 nm / 630 nm within 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	<0.120
Negative Control	<0.120
Positive Control	>1.200

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula: Cut Off = (Average Positive Control Absorbance x 0.02) + 0.080

Example:

ITEM	ABSORBANCES
Positive Control	A1 = 2.045
	A2 = 2.033
Cut Off = (Average Positive Control Absorbance x 0.02) + 0.080	((2.045 + 2.033) / 2) x 0.02 + 0.080 = 0.121

Calculate the Index by dividing the Sample's absorbance by the Cut Off value.

Example:

ITEM	ABSORBANCES
Sample	2.038
Cut Off Value	0.121
Index = Samples / Cut Off Value	2.038 / 0.121 = 16.84

Note: The data presented in the examples are for illustration only and cannot be used to calculate the results.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indices below to determine the results.

RESULTS	QUALITATITE FOR SERUM AND PLASMA SAMPLES
	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	Between 0.9 - 1.1

Note: In the case of an undetermined result, the sample must be reanalyzed. Samples that obtain repeatedly indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be collected in two weeks. The result of the last collected sample must prevail. The results provided by this kit must be interpreted by the responsible medical professional, and it is not the only criterion for determining the diagnosis and / or treatment of the patient.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, before the descriptive diagnosis of the disease.

INTERFERENTS

No interference was observed for Triglycerides up to 1500 mg/dL, Acetylsalicylic Acid up to 20 mg/dL, Ascorbic Acid up to 2 g/dL, Creatine up to 200 mg/dL, Bilirubin up to 1 g/dL, Albumin up to 2 g/dL, Hemoglobin up to 1000 mg/dL, Oxalic Acid up to 60 mg/dL, Rheumatoid Factor up to 980IU/mL, C Reactive Protein up to 41.2 mg/dL and Anti Streptolysin O up to 1023 IU/mL.

CROSS REACTIVITY

A cross-reactivity study was carried out, evaluating 124 samples negative for HCV, but positive for other infections. Among them 10 HIV positive samples, 10 HTLV positive samples, 10 HBsAg positive samples, 10 anti-HBs positive samples, 10 Syphilis positive samples, 10 Dengue positive samples, 10 Zika positive samples, 4 positive samples for Chikungunya, 10 samples positive for Chagas Disease, 10 samples positive for CMV, 10 samples positive for Rubella, 10 samples positive for Toxoplasmosis and 10 samples positive for H1N1. No cross-reactivity was observed with positive samples for HIV, HTLV, HBsAg, anti-HBs, Syphilis, Dengue, Zika, Chikungunya, Chagas Disease, CMV, Rubella, Toxoplasmosis and H1N1. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

ACCURACY

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	2
Average	1.709	1.500	0.019
Standard Deviation	0.075	0.119	0.002
Coefficient of Variation (%)	4.40	7.94	11.66

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	2
Average	1.575	1.679	0.041
Standard Deviation	0.057	0.105	0.005
Coefficient of Variation (%)	3.55	6.34	13.21

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA HCV kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA HCV kit is >99.9% and the clinical specificity is >99.9%.

CLINICAL SENSITIVITY AND SPECIFICITY	EXPECTED RESULT	BIOLISA HCV
Positive Sample	85	85
Negative Sample	152	152
Total Tested Samples	237	237

Clinical Sensitivity: >99.9% (85/85) - IC 95%: 95.8% - 100.0%

Clinical Specificity: >99.9% (152/152) - IC 95%: 97.6% - 100.0%

CLINICAL SIGNIFICANCE

Hepatitis C, formerly called non-A and non-B hepatitis, is a disease caused by the parenteral transmission of Hepatitis C Virus (HCV). HCV is an enveloped virus consisting of only one RNA molecule in its genome, which encodes 3 structural proteins (a capsid or core protein and two envelope proteins, E1 and E2) and 7 non-structural proteins (NS1, NS2, NS3, NS4A, NS4B, NS5A and NS5B). The main routes of transmission of the virus are through transfusion of blood and blood products, organ transplantation and sharing of contaminated needles and syringes. Most cases of HCV infection are asymptomatic and about 15 to 20% have nonspecific symptoms such as fatigue, loss of appetite, weight loss, fever and abdominal pain. About 50 to 80% of cases of HCV infection become chronic and can progress to liver fibrosis, cirrhosis and hepatocellular carcinoma. Serological diagnosis of Hepatitis C is the first step in identifying infected individuals, with early diagnosis being the main means of reducing morbidity and mortality caused by HCV.

BIBLIOGRAPHIC REFERENCES

1- Cancer Epidemiol Biomarkers Prev; 26(8); 1337–44.2017 AACR

2- WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and sérum samples. WHO/DIL/LAB/99.1 rev. 2, 2002: 31.

3- Warkad, S. D. et al. Developments in the HCV Screening Technologies Based on the Detection of Antigens and Antibodies. Sensors. 2019;19;4257.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 - E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 - E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA HCV kit: 10269360304

Review: September/2024

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER