BIOLISA HIV 1/2/O Ag/Ab DBS

REF K250

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Fourth generation sandwich enzyme-linked immunosorbent assay (ELISA) for the qualitative detection of the presence of p24 antigen, anti-HIV-1, anti-HIV-2, and/or subtype O in **dried blood samples collected on filter paper (DBS)**, for screening in adults and neonates, using an enzyme-linked immunosorbent assay. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic
The BIOLISA HIV 1/2/O Ag/Ab DBS kit is a fourth-generation solid-phase enzyme-linked immunosorbent assay based on the "sandwich" principle for the detection of p24 antigen, anti-HIV-1, anti-HIV-2 and/or subtype O antibodies in dried blood samples collected on filter paper. The microplate is coated with HIV-1 and HIV-2 specific antigens (HIV-1 gp41 and gp120 antigen and HIV-2 gp36 antigen) and HIV-1 anti-p24 antibody. Antibodies present in the dried blood sample are eluted and bind to the microplate coated with specific recombinant antigens, forming antigen-antibody immunocomplexes. Together with the first incubation, conjugate 1 is added, consisting of biotin-conjugated anti-HIV p24, which binds to the p24 antigens present in the sample, forming antibody-antigen-antibody immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Conjugate 2 containing HIV-1, HIV-2 and streptavidin antigens, all bound to peroxidase, is added to the microplate and then incubated. The antigens in conjugate 2 bind to the antibodies in the sample that have been captured by the antigens in the microplate, forming antigen-antibody-antigen immunocomplexes. Streptavidin linked to peroxidase, also present in conjugate 2, binds to biotin linked to the anti-HIV p24 of the antibody-antigen-antibody immunocomplexes. After these steps, another wash is carried out to remove unbound materials. The substrate is then added and incubated, producing a blue color that indicates the detection of anti-HIV antibodies and/ or p24 antigen present in the samples. The Stop Solution is added to stop the reaction with a color change from blue to yellow. The intensity of the color is measured at 450/620 nm.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8 °C. Contains: Microplate coated with HIV-1 and HIV-2 specific antigens (HIV-1 gp41 and gp120 antigen and HIV-2 gp36 antigen) and HIV-1 anti-p24 antibody.
- 2- Conjugate 1** - Store between 2 and 8 °C. Contains: Buffer solution containing biotin-linked anti-HIV p24 antibody, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Wash** - Store between 2 and 8°C. Contains: Buffer solution (phosphate < 0.5 mol/L, potassium chloride < 100 mmol/L, sodium chloride < 5 mol/L), surfactant and preservative.
- 4- Conjugate 2** - Store between 2 and 8 °C. Contains: Specific antigens for HIV-1, HIV-2 and streptavidin, all bound to peroxidase, surfactant, stabilizers, dye and preservative.
- 5- Substrate** - Store between 2 and 8 °C. Contains: Buffer solution containing urea peroxide, tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8 °C. Contains: 1 M hydrochloric acid solution.
- 7- Elution Buffer** - Store between 2 and 8 °C. Contains: Buffer solution, surfactant, stabilizer and preservative.
- 8- DBS Negative Control** - Store between 2 and 8 °C. Contains: HIV non-reactive solution impregnated on filter paper. **Potentially infectious.**
- 9- DBS Positive Control** - Store between 2 and 8 °C. Contains: HIV antibody solution impregnated on filter paper. **Potentially infectious.**
- 10- HIV-1 p24 Positive Control** - Store between 2 and 8 °C. Solution with p24 antigens, stabilizer, dye and preservative. **Potentially infectious.**

PRESENTATION

REAGENTS	1	2	3
	96 Cavities	192 Cavities	480 Cavities
1- Sensitized Plate	1 Unit	2 Units	5 Units
2- Conjugate 1	1 Vial x 6 mL	1 Vials x 12 mL	1 Vials x 30 mL
3- Concentrated Wash	1 Vial x 50 mL	1 Vials x 100 mL	1 Vials x 250 mL
4- Conjugate 2	1 Vial x 12 mL	1 Vials x 24 mL	1 Vials x 60 mL
5- Substrate	1 Vial x 12 mL	1 Vials x 24 mL	1 Vials x 60 mL
6- Stop Solution	1 Vial x 12 mL	1 Vials x 24 mL	1 Vials x 60 mL
7- Elution Buffer	1 Vial x 12 mL	1 Vials x 24 mL	1 Vials x 60 mL
8- DBS Negative Control	1 Unit	2 Units	3 Units
9- DBS Positive Control	1 Unit	2 Units	3 Units
10- HIV-1 p24 Positive Control	1 Vial x 0.5 mL	1 Vials x 1.0 mL	1 Vials x 1.0 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Materials in the kit:

- Reagents described in the above table

Required materials not contained in the kit:

- 1-** Pipettes capable of dispensing volumes of 50 to 500 µL with a coefficient of variation of less than 1.5%.
- 2-** Repipettor for repetitive pipetting of volumes of 300 µL with a coefficient of variation of less than 1.5% or multichannel pipette (optional).
- 3-** Microplate washer for filter paper technique.
- 4-** ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5-** Absorbent paper to dry the microcavities.
- 6-** Stopwatch or clock.
- 7-** Bottle to store the washing solution after dilution.
- 8-** Distilled or deionized water.
- 9-** Quality control tools.
- 10-** Incubator 37 °C ± 2 °C.
- 11-** Paper mincer (diameter 3 ± 0.2 mm).

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8 °C. Transportation at temperatures up to 30 °C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.**
- 2-** Strictly follow the proposed methodology in order to obtain accurate results.
- 3-** The sachet containing the microplate should only be opened after it has reached room temperature. Place the unused microcavity strips back in the sachet, seal and store between 2 and 8 °C.
- 4-** The water used to clean the material must be fresh and free of contaminants.

- 5-** Saturated deionizing columns release alkaline water, various ions, oxidizing and reducing agents which can significantly alter the results.
- 6-** The Stop Solution contains hydrochloric acid, which is a strong acid. Therefore, handle it with due care.
- 7-** All the product's raw materials are tested as non-reactive for HBsAg and Anti-HCV. However, these tests do not provide complete assurance of the absence of infectious agents. The handling of any product containing a biological sample is potentially capable of transmitting diseases. Biosafety precautions must therefore be taken when handling these products.
- 8-** Always pipette the reagents in the same order to minimize the difference in reaction time between the microcavities.
- 9-** As a protective measure, cover the plate during the reaction.
- 10-** Ensure that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the test.
- 11-** Do not expose reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.
- 12-** We recommend applying local, state and federal environmental protection regulations so that reagents and biological material are disposed of in accordance with current legislation.
- 13-** For information related to biosafety or in the event of accidents with the product, consult the SDS (Safety Data Sheet) available at www.bioclin.com.br or by contacting Quibasa's Customer Service Center.
- 14-** Do not use the product if the packaging is damaged.
- 15-** It is essential that the instruments and equipment used are properly calibrated and undergo periodic maintenance.

SAMPLES

Dried blood collected on filter paper (DBS) - EDTA, Puncture or Heel stick (for neonates)

The filter paper used must be specific for collecting dried blood samples, such as the S&S 903 and Ahlstrom 226 models. The laboratory must check the quality of the samples before use. Dried blood samples on filter paper can be stored for up to 35 days at room temperature, as long as they are protected from direct sunlight and in low humidity. For storage of up to 2 years, the samples should be kept refrigerated between 2 and 8 °C. Storage for longer than 2 years should be at -20 °C.⁷

PROCESS DESCRIPTION

Stability After Opening

The stability test results show that the BIOLISA HIV 1/2/O Ag/Ab DBS kit is stable after opening for at least 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the product's performance using the kit's internal controls and the technique's validation criteria.

PREPARATION OF WORKING REAGENT

Washing Solution

Dilute the contents of bottle No. 3 (Concentrated Wash) 1:20 with distilled or deionized water. In this way, dilute 50 mL of the Concentrated Wash solution in distilled or deionized water to a final volume of 1000 mL. After preparation, the solution can be stored between 2 and 30 °C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolved.

All other reagents are ready for use.

TECHNIQUE

For use in automatic equipment, consult the Customer Service Department (SAC).

Before starting the test, leave all reagents, samples and controls to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.

- 1-** Separate the wells to be used considering: p24 HIV-1 Positive Control, DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper (we recommend testing

in duplicate). Return the unused microplate strips to their original sealed packaging.

- 2-** Separate the first well for the blank (OPTIONAL).

IMPORTANT! The following steps must be carried out strictly in the sequence described below.

- 3-** Add a 3 mm disk of DBS Positive Control, DBS Negative Control and dried blood samples collected on filter paper, previously minced, to the determined wells.

- 4-** Pipette 50 µL of the p24 HIV-1 Positive Control into the previously determined well.

- 5-** Pipette 50 µL of Elution Buffer into all the wells **except** the well for the p24 HIV-1 Positive Control, including the well for the Blank.

- 6-** Pipette 50 µL of conjugate 1 into all wells, **including** the p24 HIV-1 Positive Control well and the Blank well.

- 7-** Homogenize gently for ± 30 seconds, cover the wells with a plate sealer.

- 8-** Incubate for 60 minutes ± 2 minutes in an incubator at (approximately 25 °C) and under constant agitation of 650 rpm, which can vary between 550 and 750 rpm.

Note: Stirring at the indicated speed is a crucial step for the efficient elution of dried blood samples and is necessary for the correct performance of the assay.

- 9-** After incubation, discard the contents of the wells by aspiration (washer for filter paper technique). Use approximately 300 µL of previously prepared washing solution and carry out a total of five (5) washing cycles with a three (3) second shake at the end of each cycle. To ensure that the plate is dry, at the end of the wash, tap the plate for a few seconds on absorbent paper.

Note: Poor washing/drying can cause inadequate results.

- 10-** Pipette 100 µL of Conjugate into all the wells, **including** the blank well.

- 11-** Homogenize gently for ± 30 seconds. Cover the wells with the plate sealer.

- 12-** Incubate for 60 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

- 13-** Remove the plate sealer from the wells.

- 14-** Repeat item 9.

- 15-** Pipette 100 µL of Substrate into all the wells. **Including** the blank cavity.

- 16-** Homogenize gently for ± 30 seconds. Cover the wells with the plate sealer.

- 17-** Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

- 18-** Remove the plate sealer from the wells.

- 19-** Pipette 100 µL of Stop Solution into all wells.

- 20-** Homogenize gently for ± 30 seconds.

- 21-** Read using a double filter: 450 nm / 630 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check that the results obtained for reading the Blank and the Reference Standards are compatible with the values shown below:

ITEM	ABSORBANCE
Blank	< 0.200
DBS Negative Control	< 0.200
DBS Positive Control	> 0.500
p24 HIV-1 Positive Control	> 0.500

If the values are outside the expected values, the technique should be repeated.

CALCULATIONS
QUALITATIVE

Calculate the Cut-Off according to the following formula:
Cut Off = Average Absorbance of Negative Control DBS + 0.180.

Example:

ITEM	ABSORBANCE
DBS Negative Control	0.086
	0.088
Cut Off = Average Absorbance of Negative Control DBS + 0.180	(0.086 + 0.088) / 2 + 0.180 = 0.267

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCIA
Sample	1.145
Cut-Off Value	0.267
Index = Sample / Cut-Off Value	1.145 / 0.267 = 4.288

Note: The data presented in the examples is for illustrative purposes only and cannot be used to calculate results.

Each laboratory must validate the cut-off according to the instrumentation used and the population surveyed.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	INDEX
Negative (No Reactive)	< 0.8
Indeterminate	Between 0.8 and 1.2
Positive (Reactive)	> 1.2

Non-Reactive: Sample with absorbance lower than the cut-off is considered non-reactive for p24 antigen, anti-HIV1, anti-HIV2 antibodies, and/or subtype O, and can be considered negative.

Reactive: Sample with absorbance higher than the cut-off is considered initially reactive for p24 antigen, anti-HIV1, anti-HIV2 antibodies, and/or subtype O. The sample must be re-analyzed in duplicate before the end of interpretation. A sample that is reactive in at least one of the reanalyses is presumed to be reactive and must be confirmed by another diagnostic method.

Observations: In the event of an indeterminate result, the sample must be re-analyzed in duplicate. Samples that repeatedly obtain indeterminate results should be retested using an alternative method. If the results remain indeterminate, a new sample should be taken within two weeks. The result of the last sample collected should prevail. The interpretation of a diagnostic test should not be based on a single test. Other confirmatory tests should be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information before a descriptive diagnosis of the disease can be made.

LIMITATIONS OF THE PROCESS

The interpretation of a diagnostic test should not be established on the basis of a single assay. All results must be interpreted in conjunction with other available clinical information before a definitive diagnosis can be made. The results provided by this kit must be interpreted by the medical professional responsible and are not the sole criterion for determining the diagnosis and/or treatment of the patient.

INTERFERENTS

No interference was observed for the concentrations of Triglycerides 1500 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, C-Reactive Protein 41.2 mg/dL, anti-Streptolysin O 1023 IU/mL and Rheumatoid Factor 1080 IU/mL.

CROSS-REACTIVITY

A study was carried out with 90 dried blood samples collected on filter paper that were negative for HIV but positive for other infections, in order to assess the possibility of cross-reactivity of these interferents in the BIOLISA HIV 1/2/O AG/AB DBS result. These included 10 samples positive for HTLV, 10 samples positive for Syphilis, 10 samples positive for HBsAg, 10 samples positive for HCV, 10 samples positive for Chagas Disease, 10 samples positive for Zika, 10 samples positive for CMV, 10 samples positive for Rubella and 10 samples positive for Toxoplasmosis. No cross-reactivity was observed with samples positive for HTLV, Syphilis, HBsAg, HCV, Chagas Disease, Zika, CMV, Rubella and Toxoplasmosis. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data together with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is advisable to use controls to assess the precision and accuracy of dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	SAMPLE		
	1	2	3
Average	1.598	0.756	0.098
Standard Deviation	0.121	0.104	0.013
Coefficient of Variation (%)	7.55	13.72	13.46

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	SAMPLE		
	1	2	3
Average	0.098	1.590	0.719
Standard Deviation	0.018	0.211	0.103
Coefficient of Variation (%)	18.32	13.25	14.33

Clinical Sensitivity and Specificity

The BIOLISA HIV 1/2/O Ag/Ab DBS kit was used for the analysis of clinical samples previously characterized and confirmed using another reference enzyme immunoassay method. The results show that the clinical sensitivity of the BIOLISA HIV 1/2/O Ag/Ab DBS kit is >99.9% and the clinical specificity is >99.9%.

RESULT		REFERENCE		
		Positive	Negative	Total
BIOLISA HIV 1/2/O Ag/Ab DBS	Positive	165	0	165
	Negative	0	247	247
	Total	165	247	412

Clinical Sensitivity: > 99.9% (165/165) IC 95% = 97.8 a 100%

Clinical Specificity: > 99.9% (247/247) IC 95% = 98.5 a 100%

DIAGNOSTIC SIGNIFICANCE

HIV is the etiological agent of Acquired Immunodeficiency Syndrome (AIDS). The main routes of transmission include exposure to blood and blood products, including sharing needles and syringes, sexual contact, mother-to-child transmission. The virion is surrounded by a lipid envelope that is derived from the host cell membrane. Several viral glycoproteins are on the envelope. Each virus contains two copies of genomic positive sense RNA. HIV-1 has been isolated from AIDS patients and healthy people at high risk for developing AIDS. HIV-1 infection is identified by an initial phase of antigenemia in which HIV-1 antigen (Ag) is detectable in the blood. In most cases, antigen levels are often difficult to detect, however, increasing immune system failure and rising virus levels can again stimulate detectable antigen levels. The important internal structural protein of HIV-1, the core protein P24, is one of the viral components found in the blood during antigenemia. Furthermore, HIV-1 consists of subtype M and subtype O. Highly divergent strains of HIV-1 were recognized in 1990 and provisionally grouped as subtype O, as this variation was similar to HIV-1 glycoprotein markers, but with a slight variation for the protein marker. Although rarely compared for HIV-1 and HIV-2, infections caused by Subtype O have so far been identified in Africa (Cameroon), France and Germany. HIV-2 has been isolated from AIDS patients and asymptomatic HIV-positive individuals in West Africa. HIV-1, HIV-2, and Subtype O induce an immune response. Immunological detection of anti-HIV antigens and antibodies in serum, plasma, or whole blood is the most efficient and common way to determine whether an individual has been exposed to HIV. Despite differences in their biological characteristics, serological activities and genome sequences, HIV-1, HIV-2, subtype O shows strong antigenic cross-reactivity. Antibodies against the HIV-1 P24 antigen are also included for detection of the pre-serum conversion phase of infection. Most HIV-2 positive sera can be identified through HIV-1 serological testing.

BIBLIOGRAPHIC REFERENCES

1. Chang, SY, Bowman, BH, Weiss, JB, Garcia, RE and White, TJ. The origin of HIV-1 isolate HTLV-IIIB. Nature (1993) 3;363:466-9.
2. Arya, SK, Beaver, B, Jagodzinski, L, Ensoli, B, Kanki, PJ, Albert, J, Fenyo, EM, Biberfeld, G, Zagury, JF and Laure, F. New human and simian HIV-related retroviruses possess functional transactivator (tat) gene. Nature (1987) 328:548-550.
3. Caetano JA Immunologic aspects of HIV infection. Acta Med Port (1991) 4 Suppl 1:52S-58S.
4. Janssen, RS, Satten, GA, Stramer, SL, Rawal, BD, O'Brien, TR, Weiblen, BJ, Hecht, FM, Jack, N, Cleghorn, FR, Kahn, JO, Chesney, MA and Busch MP. New testing strategy to detect early HIV-1 infection for use in incidence estimates and for clinical and prevention purposes. JAMA (1998) 280(1):42-48.
5. Travers, K, Mboup, S, Marlink, R, Gueye-Nidaye, A, Siby, T, Thior, I, Traore, I, Dieng-Sarr, A, Sankale, JL and Mullins, C. Natural protection against HIV-1 infection provided by HIV-2. Science (1995) 268:1612-161.
6. Greenberg, AE, Wiktor, SZ, DeCock, KM, Smith, P, Jaffe HW and Dondero, TJ, Jr. HIV-2 and natural protection against HIV-1 infection. Science (1996) 272:1959-1960.
7. TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
8. Bioclin - Dados de arquivos.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for BIOLISA HIV 1/2/O Ag/Ab DBS kit:
10269360458

Review: February/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER