BIOLISA HIV 1/2/O ANTÍGENO/ANTICORPO

REF K119

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Fourth-generation sandwich enzyme immunoassay (ELISA) for the qualitative detection of the presence of P24 Antigen, Anti-HIV-1, Anti-HIV-2, and/or subtype O antibodies in biological samples (human serum and plasma) through enzyme-immunoassay. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or enzyme immunoassay.
The BIOLISA HIV 1/2/O ANTIGEN/ANTIBODY kit is a fourth-generation, solid-phase, enzyme-linked immunoassay based on the "sandwich" principle for the detection of P24 Antigen, Anti-HIV-1, Anti-HIV-2, and/or subtype O in serum and plasma samples. The microplate is coated with HIV-1 and HIV-2 specific Antigens (HIV-1 Antigen gp41 and gp120 and HIV-2 Antigen gp36) and HIV-1 Anti-P24 Antibody. Antibodies present in the serum and plasma sample bind to the microplate coated with specific recombinant Antigens forming Antigen-Antibody immunocomplexes. Along with the first incubation, Conjugate 1 composed of Anti-HIV P24 conjugated to Biotin is added, which binds to P24 Antigens present in the sample, forming Antibody-Antigen-Antibody immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. Conjugate 2 containing HIV-1, HIV-2 Antigens and Streptavidin, all linked to Peroxidase, are added to the microplate and then incubated. Conjugate 2 Antigens bind to sample Antibodies that have been captured by microplate Antigens, forming Antigen-Antibody-Antigen immune complexes. Streptavidin linked to Peroxidase, also present in Conjugate 2, binds to Biotin linked to Anti-HIV P24 of the Antibody-Antigen-Antibody immune complexes. After these steps, a new wash is performed to remove unbound materials. The Substrate is then added and incubated, producing a blue color that indicates the detection of Anti-HIV Antibodies and/or P24 Antigen present in the samples. Stop Solution is added to stop the reaction with a color change from blue to yellow. Color intensity is measured at 450/620 nm.

REAGENTS

- 1- Sensitized Plate** – Store between 2 and 8 °C. Contains: Microplate coated with specific Antigens for HIV-1 and HIV-2 (HIV-1 Antigen gp41 and gp120 and HIV-2 Antigen gp36) and HIV-1 Anti-P24 Antibody.
- 2- Conjugate 1** – Store between 2 and 8 °C. Contains: Buffer Solution containing P24 Anti-HIV Antibodies linked to Biotin, surfactant, stabilizers, dye and preservative.
- 3- Concentrated Washing** – Store between 2 and 8°C. Contains: Buffer Solution (phosphate < 0.5 mol/L, potassium chloride < 100 mmol/L, sodium chloride < 5 mol/L), surfactant and preservative.
- 4- Conjugate 2** – Store between 2 and 8 °C. Contains: Specific Antigens for HIV-1, HIV-2 and Streptavidin, all linked to Peroxidase, surfactant, stabilizers, dye and preservative.
- 5- Substrate** – Store between 2 and 8 °C. Contains: Buffer Solution containing Urea Peroxide, Tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** – Store between 2 and 8 °C. Contains: 1 M Hydrochloric Acid Solution.
- 7- Negative Control** – Store between 2 and 8°C. HIV non-reactive solution, dye and preservative. **Potentially infectious.**
- 8- HIV-1 Positive Control** – Store between 2 and 8 °C. Reactive solution for HIV-1, stabilizer, dye and preservative. **Potentially infectious.**
- 9- HIV-2 Positive Control** – Store between 2 and 8 °C. HIV-2 reactive solution, stabilizer, dye and preservative. **Potentially infectious.**
- 10- HIV-1 P24 Positive Control** – Store between 2 and 8 °C. Solution with P24 Antigens, stabilizer, dye and preservative. **Potentially infectious.**
- 11- Plate Sealers.**

PRESENTATION

REAGENTS	1	2	3
	96 cavities	192 cavities	480 cavities
1- Sensitized Plate	1 Unit	2 Units	5 Units
2- Conjugate 1	1 Vial x 6 mL	1 Vial x 12 mL	1 Vial x 30 mL
3- Concentrated Washing	1 Vial x 50 mL	1 Vial x 100 mL	1 Vial x 250 mL
4- Conjugate 2	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
5- Substrate	1 Vial x 12 mL	1 Vial x 24 mL	1 Vial x 60 mL
6- Stop Solution	1 Vial x 8 mL	1 Vial x 16 mL	1 Vial x 40 mL
7- Negative Control	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 1 mL
8- HIV-1 Positive Control	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 1 mL
9- HIV-2 Positive Control	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 1 mL
10- HIV-1 P24 Positive Control	1 Vial x 1 mL	1 Vial x 1 mL	1 Vial x 1 mL
11- Plate Sealers	3 Units	1 Unit	1 Unit

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the table above.

Required materials not contained in the kit:

- 1- Pipettes capable of dispensing volumes from 50 to 500 µL with a coefficient of variation less than 1.5%.
- 2- Repeater for repetitive pipetting of volumes of 300 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3- Microplate washer for technique on filter paper.
- 4- ELISA reader with absorbance capacity at 450 and 620 nm of wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or watch.
- 7- Bottle to store the Washing Solution after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.
- 10- Incubator 37 ± 2 °C.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be between 2 and 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

1- For *in vitro* diagnostic use only.

- 2- Strictly follow the proposed methodology to obtain accurate results.
- 3- The sachet containing the microplate should only be opened after reaching room temperature. Replace the unused microwell strips in the sachet, seal and store at 2 to 8 °C.
- 4- The water used to clean the material must be recent and free of contaminants.
- 5- Saturated deionizing columns release alkaline water, various ions, and oxidizing and reducing agents that can significantly alter the results.
- 6- The Stop Solution contains Hydrochloric Acid which is a strong acid. Therefore, handle it with due care.
- 7- All raw material of the product is tested and must be non-reactive for HbsAg and Anti-HCV. However, these tests do not offer complete assurance of the absence of infectious agents. Handling any product containing serum is potentially capable of transmitting diseases. Therefore, it is necessary to take the proper biosafety precautions when handling these products.

8- Always pipette the reagents in the same order to minimize the difference in reaction time between the microwells.

9- As a protection measure, the plate must be covered during the reaction.

10- It must be ensured that the bottom of the well is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow the wells to dry out during the assay.

11- Do not expose the reagents, especially the Substrate, to strong light or hypochlorite vapors during storage or incubation steps.

12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advice Service) by Quibasa.

14- Do not use the product if the packaging is damaged.

15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin)

Hemolyzed and highly lipemic samples should not be used. The samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C. Plasma samples stored for periods longer than recommended should not be used.

PROCESS DESCRIPTION

Stability After Opening

The results of the stability test show that the BIOLISA HIV 1/2/O ANTIGEN/ANTIBODY kit is stable after opening for up to 30 days. This stability may vary according to test and environmental conditions. Therefore, it is suggested to monitor the performance of the product using the kit's internal controls and the technique validation criteria.

PREPARATION OF WORKING REAGENTS

Washing Solution

Dilute the contents of vial N° 3 (Concentrated Washing) in a 1:20 ratio of distilled or deionized water; for example: 50 mL of Concentrated Washing solution in 1000 mL of distilled or deionized water. After preparation, the solution can be stored between 2 and 30°C for up to 30 days. If crystallization occurs, heat to 37 °C until dissolution.

All other reagents are ready to use.

TECHNIQUE

For use in automatic equipment, contact SAC (Customer Advice Service).

Before starting the assay, place all reagents, samples and controls to stabilize at room temperature (15 – 30 °C) for at least 40 minutes.

1- Separate the wells to be used considering: HIV-1 Positive Control, HIV-2 Positive Control, HIV-1 P24 Positive Control, and Samples (it is recommended to test in duplicate). Return unused microplate strips to the original sealed package.

2- Separate the first cavity for the Blank (OPTIONAL).

IMPORTANT! The following steps must be carried out strictly in the sequence described below.

3- Pipette 50 µL of Conjugate 1 into all wells, **except the well for the Blank.**

4- Pipette 50 µL of HIV-1 Positive Control, HIV-2 Positive Control, HIV-1 P24 Positive Control, Negative Control and Samples into the previously determined wells.

5- Gently homogenize for ± 30 seconds, cover the wells with plate sealer.

6- Incubate for 60 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

7- After incubation, discard the contents of the wells by aspiration. Use approximately 300 µL of **previously prepared** Washing Solution and perform a total of five (5) second shake wash cycles. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.

Note: Poor washing/drying may cause poor results.

8- Pipette 100 µL of Conjugate 2 into all wells, **including the Blank well.**

9- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.

10- Incubate for 60 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

11- Remove the plate sealer from the wells.

12- Repeat item 7.

13- Pipette 100 µL of Substrate in all wells, **including the Blank well.**

14- Gently homogenize for ± 30 seconds. Cover the wells with the plate sealer.

15- Incubate for 10 minutes ± 2 minutes in an incubator at 37 °C ± 2 °C.

16- Remove the plate sealer from the wells.

17- Pipette 50 µL of Stop Solution into all wells.

18- Gently homogenize for ± 30 seconds.

19- Read using a double filter: 450 nm / 620 nm in up to 15 minutes (maximum).

TECHNIQUE VERIFICATION

Check if the results obtained for reading the Blank and the Reference Standards are compatible with the values presented below:

ITEM	ABSORBANCES
Blank	< 0.250
Negative Control	< 0.250
HIV-1 Positive Control	> 0.500
HIV-2 Positive Control	> 0.500
HIV-1 P24 Positive Control	> 0.500

If the values are outside the expected values, the technique must be repeated.

CALCULATIONS

QUALITATIVE

Calculate Cut Off according to the following formula:

Cut Off = Mean Absorbance of Negative Control + 0.130.

Example:

ITEM	ABSORBANCES
Negative Control	0.105
	0.101
Cut Off = (Average Absorbance of the Negative Control) + 0.130.	((0.105 + 0.101) / 2) + 0.130 = 0.233

Calculate the Index by dividing the absorbance of the sample by the Cut Off value.

Example:

ITEM	ABSORBANCES
Sample	2.102
Cut Off Value	0.233
Index = Sample / Cut Off Value	2.102 / 0.233 = 9.02

Note: Data shown in examples is for illustration only and cannot be used for calculation of results.

INTERPRETATION OF RESULTS

After calculating the sample index, consider the indexes below to determine the results.

RESULTS	QUALITATIVE
	INDEX
Negative (Non-Reagent)	< 0.9
Positive (Reagent)	> 1.1
Undetermined	Between 0.9 - 1.1

Non-Reactive: Sample with absorbance lower than the cut-off is considered non-reactive for P24 Antigen, Anti-HIV-1, Anti-HIV-2, and/or subtype O Antibodies, and may be considered negative.

Reagent: A sample with an absorbance greater than the cut-off is considered initially reactive for P24 Antigen, Anti-HIV1, Anti-HIV2, and/or subtype O Antibodies. The sample must be reanalyzed in duplicate before the end of the interpretation. The sample that is reactive in at least one of the reanalyses is presumed to be reactive and must be confirmed through another diagnostic method.

Notes: In case of indeterminate results, the sample must be reanalyzed in duplicate. Samples that give repeatedly indeterminate results should be retested using an alternate method. If results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected shall prevail. The interpretation of a diagnostic test should not be based on a single test. Additional confirmatory tests must be included before a specimen is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other available clinical information, prior to making a descriptive diagnosis of the disease.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be based on a single test. All results must be interpreted in conjunction with other available clinical information prior to definitive diagnosis. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.

INTERFERING

No interference was observed by Triglycerides 1200 mg/dL, Acetylsalicylic Acid 20 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Acid Oxalic 60 mg/dL, C-Reactive Protein 41.2 mg/dL and anti-Streptolysin O 1023 IU/mL. High concentrations of Rheumatoid Factor do not interfere with the result obtained from positive or negative HIV samples. Plasma samples stored for periods longer than recommended (30 days) may show fibrin and fibronectin precipitation that may interfere with the test.

CROSS REACTIVITY

A study was carried out with 81 serum and plasma samples negative for HIV, but positive for other infections, in order to evaluate the possibility of cross-reactivity of these interfering agents in the result of the BIOLISA HIV 1/2/O ANTIGEN/ ANTIBODY. Among them 10 positive samples for HBsAg, 6 positive samples for HTLV, 5 positive samples for HCV, 10 positive samples for COVID-19, 8 positive samples for Dengue, 10 positive samples for Rubella, 5 positive samples for CMV, 8 positive samples for Chagas disease, 10 positive samples for Syphilis and 9 positive samples for Toxoplasmosis. No cross-reactivity was observed with positive samples for HBsAg, HTLV, HCV, COVID-19, Dengue, Rubella, CMV, Chagas Disease, Syphilis and Toxoplasmosis. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis must consider the patient's clinical data along with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following absorbance results:

Repeatability	Sample		
	1	2	3
Average	2.356	0.552	0.044
Standard Deviation	0.039	0.030	0.004
Coefficient of Variation (%)	1.67	5.52	8.38

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 3 samples with different values, obtaining the following absorbance results:

Reproducibility	Sample		
	1	2	3
Average	2.281	0.546	0.043
Standard Deviation	0.102	0.029	0.008
Coefficient of Variation (%)	4.47	5.31	18.56

CLINICAL SENSITIVITY AND SPECIFICITY

The BIOLISA HIV 1/2/O ANTIGEN/ANTIBODY kit was used for the analysis of clinical samples previously characterized and confirmed by another reference enzyme-immunoassay method. The results show that the clinical sensitivity of the BIOLISA HIV 1/2/O ANTIGEN/ANTIBODY kit is >99.9% and the clinical specificity is >99.9%.

BIOLISA HIV 1/2/O ANTIGEN/ANTIBODY	Result Reference		
	Positive	Negative	Total
Positive	84	0	84
Negative	0	341	341
Total	84	341	425

Clinical Sensitivity: > 99.9% (84/84) IC 95% = 95.7 to 99.9%

Clinical Specificity: > 99.9% (341/341) IC 95% = 98.9 to 99.9%

DIAGNOSTIC SIGNIFICANCE

HIV is the etiologic agent of Acquired Immunodeficiency Syndrome (AIDS). The main routes of transmission include exposure to blood and blood products, including sharing needles and syringes, sexual contact, mother-to-child transmission. The virion is surrounded by a lipid envelope that is derived from the host cell membrane. Several viral glycoproteins are on the envelope. Each virus contains two copies of genomic positive sense RNA. HIV-1 has been isolated from AIDS patients and from healthy people at high potential risk for developing AIDS. HIV-1 infection is identified by an early stage of antigenemia in which HIV-1 Antigen (Ag) is detectable in the blood. In most cases, antigen levels are often difficult to detect, however, increasing immune system failure and rising virus levels can again stimulate detectable Antigen levels.

The important internal structural protein of HIV-1, the core protein P24, is one of the viral components found in the blood during antigenemia. Furthermore, HIV-1 consists of subtype M and subtype O. Highly divergent strains of HIV-1 were recognized in 1990 and provisionally grouped as subtype O, as this variation was similar to HIV-1 glycoprotein markers, but with a slight variation for the protein marker. Although rarely compared to HIV-1 and HIV-2, infections caused by Subtype O have so far been identified in Africa (Cameroon), France and Germany. HIV-2 has been isolated from AIDS patients and asymptomatic HIV-positive individuals in West Africa. HIV-1, HIV-2, and Subtype O induce an immune response. Immunological detection of anti-HIV antigens and antibodies in serum, plasma, or whole blood is the most efficient and common way to determine whether an individual has been exposed to HIV.

Despite differences in their biological characteristics, serological activities and genome sequences, HIV-1, HIV-2, subtype O show strong antigenic cross-reactivity. Antibodies against the HIV-1 P24 Antigen are also included for detection of the pre-serum conversion phase of infection. Most HIV-2 positive sera can be identified using HIV-1 serological tests.

BIBLIOGRAPHIC REFERENCES

1. Chang, SY, Bowman, BH, Weiss, JB, Garcia, RE and White, TJ. The origin of HIV-1 isolate HTLV-IIIB. Nature (1993) 3;363:466-9.
2. Arya, SK, Beaver, B, Jagodzinski, L, Ensoli, B, Kanki, PJ,Albert, J, Fenyo, EM, Biberfeld, G,Zagury, JF and Laure, F. New human and simian HIVrelated retroviruses possess functional transactivator (tat) gene. Nature (1987) 328:548-550.
3. Caetano JA Immunologic aspects of HIV infection. Acta Med Port (1991) 4 Suppl 1:52S-58S.
4. Janssen, RS, Satten, GA, Stramer, SL, Rawal, BD, O'Brien, TR, Weiblen, BJ, Hecht, FM, Jack, N, Cleghorn, FR, Kahn, JO, Chesney, MA and Busch MP. New testing strategy to detect early HIV-1 infection for use in incidence estimates and for clinical and prevention purposes. JAMA (1998) 280(1):42-48.
5. Travers, K, Mboup, S, Marlink, R, Gueye-Nidaye, A, Siby, T, Thior, I, Traore, I, Dieng-Sarr, A, Sankale, JL and Mullins, C. Natural protection against HIV-1 infection provided by HIV-2. Science (1995) 268:1612-161.
6. Greenberg, AE, Wiktor, SZ, DeCock, KM, Smith, P, Jaffe HW and Dondero, TJ, Jr. HIV-2 and natural protection against HIV-1 infection. Science (1996) 272:1959-1960.
- 7.TELELAB. Manual de Coleta de Sangue - Diagnóstico e Monitoramento das DST, Aids e Hepatites Virais.
- 8.QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiry date mentioned on the presentation packaging, provided they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of the BIOLISA HIV 1/2/O ANTIGEN/ANTIBODY kit at ANVISA: 10269360199

Review: September/2024

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER